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Trademarks: QIAGEN®, HotStarTaq®, QuantiTect® (QIAGEN Group); ABI PRISM®, Applied Biosystems® (Applera Corporation or its subsidiaries); Bio-Rad®, iCycler iQ® (Bio-Rad Laboratories, Inc.); Eppendorf® [Eppendorf AG]; Roche®, LightCycler® (Roche Group); Cepheid®, SmartCycler® (Cepheid); SYBR® (Molecular Probes, Inc.); Stratagene® [Stratagene].

The purchase price of this product (QuantiTect SYBR Green RT-PCR Kit) includes a limited, non-transferable immunity from suit under U.S. Patents Nos. 5,994,056 and 6,171,785 and corresponding patent claims outside the United States, owned by Roche Molecular Systems, Inc. or F. Hoffmann-La Roche Ltd (Roche), to use only this amount of the product for dsDNA-binding dye processes covered by said patents solely for the purchaser's own internal research. This product is also a Licensed dsDNA Dye Binding Kit for use with service sublicenses available from Applied Biosystems. No right under any other patent claims (such as apparatus or system claims in U.S. Patent No. 6,814,934) and no right to use this product for any other purpose or for commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted expressly, by implication or by estoppel. This product is for research use only. Diagnostic uses require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

NOTICE TO PURCHASER: LIMITED LICENSE

The purchase price of this product (QuantiTect SYBR Green RT-PCR Kit) includes a limited, non-transferable license under U.S. Patents Nos. 5,407,800, 5,322,770, 5,310,652 and corresponding patent claims outside the United States, owned by Roche Molecular Systems, Inc. or F. Hoffmann-La Roche Ltd (Roche), to use only this amount of product solely for the purchaser's own internal research. No right under any other patent claims (such as apparatus or system claims) and no right to use this product for any other purpose or for commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is hereby granted expressly, by implication or by estoppel. This product is for research use only. Diagnostic uses require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

NOTICE TO PURCHASER: LIMITED LICENSE

This product (QuantiTect SYBR Green RT-PCR Kit) or its use is covered by at least one claim of U.S. Pat. Nos. 5,035,996; 5,945,313; 6,287,823; or 6,518,026, owned by Invitrogen Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product, (b) its components, or (c) materials made by the employment of this product or its components to a third party or otherwise use this product or its components or materials made by the employment of this product or its components for Commercial Purposes. Commercial Purposes means any activity for which a party receives or is due to receive consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. The buyer cannot use this product or its components or materials made using this product or its components for therapeutic, diagnostic or prophylactic purposes. Further information on purchasing licenses under the above patents may be obtained by contacting the Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, CA 92008. Email: outlicensing@invitrogen.com.

NOTICE TO PURCHASER: DISCLAIMER OF LICENSE

This product (QuantiTect Primer Assays) is compatible for use in real-time PCR methodologies, including 5' nuclease processes and dsDNA-binding dye processes claimed in patents owned by Roche Molecular Systems, Inc., F. Hoffmann-La Roche Ltd, and Applera Corporation. No right to perform any of these patented processes is conveyed, expressly, by implication or by estoppel to the purchaser by the purchase of this product. A license to use these patented processes for the purchaser's own internal research accompanies the purchase of certain reagents of Applied Biosystems and other licensed suppliers, or is available from Applied Biosystems. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

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© 2001–2008 QIAGEN, all rights reserved.