

Phoenix Hot Start Taq DNA Polymerase

How stable is Phoenix™ Hot Start Taq when incubated in a PCR reaction mix at room temperature?

Phoenix Hot Start Taq remains stable at room temperature for at least 72 hours when incubated in a PCR reaction mix containing 1X of the supplied, optimized Phoenix Hot Start Taq reaction buffer (standard or GC). Additionally, endpoint PCR tests on a smaller set of amplicons demonstrated functional stability after 20 days of room temperature incubation.

How can PCR cycling conditions be optimized for Phoenix Hot Start Taq?

General Taq protocols (see recommended cycling conditions on the Phoenix Hot Start Taq Product Information Sheet) are a good starting point. To fine-tune the annealing step, vary the temperature in 2°C increments, beginning about 10°C below the lower primer T_m and increasing up to approximately 5°C above the higher primer T_m . Phoenix Hot Start Taq is also compatible with Touchdown PCR cycling protocols.

Can Phoenix Hot Start Taq use cDNA as a template for PCR?

Yes.

Is Phoenix Hot Start Taq capable of multiplex PCR?

Yes. Phoenix Hot Start Taq can amplify multiple targets in a single 50 µl PCR reaction. Increasing the amount of polymerase (up to 2.5 U) may improve multiplex results.

How can I optimize Mg^{2+} conditions for a specific amplicon when using Phoenix Hot Start Taq and the supplied reaction buffer?

The supplied Phoenix Hot Start Taq buffers provide 2 mM Mg^{2+} in the final reaction. If needed, you can adjust the final Mg^{2+} concentration by adding a concentrated Mg^{2+} solution. Keep in mind that increasing Mg^{2+} levels may reduce PCR fidelity (1).

When should I use Phoenix Hot Start Taq GC reaction buffer?

The GC reaction buffer is recommended for challenging or GC-rich amplicons ($\geq 55\%$ GC content).

Ordering Information

Product	Contents	Cat. no.
Phoenix Hot Start Taq DNA Polymerase	500 U of Phoenix Hot Start Taq DNA Polymerase (0.10mL at 5,000 U/mL), 5x Phoenix Hot Start Taq Reaction Buffer (4 x 1.5 mL), and 5x Phoenix Hot Start Taq GC Reaction Buffer (2 x 1.5 mL)	P7590L

References

1. Eckert K and Kunkel TA. (1991) *Genome Research*, **1**, 17 - 24.
2. Frey B and Suppmann B (1995) *Biochemica*, **2**, 8 - 9.
3. Ralser M, et al. (2006) *Biochemical and Biophysical Research Communications*, **347**, 747 - 751.



Click on the QR code or visit www.qiagen.com/enzymes to learn more.



Need bulk, customized or optimized products for commercial purposes? We also offer support with logistics, compliance and more. Reach out to cooperate with QIAGEN Strategic Partnerships & OEM: OEM@qiagen.com.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit instructions for use or user operator manual. QIAGEN kit instructions for use and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services (or your local distributor).

A portion of this product is licensed from Johnson & Johnson Services, Inc. under U.S. Patent Nos. 5,338,671 and 5,587,287. This product was developed, manufactured, and sold for research use only. Diagnostics rights are available, please contact us for more information.

Trademarks: QIAGEN®, Sample to Insight®, VeraSeq™ (QIAGEN Group).

Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

04/2025 QPRO-10096 © 2025 QIAGEN, all rights reserved.