

August 2023

Quick-Start Protocol

Tth DNA Ligase

Tth DNA Ligase (cat. nos. EN13-025, EN13-250) catalyzes the NAD-dependent formation of phosphodiester bonds between adjacent 3'-hydroxyl and 5'-phosphate termini in double stranded DNA. It is not active against single-stranded DNA or RNA and blunt-ended DNA. Enzyme is isolated from *E. coli* strain containing plasmid carrying the *Thermus thermophilus* DNA ligase gene.

Tth DNA Ligase is stable and active in optimum ligation temperature range of 45–65°C, which is 7–10°C higher than that of T4 DNA ligase. The final reaction ligation temperature is determined by the T_m (melting temperature) of the substrates. High ligation temperature eliminates nonspecific ligation.

The Tth DNA Ligase and its components should be shipped on dry ice, and stored at –20°C.

Further information

- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- Enzyme retains full activity after incubation for 1 week at 37°C.
- The half-life of the enzyme is about 48 hours at 65°C.
- 10x Tth Ligation Buffer is stable for 1 week at 37°C.

- Up to 20 freeze/thaw cycles will not compromise 10x Tth Ligation Buffer performance.
- The concentration of the enzyme is 5 U/ μ L.
- One unit of Tth DNA Ligase catalyzes the ligation of 50% of the cos sites present in 1 μ g of bacteriophage lambda DNA in 1 minute at 45°C.

Note: 1 U (Unit) of Tth DNA Ligase is equivalent to 15 cohesive end units (CEU)

Protocol

1. Add the reaction reagents listed below to a sterile nuclease-free tube. The reaction agents should be added in the following order:

Table 1. Reagents to be added and corresponding quantities

Component	Quantity
Nuclease-free water	Up to 25 μ L
10x Tth Ligation Buffer	2.5 μ L
Tth DNA Ligase	0.5–1 μ L
DNA	0.5–1 μ g

2. Mix gently and spin briefly.
3. Incubate for 10 min at 45–65°C for 30 min. Optimum ligation temperature range is determined by the T_m of the substrates.

Document Revision History

Date	Changes
August 2023	Initial release

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

08/2023 HB-3445-001 © 2023 QIAGEN, all rights reserved.