

QuantiNova® Internal Control RNA and Assay

The QuantiNova Internal Control RNA (QN IC RNA) in the QuantiNova Reverse Transcription Kit (cat. nos. 205411, 205413), QuantiNova Probe RT-PCR Kit (cat. nos. 208352, 208354, 208356), QuantiNova Multiplex RT-PCR Kit (cat. no. 208552, 208554, 208556), and QuantiNova SYBR® Green RT-PCR Kit (cat. nos. 208152, 208154, 208156) should be stored immediately upon receipt at -30°C to -15°C in a constant-temperature freezer.

Further information

- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- The QN IC RNA is a synthetic RNA that can be optionally used to monitor successful reverse transcription. The QN IC RNA is intended to report instrument or chemistry failures, errors in assay setup and the presence of inhibitors. Inhibitors such as phenol, ethanol, sodium dodecyl sulfate (SDS), or ethylene diaminetetraacetic acid (EDTA) may remain from the lysis and purification steps of the RNA isolation procedure.
- The primer and probe sequences for the detection of the QN IC RNA have been bioinformatically validated for non-homology against hundreds of eukaryotic and prokaryotic organisms. Additionally, they have been experimentally tested against a multitude of human, mouse, and rat RNA samples from multiple tissues and cell lines.
- The QN IC RNA is detected as a 200 bp amplicon. For SYBR® Green-based detection, use the QuantiNova IC SYBR Green Assay (cat. no. SBH1218551, ordered via GeneGlobe®, available under the name HS_QIC_2467742 QuantiNova LNA PCR Reference Assay) and the QuantiNova SYBR Green PCR and RT-PCR Kits. For probe-based detection, use the QuantiNova IC Probe Assay (cat. no. 205813) or QuantiNova IC Probe Assay Red 650 (cat. no. 205824). The QN IC RNA can be detected using the VIC®/HEX or yellow dye channel (cat. no. 205813) or the Cy5 or red dye channel (cat

no. 205824) of your real-time PCR instrument and the QuantiNova Probe PCR and RT-PCR Kits.

- The QN IC RNA (optional) is added to the experimental RNA sample. An additional no-template RNA control sample, which only contains the QN IC RNA, should also be set up. C_q shifts >2 between the template RNA+QN IC RNA compared to the QN IC RNA only samples, and between the different template RNA containing samples indicate inhibition of the RT-PCR.
- For the QuantiNova Reverse Transcription Kit, the QN IC RNA is used undiluted. Before use with the QuantiNova Probe or Multiplex RT-PCR Kit and the QuantiNova SYBR Green RT-PCR Kit, dilute the QN IC RNA 1:10 by adding 180 μL (20 μL tube) or 900 μL (100 μL tube) of RNase-free water to the RNA. Mix thoroughly by vortexing.

Procedure

1. Mix the QuantiNova Internal Control RNA thoroughly by vortexing, and dispense 1 μL into each sample as described in the quick-start protocols for the QuantiNova Reverse Transcription Kit (undiluted IC) and the QuantiNova SYBR Green or Probe RT-PCR Kits (pre-diluted IC).
2. Analyse depending on the assay used:
 - 2a. For analysis with the QuantiNova Probe or Multiplex RT-PCR Kit, add the appropriate volume of 10x QuantiNova IC Probe Assay to the sample. Signal detection is performed on the filter/channel for HEX/VIC (cat. no. 205813), or for Cy5 (cat no. 205824) of your real-time PCR instrument. For analysis with the QuantiNova SYBR Green RT-PCR Kit, use the appropriate volume of 10x QuantiNova IC SYBR Green Assay in the reaction.

Note: QuantiNova IC Assays, for the detection of QN IC RNA, need to be ordered separately as the QuantiNova IC Probe Assay (cat. no. 205813), QuantiNova IC Probe Assay Red 650 (cat no. 205824), or QuantiNova IC SYBR Green Assay (cat. no. SBH1218551, available under the name HS_QIC_2467742 QuantiNova LNA PCR Reference Assay, for SYBR[®] Green detection) at GeneGlobe (www.qiagen.com/geneglobe).

- 2b. If using the QuantiNova Reverse Transcription Kit, cDNA should be diluted (1:10–1:100) and an aliquot of the reaction should be used for subsequent amplification with the QuantiNova SYBR Green or Probe PCR Kits.

For probe-based detection, use the QuantiNova IC Probe Assays and detection in the VIC/HEX dye or yellow channel (cat. no. 205813), or in the Cy5 or red channel (cat. no. 205824) of your real-time PCR instrument. For detection with the QuantiNova SYBR Green PCR Kit, use the QuantiNova IC SYBR Green® Assay.

Data analysis and interpretation of results

1. After amplification, perform data analysis as recommended for your real-time PCR instrument.
2. Compare C_q values between the QN IC RNA only and samples containing QN IC RNA plus template RNA.
3. Consistent C_q values indicate successful RT-PCR and reliable results. A C_q difference >2 is likely to indicate inhibition or sample failure.
4. If a shifted C_q of >2 appears, indicating inhibition or failure of a specific sample, we recommend the following:
 - 4a. Check equipment for accurate performance and repeat sample/experiment to rule out pipetting or handling errors.
 - 4b. Dilute the affected template RNA using RNase-free water before repeating the experiment. This dilutes inhibitors present in the sample.
 - 4c. Consider repeating the RNA extraction and avoid contamination or carry-over of inhibitors (e.g., use an appropriate RNeasy® Kit).

Alternatively, the RNeasy MinElute Cleanup Kit (cat. no. 74204) can be used to remove potential inhibitors and concentrate the RNA template.

Expected C_q value for successful RT-PCR (e.g., no inhibition of QN IC RNA)

The C_q value for the QN IC RNA in 2-step RT-PCR using the QuantiNova Reverse Transcription Kit will vary depending on various parameters, such as:

- Detection format (SYBR Green or probe)
- Dilution factor of the cDNA after reverse transcription
- Volume of template cDNA added to the PCR
- Type of real-time PCR instrument
- C_q determination (threshold setting, Auto- C_q determination)

An example of typical results for the QN IC RNA:

- Dilution factor: 10
- Volume of cDNA added to PCR: 2 μ L
- Real-time instrument: Rotor-Gene® Q
- C_q determination: manual threshold setting
- C_q value with QuantiNova SYBR® Green PCR Kit: 21
- C_q value with QuantiNova Probe PCR Kit: 27

This example illustrates that C_q values for the QN IC RNA in 2-step RT-PCR, using the QuantiNova Reverse Transcription Kit, will not be identical under all circumstances. However, the C_q values should be consistent under identical conditions; therefore, the QN IC RNA can be used to monitor successful reverse transcription and RT-PCR. A C_q difference >2 is likely to indicate sample inhibition or failure.

Typical results for the QN IC RNA in 1-step RT-PCR:

The C_q value for the QN IC RNA in the QuantiNova Probe RT-PCR Kit depends on the real-time PCR instrument used and can be expected within a C_q range of 23–25.

The C_q value for the QN IC RNA in the QuantiNova SYBR Green RT-PCR Kit depends on the real-time PCR instrument used and can be expected within a C_q range of 17–19.

Document Revision History

Date	Changes
04/2024	Converted to the new branding template. Replaced QuantiTect Primer Assay with QuantiNova IC SYBR Green Assay and added QuantiNova IC Probe Assay Red 650 (cat no. 205824).



Scan QR code for *handbook*.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, QuantiNova™, QuantiTect®, RNeasy®, Rotor-Gene®, GeneGlobe® (QIAGEN Group); SYBR®, VIC® (Life Technologies Corporation). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

This product is provided under an intellectual property license from Life Technologies Corporation. The sale of this product is expressly conditioned on the buyer using the product or its components only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or Qulicensing@lifetech.com

1135230 04/2024 HB-1932-002 © 2024 QIAGEN, all rights reserved.