

A NEW IDH1/2 PCR ASSAY FOR ONE-STEP DETECTION OF 12 IDH1 AND IDH2 MUTATIONS IN GLIOMA

H. Girardi ¹, F. Monville ¹, S. Carpentier ¹, M. Giry ², J. Voss ³, R. Jenkins ³, B. Boisselier ², V. Frayssinet ¹,
A. Catteau ¹, K. Mokhtari ², M. Sanson ², H. Peyro-Saint-Paul ¹, C. Giannini ³

¹ QIAGEN Marseille, France - ² AP-HP, Pitié-Salpêtrière Hospital, Paris, France - ³ Mayo Clinic, Rochester, MN, USA

BACKGROUND

- Isocitrate dehydrogenase (IDH) mutational status is a strong diagnostic and prognostic marker in glioma which will probably be introduced in the next WHO classification system.
- In addition to the established 1p/19q codeletion and MGMT methylation, a series of new biomarkers such as IDH1/2, EGFR or BRAF mutations, and FGFR gene fusions, are increasingly documented to play a role as prognostic or predictive markers, and should progressively be introduced in the diagnostic and treatment decision algorithm for glioma.
- Current IDH mutations screening is performed with an IHC assay specific for IDH1 R132H, the most common mutation. Sequencing is recommended as a second-step test for IHC-negative or -equivocal cases. However, sequencing is not readily accessible in all centers, and its use generally leads to additional delay in providing a comprehensive IDH1/2 mutational status assessment. Moreover, IDH sequencing procedures can sometimes lead to inter-laboratory variability.
- Beyond IDH1 R132H, 11 other IDH mutations, 6 IDH1 and 5 IDH2, have been reported so far in literature, in large cohorts of gliomas.
- Recent data indicate that mutant specific IDH1-inhibitors impair growth of mutant-IDH1 gliomas in mice.
- A real-time PCR assay was designed to detect the 12 IDH1/2 mutations in one single step in FFPE samples, and identify the most frequent ones.

OBJECTIVES

- Establish the analytical performance of the new IDH1/2 PCR assay
- Validate the IDH1/2 PCR assay performance on FFPE glioma clinical samples by comparing PCR IDH mutational status to IHC and Sanger sequencing.

MATERIAL & METHODS

IDH1/2 one-step qPCR assay:

- PCR Clamping was used for the qualitative detection of IDH1 R132H and 11 additional IDH1/2 mutations. ARMS PCR technology was combined to selectively identify the most frequent IDH1 (R132H / R132C) and IDH2 (R172K) mutations (Table 1).

Evaluation of analytical sensitivity:

- Limit of Detection (LOD) (min. % mutant DNA detected in a WT background) was determined following CLSI/NCCLS EP17-A guidelines.
- 5 low positive samples (2-5-10-15 and 20%) obtained by mixing IDH mutant plasmid DNA with glioma IDH1/2 WT DNA were tested per mutation (n= 30 to 110 measurements per mutation and mutation percentage)

Validation of the IDH1/2 PCR assay:

- Samples**
 - 171 FFPE glioma samples: 121 samples retrospectively collected in a reverse chronological order from 2 academic centers (C1=102; C2=19) and 50 additional commercial samples. No specific selection criteria beyond tumor characteristics assessed by local pathologists.
 - Samples selection: < 10 yrs; ≥ 50mm² tissue area with ≥ 40% tumoral cells.
 - DNA extracted from 10 μm FFPE sections (QIAamp DNA FFPE Tissue Kit, Qiagen)
 - qPCRs performed on 25 ng DNA acc. to *therascreen* IDH1/2 RGQ PCR Kit IFU (Qiagen) and Rotor-Gene Q 5plex HRM instrument (Qiagen)
- Molecular biology methods**
 - IHC performed locally using the IDH1 R132H monoclonal antibody Clone H09 (Dianova)
 - Bidirectional Sequencing (centrally) using recommended primers for IDH1 codon 132 and IDH2 codon 172⁽³⁾ and newly designed primers for IDH1 codon 100.
 - Discordant PCR/Sanger Sequencing cases additionally tested by pyrosequencing +/- LNA-based sequencing (centrally).
- Rare mutations testing**
 - 22 synthetic samples (30% and 40% Mutant DNA in WT DNA) for the 11 rare mutations, processed similarly to clinical samples

IDH1/2 ASSAY DESIGN

Gene / Codon	Mutation	Base change
IDH1 / R132	R132H *	395 G>A
	R132C *	394 C>T
	R132S	394 C>A
	R132G	394 C>G
	R132L	394 G>T
	R132V	394_395 CG>GT
IDH1 / R100	R100Q	299 G>A
IDH2 / R172	R172K *	515 G>A
	R172M	515 G>T
	R172W	514 A>T
	R172S	516 G>T
	R172G	514 A>G

Table 1 - IDH1/2 mutations detected and identified* with the IDH1/2 PCR Assay

Detection of 12 mutations by PCR CLAMPING:

- 6 within IDH1 codon R132
- 5 within the homologous codon 172 of IDH2
- one within IDH1 codon 100

Identification of 3 major IDH1/2 mutations by ARMS:

- IDH1 R132H
- IDH1 R132C
- IDH2 R172K

ANALYTICAL SENSITIVITY

- Analytical sensitivity ranged from 0,6 % to 15% according to mutations (mean = 3.3 %)
- LOD was < 5% for 11/12 mutations and ≤ 3% for 9 of them.
- The identification of the 3 major IDH1/2 mutations showed very high sensitivity with LOD of 0.78% (R132H), 1.19% (R132C) and 0.61 % (R172K) respectively




Fig 1 - IDH1/2 PCR test sensitivity

CLINICAL VALIDATION COHORT

- 147 samples met the selection criteria (103 academic and 44 commercial)
- The histological distribution reflected the observed subtypes in clinical routine in an academic center, with > 40% GBMs (Fig 2)
- IDH1/2 PCR technical success was 100% on samples collected in academic centers (Fig 3)




Fig 2 - Distribution of the samples by histology (n=147)




Fig 3 - IDH1/2 kit experimental flowchart: technical success rate

RESULTS

CLINICAL PERFORMANCE

1- Comparisons of PCR to IHC and Sanger sequencing

- Overall concordance between IHC and PCR for IDH1 R132H detection was 99% (Tab 2) - The only PCR/IHC discordant case was a sample of the commercial series
- Overall concordance between Sanger sequencing and PCR was 96%
 - PCR detected 5 additional mutated cases (2 IDH1 R132H, 1 IDH1 R132C, 1 IDH1 R132, 1 IDH2 172) compared to Sequencing (Tab 3)
- Positive agreement between PCR and IHC was 98.4 % [91.3;99.7] and was 100% between PCR and sequencing [94.6;100], meeting predefined target (PA ≥95%, lower CI limit ≥90%)
- Sequencing did not identify 3 IHC-positive cases (incl. the PCR-neg case, Tab 2-4)
- Out of the IHC-negative cases (n=72), PCR identified 12 rare mutations (17%, 10 IDH1, 2 IDH2), while only 9 (12 %) were detected by sequencing (Tab 4)

IDH1/2 PCR	IHC			TOTAL
	R132H POS	R132H NEG	TOTAL	
	R132H	NON-R132H	TOTAL	
	60	0	60	
	1	72	73	
	61	72	133	

Table 2 - IHC vs PCR for detection of IDH1/2 R132H mutation

IDH1/2 PCR	SEQUENCING		TOTAL
	IDH1/2 MUT+	IDH1/2 MUT-	
	IDH1/2 MUT+	IDH1/2 MUT-	
	67	5	72
	0	61	61
	67	66	133

Table 3 - Sequencing vs PCR for IDH1/2 mutational status

	WT	IDH1 MUT				IDH2 Mut		% Mutated Cases
		R132H	R132C	R132 OTHER	R100	R172K	R172 OTHER	
IHC	72	61	0	0	0	0	0	46 %
SEQUENCING	66	58	2	6*	0	1	0	50 %
IDH1/2 PCR	61	60	3	7	0	1	1	54 %

* 3 R132S; 2 R132G; 1 R132L

Table 4 - Number and types of IDH 1/2 mutations detected by IHC, sequencing and PCR (n=133)

2- Analysis of discordant cases

CASE	IHC	IDH1/2 PCR	SANGER Seq.	PYRO Seq.	LNA-Seq.	CONCL.	COMMENT
# 1 (COMMERCIAL)	POS	WT	WT	WT	WT	WT	IHC False-Pos
# 2	POS	R132H	WT	-	-	R132H	Low Mut Allele % (15%)*
# 3	POS	R132H	WT	-	-	R132H	Low Mut Allele % (10%)*
# 4 (COMMERCIAL)	NEG	R132C	WT	R132C	-	R132C	Low Mut Allele % (14%)*
# 5 (COMMERCIAL)	NEG	R132	WT	WT	WT	WT	PCR False-Pos**
# 6	NEG	R172	WT	WT	WT	WT	PCR False-Pos**

* Theoretical Mutant Allele % acc. to PCR LOD analysis - ** PCR test result close to LOD value

Table 5 - Analysis of the IHC/PCR and PCR/Sanger sequencing discordant cases

- Cases #1, #5, #6:
Based on consistency of results obtained with the additional and highly sensitive techniques, the mutations detected by IHC (1 case) or PCR (2 cases) are likely false-positive
- Cases #2, #3, #4:
Lack of mutation detection by Sanger Sequencing likely results from a low mutant allele content in the respective samples

3- Synthetic samples

- The IDH1/2 assay correctly detected the 11 IDH1/2 rare mutations at the two tested mutation frequencies (30% and 40%).

SUMMARY AND CONCLUSION

- The newly developed IDH1/2 PCR assay showed
 - High technical success rate
 - High analytical sensitivity, with an LOD < 5% for all but one (rare IDH2) mutations, below published references for sequencing techniques
- Positive concordance with IHC (R132H) and sequencing was high (98% & 100% resp.)
- Out of the 5 PCR/Sanger sequencing discordant cases (< 4%), 3 mutations detected by the PCR assay but not by Sanger Sequencing were confirmed by other sensitive techniques highlighting a higher sensitivity of the IDH1/2 PCR assay.
- The IDH1/2 PCR assay can reliably be performed on FFPE samples of up to 10 yrs of age, which should allow the assay to be used to retrospectively analyse clinical cohorts
- This new IDH1/2 PCR assay is able to detect the major IDH1R132H mutation and 11 rare IDH1/2 mutations in **one step**. This should facilitate the implementation of a comprehensive IDH1/2 testing protocol in routine clinical practice

REFERENCES

1. Parsons DW et al. 2008; Science 321(5897):1807-12
2. Yan H et al. 2009 N Engl J Med 360(8):765-73.
3. Preusser M et al. 2011; Clin Neuropath 30:217-230

4. Hartmann C, et al. Acta Neuropathol. 2009
5. Pusch S et al. 2011; Neuropathol Appl Neurobiol. 37(4):428-30
6. Rohle D et al. 2013; Science. 340(6132):626-30

7. van den Bent MJ et al. 2013; J Neurooncol. 112(2):173-8