

Quick-Start Protocol

RNase H MBG

RNase H MBG (cat. nos. RT34-025 and RT34-125) is a 18.9 kDa recombinant endoribonuclease purified from an *Escherichia coli* strain, which over-expresses cloned RNase H gene (*rnh*). The enzyme hydrolyses specifically the phosphodiester bonds of RNA hybridized to DNA and produces 5' phosphate-terminated oligoribonucleotides and single-stranded DNA. RNase H MBG does not degrade single and double-stranded DNA or unhybridized RNA. It is a key enzyme in the removal of mRNA after first-strand cDNA synthesis. Treating cDNA with RNase H MBG prior to PCR can improve sensitivity as RNA bonded to the cDNA template may prevent binding of the amplification primers in a PCR reaction. RNase H MBG treatment is often necessary when amplifying longer, full-length cDNA targets. In addition, RNase H MBG is useful for the removal of poly(A) tails on mRNAs after hybridization with oligo(dT) and also for the site-specific enzymatic cleavage of RNA. RNase H MBG must be shipped on dry ice and stored at -20°C .

Further information

- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- The activity of RNase H MBG is inhibited by metal chelators (e.g., EDTA) and sulfhydryl SH-blocking reagents.
- Inactivate enzyme by heating at 65°C for 10 minutes.

- RNase H MBG is >90% pure as judged by SDS polyacrylamide gel. The absence of DNase, RNase, and protease activity has been confirmed using the relevant procedures.
- One unit catalyses the hydrolysis of 1 nmol of RNA in [³H]-labeled poly(A)×poly(dT) to acid-soluble ribonucleotides in a total reaction volume of 50 µL in 20 minutes at 37°C in 1x RNase H Buffer.

Procedure

1. Use 5 U of enzyme to remove RNA from a RNA:DNA duplex after reverse transcription in a 20 µL reaction. If 50 µL reaction is desired, the use of 12.5 U of enzyme is recommended.
2. The reaction mixture should be incubated at 37°C for 20 min.

Document Revision History

Date	Changes
08/2023	Initial release

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

August-2023 HB-3453-001 © 2023 QIAGEN, all rights reserved.