

QIAEX II[®] Gel Extraction Kit

The QIAEX II Gel Extraction Kit (cat. nos. 20021 and 20051) can be stored at room temperature (15–25°C) for up to 12 months if not otherwise stated on label.

Further information

- *QIAEX II Handbook*: www.qiagen.com/HB-1167
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- This protocol is for cleanup of DNA fragments of 40 bp to 50 kb.
 - The yellow color of Buffer QX1 indicates a pH ≤ 7.5 .
 - Add ethanol (96–100%) to Buffer PE concentrate before use (see bottle label for volume).
 - A heating block or water bath at 50°C is required.
 - All centrifugation steps are carried out at 17,900 $\times g$ (~13,000 rpm) in a conventional tabletop microcentrifuge at room temperature (15–25°C).
 - For purification of DNA from polyacrylamide gels or aqueous solutions, see the handbook.
1. Excise the DNA band from the agarose gel with a clean, sharp scalpel. Use a 1.5 ml microfuge tube for processing up to 250 mg agarose per tube.
 2. Weigh the gel slice in a colorless tube. Add Buffer QX1 according to DNA fragment size: 6 volumes for <100 bp; 3 volumes for 100 bp – 4 kb; 3 volumes with 2 volumes of water for >4 kb. Add 6 volumes of Buffer QX1 when using >2% or Metaphor agarose gels.
 3. Resuspend QIAEX II by vortexing for 30 s. Add QIAEX II to the sample and mix: Use 10 μ l QIAEX II for ≤ 2 μ g DNA; 30 μ l for 2–10 μ g DNA; and an additional 30 μ l for each additional 10 μ g DNA.



4. Incubate at 50°C for 10 min to solubilize the agarose and bind the DNA. Mix by vortexing* every 2 min to keep QIAEX II in suspension. Check that the color of the mixture is yellow. If the color of the mixture is orange or violet, add 10 µl 3 M sodium acetate, pH 5.0, and mix. The color should turn to yellow. The incubation should then be continued for at least 5 min.
5. Centrifuge the sample for 30 s and carefully remove supernatant with a pipet.
6. Wash the pellet with 500 µl Buffer QX1. Resuspend the pellet by vortexing.* Centrifuge the sample for 30 s and remove all traces of supernatant with a pipet. This wash step removes residual agarose contaminants.
7. Wash the pellet twice with 500 µl Buffer PE. Resuspend the pellet by vortexing.* Centrifuge the sample for 30 s and carefully remove all traces of supernatant with a pipet. This step removes residual salt contaminants.
8. Air-dry the pellet for 10–15 min or until the pellet becomes white. If 30 µl QIAEX II suspension is used, air-dry the pellet for approximately 30 min. Do not vacuum dry, as overdrying, may lead to decreased elution efficiency.
9. To elute DNA, add 20 µl of 10 mM Tris-Cl, pH 8.5, TE buffer (10 mM Tris-Cl, 1 mM EDTA, pH 8.0) or water and resuspend the pellet by vortexing.* Incubate according to the DNA fragment size: 5 min at room temperature (15–25°C) for ≤4 kb; 5 min at 50°C for 4–10 kb; or 10 min at 50°C for >10 kb.
10. Centrifuge for 30 s, and carefully pipet the supernatant into a clean tube. The supernatant now contains the purified DNA.
11. **Optional:** repeat steps 9 and 10 and combine the eluates. A second elution step will increase the yield by approximately 10–15%.

* For fragments larger than 10 kb, resuspend the pellet by inverting and flicking the tube. Vortexing can cause shearing of large DNA fragments.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. Trademarks: QIAGEN®, Sample to Insight®, QIAEX® (QIAGEN Group). 1101214 03/2016 HB-0583-002 © 2016 QIAGEN, all rights reserved.