

## Quick-Start Protocol

# miRNeasy Tissue/Cells Advanced Micro Kit

The miRNeasy Tissue/Cells Advanced Micro Kit (cat. no. 217684) is shipped at ambient temperature. Store the RNeasy® UCP MinElute® spin columns immediately at 2–8°C. Store the remaining components dry at room temperature (15–25°C). All kit components are stable for at least 9 months under these conditions if not otherwise stated on label. This protocol is for purification of total RNA, including small RNAs from animal cells and tissues.

### Further information

- *miRNeasy Tissue/Cells Advanced Micro Kit Handbook*: [www.qiagen.com/HB-2902](http://www.qiagen.com/HB-2902)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](mailto:support.qiagen.com)

### Notes before starting

- Buffers RWT, AL, and RLT contain guanidine salt and are therefore not compatible with disinfecting reagents containing bleach.
- If necessary, redissolve any precipitate in Buffer RLT or Buffer RWT by warming.
- Equilibrate buffers to room temperature.
- All steps should be performed at room temperature. Work quickly.
- Add ethanol (96–100%) to Buffer RWT and Buffer RPE concentrates before use (see bottle label for volume).

- The RNA Spike-in Kit, for RT (cat no. 339390), may be purchased separately. For instructions on preparing a working solution, please refer to the handbook.
- If purifying RNA from cell lines rich in RNases or tissue, add either 10  $\mu$ l  $\beta$ -mercaptoethanol ( $\beta$ -ME) or 20  $\mu$ l 2 M dithiothreitol (DTT) to 1 ml Buffer RLT before use. Buffer RLT containing DTT or  $\beta$ -ME can be stored at room temperature for up to 1 month.

## Procedure

1. **Cells:** Harvest a maximum of  $5 \times 10^5$  cells as a cell pellet or lysed directly in the vessel. Add 260  $\mu$ l of Buffer RLT. Vortex for 30 s or homogenize.  
**Tissues:** Disrupt the tissue ( $\leq 5$  mg \*) and homogenize the lysate in 260  $\mu$ l of Buffer RLT.
2. Add 40  $\mu$ l Buffer AL and mix thoroughly. Incubate at room temperature for 3 min. Transfer the homogenized lysate to a gDNA Eliminator spin column placed in a 2 ml collection tube (supplied).
3. Centrifuge for 30 s at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm). Discard the column and save the flow-through.  
**Optional:** Steps 4 and 5 do not need to be carried out when working with cell samples.
4. Transfer the flow-through to a new 2 ml reaction vessel (not provided). Add 20  $\mu$ l Buffer RPP. Close the tube cap and mix vigorously by vortexing for  $>20$  s. Incubate at room temperature for 3 min.
5. Centrifuge at  $12,000 \times g$  for 3 min at room temperature to pellet the precipitate.  
**Note:** Supernatant should be clear and colorless. Transfer supernatant (for cells: approx. 300  $\mu$ l; for tissue: approx. 270  $\mu$ l) to a new 2 ml reaction tube.
6. Add 1.3 volume isopropanol and mix well by pipetting. Transfer the entire sample to an RNeasy® UCP MinElute spin column. Close the lid and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
7. Repeat step 6 using the remainder of the sample.

\* The amount indicated is applicable for fresh or frozen tissue only. For stabilized tissue, use only half of the amount indicated.

- 
8. Pipet 700  $\mu$ l Buffer RWT to the RNeasy UCP MinElute spin column. Close the lid and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
  9. Pipet 500  $\mu$ l Buffer RPE onto the RNeasy UCP MinElute spin column. Close the lid and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
  10. Add 500  $\mu$ l of 80% ethanol to the RNeasy UCP MinElute spin column. Close the lid and centrifuge for 2 min at  $\geq 8000 \times g$ . Discard the flow-through and the collection tube.
  11. Place the RNeasy UCP MinElute spin column in a new 2 ml collection tube (supplied). Open the lid of the spin column, and centrifuge at full speed for 5 min to dry the membrane. Discard the flow-through and collection tube.
  12. Place the RNeasy UCP MinElute spin column in a new 1.5 ml collection tube (supplied). Add 20  $\mu$ l RNase-free water directly to the center of the spin column membrane and incubate for 1 min. Close the lid and centrifuge for 1 min at full speed to elute the RNA.

## Document Revision History

Date	Changes
04/2021	Initial release



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN® kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, MinElute®, RNeasy® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1124397 04/2021 HB-2904-001 © 2021 QIAGEN, all rights reserved.