

**Second Edition**

**June 2011**

---

# **Qproteome<sup>®</sup> Cell Compartment Handbook**

For the subcellular fractionation of proteomic samples



---

**Sample & Assay Technologies**

## **QIAGEN Sample and Assay Technologies**

QIAGEN is the leading provider of innovative sample and assay technologies, enabling the isolation and detection of contents of any biological sample. Our advanced, high-quality products and services ensure success from sample to result.

### **QIAGEN sets standards in:**

- Purification of DNA, RNA, and proteins
- Nucleic acid and protein assays
- microRNA research and RNAi
- Automation of sample and assay technologies

Our mission is to enable you to achieve outstanding success and breakthroughs. For more information, visit [www.qiagen.com](http://www.qiagen.com).

# Contents

<b>Kit Contents</b>	<b>4</b>
<b>Storage</b>	<b>4</b>
<b>Safety Information</b>	<b>4</b>
<b>Quality Control</b>	<b>5</b>
<b>Product Use Limitations</b>	<b>5</b>
<b>Product Warranty and Satisfaction Guarantee</b>	<b>5</b>
<b>Technical Assistance</b>	<b>6</b>
<b>Introduction</b>	<b>8</b>
<b>Principle and Procedure</b>	<b>8</b>
<b>Protocols</b>	
■ <b>Subcellular Fractionation of Cultured Cell Samples</b>	<b>11</b>
■ <b>Subcellular Fractionation of Tissue Samples</b>	<b>13</b>
■ <b>Acetone Precipitation of Protein Fractions</b>	<b>16</b>
<b>Troubleshooting Guide</b>	<b>17</b>
<b>Ordering Information</b>	<b>18</b>

## Kit Contents

<b>Qproteome Cell Compartment Kit</b>	
<b>Catalog no.</b>	<b>37502</b>
<b>Number of preps</b>	<b>Up to 10</b>
Extraction Buffer CE1	1 x 10 ml
Extraction Buffer CE2	1 x 10 ml
Extraction Buffer CE3	1 x 5 ml
Extraction Buffer CE4	1 x 5 ml
Benzonase® Nuclease	1 x 80 µl
Protease Inhibitor Solution (100x)	1 x 300 µl
QIAshredder	10
Handbook	1

## Storage

Benzonase® Nuclease and Extraction Buffers CE1, CE2, and CE3 should be stored at –20°C.

Protease Inhibitor Solution (100x) should be stored at 2–8°C.

Extraction Buffer CE4 should be stored at room temperature (15–25°C).

## Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs). These are available online in convenient and compact PDF format at [www.qiagen.com/ts/msds.asp](http://www.qiagen.com/ts/msds.asp) where you can find, view, and print the MSDS for each QIAGEN® kit and kit component.

### 24-hour emergency information

Emergency medical information in English, French, and German can be obtained 24 hours a day from:

Poison Information Center Mainz, Germany

Tel: +49-6131-19240

## Quality Control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of Qproteome Cell Compartment Kits is tested against predetermined specifications to ensure consistent product quality.

## Product Use Limitations

The Qproteome Cell Compartment Kit is intended for molecular biology applications. This product is not intended for the diagnosis, prevention, or treatment of a disease.

All due care and attention should be exercised in the handling of the products. We recommend all users of QIAGEN products to adhere to the NIH guidelines that have been developed for recombinant DNA experiments, or to other applicable guidelines.

## Product Warranty and Satisfaction Guarantee

QIAGEN guarantees the performance of all products in the manner described in our product literature. The purchaser must determine the suitability of the product for its particular use. Should any product fail to perform satisfactorily due to any reason other than misuse, QIAGEN will replace it free of charge or refund the purchase price. We reserve the right to change, alter, or modify any product to enhance its performance and design. If a QIAGEN product does not meet your expectations, simply call your local Technical Service Department or distributor. We will credit your account or exchange the product — as you wish.

A copy of QIAGEN terms and conditions can be obtained on request, and is also provided on the back of our invoices. If you have questions about product specifications or performance, please call QIAGEN Technical Services or your local distributor (see back cover or visit [www.qiagen.com](http://www.qiagen.com) ).

## Technical Assistance

At QIAGEN, we pride ourselves on the quality and availability of our technical support. Our Technical Service Departments are staffed by experienced scientists with extensive practical and theoretical expertise in sample and assay technologies and the use of QIAGEN products. If you have any questions or experience any difficulties regarding the Qproteome Cell Compartment Kit or QIAGEN products in general, please do not hesitate to contact us.

QIAGEN customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at QIAGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance and more information, please see our Technical Support center at [www.qiagen.com/goto/TechSupportCenter](http://www.qiagen.com/goto/TechSupportCenter) or call one of the QIAGEN Technical Service Departments or local distributors (see back cover or visit [www.qiagen.com](http://www.qiagen.com) ).

## Cell Compartment Kit Fractionation Procedure

Cell pellet (approx.  
 $5 \times 10^6$  cells)



Add Extraction  
Buffer CE1



Fraction 1  
Cytosolic proteins

Add Extraction  
Buffer CE2



Fraction 2  
Membrane proteins

Add Benzonase and  
Extraction Buffer CE3



Fraction 3  
Nuclear proteins

Add Extraction  
Buffer CE4



Fraction 4  
Cytoskeletal proteins

## Introduction

Eukaryotic cells are complex, well-ordered, and highly structured systems. The Cell Compartment Kit is designed for fast and easy subcellular fractionation of intact eukaryotic cells and tissue. By sequential addition of different extraction buffers to a cell pellet, proteins in the different cellular compartments can be selectively isolated (see Table 1, Figure 1).

## Principle and Procedure

Extraction Buffer CE1 is added to cells and selectively disrupts the plasma membrane without solubilizing it, resulting in the isolation of cytosolic proteins. Plasma membranes and compartmentalized organelles, such as nuclei, mitochondria, and the endoplasmic reticulum (ER), remain intact and are pelleted by centrifugation.

The pellet from the first step is resuspended in Extraction Buffer CE2, which solubilizes the plasma membrane as well as all organelle membranes except the nuclear membrane. After solubilization, the sample is centrifuged. The supernatant contains membrane proteins and proteins from the lumen of organelles (e.g., the ER and mitochondria). The pellet consists of nuclei.

In the next step nuclei are solubilized using Extraction Buffer CE3 in which all soluble and most membrane-bound nuclear proteins are extracted. Addition of Benzonase<sup>®</sup> Nuclease allows the release of proteins tightly bound to nucleic acids (e.g., histones).

After another centrifugation, Extraction Buffer CE4 is used to solubilize all residual — mainly cytoskeletal — proteins in the pellet.

Fractions 1 to 3 contain proteins in their native state. Extraction Buffer CE4 is strongly denaturing and not compatible with isoelectric focusing. Proteins in all fractions must be desalted (e.g., by acetone precipitation, see page 16) before further analysis using isoelectric focusing.

**Table 1. Subcellular protein fractionation from cellular components**

<b>Buffer</b>	<b>Used to isolate proteins from:</b>
Extraction Buffer CE1	Cytosol
Extraction Buffer CE2	Membranes
Extraction Buffer CE3	Nucleus
Extraction Buffer CE4	Cytoskeleton

Starting material for one fractionation procedure is  $5 \times 10^6$  cells or 20 mg tissue. The procedure has been used successfully for several different mammalian cell lines including HeLa, Jurkat, NIH-3T3, HEK293, and Cos. The tissue protocol has been tested with several rat tissues, including liver, heart, lung and brain. Tables 2 and 3 give an overview of expected protein yields using different cell lines.

Subcellular fractionation of proteins enables:

- Enrichment of low-abundance species
- Definition of the subcellular localization of enzymes, regulatory, and structural proteins
- Monitoring of compartmental redistribution of biomolecules under basal and stimulated conditions

**Table 2. Typical protein yields from cultured cells in cytosolic, membrane, nuclear, and cytoskeletal fractions\***


Cell line	Cytosol	Membranes	Nucleus	Cytoskeleton
HeLa	490 $\mu$ g	290 $\mu$ g	160 $\mu$ g	4 $\mu$ g
Jurkat	380 $\mu$ g	120 $\mu$ g	110 $\mu$ g	18 $\mu$ g
NIH-3T3	430 $\mu$ g	270 $\mu$ g	150 $\mu$ g	6 $\mu$ g
Cos	450 $\mu$ g	210 $\mu$ g	130 $\mu$ g	40 $\mu$ g

\* The yields in this table are mean values from 4 independent preparations of  $5 \times 10^6$  cells. The % CV ranged between 5 and 26% for the first 3 fractions.


**Table 3. Typical protein yields from different rat tissues in cytosolic, membrane, nuclear, and cytoskeletal fractions\***

Tissue type	Cytosol	Membranes	Nucleus	Cytoskeleton
Brain	300 $\mu$ g	194 $\mu$ g	29 $\mu$ g	0 $\mu$ g
Liver	630 $\mu$ g	360 $\mu$ g	15 $\mu$ g	0 $\mu$ g
Heart	324 $\mu$ g	122 $\mu$ g	29 $\mu$ g	28 $\mu$ g
Lung	419 $\mu$ g	42 $\mu$ g	35 $\mu$ g	0 $\mu$ g

\* Yields are mean values from 3 independent preparations of 20 mg tissue.



**Figure 1 Subcellular fractionation of NIH-3T3 cells.** Western blots of fractionated NIH-3T3 cells. Protein from fractions 1–4 (20  $\mu$ g) was separated by SDS-PAGE. After western blotting, proteins specific to each fraction were detected using **A** GAPDH and TIM23, **B** annexin, and **C** lamin antibodies, and an HRP-conjugated secondary antibody with chemiluminescent detection.



**Figure 2 Separation of subcellular fractions.** Coomassie<sup>®</sup>-stained gel showing fractionation of NIH-3T3 cells using the Qproteome Cell Compartment Kit.

## Protocol: Subcellular Fractionation of Cultured Cell Samples

This protocol is suitable for processing of  $5 \times 10^6$  cells. **Note:** The extraction buffers contain components that may interfere with protein quantification assays. A precipitation step (e.g., using acetone, see page 16) to remove interfering substances is required to accurately determine protein concentrations.

### Equipment and reagents to be supplied by user

- Ice-cold PBS (50 mM  $\text{NaH}_2\text{PO}_4$ , 150 mM NaCl, pH 7.2)
- Distilled water
- End-over-end shaker
- Optional: acetone stored at  $-20^\circ\text{C}$

### Procedure

1. **Thaw Protease Inhibitor Solution (100x) and Extraction Buffers CE1, CE2, and CE3. After thawing, mix well by vortexing and place on ice. For each fractionation procedure, prepare the volume of buffer supplemented with Protease Inhibitor Solution (100x) given in the table below.**

### Volume of Buffer Required per Fractionation Procedure

	Buffer CE1	Buffer CE2	Buffer CE3	Buffer CE4
Required volume	1 ml	1 ml	0.5 ml	0.5 ml
Protease Inhibitor Solution (100x)	10 $\mu\text{l}$	10 $\mu\text{l}$	5 $\mu\text{l}$	–

2. **Transfer a cell suspension containing  $5 \times 10^6$  cells into a 15 ml conical tube and centrifuge at  $500 \times g$  for 10 min at  $4^\circ\text{C}$ . Remove the supernatant carefully and discard it.**
3. **Resuspend the cell pellet in 2 ml ice-cold PBS by pipetting up and down with a 1 ml pipette tip and transfer the cell suspension into a microcentrifuge tube. Pellet cells by centrifuging at  $500 \times g$  for 10 min at  $4^\circ\text{C}$ . Remove the supernatant carefully and discard it.**
4. **Repeat step 3.**

- 5. Resuspend the cell pellet in 1 ml ice-cold Extraction Buffer CE1 by pipetting up and down using a 1 ml pipette tip. Incubate for 10 min at 4°C on an end-over-end shaker.**

Ensure that Protease Inhibitor Solution (100x) has been added to Extraction Buffer CE1.

- 6. Centrifuge the lysate at 1000 x g for 10 min at 4°C.**
- 7. Carefully transfer the supernatant (fraction 1) into a fresh microcentrifuge tube. Store on ice.**

This fraction primarily contains cytosolic proteins.

- 8. Resuspend the pellet in 1ml ice-cold Extraction Buffer CE2 by pipetting up and down using a 1ml pipette tip. Incubate for 30 min at 4°C on an end-over-end shaker.**

Ensure that Protease Inhibitor Solution (100x) has been added to Extraction Buffer CE2.

- 9. Centrifuge the suspension at 6000 x g for 10 min at 4°C.**
- 10. Carefully transfer the supernatant (fraction 2) into a fresh microcentrifuge tube. Store on ice.**

This fraction primarily contains membrane proteins.

- 11. Add 7  $\mu$ l Benzonase<sup>®</sup> Nuclease and 13  $\mu$ l distilled water to the pellet. Resuspend the pellet by gently flicking the bottom of the tube. Incubate for 15 min at room temperature (15–25°C).**

- 12. Pipet 500  $\mu$ l ice-cold Extraction Buffer CE3 into the tube and pipet up and down using a 1ml pipette tip. Incubate for 10 min at 4°C on an end-over-end shaker.**

Ensure that Protease Inhibitor Solution (100x) has been added to Extraction Buffer CE3.

- 13. Pellet insoluble material by centrifuging at 6800 x g for 10 min at 4°C.**

- 14. Transfer the supernatant (fraction 3) into a fresh sample tube. Store on ice.**

This fraction primarily contains nuclear proteins.

- 15. Resuspend the pellet from step 13 in 500  $\mu$ l Extraction Buffer CE4. Label the suspension fraction 4.**

This fraction primarily contains cytoskeletal proteins.

## Protocol: Subcellular Fractionation of Tissue Samples

This protocol is suitable for processing of 20 mg tissue. **Note:** The extraction buffers contain components that may interfere with protein quantification assays. A precipitation step (e.g., using acetone, see page 16) to remove interfering substances is required to accurately determine protein concentrations.

For most tissues tested by QIAGEN (heart, lung, and liver) there is sufficient buffer CE1 in the kit for a total of 10 preparations. However, some tissues (e.g., brain) require double the amount of Extraction Buffer CE1 and therefore a total of 5 preps can be performed per kit.

### Equipment and reagents to be supplied by user

- Ice-cold PBS (50 mM NaH<sub>2</sub>PO<sub>4</sub> , 150 mM NaCl, pH 7.2)
- Distilled water
- End-over-end shaker
- Disruption/homogenization device, e.g., TissueRuptor™
- Optional: acetone stored at -20°C

### Procedure

1. **Thaw Protease Inhibitor Solution (100x) and Extraction Buffers CE1, CE2, and CE3. After thawing, mix well by vortexing and place on ice. For each fractionation procedure, prepare the volume of buffer supplemented with Protease Inhibitor Solution (100x) given in the table below.**

When processing brain tissue sample prepare 2 ml Buffer CE1.

### Volume of Extraction Buffer Required per Fractionation Procedure

	Buffer CE1	Buffer CE1*	Buffer CE2	Buffer CE3	Buffer CE4
Required volume	1 ml	2 ml	1 ml	0.5 ml	0.5 ml
Protease Inhibitor Solution (100x)	10 µl	20 µl	10 µl	5 µl	–

\* When processing brain tissue sample prepare 2 ml Buffer CE1.

- 2. Using a clean scalpel, cut the tissue into 3-4 pieces and wash with 1 ml ice-cold PBS buffer.**
- 3. Place the tissue pieces into a clean 2 ml microcentrifuge tube and add 500  $\mu$ l Buffer CE1 supplemented with Protease Inhibitor Solution (prepared in step 1).**
- 4. Using the TissueRuptor, disrupt the tissue for 5 s at the lowest speed setting.**

Do not disrupt the tissue for longer than 5 s. If clumps of tissue remain they will be filtered out in the next step.
- 5. Transfer the tissue suspension to a QIAshredder homogenizer and centrifuge at 510 x g for 2 min at 4°C.**

Ensure that Protease Inhibitor Solution (100x) has been added to Extraction Buffer CE1.
- 6. Resuspend the pellet by gently pipetting up and down and transfer the suspension to a clean 2 ml tube.**
- 7. Add 500  $\mu$ l Buffer CE1 (prepared in step 1). NOTE: When processing brain tissue add 1.5 ml Buffer CE1 in this step.**
- 8. Incubate the suspension on an end-over-end shaker for 10 min at 4°C.**
- 9. Centrifuge the lysate at 4000 x g for 10 min at 4°C.**
- 10. Carefully transfer the supernatant (fraction 1) into a fresh microcentrifuge tube. Store on ice.**

This fraction primarily contains cytosolic proteins.
- 11. Resuspend the pellet from step 10 in 1 ml ice-cold Extraction Buffer CE2 by pipetting up and down using a 1 ml pipette tip. Incubate for 30 min at 4°C on an end-over-end shaker.**

Ensure that Protease Inhibitor Solution (100x) has been added to Extraction Buffer CE2.
- 12. Centrifuge the suspension at 6000 x g for 10 min at 4°C.**
- 13. Carefully transfer the supernatant (fraction 2) into a fresh microcentrifuge tube. Store on ice.**

This fraction primarily contains membrane proteins.
- 14. Add 7  $\mu$ l Benzonase<sup>®</sup> Nuclease and 13  $\mu$ l distilled water to the pellet. Resuspend the pellet by gently flicking the bottom of the tube. Incubate for 15 min at room temperature (15–25°C).**
- 15. Pipet 500  $\mu$ l ice-cold Extraction Buffer CE3 into the tube and pipet up and down using a 1ml pipette tip. Incubate for 10 min at 4°C on an end-over-end shaker.**

Ensure that Protease Inhibitor Solution (100x) has been added to Extraction Buffer CE3.

**16. Pellet insoluble material by centrifuging at 6800 x g for 10 min at 4°C.**

**17. Transfer the supernatant (fraction 3) into a fresh sample tube. Store on ice.**

This fraction primarily contains nuclear proteins.

**18. Resuspend the pellet from step 13 in 500  $\mu$ l Extraction Buffer CE4. Label the suspension fraction 4.**

This fraction primarily contains cytoskeletal proteins.

## Protocol: Acetone Precipitation of Protein Fractions

This protocol is suitable for concentrating and desalting protein samples for downstream applications.

- 1. Add four volumes of ice-cold acetone to the protein fraction and incubate for 15 min on ice.**
- 2. Centrifuge for 10 min at 12,000 x g in a pre-cooled microcentrifuge at 4°C. Discard the supernatant and air dry the pellet.**  
Do not overdry the pellet as this may make it difficult to resuspend.
- 3. Depending on the application, resuspend the pellet in the required sample buffer.**

For 2D-PAGE, an extra desalting step may be required.

- 4. Resuspend the pellet from step 2 in 100  $\mu$ l 8M urea.**
- 5. Desalt the sample using a gel filtration device (e.g., Bio-Spin<sup>®</sup> 6, Bio-Rad cat. no. 732-6227).**
- 6. Repeat steps 1 to 3.**

## Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about either the information and protocols in this handbook or sample and assay technologies (for contact information, see back cover or visit [www.qiagen.com](http://www.qiagen.com) ).

	<b>Comments and Suggestions</b>
Inconsistent results in protein quantification assays	The extraction buffers contain components that may interfere with protein quantification assays. A precipitation step (e.g., using acetone, see page 16) to remove interfering substances is required to accurately determine protein concentrations.
Marker proteins do not appear in expected fraction/appear in different fraction	<p>The amount of cells or tissue processed was too high. The protocol is suitable for processing of <math>5 \times 10^6</math> cells or 20 mg tissue.</p> <p>Marker protein may shuffle between compartments, for example, upon apoptotic stimulus cytochrome c can be transported from mitochondria (fraction 2) to the cytosol (fraction 1).</p>
A precipitate forms when storing Extraction Buffer CE4 on ice.	The precipitate can be redissolved by heating the buffer to 37°C with agitation and cooling to room temperature.
Protease Inhibitor does not thaw at room temperature	Heat the protease inhibitor solution to 37°C with agitation and cool to room temperature.

## Ordering Information

Product	Contents	Cat. no.
Qproteome Cell Compartment Kit	For up to 10 subcellular fractionations: Extraction buffers, Protease Inhibitor Solution, Benzonase® Nuclease	37502
<b>Related products</b>		
TissueRuptor	Handheld rotor–stator homogenizer, 5 TissueRuptor Disposable Probes	9001271* 9001272† 9001273‡ 9001274§
TissueRuptor Disposable Probes (25)	25 nonsterile plastic disposable probes for use with the TissueRuptor	990890
QIAshredder (50)	50 disposable cell-lysate homogenizers for use in nucleic acid minipreps, caps	79654

\* 120 V, 60 Hz (for North America and Japan); † 235 V, 50/60 Hz (for Europe, excluding UK and Ireland); ‡ 235 V, 50/60 Hz (for UK and Ireland); § 235 V, 50/60 Hz (for Australia)

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Qproteome®, (QIAGEN Group); Benzonase® (Merck KGaA, Germany); Bio-Spin® (Bio-Rad Laboratories); Coomassie® (ICI [Imperial Chemical Industries] Organics Inc.).

Benzonase® Nuclease is supplied by Merck KGaA and its Affiliates. Benzonase® is a registered trademark of Merck KGaA, Darmstadt, Germany.

#### **Limited License Agreement**

Use of this product signifies the agreement of any purchaser or user of the Qproteome Cell Compartment Kit to the following terms:

1. The Qproteome Cell Compartment Kit may be used solely in accordance with the *Qproteome Cell Compartment Handbook* and for use with components contained in the Kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this Kit with any components not included within this Kit except as described in the *Qproteome Cell Compartment Handbook* and additional protocols available at [www.qiagen.com](http://www.qiagen.com).
2. Other than expressly stated licenses, QIAGEN makes no warranty that this Kit and/or its use(s) do not infringe the rights of third-parties.
3. This Kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the Kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the Kit and/or its components.

For updated license terms, see [www.qiagen.com](http://www.qiagen.com).

© 2007–2011 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)

**Australia** ■ Orders 1-800-243-800 ■ Fax 03-9840-9888 ■ Technical 1-800-243-066

**Austria** ■ Orders 0800-28-10-10 ■ Fax 0800-28-10-19 ■ Technical 0800-28-10-11

**Belgium** ■ Orders 0800-79612 ■ Fax 0800-79611 ■ Technical 0800-79556

**Brazil** ■ Orders 0800-557779 ■ Fax 55-11-5079-4001 ■ Technical 0800-557779

**Canada** ■ Orders 800-572-9613 ■ Fax 800-713-5951 ■ Technical 800-DNA-PREP (800-362-7737)

**China** ■ Orders 86-21-3865-3865 ■ Fax 86-21-3865-3965 ■ Technical 800-988-0325

**Denmark** ■ Orders 80-885945 ■ Fax 80-885944 ■ Technical 80-885942

**Finland** ■ Orders 0800-914416 ■ Fax 0800-914415 ■ Technical 0800-914413

**France** ■ Orders 01-60-920-926 ■ Fax 01-60-920-925 ■ Technical 01-60-920-930 ■ Offers 01-60-920-928

**Germany** ■ Orders 02103-29-12000 ■ Fax 02103-29-22000 ■ Technical 02103-29-12400

**Hong Kong** ■ Orders 800 933 965 ■ Fax 800 930 439 ■ Technical 800 930 425

**Ireland** ■ Orders 1800 555 049 ■ Fax 1800 555 048 ■ Technical 1800 555 061

**Italy** ■ Orders 800-789-544 ■ Fax 02-334304-826 ■ Technical 800-787980

**Japan** ■ Telephone 03-6890-7300 ■ Fax 03-5547-0818 ■ Technical 03-6890-7300

**Korea (South)** ■ Orders 080-000-7146 ■ Fax 02-2626-5703 ■ Technical 080-000-7145

**Luxembourg** ■ Orders 8002-2076 ■ Fax 8002-2073 ■ Technical 8002-2067

**Mexico** ■ Orders 01-800-7742-639 ■ Fax 01-800-1122-330 ■ Technical 01-800-7742-436

**The Netherlands** ■ Orders 0800-0229592 ■ Fax 0800-0229593 ■ Technical 0800-0229602

**Norway** ■ Orders 800-18859 ■ Fax 800-18817 ■ Technical 800-18712

**Singapore** ■ Orders 1800-742-4362 ■ Fax 65-6854-8184 ■ Technical 1800-742-4368

**Spain** ■ Orders 91-630-7050 ■ Fax 91-630-5145 ■ Technical 91-630-7050

**Sweden** ■ Orders 020-790282 ■ Fax 020-790582 ■ Technical 020-798328

**Switzerland** ■ Orders 055-254-22-11 ■ Fax 055-254-22-13 ■ Technical 055-254-22-12

**UK** ■ Orders 01293-422-911 ■ Fax 01293-422-922 ■ Technical 01293-422-999

**USA** ■ Orders 800-426-8157 ■ Fax 800-718-2056 ■ Technical 800-DNA-PREP (800-362-7737)

