



October 2025

# QIAstat-Dx<sup>®</sup> Gastrointestinal Panel 2 Instructions for Use



Version 2



For In Vitro Diagnostic Use



For prescription use only

For use with QIAstat-Dx Analyzer 2.0



QIAGEN, GmbH, QIAGEN Strasse 1, 40724 Hilden, GERMANY



691421

R4

# Table of Contents

- Intended Use ..... 4
- Summary and Explanation ..... 6
  - Pathogen information ..... 6
  - Summary of detected organisms ..... 8
  - QIAstat-Dx Gastrointestinal Panel 2 Cartridge description ..... 14
- Materials Provided ..... 18
  - Kit contents ..... 18
- Materials Required but Not Provided ..... 19
- Warnings and Precautions ..... 20
  - Safety information ..... 20
  - Precautions ..... 22
- Reagent Storage and Handling ..... 24
  - In-use stability ..... 24
- Specimen Storage and Handling ..... 25
  - Specimen collection ..... 25
- Procedure ..... 26
  - Protocol: Stool samples in transport medium ..... 26
- Interpretation of Results ..... 37
  - Viewing results ..... 37
  - Viewing amplification curves ..... 39
  - Viewing test details ..... 41
  - Browsing results from previous tests ..... 43
  - Sample result interpretation ..... 47
- Quality Control ..... 51
  - Internal control interpretation ..... 51
  - External control information ..... 51
- Limitations ..... 52
- Performance Characteristics ..... 58
  - Analytical performance ..... 58
  - Clinical performance ..... 100

Troubleshooting Guide ..... 125

Contact Information ..... 126

Symbols ..... 127

Appendices ..... 130

    Appendix A: Installing the Assay Definition File ..... 130

    Appendix B: Glossary ..... 135

    Appendix C: Additional instructions for use ..... 136

Ordering Information ..... 138

References ..... 139

Document Revision History ..... 149

# Intended Use

The QIAstat-Dx Gastrointestinal Panel 2 is a multiplexed nucleic acid test intended for use with the QIAstat-Dx Analyzer 2.0 for the simultaneous in vitro qualitative detection and identification of nucleic acids from multiple viruses, bacteria, and parasites directly from preserved stool samples (Para-Pak® C&S or FecalSwab™) obtained from individuals with signs and/or symptoms of gastrointestinal infection. The following viruses, bacteria (including several diarrheagenic *E. coli*/*Shigella* pathotypes), and parasites are identified with the QIAstat-Dx Gastrointestinal Panel 2:

- Adenovirus F40/F41
- Astrovirus
- Norovirus GI/GII
- Rotavirus A
- *Campylobacter* (*C. jejuni*, *C. coli*, and *C. upsaliensis*)
- *Shigella*/Enteroinvasive *Escherichia coli* (EIEC)
- Enteropathogenic *Escherichia coli* (EPEC)
- Enterotoxigenic *Escherichia coli* (ETEC) lt/st
- Shiga- like toxin- producing *Escherichia coli* (STEC) *stx1/stx2* (including specific identification of *E. coli* O157 serogroup within STEC)
- *Salmonella*
- *Plesiomonas shigelloides*
- *Yersinia enterocolitica*
- *Cryptosporidium*

- *Cyclospora cayetanensis*
- *Entamoeba histolytica*
- *Giardia lamblia*\*

Concomitant culture is necessary for organism recovery and further typing of bacterial agents. The QIAstat-Dx Gastrointestinal Panel 2 is indicated as an aid in the diagnosis of specific agents of gastrointestinal illness, in conjunction with other clinical, laboratory, and epidemiological data. Positive results do not rule out co-infection with organisms not detected by the QIAstat-Dx Gastrointestinal Panel 2. The organisms detected may not be the sole or definitive cause of the disease.

Negative QIAstat-Dx Gastrointestinal Panel 2 results in the setting of clinical illness compatible with gastroenteritis may be due to infection by pathogens that are not detected by this assay test or non-infectious causes such as ulcerative colitis, irritable bowel syndrome, or Crohn's disease.

\*Also known as *Giardia intestinalis* and *Giardia duodenalis*

# Summary and Explanation

## Pathogen information

Acute gastrointestinal infections can be caused by a variety of pathogens, including parasites, bacteria, and viruses, and generally present with nearly indistinguishable clinical signs and symptoms (1). The rapid and accurate determination of the presence or absence of potential causative agent(s) help make timely decisions regarding treatment, hospital admission, infection control, and return of the patient to work and family (2,3,4). It may also greatly support improved antimicrobial stewardship and other important public health initiatives (3,5).

The QIAstat-Dx Gastrointestinal Panel 2 Cartridge allows detection and differentiation of 16 parasitic, viral, and bacterial pathogens that cause gastrointestinal symptoms. Testing requires a small sample volume and minimal hands-on time, and the results are available in approximately 78 minutes.

Pathogens that can be detected and identified with the QIAstat-Dx Gastrointestinal Panel 2 are listed in Table 1.

Table 1. Pathogens detected by the QIAstat-Dx Gastrointestinal Panel 2

Pathogen	Classification (genome type)
Adenovirus F40/F41	Adenovirus (DNA)
Astrovirus	Astrovirus (RNA)
Norovirus GI/GII	Calicivirus (RNA)
Rotavirus A	Reovirus (RNA)
Campylobacter ( <i>C. jejuni</i> , <i>C. upsaliensis</i> , <i>C. coli</i> )	Bacterium (DNA)
Enteroinvasive <i>E. coli</i> (EIEC)/ <i>Shigella</i>	Bacterium (DNA)
Enteropathogenic <i>E. coli</i> (EPEC)	Bacterium (DNA)
Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i>	Bacterium (DNA)
<i>Plesiomonas shigelloides</i>	Bacterium (DNA)
<i>Salmonella</i> spp.	Bacterium (DNA)
Shiga-like toxin-producing <i>E. coli</i> (STEC) <i>stx1/stx2</i> (including specific identification of <i>E. coli</i> O157 serogroup within STEC)	Bacterium (DNA)
<i>Yersinia enterocolitica</i>	Bacterium (DNA)
<i>Cryptosporidium</i>	Parasite (DNA)
<i>Cyclospora cayetanensis</i>	Parasite (DNA)
<i>Entamoeba histolytica</i>	Parasite (DNA)
<i>Giardia lamblia</i>	Parasite (DNA)

**Note:** Shiga-like toxin-producing *E. coli* (STEC) *stx1* and *stx2* are grouped together and reported as Shiga-like toxin-producing *E. coli* (STEC) *stx1/stx2*.

**Note:** *Giardia lamblia* is also known as *Giardia intestinalis* and *Giardia duodenalis*.

## Summary of detected organisms

### Bacteria

***Campylobacter* spp. (*C. jejuni*/*C. coli*/*C. upsaliensis*)** is a genus of gram-negative bacteria that includes more than 30 species (6). *Campylobacter jejuni* and *Campylobacter coli* are the most common *Campylobacter* species associated with diarrheal illness, being *C. jejuni* responsible for 90% of cases (7,9). The consumption of undercooked poultry or raw milk are the most common source of *Campylobacter* infections (10,11). *Campylobacter* are highly infectious, with an infectious dose as low as 500 bacteria (12); however, person-to-person spread is uncommon (10). Systemic disease, associated with significant morbidity and mortality, may occur in individuals who are immunocompromised (10,12). Infection can result in long-term consequences such as arthritis, irritable bowel syndrome, and Guillain-Barré syndrome (10,12).

***Plesiomonas shigelloides*** is a facultatively anaerobic gram-negative bacterium that can cause enteric disease in humans. The prevalence of *P. shigelloides* enteritis varies considerably, with higher rates reported from Southeast Asia and Africa and lower numbers from North America and Europe. It is unknown how many people suffer from disease caused by *P. shigelloides* each year, but mortality is rare. Infection especially occurs following the consumption of raw seafood or contaminated water (13).

***Salmonella*** is a gram-negative bacterium comprising more than 2600 serovars, including the distinct typhoidal serotypes, Typhi and Paratyphi A–C (14,15). Enteric (typhoid) fever is an invasive, life-threatening, systemic infection with predominantly non-gastrointestinal symptoms (14,16). Non-typhoidal salmonellosis is an acute, usually self-limiting gastroenteritis that is characterized by symptoms such as watery diarrhea, fever, abdominal pain, and sometimes vomiting (14,16,17). Less commonly, non-typhoidal *Salmonella* serovars cause invasive disease due to bloodstream infections that are not usually associated with diarrhea (14,16). There are 100–200 million cases of non-typhoidal salmonellosis each year, resulting in



approximately 85,000–155,000 deaths (16,18). The incidence of non-typhoidal *Salmonella* gastroenteritis is highest in the developing world but is also of considerable importance in developed countries (14).

***Yersinia enterocolitica*** is a gram-negative bacterium that has more than 70 serotypes (19); serotypes most commonly associated with infection are O:3, O:9, O:8, and O:5,27 (20). *Y. enterocolitica* infection have been reported frequently in northern Europe, particularly in Belgium, Norway, and the Netherlands; it is rarely observed in tropical countries (21). *Y. enterocolitica* is usually transmitted through the consumption of raw meats, unpasteurized dairy products, contaminated water, or via the fecal–oral route (22). Symptoms range from self-limiting enteritis with diarrhea, low-grade fever, and abdominal pain to severe disease such as terminal ileitis and mesenteric lymphadenitis, which also mimics appendicitis (23,24,25).

### Diarrheagenic *Escherichia coli*/*Shigella*

***E.coli* / *Shigella*** are gram-negative facultative anaerobic bacteria belonging to the Enterobacteriaceae family. In addition to being part of the normal intestinal microflora of mammals, *E. coli*/*Shigella* contains several pathotypes that cause a variety of diseases (26,27). There are four major pathotypes of diarrhoeagenic *E. coli*/*Shigella*, which each have unique features in their interaction with eukaryotic cells: Enteropathogenic *E. coli* (EPEC), Enterohaemorrhagic *E. coli*/Shiga-like toxin-producing *E. coli* (EHEC/STEC), Enterotoxigenic *E. coli* (ETEC), and Enteroinvasive *E. coli* (EIEC)/*Shigella* (26,27). *E. coli*/*Shigella* have a conserved core genome and a flexible gene pool containing virulence and fitness genes, which are carried on mobile genetic elements (26,27). Gene gain, via horizontal transfer, and gene loss afford the pathogenic traits to *E. coli*/*Shigella* that give rise to the different pathotypes (27).

**Enteroinvasive *E. coli* (EIEC) and *Shigella*.** EIEC is an invasive strain of *E. coli* that is very closely related in virulence and other pathogenic properties to *Shigella* (28,29). Sequencing

indicates that EIEC is more related to *Shigella* than to non-invasive *E. coli*; however, they are currently classified as distinct species (26,28,30). The virulence of this pathogen is primarily due a plasmid-encoding virulence factors that allow the adhesion and invasion to the epithelial cells (31). EIEC is under-represented in epidemiological studies due to its less severe manifestation and potential misclassification as *Shigella* (27). EIEC infection often leads only to self-limiting, mild watery diarrhea; in rare situations, it can cause symptoms of shigellosis, but complications are uncommon (27). *Shigella* is the second-leading cause of diarrhea mortality, causing approximately 13% diarrhea deaths (32). Numbers of deaths are greatest in young children and the elderly (32). It is recommended that individuals with shigellosis should not take anti-diarrheal medications such as loperamide, as these can make symptoms worse (33).

**Enteropathogenic *E. coli* (EPEC)** is primarily a disease of infants <2 years (27,34,35) and is commonly present in co-infections with other gastrointestinal pathogens (36). EPEC are classified into typical (tEPEC) and atypical (aEPEC) strains based on the presence of the *E. coli* adherence factor plasmid (pEAF). tEPEC is considered an important cause of infantile diarrhea in developing countries (37). Infections in adults, including travelers to developing countries, are rarely reported (27,35). aEPEC is frequently detected in both developing countries and industrialized countries and is suggested to be more prevalent than tEPEC (34). aEPEC is an important cause of both endemic diarrhea and outbreaks (34).

**Enterotoxigenic *E. coli* (ETEC)** is characterized by the production of heat-labile enterotoxins (LT) and heat-stable enterotoxins (ST) (38,39). ETEC is the most common diarrhea-associated *E. coli* and although infections are usually self-limiting (39), is the eighth-leading cause of diarrhea globally and accounts for >50,000 deaths every year (32). It also remains a major cause of diarrhea in travelers to low resource countries (39). ETEC is frequently antimicrobial resistant (39).

**Shiga-like toxin-producing *E. coli* (STEC) *stx1/stx2*, including *E. coli* O157**, is defined by the production of Shiga toxin 1 (Stx1) or 2 (Stx2), which show homology to Stx toxins from

*Shigella dysenteriae* (27). There are >400 serotypes of STEC, of which O157:H7 is the most common (27). Symptoms of STEC infection range from mild intestinal disease to hemorrhagic diarrhea and can lead to hemolytic uremic syndrome (HUS), end-stage renal disease, and death (27,40). Approximately 5–10% of individuals diagnosed with STEC infections develop HUS, which can be a life-threatening complication (41). The impacts of STEC are often greater in infants and children, compared to other ages (40). Antibiotics should not be used to treat STEC infections as there is currently no evidence that they aid recovery and have instead been associated with worsening of symptoms and the development of HUS (41).

## Parasites

***Cryptosporidium* spp.** are protozoan parasites that can infect humans and other animals, with *C. hominis* and *C. parvum* species responsible for the majority of human infections (42). *Cryptosporidium* spp. are found globally. In developing countries, particularly in sub-Saharan Africa, there is a greater risk of infection due to poorer water treatment and food sanitation (32,43). It is also one of the leading causes of diarrheal mortality in children <5 years of age (32,44).

***Cyclospora cayetanensis*** is a single-celled protozoa parasite, and the only known species of the genus *Cyclospora* to infect humans (45,46). It is endemic in tropical/subtropical areas. In non-endemic regions, cases and outbreaks of cyclosporiasis are usually linked to international travel and consumption of contaminated produce imported from endemic regions (45,46,47). Direct fecal-oral transmission cannot occur, because the unsporulated oocysts shed in feces sporulate in water and food environments before infecting another host (45,46,48).

***Entamoeba histolytica*** is an anaerobic, protozoan parasite (49). *E. histolytica* is common in developing countries, particularly those in the tropics and sub-tropics with poor sanitation (49,50,51). Only 10–20% of individuals infected with *E. histolytica* are symptomatic (1,2). Through destruction of the intestinal walls, trophozoites can also spread systemically to the

liver, lungs, and central nervous system (49,50,51,89). The liver is the most common extra-intestinal site of infection (49,50,51).

***Giardia lamblia*** (also referred to as ***G. duodenalis*** and ***G. intestinalis***) is a unicellular, protozoan parasite that can cause disease in humans and other mammals (52,53). *G. lamblia* has a global distribution and is common in both children and adults (54,55). Prevalence of infection is higher in developing regions of the world and in children (52,54,55). The majority (50–75%) of *G. lamblia* infections are asymptomatic (56). In immunocompetent individuals, infections are usually self-limiting, although some may become chronic (52).

## Viruses

**Adenovirus F40/41** is a double-stranded DNA, non-enveloped virus (57,58), with many distinct serotypes described and classified into 7 species (A–G) (57). Serotypes F40/41 are the most common cause of acute gastroenteritis in young children, causing 5–20% of reported cases. More than 80% of diagnosed infections occur in children aged <4 years (58). Adenoviruses have a worldwide distribution, and infections occur throughout the year without significant seasonal variability (57). Infections are usually mild and self-limiting in immunocompetent individuals but can be fatal in individuals who are immunocompromised (57,59,60).

**Astroviruses** are non-enveloped, positive-sense, single-stranded RNA viruses (61). Human astroviruses are distributed all over the world and are associated with 2–9% of cases of acute, nonbacterial diarrhea in children (61,62). It is estimated that 90% of the global population aged ≥9 years have antibodies against astrovirus type 1 (61). Astrovirus infection in healthy children and adults is mostly asymptomatic, although Astrovirus can cause severe diarrhea in children, older adults and immunocompromised patients or those with comorbidities (61,62).

**Noroviruses GI/GII** are small, non-enveloped, positive-stranded RNA viruses from the family Caliciviridae (63). They are responsible for >90% of viral gastroenteritis and around 50% of

all-cause gastroenteritis outbreaks globally (64), causing approximately 685 million cases every year (65). Approximately 200 million cases are in children aged <5 years, leading to 50,000 child deaths (65). Norovirus is commonly known as the “winter-vomiting bug”; outbreaks are more common during the winter months, but infection can occur at any time of year (65). Norovirus is infectious at very low doses and is transmitted via aerosolized droplets and by touching contaminated surfaces (65). Individuals infected with norovirus usually recover within 1–3 days. Infections in infants, older adults, and immunocompromised individuals can be severe and sometimes fatal (65). In some individuals, viral shedding can occur for many weeks/months after symptoms stop, a large contributing factor for outbreaks (66).

**Rotavirus A** is a non-enveloped, double-stranded RNA virus of the Reoviridae family, with 10 species that cause infection in humans (A–J). Rotavirus A is the most common species and causes >90% of all rotavirus infections (67,68). Rotavirus is a leading cause of diarrhea in children <5 years (67), with a seasonal infection pattern that differs across the world, particularly in middle–high income countries (69). Severe infection is most common in young children and infants. In adults, infections are often associated with milder symptoms (70). Two oral rotavirus vaccines are approved in the United States (71) and have been available in >100 countries since 2006 (71). These vaccines have substantially reduced the burden of rotavirus-associated illness (70).

# QIAstat-Dx Gastrointestinal Panel 2 Cartridge description


The QIAstat-Dx Gastrointestinal Panel 2 Cartridge (Figure 1) is a disposable plastic device that allows performance of fully automated molecular assays for the detection of gastrointestinal pathogens. Main features of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge include compatibility with a liquid sample type, hermetical containment of pre-loaded reagents necessary for testing, and walk-away operation. All sample preparation and assay testing steps are performed within the cartridge.

All reagents required for the complete execution of a test run are pre-loaded and self-contained in the QIAstat-Dx Gastrointestinal Panel 2 Cartridge. The user does not need to come in contact with or manipulate any reagents. The QIAstat-Dx Analyzer 2.0 houses air filters for both incoming and outgoing air, further safeguarding the environment. After testing, the cartridge stays hermetically closed at all times, greatly enhancing its safe disposal.

Within the cartridge, multiple steps are automatically performed in sequence using pneumatic pressure to transfer samples and fluids via the transfer chamber to their intended destinations.

## Description of the process

After the sample is manually loaded, the diagnostic tests within the QIAstat-Dx Gastrointestinal Panel 2 are performed on the QIAstat-Dx Analyzer 2.0. All of the sample preparation and analysis steps are performed automatically by the QIAstat-Dx Analyzer 2.0.



**Figure 1. Layout of the QIAstat-Dx Cartridge and its features. \*Note: The swab port is not used for the QIAstat-Dx Gastrointestinal Panel 2.**

## Sample collection and cartridge loading

The collection of samples and their subsequent loading into the QIAstat-Dx Gastrointestinal Panel 2 Cartridge should be performed by personnel trained in safe handling of biological samples.

The following steps are performed:

1. Fresh unpreserved stool specimen is collected and resuspended in Para-Pak C&S or FecalSwab transport medium as soon as possible after collection following the manufacturer's instructions. Attention should be given not to exceed the maximum fill line of the Para-Pak C&S or FecalSwab container or overfill the FecalSwab collection device.
2. The sample information is manually written on or a sample label is affixed to the top of a QIAstat-Dx Gastrointestinal Panel 2 Cartridge.
3. Liquid sample (stool resuspended in Para-Pak C&S or FecalSwab transport medium) is loaded manually into the QIAstat-Dx Gastrointestinal Panel 2 Cartridge.

**Note:** Preserved stool specimens should present a homogenous suspension (easily vortexed).

**Note:** The user must perform a visual check of the sample inspection window to confirm that the liquid sample has been loaded.

4. The sample barcode (if available) and the QIAstat-Dx Gastrointestinal Panel 2 Cartridge barcode are scanned by the QIAstat-Dx Analyzer 2.0. If the sample barcode is not available, the sample ID is manually written using the virtual keyboard of the touchscreen.
5. The QIAstat-Dx Gastrointestinal Panel 2 Cartridge is introduced into the QIAstat-Dx Analyzer 2.0.
6. The test is started on the QIAstat-Dx Analyzer 2.0.

## Sample preparation, nucleic acid amplification, and detection

The extraction, amplification, and detection of nucleic acids in the sample are performed automatically by the QIAstat-Dx Analyzer 2.0.

1. The sample is pre-treated with buffer and homogenized.
2. Resuspension of Internal Control using on-cartridge buffer and mixing with the sample.



3. Cells are lysed in the lysis chamber of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge, which includes a rotor that turns at high speed and silica beads that provide effective cell disruption.
4. Nucleic acids are purified from the lysed sample via binding to a silica membrane in the purification chamber of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge in the presence of chaotropic salts and alcohol.
5. The purified nucleic acids are eluted from the membrane in the purification chamber and are mixed with the lyophilized PCR chemistry in the dried-chemistry chamber of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge.
6. The mixture of sample and PCR reagents is dispensed into the QIAstat-Dx Gastrointestinal Panel 2 Cartridge PCR chambers, which contain air-dried, assay-specific primers and probes.
7. The QIAstat-Dx Analyzer 2.0 creates the optimal temperature profiles to carry out effective multiplex real-time RT-PCR and performs real-time fluorescence measurements to generate amplification curves.
8. The QIAstat-Dx Analyzer 2.0 Software interprets the resulting data and process controls and delivers a test report.

# Materials Provided

## Kit contents

<b>QIAstat-Dx Gastrointestinal Panel 2 Cartridge</b>	
<b>Catalog number</b>	<b>691421</b>
<b>Number of tests</b>	<b>6</b>
<hr/>	
QIAstat-Dx Gastrointestinal Panel 2 Cartridges*	6
Transfer pipettes†	6
QIAstat-Dx Gastrointestinal Panel 2 Product information Card	1
<hr/>	

\* 6 individually packaged cartridges containing all reagents needed for sample preparation and multiplex real-time RT-PCR, plus Internal Control.

† 6 individually packaged transfer pipettes for dispensing liquid sample into the QIAstat-Dx Gastrointestinal Panel 2 Cartridge.

# Materials Required but Not Provided

## Platform and software

The QIAstat-Dx Gastrointestinal Panel 2 is designed for use with the QIAstat-Dx Analyzer 2.0. Before beginning a test, make sure the following are available:

- QIAstat-Dx Analyzer 2.0 (at least one Operational Module PRO and one Analytical Module) with software version 1.6 or later
- *QIAstat-Dx Analyzer 2.0 User Manual* (for use with software version 1.6 or later)
- QIAstat-Dx latest Assay Definition File software for QIAstat-Dx Gastrointestinal Panel 2 installed in the Operational Module PRO

**Important:** Prior to use, ensure that instruments have been checked and calibrated according to the manufacturer's recommendations.

# Warnings and Precautions

- The QIAstat-Dx Gastrointestinal Panel 2 is for in vitro diagnostic use.
- For prescription use only.
- The QIAstat-Dx Gastrointestinal Panel 2 is to be used by laboratory professionals trained in the use of QIAstat-Dx Analyzer 2.0.
- False positives and false negatives can be the result of a variety of sources and causes. A trained healthcare professional should carefully interpret the results from the QIAstat-Dx Gastrointestinal Panel 2 in conjunction with a patient's signs and symptoms, results from other diagnostic tests, and relevant epidemiological information.
- Please be aware that you may be required to consult your local regulations for reporting serious incidents that have occurred in relation to the device to the manufacturer and the regulatory authority in which the user and/or the patient is established.

## Safety information

- When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient PDF at [www.qiagen.com/safety](http://www.qiagen.com/safety) where you can find, view, and print the SDS for each QIAGEN kit and kit component.
- Observe standard laboratory procedures for keeping the working area clean and contamination-free. Guidelines are outlined in publications such as the *Biosafety in Microbiological and Biomedical Laboratories* from the Centers for Disease Control and Prevention and the National Institutes of Health (72).
- Specimens and samples are potentially infectious. Discard sample and assay waste according to your local safety procedures.

- Always wear appropriate personal protective equipment and follow your institution's safety procedures for handling biological samples. Handle all samples, cartridges, and transfer pipettes as if they are capable of transmitting infectious agents.
- Always observe safety precautions as outlined in relevant guidelines, such as the Clinical and Laboratory Standards Institute® (CLSI) *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline (M29)*, or other appropriate documents provided by local authorities.
- The QIAstat-Dx Gastrointestinal Panel 2 Cartridge is a closed, single-use device that contains all reagents needed for sample preparation and multiplex real-time RT-PCR within the QIAstat-Dx Analyzer 2.0. Do not use a QIAstat-Dx Gastrointestinal Panel 2 Cartridge that is past its expiration date, appears damaged, or leaks fluid.
- Dispose of used or damaged cartridges in accordance with all national, state, and local health and safety regulations and laws.

## Emergency information

CHEMTREC

USA & Canada 1-800-424-9300

## Precautions

The following hazard and precautionary statements apply to components of the QIAstat-Dx Gastrointestinal Panel 2.



Contains: ethanol; guanidine hydrochloride; guanidine thiocyanate; iso-propanol; proteinase K; t-Octylphenoxypolyethoxyethanol. Danger! Highly flammable liquid and vapor. Harmful if swallowed or if inhaled. May be harmful in contact with skin. Causes severe skin burns and eye damage. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause drowsiness or dizziness. Harmful to aquatic life with long lasting effects. Contact with acids liberates very toxic gas. Corrosive to the respiratory tract. Keep away from heat/sparks/open flames/hot surfaces. No smoking. Avoid breathing dust/fume/gas/mist/vapors/spray. Wear protective gloves/protective clothing/eye protection/face protection. Wear respiratory protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF exposed or concerned: Immediately call a POISON CENTER or doctor. Rinse mouth. Do NOT induce vomiting. Remove person to fresh air and keep comfortable for breathing. Wash contaminated clothing before reuse. Store in a well-ventilated place. Keep container tightly closed.

To reduce the risk of contamination when handling stool samples, it is recommended that the below guidelines are applied (72).

- When handling the stool sample, a biosafety cabinet, dead air box, splash shield, or face shield should be used.
- The work area used for cartridge loading should be separate from the work area used for stool pathogen testing (i.e., stool culture, EIA) to prevent cross-contamination.
- Prior to sample handling, the work area should be thoroughly cleaned using 10% bleach or similar disinfectant.
- QIAstat-Dx Gastrointestinal Panel 2 Cartridges and samples should be processed one at a time.
- Change gloves prior to removing cartridges from shipping boxes.
- Change gloves and clean the work area between processing each sample.

- Dispose of used cartridges in a biohazard container immediately after the run is complete and avoid excessive handling.

## Precautions related to public health reporting

National, state, and local public health organizations have published guidelines for the notification of reportable diseases. While the list of reportable conditions varies by state, the Council of State and Territorial Epidemiologists (CSTE) has recommended that state health departments report cases of selected diseases to CDC's National Notifiable Diseases Surveillance System (NNDS). At the time of writing, the notifiable pathogens in the US per CDC included in the QIAstat-Dx Gastrointestinal Panel 2 are:

- *Campylobacter* spp.
- Certain *E. coli*
  - O157:H7
  - Shiga toxin-producing (STEC)
- *Cryptosporidium* spp.
- *Cyclospora cayetanensis*
- *Giardia lamblia*
- *Salmonella* spp.
- *Salmonella* enterica serotypes Paratyphi A, B [tartrate negative], and C [S. Paratyphi]
- *Salmonella* enterica serotype Typhi
- *Shigella* spp./EIEC

Pathogens are notifiable due to their outbreak potential or impact on public health. Laboratories are responsible for following their state or local regulations for submission of clinical material or isolates in positive specimens to their state public health laboratories.

# Reagent Storage and Handling

Store the QIAstat-Dx Gastrointestinal Panel 2 Cartridges in a dry, clean storage space at room temperature (15–25°C). Do not remove the QIAstat-Dx Gastrointestinal Panel 2 Cartridges or the transfer pipettes from their individual packaging until actual use. Under these conditions, QIAstat-Dx Gastrointestinal Panel 2 Cartridges can be stored until the expiration date printed on the individual packaging. The expiration date is also included in the QIAstat-Dx Gastrointestinal Panel 2 Cartridge barcode and is read by the QIAstat-Dx Analyzer 2.0 when the cartridge is inserted into the instrument to run a test. Once the cartridge is removed from the pouch, it should be protected from sunlight.

Attention should be paid to the expiration dates and storage conditions printed on the box and labels of all components. Do not use expired or incorrectly stored components.

## In-use stability

When stored under the specified storage conditions, the QIAstat-Dx Gastrointestinal Panel 2 is stable until the stated expiration date on box label.

After the cartridge package is opened, sample should be introduced into the QIAstat-Dx Gastrointestinal Panel 2 Cartridge within 30 minutes. Sample-loaded cartridges should be loaded into the QIAstat-Dx Analyzer 2.0 within 90 minutes.



# Specimen Storage and Handling

The QIAstat-Dx Gastrointestinal Panel 2 kit is for use with stool samples resuspended in transport medium (Para-Pak C&S (Meridian Bioscience) or FecalSwab (COPAN)). All samples should be treated as potentially infectious. Discard sample and assay waste according to your local safety procedures.

## Specimen collection

Stool samples should be collected and handled according to the transport medium manufacturer's recommended procedures.

Recommended storage conditions for stool resuspended in transport medium (Para-Pak C&S (Meridian Bioscience) or FecalSwab (COPAN)) specimens are listed below:

- Room temperature up to 4 days at 15–25°C
- Refrigerated up to 4 days at 2–8°C

# Procedure

## Important points before starting

- Ensure all materials required but not provided are available.
- The QIAstat-Dx Gastrointestinal Panel 2 Cartridge (cat. no 691421) is identified by a purple-colored (●) bar on the label and an icon indicating gastrointestinal tract (🦠, see "Symbols" on page 127).

## Protocol: Stool samples in transport medium


### Sample collection, transport, and storage

Collect and resuspend the stool sample in Para-Pak C&S (Meridian) or FecalSwab (COPAN) transport media according to the manufacturer's recommended procedures.

### Loading a sample into the QIAstat-Dx Gastrointestinal Panel 2 Cartridge


1. Open the package of a QIAstat-Dx Gastrointestinal Panel 2 Cartridge using the tear notches on the sides of the packaging (Figure 2).

**Important:** After the package is opened, sample should be introduced into the QIAstat-Dx Gastrointestinal Panel 2 Cartridge within 30 minutes. Sample-loaded cartridges should be loaded into the QIAstat-Dx Analyzer 2.0 within 90 minutes.



**Figure 2. Opening the QIAstat-Dx Gastrointestinal Panel 2 Cartridge.**

2. Remove the QIAstat-Dx Gastrointestinal Panel 2 Cartridge from the packaging and position it so that the barcode on the label faces you.
3. Manually write the sample information or place a sample information label on the top of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge. Make sure that the label is properly positioned and does not block the lid opening (Figure 3).




**Figure 3. Sample information placement on top of QIAstat-Dx Gastrointestinal Panel 2 Cartridge.**

4. Place the QIAstat-Dx Gastrointestinal Panel 2 Cartridge flat on the clean work surface so that the barcode on the label faces upwards. Open the sample lid of the main port on the front of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge (Figure 4).


**Important:** Do not flip the QIAstat-Dx Gastrointestinal Panel 2 Cartridge or agitate it while the main port lid is open. The main port contains silica beads used in the sample disruption. The silica beads could fall out of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge if it is agitated while the lid is open.

**Note:** The swab port is not used for the QIAstat-Dx Gastrointestinal Panel 2 assay.



**Figure 4. Opening the sample lid of main port.**

5. Thoroughly mix the stool in the Para-Pak C&S or FecalSwab transport medium, for example, by vigorously agitating the tube 3 times (Figure 5).




**Figure 5. Mixing stool sample in transport medium.**

6. Open the tube with the sample to be tested. Use the supplied transfer pipette to draw up fluid. Draw the sample to the second fill line on the pipette (i.e., 200  $\mu$ L) (Figure 6).

**Important:** Do not draw air, mucus, or particles into the pipette. If air, mucus, or particles are drawn into the pipette, carefully expel the sample fluid in the pipette back into the sample tube and draw up fluid again. In the event that the supplied transfer pipette is lost,


please use another one from the package or any other commercially available pipette with a minimum volume of 200  $\mu\text{L}$ .

**Note:** In case the test should be repeated due to previous cartridge error related to sample concentration too high, draw the sample to the first fill line on the pipette instead (100  $\mu\text{L}$ ) (See the Troubleshooting Guide section for further details on error codes and "Appendix C: Additional instructions for use" on page 136 for further instructions on repeating a sample with 100  $\mu\text{L}$ ).




**Figure 6. Drawing up sample into the supplied transfer pipette.**

7. Carefully transfer the sample into the main port of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge using the supplied single-use transfer pipette (Figure 7).



**Figure 7. Transferring sample to main port of QIAstat-Dx Gastrointestinal Panel 2 Cartridge.**


8. Firmly close the lid of the main port until it clicks (Figure 8).



**Figure 8. Closing the lid of the main port.**

9. Visually confirm that the sample has been loaded by checking the sample inspection window of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge (Figure 9). A mixture of sample and silica beads should be observed.

**Important:** After the sample is placed inside the QIAstat-Dx Gastrointestinal Panel 2 Cartridge, the cartridge must be loaded into the QIAstat-Dx Analyzer 2.0 within 90 minutes.



**Figure 9. Sample inspection window (red arrow).**

## Running a test with a QIAstat-Dx Analyzer 2.0

1. Power on the QIAstat-Dx Analyzer 2.0 using the **ON/OFF** button on the front of the instrument.

**Note:** The power switch at the back of the Analytical Module must be set in the "I" position. The QIAstat-Dx Analyzer 2.0 status indicators will turn blue.

2. Wait until the Main screen appears and the QIAstat-Dx Analyzer 2.0 status indicators turn green and stop blinking.
3. Enter your username and password for QIAstat-Dx Analyzer 2.0 to log in.


**Note:** The Login screen will appear if User Access Control is activated. If the User Access Control is disabled, username/password will not be required and the Main screen will appear.

4. If the Assay Definition File software is not installed on the QIAstat-Dx Analyzer 2.0, follow the installation instructions prior to running the test (see "Appendix A: Installing the Assay Definition File" on page 130 for additional information).
5. Press **Run Test** in the top right corner of the touchscreen of the QIAstat-Dx Analyzer 2.0.
6. When prompted, scan the sample ID barcode on the resuspended sample or scan the specimen information barcode located on the top of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge (Figure 3) using the integrated front barcode reader of QIAstat-Dx Analyzer 2.0 (Figure 10).

**Note:** It is also possible to enter the sample ID using the virtual keyboard of the touchscreen by selecting the **Sample ID** field.

**Note:** Depending on the selected system configuration, entering the patient ID may also be required at this point.


**Note:** Instructions from the QIAstat-Dx Analyzer 2.0 appear in the Instructions Bar at the bottom of the touchscreen.



**Figure 10. Scanning sample ID barcode.**

- When prompted, scan the barcode of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge to be used (Figure 11). The QIAstat-Dx Analyzer 2.0 will automatically recognize the assay to be run based on the cartridge barcode.


**Note:** The QIAstat-Dx Analyzer 2.0 will not accept QIAstat-Dx Gastrointestinal Panel 2 Cartridges with lapsed expiration dates, previously used cartridges, or cartridges for assays that have not been installed on the unit. An error message will be shown in these cases and the QIAstat-Dx Gastrointestinal Panel 2 Cartridge will be rejected. Refer to the *QIAstat-Dx Analyzer 2.0 User Manual* or "Appendix A: Installing the Assay Definition File" on page 130 for further details on how to install assays.



**Figure 11. Scanning QIAstat-Dx Gastrointestinal Panel 2 Cartridge barcode.**



8. The Confirm screen will appear. Review the entered data and make any necessary changes by selecting the relevant fields on the touchscreen and editing the information.
9. Press **Confirm** when all the displayed data are correct. If needed, select the appropriate field to edit its content, or press **Cancel** to cancel the test (Figure 12).



**Figure 12. Confirming data entry.**

10. Ensure that both sample lids of the swab port and main port of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge are firmly closed.
11. When the cartridge entrance port on the top of the QIAstat-Dx Analyzer 2.0 automatically opens, insert the QIAstat-Dx Gastrointestinal Panel 2 Cartridge with the barcode facing to the left and the reaction chambers facing down (Figure 13).

**Note:** Depending on the system configuration, the operator may be required to re-enter their user password to start the test run.


**Note:** Up to this point, it is possible to cancel the test run by pressing the **Cancel** button at the bottom right corner of the touchscreen.

12. Upon detecting the QIAstat-Dx Gastrointestinal Panel 2 Cartridge, the QIAstat-Dx Analyzer 2.0 will automatically close the lid of the cartridge entrance port and start the test run. No further action from the operator is required to start the run.

**Note:** There is no need to push the QIAstat-Dx Gastrointestinal Panel 2 Cartridge into the QIAstat-Dx Analyzer 2.0.

**Note:** The QIAstat-Dx Analyzer 2.0 will not accept a QIAstat-Dx Gastrointestinal Panel 2 Cartridge other than the one used and scanned during the test setup. If a cartridge other than the one scanned is inserted, an error will be generated, and the cartridge will be automatically ejected.

**Note:** The lid of the cartridge entrance port will close automatically after 30 seconds if a QIAstat-Dx Gastrointestinal Panel 2 Cartridge is not positioned in the port. If this occurs, repeat the procedure starting with step 5.




**Figure 13. Inserting QIAstat-Dx Gastrointestinal Panel 2 Cartridge into the QIAstat-Dx Analyzer 2.0.**


13. While the test is running, the remaining run time is displayed on the touchscreen.
14. After the test run is completed, the Eject screen will appear (Figure 14) and the Module status bar will display the test result as one of the following options:
  - TEST COMPLETED: The test was completed successfully
  - TEST FAILED: An error occurred during the test

- **TEST CANCELED:** The user canceled the test

**Important:** If the test fails, refer to the “Troubleshooting” section in the *QIAstat-Dx Analyzer 2.0 User Manual* for possible reasons and instructions on how to proceed. For additional information about specific QIAstat-Dx Gastrointestinal Panel 2 error codes and messages, please see the Troubleshooting Guide section of this document.



**Figure 14. Eject screen display.**

15. Press  **Eject** on the touchscreen to remove the QIAstat-Dx Gastrointestinal Panel 2 Cartridge and dispose of it as biohazardous waste in accordance with all national, state, and local health and safety regulations and laws. The QIAstat-Dx Gastrointestinal Panel 2 Cartridge should be removed when the cartridge entrance port opens and ejects the cartridge. If the cartridge is not removed after 30 seconds, it will automatically move back into the QIAstat-Dx Analyzer 2.0 and the cartridge entrance port lid will close. If this occurs, press **Eject** to open the lid of the cartridge entrance port again and then remove

the cartridge.

**Important:** Used QIAstat-Dx Gastrointestinal Panel 2 Cartridges must be discarded. It is not possible to re-use cartridges for tests for which the execution was started but then subsequently canceled by the operator, or for which an error was detected.

16. After the QIAstat-Dx Gastrointestinal Panel 2 Cartridge has been ejected, the results Summary screen will appear. Refer to "Interpretation of Results" on the facing page for further details. To begin the process for running another test, press **Run Test**.

**Note:** For further information on the use of the QIAstat-Dx Analyzer 2.0, refer to the *QIAstat-Dx Analyzer 2.0 User Manual*.

# Interpretation of Results

## Viewing results

The QIAstat-Dx Analyzer 2.0 automatically interprets and saves test results. After ejecting the QIAstat-Dx Gastrointestinal Panel 2 Cartridge, the results Summary screen is automatically displayed (Figure 15).







Figure 15. QIAstat-Dx Analyzer 2.0 screen.

Figure 15 shows the screen for the QIAstat-Dx Analyzer 2.0.

QIAstat-Dx Analyzer 2.0 includes an additional tab, AMR Genes, which is disabled for QIAstat-Dx Gastrointestinal Panel 2.

The main part of the screen provides the following lists and uses color-coding and symbols to indicate the results:

- The first list, under the heading “Detected”, includes all pathogens detected and identified in the sample, which are preceded by a  sign and are colored red.
- The second list, under the heading “Tested”, includes all pathogens tested in the sample. Pathogens detected and identified in the sample are preceded by a  sign and are colored red. Pathogens that were tested but not detected are preceded by a  sign and are colored green. Invalid and not applicable pathogens are also displayed in this list.

**Note:** Pathogens detected and identified in the sample are shown in both the **Detected** and **Tested** lists.

If the test failed to complete successfully, a message will indicate **Failed** followed by the specific error code.

The following Test Data are shown on the left side of the screen:

- Sample ID
- Patient ID (if available)
- Assay Type
- Sample Type

Further data about the assay is available, depending on the operator’s access rights, through the tabs at the bottom of the screen (e.g., amplification plots and test details).

A report with the assay data can be exported to an external USB storage device. Insert the USB storage device into one of the USB ports of the QIAstat-Dx Analyzer 2.0 and press **Save Report** in the bottom bar of the screen. This report can be exported later at any time by selecting the test from the “View Result” List.

The report can also be sent to the printer by pressing **Print Report** in the bottom bar of the screen.

## Viewing amplification curves

To view test amplification curves of pathogens detected, press the **Amplification Curves** tab (Figure 16).



Figure 16. Amplification Curves screen (PATHOGENS tab).

Details about the tested pathogens and controls are shown on the left and the amplification curves are shown in the center.

**Note:** If User Access Control is enabled on the QIAstat-Dx Analyzer 2.0, the **Amplification Curves** screen is only available for operators with access rights.

Press the **PATHOGENS** tab on the left side to display the plots corresponding to the tested pathogens. Press on the pathogen name to select which pathogens are shown in the amplification plot. It is possible to select single, multiple, or no pathogens. Each pathogen in the selected list will be assigned a color corresponding to the amplification curve associated with the pathogen. Unselected pathogens will be shown in gray.

The corresponding Ct and endpoint fluorescence (EP) values are shown below each pathogen name.

Press the **CONTROLS** tab on the left side to view the controls in the amplification plot. Press the circle next to the control name to select or deselect it (Figure 17).






Figure 17. Amplification Curves screen (CONTROLS tab).



The amplification plot displays the data curve for the selected pathogens or controls. To alternate between logarithmic or linear scale for the Y-axis, press the **Lin** or **Log** button at the bottom left corner of the plot.

The scale of the X-axis and Y-axis can be adjusted using the  **blue pickers** on each axis. Press and hold a blue picker and then move it to the desired location on the axis. Move a blue picker to the axis origin to return to the default values.

## Viewing test details

Press  **Test Details** in the Tab Menu bar at the bottom of the touchscreen to review the results in more detail. Scroll down to see the complete report.

The following Test Details are shown in the center of the screen (Figure 18):

- User ID
- Cartridge SN (serial number)
- Cartridge Expiration Date
- Module SN (serial number)
- Test Status (Completed, Failed, or Canceled by operator)
- Error Code (if applicable)
- Test Start Date and Time
- Test Execution Time
- Assay Name
- Test ID

- Test Result
  - **Positive** (if at least one gastrointestinal pathogen is detected/identified)
  - **Positive with warning** (if at least one pathogen is detected, but the Internal Control failed)
  - **Negative** (if no gastrointestinal pathogen is detected)
  - **Failed** (an error occurred or the test was canceled by the user)
- List of analytes tested in the assay, with Ct and endpoint fluorescence in the event of a positive signal
- Internal Control, with Ct and endpoint fluorescence




Figure 18. Example screen showing Test Data on the left panel and Test Details in the main panel.

# Browsing results from previous tests

To view results from previous tests that are stored in the results repository, press **View Results** on the Main Menu bar (Figure 19).




Figure 19. Example View Results screen.


The following information is available for every executed test (Figure 20):

- Sample ID
- Assay (name of test assay which is “GI2” for Gastrointestinal Panel 2)
- Operator ID
- Mod (Analytical Module on which the test was executed)
- Date/Time (date and time when the test was finished)

- Result (outcome of the test: positive [pos], positive with warning [pos\*], negative [neg], failed [fail], or successful [suc])

**Note:** If User Access Control is enabled on the QIAstat-Dx Analyzer 2.0, the data for which the user has no access rights will be hidden with asterisks.

Select one or more test results by pressing the **gray circle** to the left of the sample ID. A checkmark will appear next to selected results. Unselect test results by pressing this **checkmark**. The entire list of results can be selected by pressing the **checkmark circle** in the top row (Figure 20).








**Figure 20. Example of selecting Test Results in the View Results screen.**

Press anywhere in the test row to view the result for a particular test.

Press a column headline (e.g., Sample ID) to sort the list in ascending or descending order according to that parameter. The list can be sorted according to only one column at a time.

The Result column shows the outcome of each test (Table 2).

Table 2. Descriptions of the test results displayed in View Results screen

Outcome	Result	Description	Action
Positive	 pos	At least one pathogen is positive	Refer to the Summary Result Screen or Result Printout for pathogen specific results. Description of pathogen results can be found in Table 5.
Positive with warning	 pos*	At least one pathogen is positive, but the Internal Control failed	Refer to the Summary Result Screen or Result Printout for pathogen specific results. Description of pathogen results can be found in Table 5.
Negative	 neg	No pathogens were detected	Refer to the Summary Result Screen or Result Printout for pathogen specific results. Description of pathogen results can be found in Table 5.
Failed	 fail	The test failed because either an error occurred, the test was canceled by the user, or no pathogens were detected and the internal control failed.	Repeat the test using a new cartridge. Accept the results of the repeat testing. If the error persists, contact QIAGEN Technical Services for further instructions.
Successful	 Suc	The test is either positive or negative, but the user does not have the access rights to view the test results.	Login from a user profile with rights to view the results.

Press **Save Report** to save the report(s) for the selected result(s) in PDF to an external USB storage device.

Select the report type: **List of Tests** or **Test Reports**.

Press **Search** to search the test results by Sample ID, Assay, and Operator ID. Enter the search string using the virtual keyboard and press **Enter** to start the search. Only the records containing the searched text will be displayed in the search results.

If the results list has been filtered, the search will only apply to the filtered list.

Press and hold a column headline to apply a filter based on that parameter. For some parameters, such as Sample ID, the virtual keyboard will appear so the search string for the filter can be entered.

For other parameters, such as Assay, a dialog box will open with a list of assays stored in the repository. Select one or more assays to filter only the tests that were performed with the selected assays.

The  symbol to the left of a column headline indicates that the column's filter is active.

A filter can be removed by pressing **Remove Filter** in the Submenu bar.

## Exporting results to a USB drive


From any tab of the View Results screen, select **Save Report** to export and save a copy of the test results in PDF to a USB drive. The USB port is located on the front of the QIAstat-Dx Analyzer 2.0.

## Backup and regular data upload to SharePoint

The results can be exported from the instrument following these steps:

1. Press **Options > System Configuration > System Backup** (Figure 21). Insert a USB storage device into the front USB port.

2. Press **Perform Backup**. A file with the extension **.dbk** will be generated in the USB with a default file name.



**Figure 21. Perform a backup.**

Printing results

Make sure a printer is connected to the QIAstat-Dx Analyzer 2.0 and the proper driver is installed. Press **Print Report** to send a copy of the PDF test results to the printer.

Sample result interpretation

For samples in Para-Pak C&S and FecalSwab collection devices, a result for a gastrointestinal organism is interpreted as “Positive” when the corresponding PCR assay is positive except for EPEC, STEC, and *E. coli* O157. The result interpretation for EPEC, STEC and *E. coli* O157 follows the rationale explained in Table 3.

Table 3. Interpretation of EPEC, STEC, and E. coli O157 results only applicable for Para-Pak C&S samples

EPEC result	STEC <i>stx1/stx2</i> result	E. coli O157 result	Description
Negative	Negative	N/A	Enteropathogenic <i>E. coli</i> (EPEC) was not detected and Shiga-like toxin-producing <i>E. coli</i> (STEC) <i>stx1/stx2</i> is negative as both <i>stx1</i> and <i>stx2</i> have not been detected.  <i>E. coli</i> O157 result is not applicable (N/A) when STEC <i>stx1/stx2</i> is not detected due to <i>E. coli</i> O157 being a specific serotype of STEC.
Positive	Negative	N/A	EPEC was detected and STEC <i>stx1/stx2</i> is negative as both <i>stx1</i> and <i>stx2</i> have not been detected.  <i>E. coli</i> O157 result is not applicable (N/A) when STEC <i>stx1/stx2</i> is not detected due to <i>E. coli</i> O157 being a specific serotype of STEC.
N/A	Positive	Negative	EPEC result is not applicable because EPEC detection cannot be differentiated when STEC <i>stx1</i> or <i>stx2</i> is detected.  <i>E. coli</i> O157 was not detected.
N/A	Positive	Positive	EPEC result is not applicable because EPEC detection cannot be differentiated when STEC <i>stx1</i> or <i>stx2</i> is detected.  <i>E. coli</i> O157 was detected.

For both sample collection devices, internal control results are to be interpreted according to Table 4.

Table 4. Interpretation of Internal Control results

Control result	Explanation	Action
Passed	The Internal Control amplified successfully.	The run was completed with success. All results are validated and can be reported. Detected pathogens are reported as “positive” and undetected pathogens are reported as “negative”.



Table 4. Interpretation of Internal Control results (continued)

Control result	Explanation	Action
Failed	The Internal Control failed.	Positively detected pathogen(s) are reported, but all negative results (tested but not detected pathogen[s]) are invalid. Repeat the testing using a new Cartridge. Accept the results of the repeat testing. If the invalid result persists, contact QIAGEN Technical Services for further instruction.

The software provides an overall test result (Table 2) as well as a result for individual pathogens. Possible results for each organism include Detected/Positive, Not Detected/Negative, N/A, and Invalid (Table 5). If the internal control has failed and no positive signal was detected or if there is an instrument error, there will be no pathogen results provided.

Table 5. Description of pathogen results as displayed on Summary Result screen and the Result Printout






Result	Symbol	Explanation	Action
Positive/ Detected		A positive signal was detected for this pathogen. Result of the Internal Control is passed.	None. Report results.
Positive/ Detected with Warning	 !pos*	A positive signal was detected for this pathogen, but the result of the Internal Control has failed.	Report positive analyte. Repeat the test using a new cartridge. Accept the results of the repeat testing. If the invalid result persists, contact QIAGEN Technical Services for further instructions.
Negative/ Not Detected		No signal was detected for this pathogen. The Internal Control passed.	None. Report results.

Table 5. Description of pathogen results as displayed on Summary Result screen and the Result Printout (continued)

Result	Symbol	Explanation	Action
N/A (applies to <i>E. coli</i> O157 and EPEC only)		The run was successfully completed and the Internal Control passed.  For <i>E. coli</i> O157 N/A: Shiga-like toxin-producing <i>E. coli</i> (STEC) was not detected. For EPEC N/A: Shiga-like toxin producing <i>E. coli</i> (STEC) was detected.	None. Report results.
Invalid		No signal was detected for this pathogen and the Internal Control failed (but other pathogens have been detected).	Repeat the test using a new cartridge. Accept the results of the repeat testing. If the invalid result persists, contact QIAGEN Technical Services for further instructions.

# Quality Control

## Internal control interpretation

The QIAstat-Dx Gastrointestinal Panel 2 Cartridge includes a full process Internal Control, which is titrated *Schizosaccharomyces pombe*. *Schizosaccharomyces pombe* is a yeast (fungi) that is included in the cartridge in dried form and is rehydrated upon sample loading. This Internal Control material verifies all steps of the analysis process, including sample homogenization, lysis of viral and cellular structures (by means of chemical and mechanical disruption), nucleic acid purification, reverse transcription, and real-time PCR.

A passed result for the Internal Control indicates that all processing steps performed by the QIAstat-Dx Gastrointestinal Panel 2 Cartridge were successful.

A failed result of the Internal Control does not negate any positive results for detected and identified targets, but it does invalidate all negative results in the analysis. Therefore, the test should be repeated if the Internal Control signal is negative.

## External control information

All external quality control requirements and testing should be performed in accordance with local, state, and federal regulations or accreditation organizations and should follow the user's laboratory standard quality control procedures.

Blank controls are not applicable to the device because it is a single test disposable cartridge. Regular testing of negative and positive external controls is recommended by the company but controls are not provided with the QIAstat-Dx Gastrointestinal Panel 2. Use transport media as the external Negative Control and previously characterized positive samples or negative sample spiked with well characterized target organisms as external positive controls.

## Limitations

- The QIAstat-Dx Gastrointestinal Panel 2 is intended for professional use only and is not intended for self-testing. The QIAstat-Dx Gastrointestinal Panel 2 is intended for in vitro diagnostic use.
- Results from the QIAstat-Dx Gastrointestinal Panel 2 are not intended to be used as the sole basis for diagnosis, treatment, or other patient management decisions.
- All assay results should be used and interpreted by a trained healthcare professional in the context of a full clinical evaluation, laboratory, and epidemiological findings, as an aid in the diagnosis of gastrointestinal infection.
- The performance of this test has not been determined for patients without signs and symptoms of gastrointestinal illness.
- The QIAstat-Dx Gastrointestinal Panel 2 is not intended for testing of samples other than those described in this Instructions for Use (IFU). The performance of this test has only been validated with human stool preserved in transport medium (Para-Pak C&S or FecalSwab), according to the media manufacturers' instructions. It has not been validated for use with other stool transport media, rectal swabs, raw stool, vomitus, or endoscopy stool aspirates. The QIAstat-Dx Gastrointestinal Panel 2 should not be used to test samples within Para-Pak C&S or FecalSwab collection devices that have been overfilled with stool. Only stool resuspended following the collection device manufacturer's instructions should be used. The overfilling of Para-Pak C&S or FecalSwab collection devices can result in a failed test with an error indicating "Sample concentration too high".
- The detection of viral, bacterial, or parasitic sequences is dependent upon proper specimen collection, handling, transportation, storage, and preparation (including extraction). Failure to observe proper procedures in any one of these steps can lead to incorrect results. There is a risk of false negative values resulting from improperly collected, transported, or handled specimens.

- Positive results do not rule out co-infection with organisms not included in the QIAstat-Dx Gastrointestinal Panel 2. The agent detected may not be the definitive cause of the disease.
- Not all agents of acute gastrointestinal infection are detected by this assay.
- The QIAstat-Dx Gastrointestinal Panel 2 is intended to be used in conjunction with standard of care culture for organism recovery, serotyping, and/or antimicrobial susceptibility testing where applicable.
- The QIAstat-Dx Gastrointestinal Panel 2 can be used only with the QIAstat-Dx Analyzer 2.0.
- The identification of multiple diarrheagenic *E. coli* pathotypes has historically relied upon phenotypic characteristics, such as adherence patterns or toxigenicity in certain tissue culture cell lines (27). The QIAstat-Dx Gastrointestinal Panel 2 targets genetic determinants characteristic of most pathogenic strains of these organisms but may not detect all strains having phenotypic characteristics of a pathotype.
- Genetic virulence markers associated with diarrheagenic *E. coli* / *Shigella* pathotypes are often carried on mobile genetic elements (MGEs) that can be transferred horizontally between different strains (27); therefore, “Detected” results for multiple diarrheagenic *E. coli* / *Shigella* may be due to co-infection with multiple pathotypes or, less frequently, may be due to the presence of a single organism containing genes characteristic of multiple pathotypes. An example of the latter is the 2019 *E. coli* hybrid ETEC/STEC strains found in Sweden (73).
- The QIAstat-Dx Gastrointestinal Panel 2 detects Enteropathogenic *E. coli* (EPEC) through targeting of the *eae* gene, which encodes the adhesin intimin. As some Shiga-like toxin-producing *E. coli* (STEC) also carry *eae* (in particular, strains identified as enterohemorrhagic *E. coli*; EHEC) (27), the QIAstat-Dx Gastrointestinal Panel 2 cannot distinguish between STEC containing *eae* and a co-infection of EPEC and STEC. Therefore, the EPEC result is not applicable (N/A) and not reported for specimens in which STEC has

also been detected. In rare cases, STEC may be reported as EPEC when a STEC carrying *eae* (EHEC) is present in a specimen below the LoD of the STEC oligonucleotide design(s). Rare instances of other organisms carrying *eae* have been documented (e.g., *Escherichia albertii*, and *Shigella boydii*)(74).

- *Shigella dysenteriae* serotype 1 possess a shiga toxin gene (*stx*) that is identical to the *stx1* gene of STEC (27). *Stx* genes have been more recently found in other *Shigella* species (e.g., *S. sonnei* and *S. flexneri*) (75,76). The detection of both *Shigella*/Enteroinvasive *E. coli* (EIEC) and STEC *stx1/stx2* analytes in the same specimen may indicate the presence of *Shigella* species such as *S. dysenteriae*. Rare instances of the detection of Shiga-like toxin genes in other genera/species have been reported (e.g., *Acinetobacter haemolyticus*, *Enterobacter cloacae*, and *Citrobacter freundii* (77,78,79).
- *E. coli* O157 result is only reported as specific serogroup identification in association with STEC *stx1/stx2*. While non-STEC O157 strains have been detected in human stool (80), their role in disease has not been established (81). Serotype O157 EPEC have been identified and will be detected by the QIAstat-Dx Gastrointestinal Panel 2 (by the EPEC oligonucleotides design) due to their carriage of the *eae* gene.
- The QIAstat-Dx Gastrointestinal Panel 2 cannot distinguish between infections with a single toxigenic STEC O157 or rare co-infections of STEC (non-O157) with a *stx1/stx2*-negative *E. coli* O157.
- This test only detects *Campylobacter jejuni*, *C. coli*, and *C. upsaliensis*, and does not differentiate between these three species of *Campylobacter*. Additional testing is required to differentiate between these species and to detect other *Campylobacter* species that may be present in stool specimens. In particular, the *Campylobacter upsaliensis* oligonucleotides design may cross-react with the *Campylobacter* species, *C. lari* and *C. helveticus* organisms.
- Negative results do not exclude the possibility of gastrointestinal infection. Negative test results may occur from sequence variants in the region targeted by the assay, the presence

of inhibitors, technical errors, sample mix-ups, or an infection caused by an organism not detected by the panel. Test results may also be affected by use of certain medications (e.g., calcium carbonate), concurrent antimicrobial therapy, or levels of organism in the sample that are below the limit of detection for the test. Sensitivity in some clinical settings may differ from that described in the Instructions for Use. Negative results should not be used as the sole basis for diagnosis, treatment, or other management decisions.

- Organism and amplicon contamination may produce erroneous results for this test. Particular attention should be given to the Laboratory Precautions noted under the Laboratory Precautions section.
- There is a risk of false-positive values resulting from cross-contamination by target organisms, their nucleic acids or the amplified product, or from non-specific signals in the assay.
- There is a risk of false negative results due to the presence of strains with sequence variability in the target regions of the oligonucleotides design. Refer to the Inclusivity testing section of this document for additional information.
- Not all *Salmonella* serotypes were tested in validation studies; however, representatives of the 20 most prevalent serotypes recently circulating in the US (CDC National *Salmonella* Surveillance Annual Summary 2016) were evaluated during analytical reactivity studies. In silico sequence analysis supports detection of all subspecies and serotypes of *Salmonella*.
- The performance of the QIAstat-Dx Gastrointestinal Panel 2 has not been established in individuals who received Rotavirus A vaccine. Recent oral administration of a Rotavirus A vaccine may cause positive results for Rotavirus A if the virus is passed in the stool.
- The performance of this test has not been evaluated for immunocompromised individuals.
- Underlying polymorphisms in primer-binding regions can affect the targets being detected and subsequently the test results returned.

- Positive and negative predictive values are highly dependent on prevalence. False negative test results are more likely when prevalence of disease is high. False positive test results are more likely when prevalence is low.
- The effect of interfering substances has only been evaluated for those listed in the labeling at its indicated amount or concentration. Interference by substances other than those described in the "Interfering Substances" section of the Instructions for Use can lead to erroneous results.
- Cross-reactivity with gastrointestinal tract organisms other than those listed in the "Analytical Specificity" section of the Instructions for Use may lead to erroneous results.
- This test is a qualitative test and does not provide the quantitative value of detected organism present.
- The assay sensitivity to detect *Cyclospora cayetanensis*, Adenovirus F40/F41, *Entamoeba histolytica* and the Shiga-like toxin- producing *Escherichia coli* (STEC) might be reduced up to 3.16-fold when using half-input sample volume (100 µL) workflow detailed in "Appendix C: Additional instructions for use" on page 136.
- Due to the small number of positive specimens collected for certain analytes during the prospective clinical study, performance characteristics for Adenovirus 40/41, ETEC, *Plesiomonas shigelloides*, *Shigella*/EIEC, STEC, *E. coli* 0157, *Yersinia enterocolitica*, *Cryptosporidium*, and *Giardia lamblia* were established additionally with retrospective clinical specimens. Performance characteristics for Astrovirus and *Entamoeba histolytica* were established primarily with contrived clinical specimens.
- If four or more distinct organisms are detected in a specimen, retesting is recommended to confirm the polymicrobial result.
- Virus, bacteria, and parasite nucleic acid may persist in vivo independently of organism viability. Additionally, some organisms may be carried asymptotically. Detection of organism targets does not imply that the corresponding organisms are infectious or are the



causative agents for clinical symptoms.

- The performance of this test has not been established for monitoring treatment of infection with any of the panel organisms.
- For ETEC, the assay is not predicted to detect bacteria carriers of Heat-labile enterotoxin gene subtype LT-II and/or of Heat-stable enterotoxin gene variant Stb e.
- The assay is not predicted to detect Human Astrovirus types MLB1-3 and VA1-5.
- The potential for competitive inhibition at high concentrations between on-panel analytes was evaluated for a limited number of pathogens (See Table 14). The potential for competitive inhibition between other on-panel analytes is unknown.

# Performance Characteristics

The analytical and clinical performance shown below was demonstrated using QIAstat-Dx Analyzer 1.0. The QIAstat-Dx Gastrointestinal Panel 2 is no longer commercialized for use with the QIAstat-Dx Analyzer 1.0 and it can only be used with the QIAstat-Dx Analyzer 2.0. However, the performance on QIAstat-Dx Analyzer 1.0 is applicable for QIAstat-Dx Gastrointestinal Panel 2 and remains in the Instructions for Use. The QIAstat-Dx Analyzer 2.0 uses the same Analytical Module as QIAstat-Dx Analyzer 1.0; therefore, the performance is not impacted by QIAstat-Dx Analyzer 2.0.

## Analytical performance

### Limit of detection

The Limit of Detection (LoD) is defined as the lowest concentration at which  $\geq 95\%$  of the tested samples generate a positive call.

The LoD for each of the QIAstat-Dx Gastrointestinal Panel 2 target pathogenic organisms was assessed, using in total 36 pathogen strains, by analyzing serial dilutions of analytical samples prepared from culture isolates from commercial suppliers (e.g., ZeptoMetrix® and ATCC®), confirmed clinical isolates, or artificial samples for target analytes commercially unavailable. Each sample tested was prepared in human stool matrix, which consists of a pool of previously tested negative clinical stool specimens resuspended in Para- Pak C&S transport medium.

Each of the 36 strains was tested in human stool matrix prepared following the manufacturer's instructions for the Para-Pak C&S collection device. The confirmed LoD was established by testing 20 replicates at the concentration determined from the preliminary LoD for each strain. The LoD for each strain was confirmed if  $\geq 95\%$  of the replicates were positive. To further

confirm the LoD, at least one dilution below the LoD was tested for each strain and was also tested in 20 replicates and was required to result in less than 95% positivity. A transport media equivalency study between Para-Pak C&S and FecalSwab transport media was conducted to support the conclusions in the section.

Individual LoD values for each QIAstat-Dx Gastrointestinal Panel 2 target are shown in Table 6.

**Table 6. LoD values obtained for the different gastrointestinal target strains tested with the QIAstat-Dx Gastrointestinal Panel 2**

Pathogen	Strain	Source	Concentration (molecular units)* (copies/mL)	Concentration (microbiological units)	Detection rate
<i>Campylobacter</i>	<i>Campylobacter coli</i> 76-GA2 [LMG 21266]	ATCC 43478	5802	1.2 CFU/mL	20/20
	<i>Campylobacter coli</i> CIP 7080	ATCC 33559	8941	0.6 CFU/mL	20/20
	<i>Campylobacter jejuni</i> Z086	ZeptoMetrix 0801650	14,491	1660 CFU/mL	20/20
	<i>Campylobacter jejuni</i> subsp. <i>Jejuni</i> RM3193	ATCC BAA-1234	7210	110 CFU/mL	19/20
	<i>Campylobacter upsaliensis</i> NCTC 11541	ZeptoMetrix 0801999	56,165	2259.4 CFU/mL	20/20
	<i>Campylobacter upsaliensis</i> RM3195	ATCC BAA-1059	7631	35 CFU/mL	19/20
<i>Plesiomonas shigelloides</i>	Z130	ZeptoMetrix 0801899	481	2291 CFU/mL	20/20
	Bader	ATCC 14029	116	2.7 CFU/mL	19/20

**Table 6. LoD values obtained for the different gastrointestinal target strains tested with the QIAstat-Dx Gastrointestinal Panel 2 (continued)**

Pathogen	Strain	Source	Concentration (molecular units)* (copies/mL)	Concentration (microbiological units)	Detection rate
<i>Salmonella</i>	<i>Salmonella enterica</i> Serovar choleraesuis	ATCC 13312	647	91.6 CFU/mL	20/20
	<i>Salmonella enterica</i> Serovar Typhimurium Z005	ZeptoMetrix 0801437	1441	4518.8 CFU/mL	20/20
<i>Yersinia enterocolitica</i>	Z036	ZeptoMetrix 0801734	719	2070 CFU/mL	20/20
	subsp. <i>enterocolitica</i> NTCC 11175, Biotype 4, serotype 3	ATCC 700822	2496	120.1 CFU/mL	20/20
Enteroinvasive <i>E. coli</i> (EIEC)/ <i>Shigella</i>	<i>Shigella sonnei</i> NCDC 1120-66	ATCC 25931	488	0.2 CFU/mL	20/20
	<i>Escherichia coli</i> CDC EDL 1282, O29:NM	ATCC 43892	1431	41.3 CFU/mL	20/20
Enteropathogenic <i>E. coli</i> (EPEC)	<i>Escherichia coli</i> O111:NM (EPEC)	ZeptoMetrix 0801747	1817	2581.7 CFU/mL	20/20
	<i>Escherichia coli</i> 7.1493; EPEC; O84:H28	Zeptomatrix 0801938	29,021	1190 CFU/mL	20/20
Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i>	<i>Escherichia coli</i> H10407, O78:H11	ATCC 35401	367	10.1 CFU/mL	19/20
	<i>Escherichia coli</i> ETEC; ST+, LT+	ZeptoMetrix 0801624	855	567 CFU/mL	20/20
Shiga-like toxin-producing <i>E. coli</i> (STEC) <i>stx1/stx2</i>	<i>Escherichia coli</i> O26:H4	ZeptoMetrix 0801748	2012	726.8 CFU/mL	20/20

**Table 6. LoD values obtained for the different gastrointestinal target strains tested with the QIAstat-Dx Gastrointestinal Panel 2 (continued)**

Pathogen	Strain	Source	Concentration (molecular units)* (copies/mL)	Concentration (microbiological units)	Detection rate
Shiga-like toxin-producing <i>E. coli</i> (STEC) <i>E. coli</i> O157	<i>Escherichia coli</i> O157:H7; EDL933	ZeptoMetrix 0801622	1217	2281.5 CFU/mL	STEC <i>stx</i> 1: 19/20 STEC <i>stx</i> 2: 19/20 O157: 19/20
<i>Cryptosporidium</i>	<i>Cryptosporidium hominis</i>	Public Health Wales UKM 84	357	N/A	20/20
	<i>Cryptosporidium parvum</i> – Iowa isolate	Waterborne® P102C	661	N/A	20/20
<i>Cyclospora cayetanensis</i>	N/A	LACNY-Clinical sample LAC2825	53	N/A	19/20
	N/A	LACNY Clinical sample LAC2827	137	N/A	20/20
<i>Entamoeba histolytica</i>	HM-1:IMSS (Mexico City 1967)	ATCC 30459	7	0.2 cells/mL	20/20
	HK-9 (Korea)	ATCC 30015	1	0.13 cells/mL	19/20
<i>Giardia lamblia</i>	WB (Bethesda)	ATCC 30957	11,850	790 cells/mL	19/20
	Portland-1	ATCC 30888	14,500	635 cells/mL	20/20
Adenovirus F40/F41	Type 40 (Dugan)	ZeptoMetrix 0810084CF	11,726	0.1 TCID <sub>50</sub> /mL	20/20
	Type 41 (Tak)	ZeptoMetrix 0810085CF	979	0.05 TCID <sub>50</sub> /mL	19/20

**Table 6. LoD values obtained for the different gastrointestinal target strains tested with the QIAstat-Dx Gastrointestinal Panel 2 (continued)**

Pathogen	Strain	Source	Concentration (molecular units)* (copies/mL)	Concentration (microbiological units)	Detection rate
Astrovirus	ERE IID 2371 (type 8)	Zeptomatrix 0810277CF	1,158,6371	11.7 TCID <sub>50</sub> /mL	20/20
	ERE IID 2868 (type 4)	Zeptomatrix 0810276CF	52,184	1.3 TCID <sub>50</sub> /mL	19/20
Norovirus GI/GII	GI.1 (recombinant)	ZeptoMatrix 0810086CF	24,629	891.1 TCID <sub>50</sub> /mL	19/20
	GII.4 (recombinant)	ZeptoMatrix 0810087CF	8998	10.5 TCID <sub>50</sub> /mL	20/20
Rotavirus A	69M	ZeptoMatrix 0810280CF	5787	436.1 TCID <sub>50</sub> /mL	19/20
	Wa	ZeptoMatrix 0810041CF	5201	14.1 TCID <sub>50</sub> /mL	19/20

\* Molecular unit titers were determined using in-house developed and validated qPCR assays.

**Exclusivity (analytical specificity)**

The Exclusivity (Analytical Specificity) study was carried out by in vitro testing and *in silico* analysis to assess the potential cross-reactivity and exclusivity of the QIAstat-Dx Gastrointestinal Panel 2. On-panel organisms were tested to assess the potential for intra-panel cross-reactivity and Off-panel organisms were tested to evaluate cross-reactivity with organisms not covered by the panel content. The On-panel and Off-panel organisms tested are shown in Table 7 and Table 8, respectively.

Samples were prepared by single spiking organisms into negative stool resuspended in Para-Pak C&S media at the highest concentration possible based on the organism stock, preferably at 10<sup>5</sup> TCID<sub>50</sub>/mL for viral, 10<sup>5</sup> cells/mL for parasite targets, and 10<sup>6</sup> CFU/mL for bacterial

targets. The pathogens were tested in 3 replicates. There was no intra-panel or Off-panel cross-reactivity for all pathogens tested in vitro, except for two non-targeted *Campylobacter* species (*C. helveticus* and *C. lari*) that cross-reacted with the *Campylobacter* assay oligonucleotides included in the QIAstat-Dx Gastrointestinal Panel 2.

**Table 7. List of analytical specificity On-panel pathogens tested**

Type	Pathogen	
Bacteria	<i>Campylobacter coli</i>	<i>Plesiomonas shigelloides</i>
	<i>Campylobacter jejuni</i>	<i>Salmonella enterica</i>
	<i>Campylobacter upsaliensis</i>	<i>Shigella sonnei</i>
	<i>Escherichia coli</i> (EPEC)*	<i>Yersinia enterocolitica</i>
	<i>Escherichia coli</i> (ETEC)	
	<i>Escherichia coli</i> (STEC)*	
Parasites	<i>Cryptosporidium parvum</i>	<i>Entamoeba histolytica</i>
	<i>Cyclospora cayetanensis</i>	<i>Giardia lamblia</i>
Viruses	Adenovirus F41	Norovirus GII
	Astrovirus	Rotavirus A
	Norovirus GI	

\*Only applicable for samples with ParaPak C&S collection device.

Table 8. List of analytical specificity Off-panel pathogens tested

Type	Pathogen (potential cross-reactant)	
Bacteria	<i>Abiotrophia defectiva</i>	<i>Enterobacter cloacae</i>
	<i>Acinetobacter baumannii</i>	<i>Enterococcus faecalis</i>
	<i>Aeromonas hydrophila</i>	<i>Enterococcus faecium</i>
	<i>Arcobacter cryaerophilus</i>	<i>Escherichia fergusonii</i>
	<i>Bacillus subtilis</i>	<i>Escherichia hermannii</i>
	<i>Bifidobacterium bifidum</i>	<i>Escherichia vulneris</i>
	<i>Campylobacter fetus</i>	<i>Faecalibacterium prausnitzii</i>
	<i>Campylobacter gracilis</i>	<i>Gardnerella vaginalis</i>
	<i>Campylobacter helveticus</i> (Cross-reactive for <i>Campylobacter</i> target)	<i>Haemophilus influenzae</i>
	<i>Campylobacter hominis</i>	<i>Helicobacter pylori</i>
	<i>Campylobacter lari</i> (Cross-reactive for <i>Campylobacter</i> target)	<i>Klebsiella pneumoniae</i>
	<i>Campylobacter mucosalis</i>	<i>Lactobacillus casei</i>
	<i>Campylobacter rectus</i>	<i>Listeria monocytogenes</i>
	<i>Chlamydia trachomatis</i>	<i>Proteus mirabilis</i>
	<i>Citrobacter freundii</i>	<i>Proteus vulgaris</i>
	<i>Clostridium difficile</i> non-toxigenic	<i>Pseudomonas aeruginosa</i>
	<i>Clostridium perfringens</i>	<i>Staphylococcus aureus</i>
	<i>Clostridium septicum</i>	<i>Staphylococcus aureus</i> subsp. <i>Aureus</i>
	<i>Clostridium tetani</i>	<i>Staphylococcus epidermidis</i>
	<i>Corynebacterium genitalium</i>	<i>Streptococcus agalactiae</i>
	<i>Enterobacter aerogenes</i>	<i>Streptococcus pyogenes</i>
Fungi	<i>Aspergillus fumigatus</i>	<i>Saccharomyces boulardii</i>
	<i>Candida albicans</i>	<i>Saccharomyces cerevisiae</i>
Parasites	<i>Babesia microti</i>	<i>Toxoplasma gondii</i>
	<i>Blastocystis hominis</i>	<i>Trichomonas tenax</i>
	<i>Giardia muris</i>	



Table 8. List of analytical specificity Off-panel pathogens tested (continued)

Type	Pathogen (potential cross-reactant)
Viruses	Adenovirus C:2
	Adenovirus B:34
	Adenovirus B3
	Adenovirus E:4a
	Adenovirus serotype 1
	Adenovirus serotype 5
	Adenovirus serotype 8
	Bocavirus Type 1
	Coronavirus 229E
	Coxsackievirus B3
	Cytomegalovirus
	Enterovirus 6 (Echovirus)
	Enterovirus 68
	Herpes Simplex Virus Type 2
	Rhinovirus 1A

In silico predictions of potential cross-reactions showed that the following cross-reactions may occur when testing stool samples with the QIAstat-Dx Gastrointestinal Panel 2 (Table 9).

Table 9. Potential cross-reactions based on in silico analysis

QIAstat-Dx Gastrointestinal Panel 2 target	Potential cross-reactive organisms
Enteropathogenic <i>E. coli</i> (EPEC) *	<i>Shigella boydii</i> †§ <i>Escherichia albertii</i> ††
<i>Campylobacter</i> spp.	<i>Campylobacter lari</i> § <i>Campylobacter helveticus</i> §
Shiga-like toxin-producing <i>E. coli</i> (STEC) <i>stx1/stx2</i>	<i>Shigella sonnei</i> †§ <i>Shigella dysenteriae</i> †§ <i>Acinetobacter haemolyticus</i> †** <i>Citrobacter freundii</i> †** <i>Enterobacter cloacae</i> †** <i>Aeromonas caviae</i> †** <i>Escherichia albertii</i> †¶

Table 9. Potential cross-reactions based on in silico analysis (continued)

QIAstat-Dx Gastrointestinal Panel 2 target	Potential cross-reactive organisms
<i>E. coli</i> O157*	Non-STEC <i>E.coli</i> O157 strains††
† Note that these predicted cross-reactivity identified by in silico analysis reflects sequences which can be acquired between species by horizontal gene transfer(27,82)	
‡ Rare or less common eae intimin carrier organisms (74).	
§ On-panel target.	
¶ Rare or less common Stx toxins producers (77,83,84,85,86,87)	
** In vitro testing of <i>Campylobacter lari</i> and <i>Campylobacter helveticus</i> strains at high concentration confirmed potential cross-reactivity of these <i>Campylobacter</i> species with the QIAstat-Dx Gastrointestinal Panel 2 assay.	
†† <i>E. coli</i> O157 will only be reported by the QIAstat-Dx Gastrointestinal Panel 2 when there is a positive amplification for the <i>E. coli</i> (STEC) design according to the calling algorithm. An infrequent case of an <i>E. coli</i> (STEC) and an <i>E. coli</i> O157 co-infection will not be differentiated from a single infection caused by an STEC O157:H7 strain.	

Inclusivity (analytical reactivity)

Inclusivity (Analytical Reactivity) was evaluated with gastrointestinal pathogen isolates/strains that were selected based on clinical relevance and genetic, temporal, and geographical diversity. Samples were prepared by spiking organisms into negative stool matrix resuspended in Para-Pak C&S transport media. Based on in vitro (wet) testing and in silico analysis, the QIAstat-Dx Gastrointestinal Panel 2 primers and probes are specific and inclusive for clinically prevalent and relevant strains for each pathogen tested.

In vitro (wet) testing

QIAstat-Dx Gastrointestinal Panel 2 is inclusive for 100% (114 out of 114) of the pathogen strains tested in vitro. Most pathogen strains evaluated in wet testing were detected at ≤3-fold (104/114) of the corresponding LoD reference strain. Less than 100% detection was observed for one strain each of ETEC, EIEC/*Shigella*, and Rotavirus and two strains each of STEC (one STEC O157), Adenovirus, and Norovirus at 3x LoD. Testing of these strains at 10x LoD generated the expected positive result for all replicates (Table 10).

**Table 10. Inclusivity test results for all the pathogens tested with the QIAstat-Dx Gastrointestinal Panel 2 Assay. LoD reference strain for every pathogen is written in bold.**

**Table 10a. Inclusivity test results for *Campylobacter* strains**

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
<i>Campylobacter</i>	<b><i>Campylobacter coli</i></b>	<b>76-GA2 [LMG 21266]*</b>	ATCC	43478	1x LoD
	<i>Campylobacter coli</i>	Z293	ZeptoMetrix	804272	1 x LoD
	<i>Campylobacter coli</i>	CIP 7080 [1407, CIP 70.80]	ATCC	33559*	3x LoD
	<b><i>Campylobacter jejuni</i></b>	<b>Z086</b>	<b>ZeptoMetrix</b>	<b>0801650*</b>	<b>1x LoD</b>
	<i>Campylobacter jejuni</i>	subsp. <i>jejuni</i> RM3193	ATCC	BAA-1234*	0.1x LoD
	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	O:19 HL7; D3180	ATCC	BAA-218	0.1x LoD
	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	AS-83-79	ATCC	33291	0.1x LoD
	<i>Campylobacter jejuni</i> subsp. <i>doylei</i>	NCTC 11951	ATCC	49349	0.1x LoD
	<b><i>Campylobacter upsaliensis</i></b>	<b>NCTC 11541</b>	<b>ZeptoMetrix</b>	<b>0801999*</b>	<b>1x LoD</b>
	<i>Campylobacter upsaliensis</i>	RM 3195 (1994)	ATCC	BAA-1059*	0.3x LoD
	<i>Campylobacter upsaliensis</i>	NCTC 11541 [C231]	ATCC	43954	1x LoD

\* LoD reference strain for every pathogen is written in bold.

† Strain tested during LoD verification study.

Table 10b. Inclusivity test results for *Plesiomonas shigelloides* strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
<i>Plesiomonas shigelloides</i>	<b><i>Plesiomonas shigelloides</i></b>	<b>Z130</b>	<b>ZeptoMetrix</b>	<b>0801899*</b>	<b>1 x LoD</b>
	<i>Plesiomonas shigelloides</i>	GNI 14	ATCC	51903	1 x LoD
	<i>Plesiomonas shigelloides</i>	CDC 3085-55 [Bader M51, NCIB 9242, NCTC 10360, RH 798]	ATCC	14029*	0.3x LoD

\*Strain tested during LoD verification study.

Table 10c. Inclusivity test results for *Salmonella* strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
Salmonella	<b><i>Salmonella enterica</i></b>	<b>Serovar Typhimurium Z005*</b>	<b>ZeptoMetrix</b>	<b>0801437*</b>	<b>1 x LoD</b>
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Bareilly	NCTC	NC05745	1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar typhi, Z152	ZeptoMetrix	0801933	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Enteridis, CDC K-1891 [ATCC 25928]	ATCC	13076	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Infantis, MZ1479 [SARB27]	ATCC	BAA-1675	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Montevideo, G4639	ATCC	BAA-710	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Javiana	NCTC	NC06495	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Thompson	NCTC	NC08496	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Saintpaul	ATCC	9712	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Berta	NCTC	NC05770	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Salame, II NCTC 10310 [JT945, SS140/61]	ATCC	700151	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. diarizonae IIIb, 62	ATCC	29934	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. houtenae IV, CIP 82.32 [264.66]	ATCC	43974	0.1 x LoD

**Table 10c. Inclusivity test results for *Salmonella* strains (continued)**

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
	<i>Salmonella enterica</i>	Subsp. Indica VI, CIP 102501 [F. Kauffmann 1240]	ATCC	43976	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Agona, CDC 873 [CDC 1111-61]	ATCC	51957	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Muenchen, 54	ATCC	8388	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Oranienburg, E1093	ATCC	9239	0.1 x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Paratyphi B var. Java, CDC 5	ATCC	51962	0.1 x LoD
	<i>Salmonella bongori</i>	CIP 82.33 [1224.72]	ATCC	43975	0.3x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Choleraesuis, NCTC 5735 [1348, K.34]	ATCC	13312 <sup>†</sup>	0.3x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Newport, C487-69	ATCC	27869	0.3x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, 4, 5, 12:7-:, serovar Typhimurium	NCTC	NC13952	0.3x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Braenderup	ATCC	700136	0.3x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Anatum	NCTC	NC05779	0.3x LoD
	<i>Salmonella enterica</i>	Subsp. arizonae IIIa, NCTC 7311 [CDAI 426]	ATCC	700156	0.3x LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Heidelberg, [16]	ATCC	8326	0.3x LoD

Table 10c. Inclusivity test results for *Salmonella* strains (continued)

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
	<i>Salmonella enterica</i>	Subsp. Enterica, serovar Mississippi, CDC 2012K-0487	ATCC	BAA-2739	0.3x LoD

\* Strain tested during LoD verification study.

Table 10d. Inclusivity test results for *Yersinia enterocolitica* strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
<i>Yersinia enterocolitica</i>	<b><i>Yersinia enterocolitica</i></b>	<b>Z036</b>	<b>ZeptoMetrix</b>	<b>0801734*</b>	<b>1x LoD</b>
	<i>Yersinia enterocolitica</i>	NTCC 11175, Biotype 4, serotype 3 (O:3)	ATCC	700822*	1x LoD
	<i>Yersinia enterocolitica</i>	33114 [CCUG 11291, CCUG 12369, CIP 80.27, DSM 4780, LMG 7899, NCTC 12982], Biovar 1, O:8	ATCC	9610	1x LoD
	<i>Yersinia enterocolitica</i>	O:9	ATCC	55075	3x LoD

\* Strain tested during LoD verification study.

Table 10e. Inclusivity test results for Enteropathogenic *E. coli* (EPEC) strains. Only applicable for Para-Pak C&S samples.

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
Enteropathogenic <i>E. coli</i> (EPEC)	<b>Enteropathogenic <i>E. coli</i> (EPEC)</b>	<b>O111:NM*</b>	<b>ZeptoMetrix</b>	<b>0801747*</b>	<b>1x LoD</b>
	Enteropathogenic <i>E. coli</i> (EPEC)	7.1493,O84:H28	ZeptoMetrix	0801938*	1x LoD
	Enteropathogenic <i>E. coli</i> (EPEC)	Stoke W,O111:K58[B4]:H-	ATCC	33780	1x LoD

\* Strain tested during LoD verification study.

Table 10f. Inclusivity test results for Enterotoxigenic *E. coli* (EPEC) strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
Enterotoxigenic <i>E. coli</i> (EPEC) lt/st	<b>Enterotoxigenic <i>E. coli</i> (EPEC) lt/st</b>	<b>ST+, LT+</b>	<b>ZeptoMetrix</b>	<b>0801624*</b>	<b>1x LoD</b>
	Enterotoxigenic <i>E. coli</i> (EPEC) lt/st	H10407,O78:H11,LT(+)/ctx A11(+)	ATCC	35401*	0.3x LoD
	Enterotoxigenic <i>E. coli</i> (EPEC) lt/st	O27:H7,ST (+)/ LT (-)	SSI Diagnostica	82173	0.1x LoD
	Enterotoxigenic <i>E. coli</i> (EPEC) lt/st	O115:H15,ST (+)/ LT (-)	SSI Diagnostica	82174	3x LoD
	Enterotoxigenic <i>E. coli</i> (EPEC) lt/st	O169:H,ST (-)/LT (+)	SSI Diagnostica	82172	10x LoD†

\* Strain tested during LoD verification study.

† Testing at a lower concentration resulted in a detection rate of <100%.



Table 10g. Inclusivity test results for Enteroinvasive *E. coli* (EIEC)/*Shigella* strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
Enteroinvasive <i>E. coli</i> (EIEC)/ <i>Shigella</i>	Enteroinvasive <i>E. coli</i> (EIEC)	CDC EDL 1282, O29:NM	ATCC	43892*	1x LoD
	Enteroinvasive <i>E. coli</i> (EIEC)	O172:H-	SSI Diagnostica	82171	3x LoD
	<i>Shigella sonnei</i>	NCDC 1120-66	ATCC	25931*	1x LoD
	<i>Shigella boydii</i> (Serogroup C)	Z131	ZeptoMetrix	0801900	1x LoD
	<i>Shigella flexneri</i> (Serogroup B)	AMC 43-G-68 [EVL 82, M134]	ATCC	9199	1x LoD
	<i>Shigella flexneri</i> (Serogroup B)	Z046	ZeptoMetrix	0801757	1x LoD
	<i>Shigella sonnei</i> (Serogroup D)	WRAIR I virulent	ATCC	29930	1x LoD
	<i>Shigella sonnei</i> (Serogroup D)	Z004	ZeptoMetrix	0801627	3x LoD
	<i>Shigella boydii</i> (Serogroup C)	AMC 43-G-58 [M44 (Type 170)]	ATCC	9207	10x LoD

\* Strain tested during LoD verification study.

† Testing at a lower concentration resulted in a detection rate of <100%.

**Table 10h. Inclusivity test results for Shiga-like toxin *E. coli* (STEC)(*stx1/stx2*-carrier strains)**

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx1/stx2</i>	<b>Shiga-like toxin producing <i>E. coli</i> (STEC) <i>stx1/stx2</i></b>	<b>O157:H7; EDL933</b>	<b>ZeptoMetrix</b>	<b>0801622*</b>	<b>1x LoD</b>
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx1</i>	O26:H4, <i>stx1</i> (+)	ZeptoMetrix	0801748*	1x LoD
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx1/stx2</i>	Reference ATCC® 35150 (EDL 931),O157:H7, <i>stx1</i> (+), <i>stx2</i> (+)	Microbiologics	617	3x LoD
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx1</i>	Reference CDC 00-3039,O45:H2,unknown	Microbiologics	1098	1x LoD
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx1</i>	O103:H2, <i>stx1</i> (+)	SSI Diagnostica	82170	3x LoD
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx1/stx2</i>	O22:H8, <i>stx1c</i> (+), <i>stx2b</i> (+)	SSI Diagnostica	91350	1x LoD
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx2</i>	O92,O107:K+:H48, <i>stx2d</i> (+)	SSI Diagnostica	91352	10x LoD
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx2</i>	O101:K32:H-, <i>stx2e</i> (+)	SSI Diagnostica	91354	0.3x LoD
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx1/stx2</i>	O128ac:H-, <i>stx2f</i> (+)	SSI Diagnostica	91355	10x LoD
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx2</i>	O26:H11, <i>stx2a</i> (+)	SSI Diagnostica	95211	1x LoD

Table 10h. Inclusivity test results for Shiga-like toxin *E. coli* (STEC)(*stx1/stx2*-carrier strains) (continued)

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
	Shiga-like toxin producing <i>E. coli</i> (STEC) - <i>stx1</i>	O8 , <i>stx1d</i> (+)	SSI Diagnostica	91349	1x LoD

\* Strain tested during LoD verification study

†Testing at a lower concentration resulted in a detection rate of <100%.

Table 10i. Inclusivity test results for Shiga-like toxin producing *E. coli* (STEC) *stx1/stx2* O157 strains.

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalogue ID	Times LoD
Shiga-like toxin producing <i>E. coli</i> (STEC) O157	<b>Shiga-like toxin producing <i>E. coli</i> (STEC) - O157</b>	<b>O157:H7; EDL933</b>	<b>ZeptoMetrix</b>	<b>0801622*</b>	<b>1x LoD</b>
	Shiga-like toxin producing <i>E. coli</i> (STEC) O157	O128ac:H-, <i>stx2f</i> (+)	SSI Diagnostica	91355†	1x LoD 10x LoD†
	Shiga-like toxin producing <i>E. coli</i> (STEC) O157	Reference ATCC 35150 (EDL 931), O157:H7, <i>stx1</i> (+), <i>stx2</i> (+)	Microbiologics	617	1x LoD

\* Strain tested during LoD verification study.

†The *E. coli* strain 91355 from SSI Diagnostica is reported as following in its catalog: vtx2f+, eae+. However, it was found to amplify for *E. coli* O157 in both QIAstat-Dx and an FDA-cleared test.

‡ Testing at a lower concentration resulted in a detection rate of <100%.

Table 10j. Inclusivity test results for *Cryptosporidium* strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
<i>Cryptosporidium</i>	<i>Cryptosporidium parvum</i>	Iowa isolate	Waterborne	P102C*	1x LoD
	<i>Cryptosporidium hominis</i>	n/a	Public Health Wales	Clinical sample; UKM 84*	0.01x LoD
	<i>Cryptosporidium parvum</i>	–	ATCC	PRA-67DQ (isolated genomic DNA)	<0.01 LoD
	<i>Cryptosporidium meleagridis</i>	–	Public Health Wales	Clinical sample; UKMEL 14	<0.01 LoD

\* Strain tested during LoD verification study.

Table 10k. Inclusivity test results for *Cyclospora cayetanensis* strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalogue ID	Times LoD
<i>Cyclospora cayetanensis</i>	<i>Cyclospora cayetanensis</i>	n/a	Clinical sample	LAC2825*	1x LoD
	<i>Cyclospora cayetanensis</i>	n/a	Clinical sample	LAC2827*	1x LoD
	<i>Cyclospora cayetanensis</i>	–	ATCC	PRA-3000SD	1x LoD

\* Strain tested during LoD verification study.

Table 10l. Inclusivity test results for *Entamoeba histolytica* strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalogue ID	Times LoD
<i>Entamoeba histolytica</i>	<i>Entamoeba histolytica</i>	HM-1:IMSS (Mexico City 1967)	ATCC	30459*	1x LoD
	<i>Entamoeba histolytica</i>	HK-9 (Korea)	ATCC	30015*	1x LoD
	<i>Entamoeba histolytica</i>	–	Vall d'Hebrón	Clinical sample; 1	1x LoD

\* Strain tested during LoD verification study.

Table 10m. Inclusivity test results for *Giardia lamblia* strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalogue ID	Times LoD
<i>Giardia lamblia</i>	<i>Giardia lamblia</i>	Portland -1 (Portland, OR, 1971)*	ATCC	30888*	1x LoD
	<i>Giardia lamblia</i>	WB (Bethesda, MD, 1979)	ATCC	30957*	1x LoD
	<i>Giardia intestinalis</i>	H3 isolate	Waterborne	P101	1x LoD

\* Strain tested during LoD verification study.

Table 10n. Inclusivity test results for Adenovirus F40/F41 targets

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
Adenovirus F40/F41	Human Adenovirus F41	Tak	ZeptoMetrix	0810085CF*	1x LoD
	Human Adenovirus F41	Tak (73-3544)	ATCC	VR-930	10x LoD†
	Human Adenovirus F40	Dugan [79-18025]	ATCC	VR-931	10x LoD†
	Human Adenovirus Type 40	Dugan	ZeptoMetrix	0810084CF*	3x LoD

\* Strain tested during LoD verification study.

† Testing at a lower concentration resulted in a detection rate of <100%.

Table 10o. Inclusivity test results for Astrovirus strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
Astrovirus	Human Astro-virus	ERE IID 2371 (type 8)	ZeptoMetrix	0810277CF*	1x LoD
	Human Astrovirus	HAsV-1	Universitat de Barcelona	Clinical sample; 160521599	1x LoD
	Human Astrovirus	ERE IID 2868 (type 4)	ZeptoMetrix	0810276CF*	1x LoD
	Human Astrovirus	HAsV-3	Universitat de Barcelona	Clinical sample; 151601306	1x LoD

\* Strain tested during LoD verification study.

Table 10p. Inclusivity test results for Norovirus GI/GII strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
Norovirus GI/GII	<b>Human Norovirus Genogroup 1</b>	<b>Recombinant GI.1</b>	<b>ZeptoMetrix</b>	<b>0810086CF*</b>	<b>1 x LoD</b>
	Human Norovirus Genogroup 1	–	Indiana University Health	Clinical sample; IU3156	1 x LoD
	Human Norovirus Genogroup 1	–	Indiana University Health	Clinical sample; IU3220	1 x LoD
	Human Norovirus Genogroup 1	–	TriCore Reference Laboratories	Clinical sample; TC4274	3x LoD
	<b>Human Norovirus Genogroup 2</b>	<b>Recombinant GII.4</b>	<b>ZeptoMetrix</b>	<b>0810087CF*</b>	<b>1 x LoD</b>
	Human Norovirus Genogroup 2	GI.2	Vall d'Hebrón	Clinical sample; 198058327	1 x LoD
	Human Norovirus Genogroup 2	GI.4	Universitat de Barcelona	Clinical sample; N26.2TA	1 x LoD
	Human Norovirus Genogroup 2	–	Lacny Hospital	Clinical sample; LAC2019	1 x LoD
	Human Norovirus Genogroup 2	–	Nationwide Children's Hospital	Clinical sample; NWC6063	1 x LoD
	Human Norovirus Genogroup 2	GI.6	QIAGEN Barcelona (STAT-Dx)	Clinical sample; GI 12	3x LoD

Table 10p. Inclusivity test results for Norovirus GI/GII strains (continued)

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
	Human Norovirus Genogroup 2	–	Lacny Hospital	Clinical sample; LAC2133	10x LoD
	Human Norovirus Genogroup 2	–	Lacny Hospital	Clinical sample; LAC2074	10x LoD <sup>†</sup>

\* Strain tested during LoD verification study.

<sup>†</sup> Testing at a lower concentration resulted in a detection rate of <100%.

Table 10q. Inclusivity test results for Rotavirus A strains

QIAstat-Dx target	Pathogen	Strain	Supplier	Catalog ID	Times LoD
Rotavirus A	Human Rotavirus A	69M	ZeptoMetrix	0810280CF*	1x LoD
	Human Rotavirus A	Wa, G1P1A[8]	ZeptoMetrix	0810041CF*	1x LoD
	Human Rotavirus A	DS-1, G2P1B[4]	ATCC	VR-2550	1x LoD
	Human Rotavirus A	Va70	ZeptoMetrix	0810281CF	1x LoD
	Human Rotavirus A	RRV	ZeptoMetrix	0810530CF	10x LoD <sup>†</sup>

\* Strain tested during LoD verification study.

<sup>†</sup> Testing at a lower concentration resulted in a detection rate of <100%.

In silico analysis

In silico analysis of potential reactivity showed that the following organisms (including species, subspecies, subtypes, serotypes, or serovars) are predicted to be detected with the QIAstat-Dx Gastrointestinal Panel 2 (Table 11).



Table 11. Organisms with predicted reactivity based on in silico analysis

QIAstat-Dx Gastrointestinal Panel 2  
target

Organisms with predicted reactivity

Bacteria

<i>Campylobacter</i>	<i>Campylobacter coli</i> *, <i>Campylobacter jejuni</i> , <i>Campylobacter jejuni</i> subsp. <i>jejuni</i> , <i>Campylobacter jejuni</i> subsp. <i>doylei</i> , <i>Campylobacter upsaliensis</i>
<i>Salmonella</i>	<i>Salmonella bongori</i> *, <i>Salmonella enterica</i> subsp. <i>salamae</i> II (e.g. serovar 55:k:z39), <i>Salmonella enterica</i> subsp. <i>arizonae</i> IIIa (e.g., serovar 63:g:z51), <i>Salmonella enterica</i> subsp. <i>diarizonae</i> IIIb (e.g., serovar 47:l,v:z), <i>Salmonella enterica</i> subsp. <i>houtenae</i> IV (e.g., serovar 43:z4), <i>Salmonella enterica</i> subsp. <i>indica</i> VI.  <i>Salmonella enterica</i> subsp. <i>enterica</i> (up to 92 different serovars including Agona, Anatum, Bareilly, Choleraesuis, Enteritidis, Heidelberg, Infantis, Kentucky, Montevideo, Newport, Paratyphi A*, Senftenberg, Tennessee, Thompson, Typhi, Typhimurium, Weltevreden*)
<i>Plesiomonas shigelloides</i>	<i>Plesiomonas shigelloides</i> (e.g. strains NCTC10360, ATCC 14029T, R4605035)
<i>Yersinia enterocolitica</i>	<i>Yersinia enterocolitica</i> , <i>Yersinia enterocolitica</i> subsp. <i>paleartica</i> , <i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i>
Enteroinvasive <i>E. coli</i> (EIEC)/ <i>Shigella</i>	Enteroinvasive <i>E. coli</i> (EIEC), <i>Escherichia coli</i> sp., <i>Shigella flexneri</i> , <i>Shigella dysenteriae</i> , <i>Shigella boydii</i> , <i>Shigella sonnei</i> .
Enteropathogenic <i>E. coli</i> (EPEC)	Enteropathogenic <i>E. coli</i> (EPEC) (e.g. including serotypes OUT: HND, OUT:H6, OUT:H34, OUT:H21, O55:H7, O119:HNM, O117)
Enterotoxigenic <i>E. coli</i> (ETEC) <sup>†</sup>	Enterotoxigenic <i>E. coli</i> (ETEC) (including H10407 and E24377A strains and serotypes O169:H41, O25:H42, O148:H28, O6:H16) carrier of: Heat-labile enterotoxin gene subtype LT-I and Heat-stable enterotoxin gene variant Sta, subtypes STp and STh

Table 11. Organisms with predicted reactivity based on in silico analysis (continued)

**QIAstat-Dx Gastrointestinal Panel 2 target**

**Organisms with predicted reactivity**

Shiga-like toxin-producing <i>E. coli</i> (STEC) - <i>stx1/stx2</i>	Shiga-like toxin-producing <i>E. coli</i> (STEC) including O157:H7 and O157:NM serotype and non-O157 serotypes (O111:NM, O111:H-, O26:H11, O145:NM, O145:H28, O45:H2, O26:H11, ONT:NM, O104:H4, O121:H19, O145:H34, O113:H21, ONT:H-, O128:H2, OUT:HNM, O124:HNM  <i>E. coli</i> strains carrier of:  <i>stx1a</i> , <i>stx1c</i> , <i>stx1d</i> , <i>stx2a</i> , <i>stx2b</i> , <i>stx2c</i> , <i>stx2d</i> , <i>stx2e</i> , <i>stx2f</i> , <i>stx2g</i> , <i>stx2h</i> , <i>stx2i</i> , <i>stx2j</i> , <i>stx2k</i> , and <i>stx2l</i>  Other <i>stx</i> -carrying bacteria: <i>Shigella sonnei</i> , <i>Shigella dysenteriae</i>
Shiga-like toxin-producing <i>E. coli</i> (STEC) O157	<i>Escherichia coli</i> O157 including: STEC O157:H7 strains (e.g., EDL933) and <i>E. coli</i> O157: non-H7 groups including non-Shiga-toxigenic <i>E. coli</i> O157 bacteria (e.g., serotype O157:H45)

**Parasites**

<i>Cryptosporidium</i> <sup>†</sup>	Common <i>Cryptosporidium</i> species involved in human disease <i>C. parvum</i> , <i>C. hominis</i> .  Less common <i>Cryptosporidium</i> species involved in human infections: <i>C. meleagridis</i> , <i>C. felis</i> , <i>C. bovis</i> , <i>C. viatorum</i> , <i>C. ubiquitum</i> , <i>C. tyzzeri</i> , <i>C. cuniculus</i> , <i>Cryptosporidium</i> sp. Chipmunk genotype I, <i>C. canis</i> <sup>*</sup>  Rare or non-human species: <i>Cryptosporidium wrairi</i>
<i>Cyclospora cayetanensis</i>	<i>Cyclospora cayetanensis</i> (including strains LG, CY9, NP20, and NP21) <sup>*</sup>
<i>Entamoeba histolytica</i>	<i>Entamoeba histolytica</i> (e.g. strains HM-1: IMSS, EHMfas1, and HK-9) <sup>*</sup>
<i>Giardia lamblia</i>	<i>Giardia lamblia</i> (a.k.a. <i>Giardia duodenalis</i> , <i>Giardia intestinalis</i> ) <sup>*</sup>

**Viruses**

Adenovirus	Human Adenovirus F40/41
Astrovirus <sup>§</sup>	Human Astrovirus (including types 1, 2, 3, 4, 5, 6, 7, 8)

Table 11. Organisms with predicted reactivity based on in silico analysis (continued)

QIAstat-Dx Gastrointestinal Panel 2  
target

Organisms with predicted reactivity

Norovirus GI/GII	Norovirus genogroup II genotypes: GII.1, GII.2, GII.3*, GII.4*, GII.5, GII.6, GII.7, GII.8, GII.9, GII.10, GII.12, GII.13, GII.14, GII.16, GII.17, GII.20, GII.21, GII.22, GII.23, GII.24*, GII.25, GII.26, GII.27, GII.NA1, and GII.NA2*  Norovirus genogroup I genotypes: GI.1, GI2, GI.3*, GI.4*, GI.5, GI.6*, GI.7*, GI.8, and GI.9
Rotavirus	Rotavirus A including genotypes: G1P[8]*, G2P[4]*, G3P[8]*, G4P[8]*, G9P[6], G9P[8]*, G12P[6]*, and G12P[8]*

\*Certain sequences are predicted to be detected with reduced sensitivity due to the presence of a reduced number of mismatches at critical positions of the primer-probe design.

†The assay is not predicted to detect bacteria carrier of Heat-labile enterotoxin gene subtype LT-II and/or of Heat-stable enterotoxin gene variant Stb e.

‡The assay is not predicted to detect other *Cryptosporidium* spp. less involved in human disease: *C. andersoni* and *C. muris* (88).

§The assay is not predicted to detect Human Astrovirus types MLB1-3 and VA1-5.

Interfering substances

The effect of potentially interfering substances on the detectability of the QIAstat-Dx Gastrointestinal Panel 2 organisms was evaluated. Thirty- four (34) potentially interfering substances were spiked into the sample mixes at a level predicted to be above the concentration of the substance likely to be found in stool specimens. Endogenous substances such as human whole blood, human genomic DNA, and several pathogens were tested alongside exogenous substances like antibiotics, other gastrointestinal-related medications, and different technique-specific substances.

Testing included samples containing negative clinical stool matrix in Para-Pak C&S media with and without addition of each potentially interfering substance. Samples containing organism mixes with one strain for each targeted pathogen were tested at 3x LoD. Testing was

performed in triplicate. Additionally, for endogenous substances, negative specimens (stool matrix in Para-Pak C&S media matrix with no organism mix) were spiked with only the test substance to evaluate the potential for false positive results due to the test substance itself.

For the vast majority of substances tested, no interference was observed, with the exceptions of mucin, calcium carbonate, nonoxynol-9 and Rotavirus reassortants, that demonstrated interference at high concentration.

Mucin at 5% w/v was found to generate false positives results for the *Yersinia* target. These signals were investigated by testing the interfering substance with an FDA-cleared method and they were confirmed to be present in the endogenous substance.

Calcium carbonate at concentrations above 0.5% w/v was found to generate false negative results for all the QIAstat-Dx Gastrointestinal Panel 2 targets and the internal control.

Nonoxynol-9 at concentrations above 0.02% v/v was found to generate false negative results for detection of *Entamoeba*.

As predicted, Rotavirus reassortants WC3:2-5, R574(9) and WI79-4,9 used in Rotavirus A vaccines generated positive results for Rotavirus A in the QIAstat-Dx Gastrointestinal Panel 2. Final concentrations without interference (i.e., no false positive results for Rotavirus) for WC3:2-5, R574(9) and WI79-4,9 were  $8.89 \times 10^{-5}$  TCID<sub>50</sub>/mL and 1.10 PFU/mL, respectively; refer to (Table 12) for other concentrations tested.

Results from the 34 interfering substances that could be present or introduced in a stool specimen are provided in Table 12.

**Table 12. Final highest concentration without observable inhibitory effect**

Substance tested	Concentration tested	Result
<b>Endogenous substances</b>		
Bovine and ovine bile	12% w/v	No Interference
Cholesterol	1.5% w/v	No Interference
Fatty acids (palmitic acid)	0.2% w/v	No Interference
Fatty acids (stearic acid)	0.4% w/v	No Interference
Human genomic DNA	20 µg/mL	No Interference
Human stool (overfill of Cary-Blair vial)	300 mg/mL	No Interference
Human urine	50% v/v	No Interference
Human whole blood with Na Citrate	40% v/v	No Interference
Mucin from bovine submaxillary	5% w/v	Interference†
	2.5% w/v	No Interference
Triglycerides	5% w/v	No Interference
<b>Exogenous substances</b>		
Bacitracin	250 U/mL	No Interference
Bisacodyl	0.3% w/v	No Interference
Bismuth subsalicylate	0.35% w/v	No Interference
Calcium carbonate (TUMS® Extra Strength 750)	5% w/v	Interference
	0.5% w/v	No Interference
Docusate sodium	2.5% w/v	No Interference
Doxycycline hydrochloride	0.05% w/v	No Interference
Glycerin	50% v/v	No Interference
Hydrocortisone	0.5% w/v	No Interference

Table 12. Final highest concentration without observable inhibitory effect (continued)

Substance tested	Concentration tested	Result
Loperamide hydrochloride	0.078% w/v	No Interference
Magnesium hydroxide	0.1% w/v	No Interference
Metronidazole	1.5% w/v	No Interference
Mineral oil	50% v/v	No Interference
Naproxen sodium	0.7% w/v	No Interference
Nonoxynol-9	1.2% v/v	No Interference
	0.6% v/v	Interference
	0.3% v/v	Interference
	0.15% v/v	Interference
	0.075% v/v	Interference
	0.02% v/v	No Interference
Nystatin	10,000 USP units/mL	No Interference
Phenylephrine hydrochloride	0.075% w/v	No Interference
Sodium phosphate	5% w/v	No Interference
Vaccine components		
Rotavirus reassortant WC3:2-5, R574 (9) - VR 2195	$8.89 \times 10^{-3}$ TCID <sub>50</sub> /mL	Interference
	$8.89 \times 10^{-4}$ TCID <sub>50</sub> /mL	Interference
	$8.89 \times 10^{-5}$ TCID <sub>50</sub> /mL	No interference
Rotavirus reassortant WI79-4,9 - VR 2415	$1.10 \times 10^2$ pfu/mL	Interference
	$1.10 \times 10^1$ pfu/mL	Interference
	1.10 pfu/mL	No interference
Technique-specific Substances, Transport Media		
Bleach	0.5% v/v	No Interference
Ethanol	0.2% v/v	No Interference

Table 12. Final highest concentration without observable inhibitory effect (continued)

Substance tested	Concentration tested	Result
Puritan Fecal Opti-Swab Collection & Transport System with Cary-Blair Medium*	100%	No Interference
Puritan PurSafe® DNA/RNA Preservative*	100%	No Interference
Sigma Fecal Transwab*	1 swab/2 mL Cary-Blair	No Interference

\* Performance not established for this transport media.

† This substance was tested by another FDA-cleared test that also detected *Yersinia* positive signals.

Microbial interference

A microbial interference study was conducted to assess the inhibitory effects of select non-target organisms on the ability to detect the QIAstat-Dx Gastrointestinal Panel 2 targets. Clinically relevant and challenging concentrations of non-target organisms ( $1 \times 10^6$  CFU/mL for bacteria,  $1 \times 10^5$  cells/mL for yeast, and  $1 \times 10^5$  TCID<sub>50</sub>/mL for viruses) were individually mixed with negative clinical stool matrix with spiked targeted pathogens at 3x LoD. Testing was performed in triplicate. All combinations and replicates successfully detected all the QIAstat-Dx Gastrointestinal Panel 2 targets. See Table 13 for a list of the non-target organisms tested and the result summary.

Table 13. Final highest concentration without observable inhibitory effect

Substance tested	Concentration tested	Result
Non-target microorganisms		
<i>Aeromonas hydrophila</i>	$1 \times 10^6$ units/mL	No Interference
<i>Bacteroides vulgatus</i>	$1 \times 10^6$ units/mL	No Interference
<i>Bifidobacterium bifidum</i>	$1 \times 10^6$ units/mL	No Interference
Enterovirus Species D, Serotype EV-D68	$1 \times 10^5$ units/mL	No Interference

Table 13. Final highest concentration without observable inhibitory effect (continued)

Substance tested	Concentration tested	Result
Non-pathogenic <i>E. coli</i>	1 x 10 <sup>6</sup> units/mL	No Interference
<i>Helicobacter pylori</i>	1 x 10 <sup>6</sup> units/mL	No Interference
<i>Saccharomyces cerevisiae</i> (deposited as <i>S. boulardii</i> )	1 x 10 <sup>5</sup> units/mL	No Interference

Competitive interference

Competitive interference was tested in a subset of pathogens. No interference was observed when evaluating competitive interference by target pathogens when two QIAstat-Dx Gastrointestinal Panel target pathogens were tested by spiking samples with one pathogen target at 3x LoD and one at 50x LoD. Results from the pathogen targets tested are provided in Table 14.

Table 14. QIAstat-Dx Gastrointestinal Panel 2 results for competitive interference

Sample Mix	Target	Final concentration (molecular units)*	Final concentration tested x LoD	Co-infection detected
Norovirus 50x - Rotavirus 3x	Norovirus GI/GII	4.5E+05 copies/mL	50x	Yes
	Rotavirus A	1.7E+04 copies/mL	3x	
Norovirus 3x - Rotavirus 50x	Norovirus GI/GII	2.7E+04 copies/mL	3x	Yes
	Rotavirus A	2.9E+05 copies/mL	50x	
Giardia 50x - Adenovirus 3x	<i>Giardia lamblia</i>	7.2 E+05 copies/mL	50x	Yes
	Adenovirus F40/F41	2.9E+03 copies/mL	3x	

\*Molecular unit titers were determined using in-house developed and validated qPCR assays.



## Carryover

A carryover study was performed to evaluate the potential occurrence of cross-contamination between consecutive runs when using the QIAstat-Dx Gastrointestinal Panel 2 on the QIAstat-Dx Analyzer 1.0.

Pathogen samples of stool sample matrix in Para-Pak C&S transport media, with alternating high-positive ( $10^6$  CFU/mL for bacteria,  $10^5$  TCID<sub>50</sub> or organism/mL for viruses and parasites) and negative samples, were conducted on two QIAstat-Dx Analyzer 1.0 instruments.

No carryover between samples was observed in the QIAstat-Dx Gastrointestinal Panel 2, demonstrating that the system design and recommended sample handling and testing practices are effective in preventing false-positive results due to carryover or cross-contamination between samples.

## Reproducibility

Reproducibility testing of contrived samples was performed at three test sites including one internal site (Site A) and two external sites (Site B and Site C). The study incorporated a range of potential variations introduced by sites, days, replicates, cartridge lots, operators, and QIAstat-Dx analyzers. For each site, testing was performed across 5 non-consecutive days with 6 replicates per day (leading to a total of 30 replicates per target, concentration, and site), 4 QIAstat-Dx Analyzers (2 analyzers per operator and per site), and at least 2 operators on each testing day. A total of 5 sample mixes (two combined samples at 1x LoD and 3x LoD plus one negative sample) were prepared. For each mix, 6 replicates were tested and evaluated.

Table 15 shows the detection rate per target and concentration for each site of the Reproducibility study. In addition, data obtained at all three sites have been compiled to calculate the exact 2-sided 95% Confidence Interval by target and concentration. During the Reproducibility study, potential variations introduced by sites, days, replicates, cartridge lots,

operators, and QIAstat-Dx analyzers were analyzed showing no significant contribution to variability (Standard Deviation and Coefficient of Variation values below 1% and 5%, respectively) caused by any of the assessed variables.

Table 15. Detection rate per target and concentration for each site of the Reproducibility study and exact 2-sided 95% Confidence Interval by target and concentration

Pathogen tested	Concentration tested	Expected result	% Agreement with expected result			
			Site A	Site B	Site C	All sites (95% Confidence Interval)
Adenovirus F41 ZeptoMetrix 0810085CF	3x LoD	Detected	30/30 100%	30/30 100%	30/30 100%	90/90 100% (95.98 – 100.00%)
	1x LoD	Detected	30/30 100%	30/30 100%	30/30 100%	90/90 100% (95.98 – 100.00%)
	None	Not Detected	29/30 96.67%		29/30 96.67%	87/90* 96.7% (90.98 – 98.9%)
Campylobacter ZeptoMetrix 801650	3x LoD	Detected	30/30 100%	30/30 100%	30/30 100%	90/90 100% (95.98 – 100.00%)
	1x LoD	Detected	30/30 100%	30/30 100%	30/30 100%	90/90 100% (95.98 – 100.00%)
	None	Not Detected	30/30 100%	30/30 100%	30/30 100%	90/90 100% (95.98 – 100.00%)

Table 15. Detection rate per target and concentration for each site of the Reproducibility study and exact 2-sided 95% Confidence Interval by target and concentration (continued)

Pathogen tested	Concentration tested	Expected result	% Agreement with expected result			All sites (95% Confidence Interval)
			Site A	Site B	Site C	
<i>Escherichia coli</i> (EPEC) ZeptoMetrix 801747	3x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100% (95.98 – 100.00%)
	1x LoD	Detected	30/30	29/30	30/30	89/90
			100%	96.67%	100%	98.89% (93.96 – 99.97%)
	None	Not Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100% (95.98– 100.00%)
<i>Entamoeba histolytica</i> ATCC 30459	3x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100% (95.98 – 100.00%)
	1x LoD	Detected	30/30	30/30	29/30	89/90
			100%	100%	96.67%	98.89% (93.96 – 99.97%)
	None	Not Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100% (95.98 – 100.00%)

Table 15. Detection rate per target and concentration for each site of the Reproducibility study and exact 2-sided 95% Confidence Interval by target and concentration (continued)

Pathogen tested	Concentration tested	Expected result	% Agreement with expected result			All sites (95% Confidence Interval)
			Site A	Site B	Site C	
<i>Giardia lamblia</i> † ATCC 30888	3x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100%
						(95.98 – 100.00%)
	1x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100%
						(95.98 – 100.00%)
	None	Not Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100%
						(95.98 – 100.00%)
Norovirus GII ZeptoMetrix 0810087CF	3x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100%
						(95.98 – 100.00%)
	1x LoD	Detected	29/30	30/30	30/30	89/90
			96.67%	100%	100%	98.89%
						(93.96 – 99.97%)
	None	Not Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100%
						(95.98 – 100.00%)

Table 15. Detection rate per target and concentration for each site of the Reproducibility study and exact 2-sided 95% Confidence Interval by target and concentration (continued)

Pathogen tested	Concentration tested	Expected result	% Agreement with expected result			All sites (95% Confidence Interval)
			Site A	Site B	Site C	
Rotavirus A†  ZeptoMetrix 0810280CF	3x LoD	Detected	29/30 96.67%	29/30 96.67%	30/30 100%	88/90 97.8% (92.20 – 99.73%)
			23/30 76.67%	26/30 86.67%	12/12 100%	61/72 84.7% (74.31 – 92.12%)
	None	Not Detected	30/30 100%	30/30 100%	30/30 100%	90/90 100% (95.98 – 100.00%)
	3x LoD	Detected	30/30 100%	30/30 100%	30/30 100%	90/90 100% (95.98 – 100.00%)
Escherichia coli (STEC) O157:H7§  ZeptoMetrix 0801622	3x LoD	Detected	30/30 100%	30/30 100%	30/30 100%	90/90 100% (95.98 – 100.00%)
	1x LoD	Detected	30/30 100%	30/30 100%	29/30 96.67%	89/90 98.89% (93.96 – 99.97%)
	None	Not Detected	30/30 100%	30/30 100%	30/30 100%	90/90 100% (95.98 – 100.00%)

Table 15. Detection rate per target and concentration for each site of the Reproducibility study and exact 2-sided 95% Confidence Interval by target and concentration (continued)

Pathogen tested	Concentration tested	Expected result	% Agreement with expected result			All sites (95% Confidence Interval)
			Site A	Site B	Site C	
<i>Escherichia coli</i> (STEC) <i>stx1/stx2</i> ZeptoMetrix 0801622	3x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100% (95.98 – 100.00%)
	1x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100% (95.98 – 100.00%)
	None	Not Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100% (95.98 – 100.00%)
<i>Salmonella enterica</i> ZeptoMetrix 0801437	3x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100% (95.98 – 100.00%)
	1x LoD	Detected	30/30	29/30	29/30	88/90
			100%	96.67%	96.67%	97.78% (92.20 – 99.73%)
	None	Not Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100% (95.98 – 100.00%)

Table 15. Detection rate per target and concentration for each site of the Reproducibility study and exact 2-sided 95% Confidence Interval by target and concentration (continued)

Pathogen tested	Concentration tested	Expected result	% Agreement with expected result			All sites (95% Confidence Interval)
			Site A	Site B	Site C	
<i>Yersinia enterocolitica</i> Zeptomatrix 0801734	3x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100%
	1x LoD	Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100%
	None	Not Detected	30/30	30/30	30/30	90/90
			100%	100%	100%	100%

\* Three (3) Adenovirus F40/41 false positives were observed when testing negative sample. Retesting of the 3 samples resulted in the expected negative results.

† One (1) *Giardia lamblia* false positive was observed when testing a positive sample not containing the pathogen. Repeat testing of this sample resulted in the expected negative result.

‡ The Reproducibility study was fully re-tested for Rotavirus A with a new sample set due to an unexpected number of false negatives for Rotavirus A at the 1x LoD concentration. This was observed during an interim data evaluation (61/72, 84.7%) and that was attributed to the sample manufacture and while unrelated to the study workflow variables (operator, lot, day, instrument, and site). Test runs derived from Rotavirus A new sample set resulted in 90/90 (100%; 95.98-100% CI) for the 3x LoD and 89/90 (98.89%; 93.96-99.97% CI) for the 1x LoD. During this testing, one (1) *Campylobacter* false positive was observed. Retesting of this sample resulted in the expected negative result.

Repeatability

A Repeatability study was conducted on the QIAstat-Dx Analyzer 1.0 instruments using a set of samples composed of low-concentrated analytes spiked into stool matrix (3x LoD and 1x LoD) and negative stool samples. QIAstat-Dx Gastrointestinal Panel 2 detected pathogens included in the positive samples were Adenovirus, *Campylobacter*, *Entamoeba histolytica*,

*Giardia lamblia*, Norovirus GII, Rotavirus, *Salmonella enterica*, *Yersinia enterocolitica*, Enteropathogenic *E. coli* (EPEC), STEC *stx1/stx2*, and *E. coli* O157. Each sample was tested with the same instrument over 12 days. In total, 60 replicates of 1x LoD and 60 replicates of 3x LoD per each of the tested targets and 60 replicates of negative samples were run. Overall results showed a 93.33–100.00% and 95.00–100.00% detection rate for 1x LoD and 3x LoD samples, respectively. Negative samples showed 100% of negative calls for all panel analytes.

Expected values

The number and percentage of positive results as determined by the QIAstat-Dx Gastrointestinal Panel 2 in the prospective clinical evaluation, stratified by age group, are presented in Table 16. Overall, the QIAstat-Dx Gastrointestinal Panel 2 detected at least 1 organism 17.4% (213/1222) and 23.8% (171/717) of the prospectively collected stool specimens in FecalSwab and Para-Pak C&S, respectively.

Table 16. Expected values summary by age group for the prospective clinical study as determined by the QIAstat-Dx Gastrointestinal Panel 2

Pathogen	Medium brand	Overall	0-5 years	6-21 years	22-49 years	50+ years	Not reported
Viruses							
Adenovirus F40/F41	FecalSwab	5 (0.4%)	3 (1.7%)	2 (1.7%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
	Para-Pak C&S	2 (0.3%)	1 (3.2%)	0 (0.0%)	0 (0.0%)	1 (0.2%)	0 (0.0%)
Astrovirus	FecalSwab	3 (0.2%)	3 (1.6%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
	Para-Pak C&S	6 (0.8%)	2 (6.5%)	0 (0.0%)	3 (1.4%)	0 (0.0%)	0 (0.0%)



Table 16. Expected values summary by age group for the prospective clinical study as determined by the QIAstat-Dx Gastrointestinal Panel 2 (continued)

Pathogen	Medium brand	Overall	0-5 years	6-21 years	22-49 years	50+ years	Not reported
Norovirus GI/GII	FecalSwab	43 (3.5%)	22 (12.1%)	1 (0.8%)	14 (4.8%)	6 (1.0%)	0 (0.0%)
	Para-Pak C&S	16 (2.3%)	3 (9.7%)	1 (2.8%)	3 (1.4%)	9 (2.2%)	0 (0.0%)
Rotavirus A	FecalSwab	23 (1.9%)	13 (7.1%)	2 (1.7%)	7 (2.4%)	1 (0.2%)	0 (0.0%)
	Para-Pak C&S	4 (0.6%)	2 (6.5%)	0 (0.0%)	0 (0.0%)	2 (0.5%)	0 (0.0%)
Bacteria							
Campylobacter	FecalSwab	69 (5.6%)	25 (13.7%)	7 (5.8%)	17 (5.9%)	20 (3.2%)	0 (0.0%)
	Para-Pak C&S	30 (4.2%)	2 (6.5%)	0 (0.0%)	10 (4.7%)	18 (4.3%)	0 (0.0%)
Plesiomonas shigelloides	FecalSwab	2 (0.2%)	0 (0.0%)	0 (0.0%)	2 (0.7%)	0 (0.0%)	0 (0.0%)
	Para-Pak C&S	7 (1.0%)	1 (3.2%)	0 (0.0%)	4 (1.9%)	2 (0.5%)	0 (0.0%)
Salmonella	FecalSwab	14 (1.1%)	5 (2.7%)	4 (3.3%)	3 (1.0%)	2 (0.3%)	0 (0.0%)
	Para-Pak C&S	17 (2.4%)	4 (12.9%)	0 (0.0%)	3 (1.4%)	10 (2.4%)	0 (0.0%)
Yersinia enterocolitica	FecalSwab	22 (1.8%)	3 (1.6%)	2 (1.7%)	9 (3.1%)	8 (1.3%)	0 (0.0%)
	Para-Pak C&S	8 (1.1%)	0 (0.0%)	0 (0.0%)	4 (1.9%)	4 (1.0%)	0 (0.0%)

**Table 16. Expected values summary by age group for the prospective clinical study as determined by the QIAstat-Dx Gastrointestinal Panel 2 (continued)**

Pathogen	Medium brand	Overall	0-5 years	6-21 years	22-49 years	50+ years	Not reported
<b>Diarrheagenic <i>E. coli</i>/Shigella</b>							
<b>Enteropathogenic <i>E. coli</i> (EPEC)</b>	FecalSwab	132 (10.8%)	47 (25.8%)	12 (9.9%)	34 (11.7%)	39 (6.2%)	0 (0.0%)
	Para-Pak C&S	56 (7.9%)	9 (29.0%)	2 (5.6%)	18 (8.4%)	27 (6.5%)	0 (0.0%)
<b>Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i></b>	FecalSwab	18 (1.5%)	2 (1.1%)	2 (1.7%)	11 (3.8%)	3 (0.5%)	0 (0.0%)
	Para-Pak C&S	17 (2.4%)	1 (3.2%)	0 (0.0%)	7 (3.3%)	9 (2.2%)	0 (0.0%)
<b>Shiga-like toxin <i>E. coli</i> (STEC) <i>stx1/stx2</i></b>	FecalSwab	15 (1.2%)	9 (4.9%)	1 (0.8%)	2 (0.7%)	3 (0.5%)	0 (0.0%)
	Para-Pak C&S	9 (1.3%)	0 (0.0%)	0 (0.0%)	6 (2.8%)	3 (0.7%)	0 (0.0%)
<b><i>E. coli</i> O157</b>	FecalSwab	3 (0.2%)	3 (1.6%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
	Para-Pak C&S	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
<b>Shigella/Enteroinvasive <i>E. coli</i> (EIEC)</b>	FecalSwab	10 (0.8%)	1 (0.5%)	0 (0.0%)	6 (2.1%)	3 (0.5%)	0 (0.0%)
	Para-Pak C&S	3 (0.4%)	0 (0.0%)	0 (0.0%)	1 (0.5%)	2 (0.5%)	0 (0.0%)
<b>Parasites</b>							

Table 16. Expected values summary by age group for the prospective clinical study as determined by the QIAstat-Dx Gastrointestinal Panel 2 (continued)

Pathogen	Medium brand	Overall	0-5 years	6-21 years	22-49 years	50+ years	Not reported
<i>Cryptosporidium</i>	FecalSwab	2 (0.2%)	0 (0.0%)	1 (0.8%)	1 (0.3%)	0 (0.0%)	0 (0.0%)
	Para-Pak C&S	7 (1.0%)	0 (0.0%)	1 (2.8%)	4 (1.9%)	2 (0.5%)	0 (0.0%)
<i>Cyclospora cayentanensis</i>	FecalSwab	3 (0.2%)	0 (0.0%)	1 (0.8%)	2 (0.7%)	0 (0.0%)	0 (0.0%)
	Para-Pak C&S	18 (2.5%)	0 (0.0%)	0 (0.0%)	6 (2.8%)	12 (2.9%)	0 (0.0%)
<i>Giardia lamblia</i>	FecalSwab	15 (1.2%)	3 (1.6%)	1 (0.8%)	7 (2.4%)	4 (0.6%)	0 (0.0%)
	Para-Pak C&S	1 (0.1%)	1 (3.2%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
<i>Entamoeba histolytica</i>	FecalSwab	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
	Para-Pak C&S	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)

# Clinical performance

The clinical performance shown below was demonstrated using QIAstat-Dx Analyzer 1.0. The QIAstat-Dx Analyzer 2.0 uses the same Analytical Modules as QIAstat-Dx Analyzer 1.0; therefore, the performance is not impacted by the QIAstat-Dx Analyzer 2.0. The clinical performance of QIAstat-Dx Gastrointestinal Panel 2 was established during a multi-center international prospective study conducted at thirteen clinical settings representatives of different geographical areas within US and Europe (9 sites in US and 4 sites in Europe) between May and July 2021. All study sites were hospital-associated or independent clinical diagnostics laboratories that perform routine diagnostics of GI infections. A total of 1939 prospectively collected stool specimens (stool preserved in Para-Pak C&S (Meridian Bioscience) or FecalSwab (COPAN)) were obtained from patients with clinical indications of diarrhea caused by gastrointestinal infection. Table 17 provides a summary of prospective specimen distribution across all study sites.

**Table 17. Prospective specimens distribution across the study sites**

Site/Country	FecalSwab	Para-Pak C&S	Total
Germany	293	46	339
Denmark	293	0	293
Spain	247	0	247
France	63	0	63
USA site 1	0	186	186
USA site 2	0	43	43
USA site 3	282	0	282
USA site 4	0	177	177
USA site 5	44	0	44
USA site 6	0	39	39

Table 17. Prospective specimens distribution across the study sites (continued)

Site/Country	FecalSwab	Para-Pak C&S	Total
USA site 7	0	0	0*
USA site 8	0	131	131
USA site 9	0	95	95
Total	1222	717	1939

\* The specimens from this site were excluded from the analysis because they were collected with another device different to Para-Pak C&S or FecalSwab.

The demographic information for the 1939 specimens evaluated in the prospective study is summarized in Table 18.

Table 18. Demographic data for prospective evaluated specimens

Demographic data	FecalSwab		Para-Pak C&S	
	N	%	N	%
<b>Gender</b>				
Female	628	32.4	442	22.8
Male	594	30.6	275	14.2
<b>Age Group</b>				
0–5 years	182	9.4	31	1.6
6–21 years	121	6.2	38	2.0
22–49 years	290	15.0	215	11.1
50+ years	629	32.4	426	22.0
Not Reported	0	0.0	7	0.4
<b>Patient population</b>				
Emergency room	46	2.4	29	1.5
Hospitalized	342	17.6	143	7.4
Immunocompromised	3	0.2	0	0.0
Outpatient	491	25.3	325	16.8
No information available	340	17.5	220	11.3
<b>No. of days between symptom onset and QIAstat-Dx testing</b>				
>7 days	89	4.6	0	0.0
≤7 days	146	7.5	16	0.8
Not Reported	987	50.9	701	36.2

The performance of the QIAstat-Dx Gastrointestinal Panel 2 was evaluated for each panel test result using one FDA-cleared test as comparator for the most analytes. A composite comparator consisting of either three independent FDA-cleared test methods or two independent FDA-cleared tests methods and two validated PCR assays followed by bi-directional sequencing was used for Norovirus GI/GII, ETEC, STEC and *Giardia lamblia* (Table 19).

**Table 19. Comparator methods for the clinical evaluation of QIAstat-Dx Gastrointestinal Panel 2**

QIAstat-Dx Gastrointestinal Panel 2 test result	Comparator method
Adenovirus F40/F41	One FDA-cleared test method
Astrovirus	
Rotavirus A	
<i>Campylobacter</i>	
<i>Plesiomonas shigelloides</i>	
<i>Salmonella</i>	
<i>Yersinia enterocolitica</i>	
<i>E. coli</i> O157	
Enteropathogenic <i>E. coli</i> (EPEC)	
<i>Shigella</i> /Enteroinvasive <i>E. coli</i> (EIEC)	
<i>Cryptosporidium</i>	
<i>Cyclospora cayetanensis</i>	
<i>Entamoeba histolytica</i>	
Norovirus GI/GII	Composite of three FDA-cleared test methods
Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i>	Composite of three FDA-cleared test methods
Shiga-like toxin <i>E. coli</i> (STEC) <i>stx1/stx2</i>	Composite of three FDA-cleared test methods

Table 19. Comparator methods for the clinical evaluation of QIAstat-Dx Gastrointestinal Panel 2 (continued)

QIAstat-Dx Gastrointestinal Panel 2 test result	Comparator method
<i>Giardia lamblia</i>	Composite of two FDA-cleared test methods and two validated PCR tests followed by bi-directional sequencing†

† Each of the two PCR assays used were well-characterized and validated nucleic acid amplification tests (NAAT) followed by bi-directional sequencing analysis. Both assays were designed to amplify different sequences than those targeted by the QIAstat-Dx Gastrointestinal Panel 2. Positive results required to generate sequences from bi-directional sequencing with at least 200 bases of adequate quality that by BLAST analyses matched a sequence of the expected organism or gene from NCBI GenBank database with at least 95% query coverage and at least 95% identity compared to the reference.

Additional prospective archived samples were collected for Norovirus GI/GII (81 samples) and STEC (18 samples). These were prospectively collected samples from four different collection sites (3 US and 1 EU), where only those positive for the pathogen by standard of care method were archived for analysis alongside 20 negative specimens. A second collection of 75 prospective archived samples positive for STEC preserved in FecalSwab from three different collection sites in the US and 17 negative specimens were also analyzed.

In addition, to supplement the results of the prospective clinical studies, a total of 750 preselected archived frozen (retrospective) specimens were also evaluated. These specimens served to increase the sample size for analytes that showed low prevalence in the clinical prospective study or that were less represented in a particular sample type (Para-Pak C&S or FecalSwab). The same Comparator Methods detailed in Table 19 were used as confirmatory testing for the presence of the nucleic acids from the expected analytes. In total, 2900 specimens (1939 prospective, 211 prospective archived, and 750 retrospective) were evaluated in the clinical study. These specimens were collected using Para-Pak C&S (1217) or FecalSwab (1683).

The positive percentage agreement (PPA) and the negative percentage agreement (NPA) were calculated for the prospective and retrospective studies and for each sample type (Para-Pak C&S and FecalSwab) separately.



The PPA was calculated as  $100\% \times (TP / (TP + FN))$ . True positive (TP) indicates that both the QIAstat-Dx Gastrointestinal Panel 2 and comparator method showed a positive result for this specific target, and false negative (FN) indicates that the QIAstat-Dx Gastrointestinal Panel 2 result was negative while the comparator method result was positive. The NPA was calculated as  $100\% \times (TN / (TN + FP))$ . True negative (TN) indicates that both the QIAstat-Dx Gastrointestinal Panel 2 and the comparator method showed negative results, and a false positive (FP) indicates that the QIAstat-Dx Gastrointestinal Panel 2 result was positive, but the comparator method result was negative. The exact binomial two-sided 95% Confidence Intervals for PPA and NPA were calculated.

Where a composite comparator was used (Table 19), the result was determined as the majority of the three individual test results (i.e., a positive composite comparator result is based on positive results for at least two comparator tests and a negative composite comparator result is based on negative results for at least two comparator tests). If insufficient pathogen positive sample was available to obtain a majority test result, a worst-case model was applied in the PPA calculation. In this model, the PPA was calculated including all observed true positive and false negative samples between QIAstat-Dx and the composite comparator, while for the samples where it was not possible to conduct testing with the complete comparator due to insufficient sample volume, the following was done:

- Samples that were negative in QIAstat-Dx and positive for one comparator assay, negative (or insufficient volume) for a second comparator and insufficient volume for a third comparator were included in the calculations as worst-case false negatives;
- Samples that were positive in QIAstat-Dx and positive in one comparator test, negative (or insufficient volume) for a second comparator and insufficient volume for the third comparator, were considered as worst-case false positives and, therefore, excluded in the PPA calculations.

The results of the clinical performance of the prospective, prospective archived, and retrospective studies are summarized in Table 20, Table 21, and Table 22, respectively.

Discrepancies between the QIAstat-Dx Gastrointestinal Panel 2 and the comparator methods were investigated for the analytes that the QIAstat-Dx Gastrointestinal Panel 2 test result was compared to one FDA-cleared method. Discrepancies analyses are footnoted on each clinical performance summary (Table 20 and Table 22).

Table 20. Clinical Performance in the Prospective study

Analyte	Medium brand	Positive Percent Agreement		Negative Percent Agreement			
		TP/TP+FN	%	95% CI	TN/TN+FP	%	95% CI
Virus							
Adenovirus F40/F41	FecalSwab	5/6 <sup>a</sup>	83.3	43.7–97.0	1216/1216	100.0	99.7–100.0
	Para-Pak C&S	1/2 <sup>b</sup>	50.0	9.5–90.6	703/704 <sup>b</sup>	99.9	99.2–100.0
Astrovirus	FecalSwab	3/3	100.0	43.9–100.0	1219/1219	100.0	99.7–100.0
	Para-Pak C&S	6/6	100.0	61.0–100.0	700/700	100.0	99.5–100.0
Norovirus GI/GII	FecalSwab	31/33 <sup>c</sup>	93.9	80.4–98.3	493/495 <sup>c</sup>	99.6	98.6–100.0
	Para-Pak C&S	14/18 <sup>d</sup>	77.8	54.8–98.3	399/399 <sup>d</sup>	100.0	99.1–100.0
Rotavirus A	FecalSwab	21/23 <sup>e</sup>	91.3	73.2–97.6	1197/1199 <sup>e</sup>	99.8	99.4–100.0
	Para-Pak C&S	3/3 <sup>f</sup>	100.0	43.9–100.0	702/703 <sup>f</sup>	99.9	99.2–100.0
Bacteria							
Campylobacter	FecalSwab	65/67 <sup>g</sup>	97.0	89.8–99.2	1151/1155 <sup>g</sup>	99.7	99.1–99.9
	Para-Pak C&S	30/31 <sup>h</sup>	96.8	83.8–99.4	675/677 <sup>h</sup>	99.7	98.9–99.9
Plesiomonas shigelloides	FecalSwab	0/0	N/A	N/A	1220/1222 <sup>i</sup>	99.8	99.4–100.0
	Para-Pak C&S	5/6 <sup>i</sup>	83.3	43.7–97.0	698/700 <sup>i</sup>	99.7	99.0–99.9
Salmonella	FecalSwab	14/16 <sup>k</sup>	87.5	64.0–96.5	1206/1206	100.0	99.7–100.0
	Para-Pak C&S	19/20 <sup>j</sup>	95.0	76.4–99.1	688/688	100.0	99.4–100.0

Table 20. Clinical Performance in the Prospective study (continued)

Analyte	Medium brand	Positive Percent Agreement		Negative Percent Agreement		
		TP/TP+FN	%	95% CI	TN/TN+FP	%
<i>Yersinia enterocolitica</i>	FecalSwab	1.5 / 16 <sup>m</sup>	93.8	71.7–99.1	1199 / 1206 <sup>m</sup>	99.4
	Para-Pak C&S	3/3	100.0	43.9–100.0	698/703 <sup>n</sup>	99.3
Diarrheagenic <i>E. coli</i> / <i>Shigella</i>						
Enteropathogenic <i>E. coli</i> (EPEC)	FecalSwab	126/145	86.9	80.4–91.5	1059/1063	99.6
	Para-Pak C&S	57/65	87.7	77.6–93.6	632/632	100.0
Enterotoxigenic <i>E. coli</i> (ETEC) lt/st	FecalSwab	9/10 <sup>o</sup>	90.0	59.6–99.2	427/430 <sup>o</sup>	99.3
	Para-Pak C&S	9/10 <sup>o</sup>	90.0	59.6–99.2	390/395 <sup>p</sup>	98.7
Shiga-like toxin <i>E. coli</i> (STEC) <i>stx1/stx2</i>	FecalSwab	3 / 5 <sup>q</sup>	60.0	23.1–88.2	434–438 <sup>q</sup>	99.1
	Para-Pak C&S	5/6 <sup>q</sup>	83.3	43.7–97.0	397/400 <sup>q</sup>	99.3
<i>E. coli</i> O157	FecalSwab	0/0 <sup>s</sup>	N/A	N/A	3 / 3 <sup>s</sup>	100.0
	Para-Pak C&S	0/0	N/A	N/A	5/5	100.0
<i>Shigella</i> /Enteroinvasive <i>E. coli</i> (EIEC)	FecalSwab	10/10	100.0	72.3–100.0	1212/1212	100.0
	Para-Pak C&S	2/2	100.0	34.2–100.0	703/704 <sup>t</sup>	99.9
<i>Cryptosporidium</i>	FecalSwab	2/4	50.0	15.0–85.0	1218/1218	100.0
	Para-Pak C&S	6/6	100.0	61.0–100.0	699/700 <sup>u</sup>	99.9

Table 20. Clinical Performance in the Prospective study (continued)

Analyte	Medium brand	Positive Percent Agreement		Negative Percent Agreement		
		TP / TP+FN	%	95% CI	TN / TN+FP	% 95% CI
<i>Cytophora cayetanensis</i>	FecalSwab	3/3	100.00	43.9–100.0	1219/1219	100.0 99.7–100.0
	Para-Pak C&S	18/19 <sup>a</sup>	94.7	75.4–99.1	687/687	100.0 99.4–100.0
<i>Entamoeba histolytica</i>	FecalSwab	0/0	N/A	N/A	1222/1222	100.0 99.7–100.0
	Para-Pak C&S	0/0	N/A	N/A	706/706	100.0 99.5–100.0
<i>Giardia lamblia</i>	FecalSwab	6/8 <sup>x</sup>	75.0	40.9–92.9	434/441 <sup>w</sup>	98.4 96.8–99.2
	Para-Pak C&S	1/1	100.0	20.7–100.0	406/406 <sup>y</sup>	100.0 99.1–100.0

<sup>a</sup> Adenovirus F40/41 was not detected in the single false negative specimen (0/1) in FecalSwab using a different FDA-cleared test method.

<sup>b</sup> Adenovirus F40/41 was not detected in the single false negative specimen (0/1) and in the single false positive specimen (0/1) in Para-Pak C&S using a different FDA-cleared test method.

<sup>c</sup> Ten (10) FecalSwab samples positive for Norovirus GI/GII in both QIAstat-Dx and one FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. The sample size for NPA is smaller for Norovirus GI/GII as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the prospective study.

<sup>d</sup> Two (2) Para-Pak C&S samples positive for Norovirus GI/GII in both QIAstat-Dx and one FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. One (1) Para-Pak C&S sample negative in QIAstat-Dx and positive with one FDA-cleared comparator with insufficient volume for complete composite comparator testing were classed as false negative in the PPA calculations. The sample size for NPA is smaller for Norovirus GI/GII as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the prospective study.

Table 20. Clinical Performance in the Prospective study (continued)

Analyte	Medium brand	Positive Percent Agreement		Negative Percent Agreement	
		TP/TP+FN	%	TN/TN+FP	%
95% CI					
95% CI					
e Rotavirus A was detected in one of the two false negative specimens (1/2) and was not detected in the two false positive specimens (0/2) in FecalSwab using a different FDA-cleared test method.					
f Rotavirus A was not detected in the single false positive specimen (0/1) in Para-Pak C&S using a different FDA-cleared test method.					
g <i>Campylobacter</i> was not detected in the two false negative specimens (0/2) and was detected in three of the four false positive specimens (3/4) in FecalSwab using a different FDA-cleared test method.					
h <i>Campylobacter</i> was not detected in the single false negative specimens (0/1) and was detected in one of the two false positive specimens (1/2) in Para-Pak C&S using a different FDA-cleared test method.					
i <i>Plesiomonas shigelloides</i> was not detected in the two false positive specimens (0/2) in FecalSwab using a different FDA-cleared test method.					
j <i>Plesiomonas shigelloides</i> was not detected in the single false negative specimen (0/1) and was not detected in the two false positive specimens (0/2) in Para-Pak C&S using a different FDA-cleared test method.					
k <i>Salmonella</i> was not detected in the two false negative specimens (0/2) in FecalSwab using a different FDA-cleared test method.					
l <i>Salmonella</i> was not detected in the single false negative specimen (0/1) in Para-Pak C&S using a different FDA-cleared test method.					
m <i>Yersinia enterocolitica</i> was not detected in the single false negative specimen (0/1) and was not detected in the seven false positive specimens (0/7) in FecalSwab using a different FDA-cleared test method.					
n <i>Yersinia enterocolitica</i> was not detected in the five false positive specimens (0/5) using a different FDA-cleared test method.					
o Six (6) FecalSwab samples positive for ETEC in both QIAstat-Dx and the primary FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. The sample size for NPA is smaller for ETEC as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the prospective study.					
p Three (3) Para-Pak C&S samples positive for ETEC in both QIAstat-Dx and one FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. The sample size for NPA is smaller for ETEC as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the prospective study.					

Table 20. Clinical Performance in the Prospective study (continued)

Analyte	Medium brand	Positive Percent Agreement		Negative Percent Agreement	
		TP/TP+FN	%	TN/TN+FP	%

- <sup>a</sup> Eight (8) FecalSwab sample positive for STEC in both QIAstat-Dx and one FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. The sample size for NPA is smaller for STEC as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the prospective study.
- <sup>r</sup> One (1) Para-Pak C&S sample positive for STEC in both QIAstat-Dx and one FDA-cleared comparator was excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. The sample size for NPA is smaller for STEC as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the prospective study.
- <sup>s</sup> Three (3) positive and nine (9) negative samples for E.coli O157 by QIAstat-Dx were excluded from the PPA/NPA calculations because reporting of the E.coli O157 result is dependent on the preceding STEC result (E. coli O157 is subtype within STEC) and the STEC result for all twelve (12) samples is negative, or not available or unconfirmed with the (composite) reference test.
- <sup>t</sup> *Shigella*/ EIEC was detected in the single false positive specimen (1/1) in Para-Pak C&S using a different FDA-cleared test method.
- <sup>u</sup> *Cryptosporidium* was not detected in the single false positive specimen (0/1) in Para-Pak C&S using PCR followed by bi-directional sequence analysis.
- <sup>v</sup> For *Cyclospora cayentensis*, there was one (1) false negative specimen in Para-Pak C&S that was not further investigated by discrepant analyses.
- <sup>w</sup> Two (2) FecalSwab samples positive for *Giardia lamblia* in both QIAstat-Dx and one FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. Two (2) FecalSwab samples negative in QIAstat-Dx and positive with one FDA-cleared comparator with insufficient volume for complete composite comparator testing were classed as false negative in the PPA calculations. The sample size for NPA is smaller for *Giardia lamblia* as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the prospective study.
- <sup>x</sup> The sample size for NPA is smaller for *Giardia lamblia* as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the prospective study.

Table 21. Clinical performance in the prospective archived study

Analyte	Medium brand	Positive Percent Agreement			Negative Percent Agreement		
		TP/TP+FN	%	95% CI	TN/TN+FP	%	95% CI
Norovirus GI/GII	FecalSwab	48/49*	98.0	89.3–99.6	2/4*	50.0	15.0–85.0
	Para-Pak C&S	26/28*†	92.9	77.4–98.0	37/38*	97.4	86.5–99.5
Shiga-like toxin <i>E. Coli</i> (STEC) <i>stx1/stx2</i>	FecalSwab	24 / 24 <sup>c</sup>	100.0	86.2–100.0	67 / 68 <sup>c</sup>	98.5	92.1–99.7
	Para-Pak C&S	12/13†§	92.3	66.7–98.6	51/52†	98.1	89.9–99.7

\* For Norovirus GI/GII, four out of the eighty-one (4/81) prospectively archived samples (positive by standard of care) were negative by the composite comparator and therefore included as negative samples in the NPA calculations.

† One (1) Para-Pak C&S sample negative in QIAstat-Dx and positive for Norovirus GI/GII with one FDA-cleared comparator with insufficient volume for complete composite comparator testing was classed as false negative in the PPA calculations.

<sup>c</sup> For STEC fifty-one out of the seventy-five (51/75) prospectively archived FecalSwab samples (positive by standard of care) were negative by the composite comparator and therefore included as negative samples in the NPA calculations.

† For STEC five out of the eighteen (5/18) prospectively archived samples (positive by standard of care) were negative by the composite comparator and therefore included as negative samples in the NPA calculations.

† One (1) Para-Pak C&S sample positive for STEC in both QIAstat-Dx and one FDA-cleared comparator was excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing.



Table 22. Clinical performance in the retrospective study

Analyte	Medium brand	Positive Percent Agreement		Negative Percent Agreement			
		TP/TP+FN	%	95% CI	TN/TN+FP	%	95% CI
Adenovirus F40/F41	Viruses						
	FecalSwab	23/26 <sup>a</sup>	88.5	71.0–96.0	203/203	100.0	98.1–100.0
	Para-Pak C&S	29/29	100.0	88.3–100.0	39/39	100.0	91.0–100.0
	FecalSwab	2/3 <sup>b</sup>	66.7	20.8–93.9	191/191	100.0	98.0–100.0
Astrovirus	Para-Pak C&S	0/0	N/A	N/A	14/14	100.0	78.5–100.0
	FecalSwab	28/32 <sup>c</sup>	87.5	71.9–95.0	74/75 <sup>c</sup>	98.7	92.8–99.8
Norovirus GI/GII	Para-Pak C&S	27/29	93.1	78.0–98.1	86/86 <sup>d</sup>	100.0	95.7–100.0
	FecalSwab	8/9 <sup>e</sup>	88.9	56.5–98.0	185/185	100.0	98.0–100.0
Rotavirus A	Para-Pak C&S	2/2	100.0	34.2–100.0	12/12	100.0	75.8–100.0
	Bacteria						
Campylobacter	FecalSwab	31/31	100.0	89.0–100.0	161/163 <sup>a</sup>	98.8	95.6–99.7
	Para-Pak C&S	3/3	100.0	43.9–100.0	11/11	100.0	74.1–100.0
Plesiomonas shigelloides	FecalSwab	2/2	100.0	34.2–100.0	192/192	100.0	98.0–100.0
	Para-Pak C&S	33/36 <sup>g</sup>	91.7	78.2–97.1	117/117	100.0	96.8–100.0
Salmonella	FecalSwab	30/31 <sup>h</sup>	96.8	83.8–99.4	161/163 <sup>h</sup>	98.8	95.6–99.7
	Para-Pak C&S	1/1	100.0	20.7–100.0	13/13	100.0	77.2–100.0

Table 22. Clinical performance in the retrospective study (continued)

Analyte	Medium brand	Positive Percent Agreement		Negative Percent Agreement			
		TP/TP+FN	%	95% CI	TN/TN+FP	%	95% CI
<i>Yersinia enterocolitica</i>	FecalSwab	32 / 34 <sup>c</sup>	94.1	80.9–98.4	160 / 160	100.0	97.7–100.0
	Para-Pak C&S	1 / 1	100.0	20.7–100.0	14 / 14	100.0	78.5–100.0
Diarrheagenic <i>E. coli</i> / <i>Shigella</i>							
Enteropathogenic <i>E. coli</i> (EPEC)	FecalSwab	46 / 48	95.8	86.0–98.9	164/164	100.0	97.7–100.0
	Para-Pak C&S	60/65 <sup>l</sup>	92.3 <sup>o</sup>	83.2–96.7	42/42	100.0	91.6–100.0
Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt</i> / <i>st</i>	FecalSwab	22/24 <sup>k</sup>	91.7	74.2–97.7	85/86 <sup>k</sup>	98.8	93.7–99.8
	Para-Pak C&S	23/24	95.8	79.8–99.3	61/61 <sup>l</sup>	100.0	94.1–100.0
Shiga-like toxin <i>E. coli</i> (STEC) <i>stx1</i> / <i>stx2</i>	FecalSwab	2 / 3 <sup>m</sup>	66.7	20.8–93.9	62/63 <sup>m</sup>	98.4	91.5–99.7
	Para-Pak C&S	60/64	93.8	85.0–97.5	44/44 <sup>m</sup>	100.0	92.0–100.0
<i>E. coli</i> O157	FecalSwab	0/0 <sup>o</sup>	N/A	N/A	2/2	100.0	34.2–100.0
	Para-Pak C&S	39/42 <sup>p</sup>	92.9%	80.1–99.4	16/16	100.0	80.6–100.0
<i>Shigella</i> / <i>Enteroinvasive E. coli</i> (EIEC)	FecalSwab	22/24 <sup>q</sup>	91.7	74.2–97.7	170/170	100.0	97.8–100.0
	Para-Pak C&S	0/0	N/A	N/A	14/14	100.0	78.5–100.0
Parasites							
<i>Cryptosporidium</i>	FecalSwab	6/6	100.0	61–100.0	186/188 <sup>r</sup>	98.9	96.2–99.7
	Para-Pak C&S	26/26	100.0	87.1–100.0	117/117	100.0	96.8–100.0

Table 22. Clinical performance in the retrospective study (continued)

Analyte	Medium brand	Positive Percent Agreement			Negative Percent Agreement		
		TP/TP+FN	%	95% CI	TN/TN+FP	%	95% CI
<i>Cyclospora cayetanensis</i>	FecalSwab	1/1	100.0	20.7–100.0	193/193	100.0	98.1–100.0
	Para-Pak C&S	1/1	100.0	20.7–100.0	13/13	100.0	77.2–100.0
<i>Entamoeba histolytica</i>	FecalSwab	0/0	N/A	N/A	194/194	100.0	98.1–100.0
	Para-Pak C&S	0/0	N/A	N/A	14/14	100.0	76.5–100.0
<i>Giardia lamblia</i>	FecalSwab	29/31 <sup>s</sup>	93.6	79.3–98.2	46/48 <sup>s</sup>	95.8	86.0–98.9
	Para-Pak C&S	27/28 <sup>t</sup>	96.4	82.3–99.4	1	100.0	96.0–100.0

<sup>a</sup> Adenovirus F40/41 was detected in one of the three false negatives (1/3) in FecalSwab using a different FDA-cleared test method.

<sup>b</sup> Astrovirus was detected in the single false negative specimen (1/1) in FecalSwab using a different FDA-cleared test method.

<sup>c</sup> Two (2) FecalSwab samples positive for Norovirus GI/GII in both QIAstat-Dx and one FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. The sample size for NPA is smaller for Norovirus GI/GII as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the retrospective study.

<sup>d</sup> The sample size for NPA is smaller for Norovirus GI/GII in Para-Pak C&S as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the retrospective study.

<sup>e</sup> Rotavirus A was detected in the single false negative specimen (1/1) in FecalSwab using a different FDA-cleared test method.

Table 22. Clinical performance in the retrospective study (continued)

Analyte	Medium brand	Positive Percent Agreement		Negative Percent Agreement	
		TP/TP+FN	%	TN/TN+FP	%

<sup>f</sup> <i>Campylobacter</i> was detected in one of the two false positive specimens (1/2) in FecalSwab using a different FDA-cleared test method.			95% CI		95% CI
<sup>g</sup> <i>Plesiomonas shigelloides</i> was detected in one of the three false negative specimens (1/3) in Para-Pak C&S using a different FDA-cleared test method.					
<sup>h</sup> <i>Salmonella</i> was not detected in the single false negative specimen (0/1) and was not detected in the two false positive specimens (0/2) in FecalSwab using a different FDA-cleared test method.					
<sup>i</sup> <i>Yersinia enterocolitica</i> was not detected in the two false negative specimens (0/2) in FecalSwab using a different FDA-cleared test method.					
<sup>j</sup> Enteropathogenic <i>E. coli</i> (EPEC) was detected in all three false negative specimens (3/3) in Para-Pak C&S using PCR followed by bi-directional sequence analysis. There were two (2) other false negative specimens that were not further investigated by discrepant analyses.					
<sup>k</sup> Ten (10) FecalSwab samples positive for ETEC in both QIAstat-Dx and one FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. One (1) FecalSwab sample negative in QIAstat-Dx and positive with one FDA-cleared comparator with insufficient volume for complete composite comparator testing was classed as false negative in the PPA calculations. The sample size for NPA is smaller for ETEC as only a portion of the samples with a negative result in QIAstat-Dx and one FDA-cleared comparator was tested with the complete composite comparator in the retrospective study.					
<sup>l</sup> The sample size for NPA is smaller for ETEC in Para-Pak CS&S as only a portion of the samples with a negative result in QIAstat-Dx and one FDA-cleared comparator was tested with the complete composite comparator in the retrospective study).					
<sup>m</sup> Fifteen (15) FecalSwab sample positive for STEC in both QIAstat-Dx and one FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. One (1) FecalSwab sample negative in QIAstat-Dx and positive with one (1) FDA-cleared comparator with insufficient volume for complete composite comparator testing was classed as false negative in the PPA calculations. The sample size for NPA is smaller for STEC in FecalSwab as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the retrospective study.					
<sup>n</sup> The sample size for NPA is smaller for STEC in Para-Pak C&S as only a portion of the samples with a negative result in QIAstat-Dx and in one FDA-cleared comparator was tested with the complete composite comparator in the retrospective study.					
<sup>o</sup> One (1) positive sample for <i>E. coli</i> O157 by QIAstat-Dx was excluded from the PPA calculation because reporting of the <i>E. coli</i> O157 result is dependent on the preceding STEC result [ <i>E. coli</i> O157 is subtype within STEC] and the STEC result for that sample is unconfirmed.					

Table 22. Clinical performance in the retrospective study (continued)

Analyte	Positive Percent Agreement		Negative Percent Agreement	
	Medium brand	TP/TP+FN %	TN/TN+FP %	95% CI
<sup>p</sup> <i>E. coli</i> O157 was not detected in two false negative specimens (0/2) in Para-Pak C&S using a different FDA-cleared test method. There was one (1) false negative specimen in Para-Pak C&S that was not further investigated by discrepant analyses.				
<sup>q</sup> <i>Shigella</i> /Enteroinvasive <i>E. coli</i> (EIEC) was detected in one of the two false negative specimens (1/2) in FecalSwab using a different FDA-cleared test method.				
<sup>r</sup> <i>Cryptosporidium</i> was not detected in the two false positive specimens (0/2) in FecalSwab using PCR followed by bi-directional sequence analysis.				
<sup>s</sup> Four (4) samples positive for <i>Giardia lamblia</i> in both QIAstatDx and one FDA-cleared comparator were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. Two (2) FecalSwab samples negative in QIAstatDx and positive with one FDA-cleared comparator with insufficient volume for complete composite comparator testing was classed as false negative in the PPA calculations. The sample size for NPA is smaller for <i>Giardia lamblia</i> as only a portion of the samples with a negative result in QIAstatDx and in one FDA-cleared comparator was tested with the complete composite comparator in the retrospective study.				

<sup>t</sup> One (1) Para-Pak C&S samples positive for *Giardia lamblia* in both QIAstatDx and primary FDA-cleared comparator (were excluded from the PPA calculations because the samples did not have sufficient volume for complete composite comparator testing. One (1) Para-Pak C&S sample negative in QIAstatDx and positive with one FDA-cleared comparator with insufficient volume for complete composite comparator testing was classed as false negative in the PPA calculations. The sample size for NPA is smaller for *Giardia lamblia* as only a portion of the samples with a negative result in QIAstatDx and in one FDA-cleared comparator was tested with the complete composite comparator in the retrospective study.

The proportion of failed runs on initial attempt, and following repeats are summarized in Table 23. The error breakdown due to instrument, invalid results, “sample too concentrated” failures, and other run failures are summarized in Table 24.

Table 23. Failure rates summary

Transport media	Study	Initial runs			Final runs		
		N/Total	%	95% CI	N/Total	%	95% CI
FecalSwab	Prospective	16/1227	1.3	0.8 – 2.1	3/1227	0.2	0.1 – 0.7
	Prospective Archived	0/145	0.0	0.0 – 2.6	0/145	0.0	0.0 – 2.6
	Retrospective	11/366	3.0	1.5 – 5.3	5/366	1.4	0.4 – 3.2
	Total	27/1738	1.6	1.0 – 2.3	8/1738	0.5	0.2 – 0.9
Para-Pak C&S	Prospective	66/740	8.9	7.0 – 11.2	21/740	2.8	1.8 – 4.3
	Prospective Archived	3/66	4.5	0.9 – 12.7	0/66	0.0	0.0 – 5.4
	Retrospective	46/454	10.1	7.5 – 13.3	25/454	5.5	3.6 – 8.0
	Total	115/1260	9.1	7.6 – 10.9	46/1260	3.7	2.7 – 4.8

Table 24. Failure types breakdown

Transport media	Study	Failure reason	Initial runs		Final runs	
			N/Total	%	N/Total	%
FecalSwab	Prospective	Instrument	0/1227	0.0	0/1227	0.0
		Invalid*	0/1227	0.0	0/1227	0.0
		Sample too Concentrated†	5/1227	0.4	0/1227	0.0
		Other‡	11/1227	0.9	3/1227	0.2
	Prospective Archived	Instrument	0/145	0.0	0/145	0.0
		Invalid	0/145	0.0	0/145	0.0
		Sample too Concentrated	0/145	0.0	0/145	0.0
		Other	0/145	0.0	0/145	0.0
	Retrospective	Instrument	1/366	0.3	0/366	0.0
		Invalid	1/366	0.3	0/366	0.0
		Sample too Concentrated	0/366	0.0	0/366	0.0
		Other	9/366	2.5	5/366	1.4

Table 24. Failure types breakdown (continued)

Transport media	Study	Failure reason	Initial runs		Final runs	
			N/Total	%	N/Total	%
Para-Pak C&S	Prospective	Instrument	9/740	1.2	2/740	0.3
		Invalid	5/740	0.7	5/740	0.7
		Sample too Concentrated	35/740	4.7	7/740	0.9
		Other	17/740	2.3	7/740	0.9
	Prospective Archived	Instrument	0/66	0.0	0/66	0.0
		Invalid	1/66	1.5	0/66	0.0
		Sample too Concentrated	1/66	1.5	0/66	0.0
		Other	1/66	1.5	0/66	0.0
	Retrospective	Instrument	1/454	0.2	0/454	0.0
		Invalid	10/454	2.2	6/454	1.3
		Sample too Concentrated	10/454	2.2	2/454	0.4
		Other	25/454	5.5	17/454	3.7

\* Internal Control failures with at least one analyte detected and the other analytes reported as "invalid".

† Run failures related to "sample concentration too high". These specimens were repeated with 100 microliters as detailed in "Appendix C: Additional instructions for use" on page 1.

‡ Run failures related to workflow checkpoints.

Co-infections

The QIAstat-Dx Gastrointestinal Panel 2 reported multiple organism detections (i.e., co-infections) for a total of 58 and 29 prospective specimens in FecalSwab and Para-Pak C&S, respectively. This represents 18.6% of positive specimens (58/312) in FecalSwab and 17.0% of positive specimens (29/171) in Para-Pak C&S. Most multiple detections in FecalSwab specimens (51/58; 87.9%) contained two organisms, while 8.6% (5/58) contained three



organisms and 3.4% (2/58) contained four organisms. In Para-Pak C&S specimens, most multiple detections (22/29; 75.9%) contained two organisms, while 24.1% (7/29) contained three organisms. The most common multiple infections are shown in Table 25 and Table 26.

**Table 25. Most prevalent multiple detection combinations (≥2 instances) as determined by the QIAstat-Dx Gastrointestinal Panel 2 in the prospective clinical study in FecalSwab specimens**

Multiple detection combination	Number of FecalSwab specimens
Adenovirus F40/F41 + Enteropathogenic <i>E. coli</i> (EPEC)	2
<i>Campylobacter</i> + Enteropathogenic <i>E. coli</i> (EPEC) + Rotavirus A	2
<i>Campylobacter</i> + Rotavirus A	2
<i>E. coli</i> O157 + Shiga-like toxin <i>E. coli</i> (STEC) <i>stx1/stx2</i>	2
Enteropathogenic <i>E. coli</i> (EPEC) + Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i> + Norovirus GI/GII	2
Enteropathogenic <i>E. coli</i> (EPEC) + <i>Giardia lamblia</i>	2
Enteropathogenic <i>E. coli</i> (EPEC) + Rotavirus A	2
Enteropathogenic <i>E. coli</i> (EPEC) + <i>Yersinia enterocolitica</i>	2
Enteropathogenic <i>E. coli</i> (EPEC) + Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i>	4
Enteropathogenic <i>E. coli</i> (EPEC) + Norovirus GI/GII	12
<i>Campylobacter</i> + Enteropathogenic <i>E. coli</i> (EPEC)	13

**Table 26. Most prevalent multiple detection combinations (≥2 instances) as determined by the QIAstat-Dx Gastrointestinal Panel 2 in the prospective clinical study in Para-Pak C&S specimens**

Multiple detection combination	Number of Para-Pak C&S specimens
<i>Campylobacter</i> + Enteropathogenic <i>E. coli</i> (EPEC)	3
Enteropathogenic <i>E. coli</i> (EPEC) + <i>Salmonella</i>	3
Enteropathogenic <i>E. coli</i> (EPEC) + Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i>	4

The analytes most commonly found in mixed infections in the FecalSwab specimens were EPEC (45), *Campylobacter* (22), Norovirus GI/GII (20), Rotavirus (6) and ETEC (7) as shown in Table 27, while the analytes most commonly found in mixed infections in the Para-Pak C&S specimens were EPEC (17), ETEC (8), *Campylobacter* (7), Norovirus GI/GII (5), Rotavirus (4) and STEC (5) as shown in Table 28.

**Table 27. Prevalence of analytes in mixed infections in FecalSwab specimens as determined by the QIAstat-Dx Gastrointestinal Panel 2**

Analyte	N	%
Adenovirus F40/F41	3	2.4
Astrovirus	2	1.6
<i>Campylobacter</i>	22	17.6
<i>E. coli</i> O157	3	2.4
Enteropathogenic <i>E. coli</i> (EPEC)	45	36.0
Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i>	7	5.6
<i>Giardia lamblia</i>	4	3.2
Norovirus GI/GII	20	16.0
<i>Plesiomonas shigelloides</i>	1	0.8
Rotavirus A	6	4.8
<i>Salmonella</i>	1	0.8
Shiga-like toxin <i>E. coli</i> (STEC) <i>stx1/stx2</i>	5	4.0
<i>Shigella</i> /Enteroinvasive <i>E. coli</i> (EIEC)	3	2.4
<i>Yersinia enterocolitica</i>	3	2.4

**Table 28. Prevalence of analytes in mixed infections in Para-Pak C&S specimens as determined by the QIAstat-Dx Gastrointestinal Panel 2**

Analyte	N	%
Adenovirus F40/F41	1	1.5
Astrovirus	1	1.5
<i>Campylobacter</i>	7	10.8
<i>Cryptosporidium</i>	2	3.1
<i>Cyclospora cayetanensis</i>	2	3.1
Enteropathogenic <i>E. coli</i> (EPEC)	17	26.2
Enterotoxigenic <i>E. coli</i> (ETEC) <i>lt/st</i>	8	12.3
<i>Giardia lamblia</i>	1	1.5
Norovirus GI/GII	5	7.7
<i>Plesiomonas shigelloides</i>	7	10.8
Rotavirus A	1	1.5
<i>Salmonella</i>	4	6.2
Shiga-like toxin <i>E. coli</i> (STEC) <i>stx1/stx2</i>	5	7.7
<i>Shigella</i> /Enteroinvasive <i>E. coli</i> (EIEC)	3	4.6
<i>Yersinia enterocolitica</i>	1	1.5

## Contrived specimens testing

Several analytes, such as *Entamoeba histolytica* are so rare that both prospective and retrospective testing efforts were insufficient to demonstrate system performance. To supplement the prospective and retrospective specimens' test results, an evaluation of contrived specimens was performed. Contrived specimens were prepared using negative residual specimens that had previously tested negative by QIAstat-Dx Gastrointestinal Panel 2 and comparator methods. At least 50% of these specimens were spiked at concentrations

slightly above the Limit of Detection (2x LoD) and the rest at 5x and 10x LoD, using quantified strains for each pathogen. A minimum of 50 contrived specimens were tested for each evaluated analyte. The analyte status of each contrived specimen was blinded to the users analyzing the specimens. Results are summarized in Table 29.

Table 29. Test results summary for contrived specimens.

QIAstat-Dx GI2 Target	Positive Percent Agreement (PPA)			
	Medium brand	Fraction	Percentage (%)	95% CI (%)
Astrovirus	FecalSwab	33/34	97.1	85.1–99.5
	Para-Pak C&S	34/34	100.0	89.8–100.0
Rotavirus A	FecalSwab	35/35	100.0	90.1–100.0
	Para-Pak C&S	34/35	97.1	85.5–99.5
<i>Plesiomonas shigelloides</i>	FecalSwab	33/33	100.0	89.6–100.0
	Para-Pak C&S	34/35	97.1	85.5–99.5
<i>Yersinia enterocolitica</i>	FecalSwab	34/34	100.0	89.8–100.0
	Para-Pak C&S	34/35	97.1	85.5–99.5
<i>E. coli</i> O157	FecalSwab	35 / 35	100.0	90.1–100.0
	Para-Pak C&S	34 / 34	100.0	89.9–100.0
Shigella/EIEC	FecalSwab	35/35	100.0	90.1–100.0
	Para-Pak C&S	34/34	100.0	89.8–100.0
<i>Cryptosporidium</i>	FecalSwab	27/27	100.0	87.5–100.0
	Para-Pak C&S	31/31	100.0	89.0–100.0
<i>Cyclospora cayetanensis</i>	FecalSwab	26/26	100.0	87.1–100.0
	Para-Pak C&S	30/30	100.0	88.6–100.0
<i>Entamoeba histolytica</i>	FecalSwab	35/35	100.0	90.1–100.0
	Para-Pak C&S	34/35	97.1	85.5–99.5

# Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. For more information, see also the Frequently Asked Questions page at our Technical Support Center: [www.qiagen.com/FAQ/FAQList.aspx](http://www.qiagen.com/FAQ/FAQList.aspx) (for contact information, visit [www.qiagen.com](http://www.qiagen.com)).

Additional information about specific QIAstat-Dx Gastrointestinal Panel 2 error codes and messages can be found in Table 30:

**Table 30. Information about specific QIAstat-Dx Gastrointestinal Panel 2 error codes and messages**

Error Code	Error message displayed
0x02C9	Cartridge execution failure: Sample concentration too high.
0x032D	Please repeat by loading 100 microliters of the sample in a new cartridge (as per Appendix C
0x0459	explanation)
0x045A	
0x04BF	
0x0524	
0x058B	
0x05E9	
0x0778	
0x077D	
0x14023	











When the sample concentration is too high and the test must be repeated by loading 100 µL, follow the workflow detailed in "Appendix C: Additional instructions for use" on page 136 of this document.










# Contact Information

For technical assistance and more information, please see our Technical Support Center at [www.qiagen.com/Support](http://www.qiagen.com/Support), call 800-426-8157, or contact one of the QIAGEN Technical Service Departments or local distributors (see back cover or visit [www.qiagen.com](http://www.qiagen.com)).

# Symbols

The following symbols appear in the instructions for use or on the packaging and labeling:

Symbol	Symbol definition
 <N>	Contains reagents sufficient for <N> reactions
 Rx Only	Prescription Use only
	Use by
	For in vitro diagnostic use
	Catalog number
	Lot number
	Material number (i.e., component labeling)
	Unique Device Identifier
	Contains
	Component

Symbol	Symbol definition
	Number
	Gastrointestinal application
<b>Rn</b>	R is for revision of the Instructions for Use and n is the revision number
	Temperature limitation
	Manufacturer
	Consult instructions for use
	Protect from light
	Do not reuse
	Caution, consult accompanying documents
	Serial number



## Symbol

## Symbol definition



Do not use if package is damaged



Flammable, risk of fire



Corrosive, risk of chemical burn



Health Hazard, risk of sensitization, carcinogenicity



Risk of harm

# Appendices

## Appendix A: Installing the Assay Definition File

The Assay Definition File of the QIAstat-Dx Gastrointestinal Panel 2 must be installed on the QIAstat-Dx Analyzer 2.0 prior to testing with QIAstat-Dx Gastrointestinal Panel 2 Cartridges.

**Note:** Whenever a new version of the QIAstat-Dx Gastrointestinal Panel 2 assay is released, the new QIAstat-Dx Gastrointestinal Panel 2 Assay Definition File must be installed prior to testing.


The Assay Definition File (.asy file type) is available at [www.qiagen.com](http://www.qiagen.com).

The Assay Definition File (.asy file type) must be saved onto a USB drive prior to installation on the QIAstat-Dx Analyzer 2.0. This USB drive must be formatted with a FAT32 file system.

To import an ADF from the USB to the QIAstat-Dx Analyzer 2.0, proceed with the following steps:


1. Insert the USB drive containing the Assay Definition File into one of the USB ports on the QIAstat-Dx Analyzer 2.0.
2. Press **Options > Assay Management**.

The Assay Management screen appears in the Content area of the display (Figure 22).




**Figure 22. Assay Management screen.**

3. Press the **Import** icon in the bottom left of the screen.
4. Select the file corresponding to the assay to be imported from the USB drive.  
A dialog box will appear to confirm the upload of the file.
5. A dialog box may appear to overwrite the current version by a new one. Press **Yes** to overwrite (Figure 23).



**Figure 23.** Dialog box that appears when upgrading the ADF version.

The assay becomes active by selecting **Assay Active** (Figure 24).



**Figure 24. Activating the assay.**

6. To assign the active assay to a user, perform the following steps:

- a. Go to **Options > User Management**.
- b. Select the user who should be allowed to run the assay.

**Note:** If needed, this step can be repeated for every user created in the system.

- c. Select **Assign Assays** from the User Options tab.
- d. Enable the assay and press **Save** (Figure 25).

technician

User Management

08:24 2025-10-16

1 Available

2 Not installed

3 Not installed

4 Not installed

USER

administrator

ADMINISTRATOR

labuser

LABORATORY SUPERVISOR

technician

SERVICE TECHNICIAN

USER OPTIONS

User Name

labuser

Password

\*\*\*\*\*

User Active

☒

Assign user profile

>

Assign Assays

>

Assay Statistics

>

ASSAYS

GI2\_691421

☒

Run Test

View Results

Options

Log Out

Add User

Save

Cancel

**Figure 25. Assigning the active assay.**

## Appendix B: Glossary

**Amplification curve:** Graphical representation of the multiplex real-time RT-PCR amplification data.

**Analytical Module (AM):** The main QIAstat-Dx Analyzer 2.0 hardware module, in charge of executing tests on QIAstat-Dx Gastrointestinal Panel 2 Cartridges. It is controlled by the Operational Module. Several Analytical Modules can be connected to one Operational Module.

**IFU:** Instructions For Use.

**Main port:** In the QIAstat-Dx Gastrointestinal Panel 2 Cartridge, inlet for transport medium liquid samples.

**Nucleic acids:** Biopolymers, or small biomolecules composed of nucleotides, which are monomers made of three components: a 5-carbon sugar, a phosphate group and a nitrogenous base.

**Operational Module (OM):** The dedicated QIAstat-Dx Analyzer 1.0 hardware that provides the user interface for 1–4 Analytical Modules (AMs).

**Operational Module PRO (OM PRO):** The dedicated QIAstat-Dx Analyzer 2.0 hardware that provides the user interface for 1–4 Analytical Modules (AMs).

**PCR:** Polymerase Chain Reaction.

**QIAstat-Dx Analyzer 1.0:** The QIAstat-Dx Analyzer 1.0 consists of an Operational Module and an Analytical Module. The Operational Module includes elements that provide connectivity to the Analytical Module and enables user interaction with the QIAstat-Dx

Analyzer 1.0. The Analytical Module contains the hardware and software for sample testing and analysis.

**QIAstat-Dx Analyzer 2.0:** The QIAstat-Dx Analyzer 2.0 consists of an Operational Module PRO and Analytical Module. The Operational Module PRO includes elements that provide connectivity to the Analytical Module and enables user interaction with the QIAstat-Dx Analyzer 2.0. The Analytical Module contains the hardware and software for sample testing and analysis.

**QIAstat-Dx Gastrointestinal Panel 2 Cartridge:** A self-contained disposable plastic device with all pre-loaded reagents required for the complete execution of fully automated molecular assays for the detection of gastrointestinal pathogens.

**RT:** Reverse Transcription.

**Swab port:** In the QIAstat-Dx Gastrointestinal Panel 2 Cartridge, inlet for dry swabs. The swab port is not used for the QIAstat-Dx .

**User:** A person who operates the QIAstat-Dx Analyzer 2.0 / QIAstat-Dx Gastrointestinal Panel 2 Cartridge in the intended way.

## Appendix C: Additional instructions for use

In case cartridge execution failures corresponding to error codes (0x02C9, 0x032D, 0x0459, 0x045A, 0x04BF, 0x0524, 0x058B, 0x05E9, 0x0778, 0x077D, 0x14023) occur during the testing, the following error message will be displayed in the QIAstat-Dx Analyzer 2.0 screen after the run has finalized.

"Cartridge execution failure: Sample concentration too high. Please repeat by loading 100 microliters of the sample in a new cartridge (as per IFU explanation)."



In this case, the test should be repeated using 100 µL of the same sample following equivalent testing procedures detailed in the 'Procedure' section of the IFU adapted to 100 µL sample input volume:

1. Open the package of a new QIAstat-Dx Gastrointestinal Panel 2 Cartridge using the tear notches on the sides of the packaging.
2. Remove the QIAstat-Dx Gastrointestinal Panel 2 Cartridge from the packaging.
3. Manually write the sample information, or place a sample information label, on the top of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge. Make sure that the label is properly positioned and does not block the lid opening.
4. Place the QIAstat-Dx Gastrointestinal Panel 2 Cartridge flat on the clean work surface so that the bar code on the label faces upwards. Open the sample lid of the main port on the front of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge.
5. Thoroughly mix the stool in the transport medium, for example, by vigorously agitating the tube 3 times.
6. Open the tube with the sample to be tested. Use the supplied transfer pipette to draw up fluid. Draw the sample to the first fill line on the pipette (i.e., 100 µL).

**Important:** Do not draw air, mucus, or particles into the pipette. If air, mucus, or particles are drawn into the pipette, carefully expel the sample fluid in the pipette back into the sample tube and draw up fluid again.

7. Carefully transfer the sample into the main port of the QIAstat-Dx Gastrointestinal Panel 2 Cartridge using the supplied single-use transfer pipette.
8. Firmly close the lid of the main port until it clicks.
9. From this point, proceed following the instructions described in the IFU.

# Ordering Information

Product	Contents	Cat. no.
QIAstat-Dx Gastrointestinal Panel 2	For 6 tests: 6 individually packaged QIAstat-Dx Gastrointestinal Panel 2 Cartridges and 6 individually packaged transfer pipettes	691421
<b>Related Products</b>		
QIAstat-Dx Analyzer 2.0	1 QIAstat-Dx Analytical Module, 1 QIAstat-Dx Operational Module PRO, and related hardware and software to run molecular diagnostic QIAstat-Dx assay cartridges	9002828

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit Instructions for Use or User Manual. QIAGEN kit Instructions for Use and User Manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

# References

1. Hata DJ, Powell EA, Starolis MW. Utility and Recommendations for the Use of Multiplex Molecular Gastrointestinal Pathogen Panels. *J Appl Lab Med*. 2023 Nov 2; 8(6): p. 1148-1159.
2. Beal SG, Tremblay EE, Toffel S, Velez L, Rand KH. A Gastrointestinal PCR Panel Improves Clinical Management and Lowers Health Care Costs. *J Clin Microbiol*. 2017 Dec 26; 56(1): p. e01457-17.
3. Cybulski Jr RJ, Bateman AC, Bourassa L, Bryan A, Beail B, Matsumoto J, et al. Clinical Impact of a Multiplex Gastrointestinal Polymerase Chain Reaction Panel in Patients With Acute Gastroenteritis. *Clin Infect Dis*. 2018 Nov 13; 67(11): p. 1688-1696.
4. Maldonado-Garza HJ, Garza-González E, Bocanegra-Ibarias P, Flores-Treviño S. Diagnostic syndromic multiplex approaches for gastrointestinal infections. *Expert Rev Gastroenterol Hepatol*. 2021 Jul; 15(7): p. 743-757.
5. Castany-Feixas M, Simo, Garcia-Garcia S, Fernandez de Sevilla M, Launes C, Kalkgruber M, et al. Rapid molecular syndromic testing for aetiological diagnosis of gastrointestinal infections and targeted antimicrobial prescription: experience from a reference paediatric hospital in Spain. *Eur J Clin Microbiol Infect Dis*. 2021 May 8; 40(10): p. 2153-2160.
6. Costa D, Iraola G. Pathogenomics of Emerging *Campylobacter* Species. *Clin Microbiol Rev*. 2019 Jul 3; 32(4): p. e00072-18.
7. Mishu Allos B, Iovine NM, Blaser MJ. *Campylobacter jejuni* and Related Species. In Bennett JE, Dolin R, Blaser MJ, editors. *Principles and Practice of Infectious Diseases*. 8th ed. Philadelphia: ELSEIVER SAUNDERS; 2015. p. 2485-2493.e4.
8. Humphries RM, Linscott AJ. Practical Guidance for Clinical Microbiology Laboratories: Diagnosis of Bacterial Gastroenteritis. *Clin Microbiol Rev*. 2015 Jan 1; 28(1): p. 3-31.

9. Fitzgerald, Nachamkin. *Manual of Clinical Microbiology*. 11th ed. Jorgensen JH, Carroll KC, Funke G, Pfaller MA, Landry ML, Richter SS, et al., editors.: Wiley; 2015.
10. CDC. Centers for Disease Control and Prevention. [Online].; 2021 [cited 2024 Jan 10. Available from: <https://www.cdc.gov/campylobacter/index.html>.
11. WHO. World Health Organization. [Online].; 2020 [cited 2024 Jan 10. Available from: <https://www.who.int/news-room/fact-sheets/xdetail/campylobacter>.
12. Janssen R, Krogfelt KA, Cawthraw SA, van Pelt W, Wagenaar JA, Owen RJ. Host-Pathogen Interactions in *Campylobacter* Infections: the. *Clin Microb Rev*. 2008 Jul 1; 21 (3): p. 505–518.
13. Janda JM, Abbott SL. The Changing Face of the Family Enterobacteriaceae (Order: "Enterobacterales"): New Members, Taxonomic Issues, Geographic Expansion, and New Diseases and Disease Syndromes. *Clin Microbiol Rev*. 2021 Feb 24; 34(2).
14. Kurtz JR, Goggins JA, McLachlan JB. *Salmonella* infection: interplay between the bacteria and host immune system. *Immunol Lett*. 2017 Oct 1; 190(n/a): p. 42-50.
15. Gal-Mor, C Boyle , A Grassl. Same species, different diseases: how and why typhoidal and non-typhoidal *Salmonella enterica* serovars differ. *Front Microbiol*. 2014 Aug 04; 5 (391).
16. Collaborators G2DD. Estimates of the global, regional, and national morbidity, mortality, and aetiologies of diarrhoea in 195 countries: a systematic analysis for the Global Burden of Disease Study 2016. *Lancet Infect Dis*. 2018 Nov; 18(11): p. 1211-1228.
17. Achtman M, Wain J, Weill FX, Nair S, Zhou Z, Sangal V, et al. Multilocus Sequence Typing as a Replacement for. *PLOS Path*. 2012 Jun 21; 8(6).
18. Majowicz SE, Musto , Scallan , Angulo FJ, Kirk , O'Brien SJ, et al. The Global Burden of Nontyphoidal *Salmonella*. *Clin Infect Dis*. 2010 Mar 15; 50(6): p. 882-889.

19. Bottone EJ. *Yersinia enterocolitica*: the charisma continues. Clin Microbiol Rev. 1997 Apr 1; 10(2).
20. Garzetti D, Susen R, Fruth , Tietze , Heesemann , Rakin. A molecular scheme for *Yersinia enterocolitica* patho-serotyping derived from genome-wide analysis. Int J Med Microbiol. 2014 May; 304(3-4): p. 275-283.
21. Ostroff S. *Yersinia* as an emerging infection: epidemiologic aspects of Yersiniosis. Contrib Microbiol Immunol. 1995; 13: p. 5-10.
22. CDC. Centers for Disease Control and Prevention. [Online].; 2016 [cited 2024 Jan 14]. Available from: <https://www.cdc.gov/yersinia/index.html>.
23. Ong KL, Gould LH, Chen DL, Jones TF, Scheftel , Webb TH, et al. Changing Epidemiology of *Yersinia enterocolitica* Infections: Markedly Decreased Rates in Young Black Children, Foodborne Diseases Active Surveillance Network (FoodNet), 1996–2009. Clin Infect Dis. 2012 Jun 1; 54(5): p. 385-390.
24. Fredriksson-Ahomaa, Cernela , Hächler , Stephan. *Yersinia enterocolitica* strains associated with human infections in Switzerland 2001-2010. Eur J Clin Microbiol Infect Dis. 2012 Jul; 31(7): p. 1543-1550.
25. Rosner BM, Werber , Höhle , Stark K. Clinical aspects and self-reported symptoms of sequelae of *Yersinia enterocolitica* infections in a population-based study, Germany 2009-2010. BMC Infect Dis. 2013 May 23; 13(236).
26. Kaper JB, Nataro JP, Mobley HL. Pathogenic *Escherichia coli*. Nat Rev Microbiol. 2004 Feb; 2(2): p. 123-140.
27. Croxen MA, Law RJ, Scholz R, Keeney KM, Wlodarska M, Finlay BB. Recent advances in understanding enteric pathogenic *Escherichia coli*. Clin Microbiol Rev. 2013 Oct; 26(4): p. 822-80.

28. Pasqua M, Michelacci V, Di Martino ML, Tozzoli R, Grossi M, Colonna B, et al. The Intriguing Evolutionary Journey of Enteroinvasive E. coli (EIEC) toward Pathogenicity. *Front. Microbiol.* 2017 Dec 5; 8(2390).
29. Clements, Young JC, Const. Infection strategies of enteric pathogenic *Escherichia coli*. *Gut Microbes.* 2012 Mar; 3(2): p. 71-87.
30. Devanga Ragupathi NK, Muthuirulandi Sethuvel DP, Inbanathan FY, Veeraraghavan. Accurate differentiation of *Escherichia coli* and *Shigella* serogroups: challenges and strategies. *New Microbes New Infect.* 2017 Sep 23; 21(n/a): p. 58-62.
31. Pakbin, Brück WM, Rossen JWA. Virulence Factors of Enteric Pathogenic *Escherichia coli*: A Review. *Int. J. Mol. Sci.* 2021 Aug; 22(18).
32. Khalil IA, Troeger , Blacker BF, Rao PC, Brown A, Atherly DE, et al. Morbidity and mortality due to *Shigella* and enterotoxigenic *Escherichia coli* diarrhoea: the Global Burden of Disease Study 1990-2016. *Lancet Infect Dis.* 2018 Nov; 18(11): p. 1229-1240.
33. CDC. Centers for Disease Control and Prevention. [Online].; 2023 [cited 2024 Jan 11. Available from: <https://www.cdc.gov/Shigella/index.html>.
34. Ochoa TJ, Contreras CA. Enteropathogenic *Escherichia coli* infection in children. *Curr Opin Infect Dis.* 2011 Oct; 24(5).
35. Nataro JP, Kaper JB. Diarrheagenic *Escherichia coli*. *Clin Microbiol Rev.* 1998 Jan; 11(1): p. 142-201.
36. Tobias, Kassem E, Rubinstein, Bialik, Vutukuru SR, Navaro , et al. Involvement of main diarrheagenic *Escherichia coli*, with emphasis on enteroaggregative E. coli, in severe non-epidemic pediatric diarrhea in a high-income country. *BMC Infectious Diseases.* 2015 Feb 21; 15(79).

37. Hu J, Torres AG. Enteropathogenic *Escherichia coli*: foe or innocent bystander? Clin Microbiol Infect. 2015 Jan 28; 21(8): p. 729-734.
38. Qadri, Svennerholm AM, Faruque ASG, Sack RB. Enterotoxigenic *Escherichia coli* in developing countries: epidemiology, microbiology, clinical features, treatment, and prevention. Clin Microbiol Rev. 2005 Jul; 18(3): p. 465-483.
39. CDC. Centers for Disease Control and Prevention. [Online].; 2014 [cited 2024 Jan 11]. Available from: <https://www.cdc.gov/ecoli/etec.html>.
40. Majowicz SE, Scallan , Jones-Bitton , Sargeant JM, Stapleton J, Angulo FJ, et al. Global incidence of human Shiga toxin-producing *Escherichia coli* infections and deaths: a systematic review and knowledge synthesis. Foodborne Pathog Dis. 2014 Apr 21; 11(6): p. 447-455.
41. CDC. Centers for Disease Control and Prevention. [Online].; 2014 [cited 15 Jan 2024]. Available from: <https://www.cdc.gov/ecoli/general/index.html>.
42. Leitch GJ, He. Cryptosporidiosis-an overview. J Biomed Res. 2012 Jan; 25(1): p. 1-16.
43. CDC. Centers for Disease Control and Prevention. [Online].; 2019 [cited 2024 Jan 10]. Available from: <https://www.cdc.gov/parasites/crypto/>.
44. Collaborators GDD. Estimates of global, regional, and national morbidity, mortality, and aetiologies of diarrhoeal diseases: a systematic analysis for the Global Burden of Disease Study 2015. Lancet Infect Dis. 2017 Jun 1; 17(9): p. 909-948.
45. Almeria, Cinar HN, Dubey JP. Cyclospora cayetanensis and Cyclosporiasis: An Update. Microorganisms. 2019 Sep 4; 7(9).
46. Ortega YR, Sanchez. Update on Cyclospora cayetanensis, a Food-Borne and Waterborne Parasite. Clin Microbiol Rev. 2010 Jan; 23(1): p. 218-234.
47. Chacin-Bonilla. CYCLOSPORA CAYETANENSIS. In PROJECT GWP. SPECIFIC EXCRETED PATHOGENS: ENVIRONMENTAL AND. Maracaibo; 2017. p. 36.

48. CDC. Centers for Disease Control and Prevention. [Online].; 2023 [cited 2024 Jan 10]. Available from: <https://www.cdc.gov/parasites/cyclosporiasis/index.html>.
49. Shirley DAT, Farr , Watanabe , Moonah. A Review of the Global Burden, New Diagnostics, and Current Therapeutics for Amebiasis. Open Forum Infect Dis. 2018 Jul 5; 5(7).
50. Kantor, Abrantes, Estevez, Schiller, Torrent , Gascon , et al. Entamoeba Histolytica: Updates in Clinical Manifestation, Pathogenesis, and Vaccine Development. Can J Gastroenterol Hepatol. 2018 Dec 2; 2018(4601420).
51. Ben Ayed, Sabbahi. ENTAMOEBA HISTOLYTICA. In PROJECT GWP. SPECIFIC EXCRETED PATHOGENS: ENVIRONMENTAL AND. Michigan; 2017. p. 35.
52. Rumsey, Waseem. Giardia lamblia Enteritis. StatPearls. 2023 Jul 4.
53. CDC. Centers for Disease for Control and Prevention. [Online].; 2022 [cited 2024 Jan 11]. Available from: <https://www.cdc.gov/parasites/giardia/index.html>.
54. Boarato-David, Guimarães , Cacciò. Global Water Pathogen Project. [Online].; 207 [cited 2024 Jan 11]. Available from: <https://www.waterpathogens.org/book/giardia-duodenalis>.
55. Hooshyar, Rostamkhani , Arbabi , Delavari. Giardia lamblia infection: review of current diagnostic strategies. Gastroenterol Hepatol Bed Bench. 2019 Winter; 12(1): p. 3-12.
56. Leung AKC, Leung AAM, Wong AHC, Hon KL. Human Ascariasis: An Updated Review. Recent Pat Inflamm Allergy Drug Discov. 2020; 14(2): p. 133-145.
57. Ison MG, Hayden RT. Adenovirus. Microbiol Spectr. 2016 Aug; 4(4).
58. Lynch 3rd JP, Kajon AE. Adenovirus: Epidemiology, Global Spread of Novel Serotypes, and Advances in Treatment and Prevention. Semin Respir Crit Care Med. 2016 Aug; 37(4): p. 586-602.



59. Lion T. Adenovirus infections in immunocompetent and immunocompromised patients. Clin Microbiol Rev. 2014 Jul; 27(3): p. 441-462.
60. CDC. Centers for Disease Controls and Prevention. [Online].; 2024 [cited 2024 Jan 10]. Available from: <https://www.cdc.gov/adenovirus/>.
61. Johnson, Hargest , Cortez , Meliopoulos VA, Schultz-Cherry S. Astrovirus Pathogenesis. Viruses. 2017 Jan 22; 9(1).
62. Bosch A, Pintó RM, Guix. Human astroviruses. Clin Microbiol Rev. 2014 Oct; 27(4): p. 1048-1074.
63. Robilotti, Deresinski, Pinsky BA. Norovirus. Clin Microbiol Rev. 2015 Jan; 28(1): p. 134-164.
64. Patel MM, Hall AJ, Vinjé, Parashar UD. Noroviruses: a comprehensive review. J Clin Virol. 2009 Jan; 44(1): p. 1-8.
65. CDC. Centers for Disease Control and Prevention. [Online].; 2024 [cited 2024 Jan 24]. Available from: <https://www.cdc.gov/norovirus/index.html>.
66. Milbrath MO, Spicknall IH, Zelner JL, Moe CL, Eisenberg JNS. Heterogeneity in norovirus shedding duration affects community risk. 2013 August; 141(8): p. 1572–1584.
67. Crawford SE, Ramani S, Tate JE, Parashar UD, Svensson L, Hagbom , et al. Rotavirus infection. Nat Rev Dis Primers. 2017 Nov 09; 3(17083).
68. Desselberger. Rotaviruses. Virus Res. 2014 Sep 22; 190: p. 75-96.
69. Patel MM, E Pitzer, Alonso WJ, Vera, Lopman, Tate , et al. Global seasonality of rotavirus disease. Pediatr Infect Dis J. 2013 Apr; 32(4): p. e134-e147.
70. CDC. Center for Disease Control and Prevention. [Online].; 2021 [cited 2024 Jan 15]. Available from: <https://www.cdc.gov/rotavirus/index.html>.

71. Burnett, Parashar , Tate. Rotavirus Vaccines: Effectiveness, Safety, and Future Directions. *Paediatr Drugs*. 2018 Jun; 20(3): p. 223-233.
72. CDC. Biosafety in Microbiological and Biomedical Laboratories (BMBL) 6th Edition. [Online].; 2023 [cited 2024 Apr 23. Available from: <https://www.cdc.gov/labs/BMBL.html>.
73. Bai X, Zhang J, Ambikan A, Jernberg C, Ehricht R, Scheutz F, et al. Molecular Characterization and Comparative Genomics of Clinical Hybrid Shiga Toxin-Producing and Enterotoxigenic *Escherichia coli* (STEC/ETEC) Strains in Sweden. *Sci Rep*. 2019 Apr 4; 9(1): p. 5619.
74. Hyma KE, Lacher DW, Nelson AM, Bumbaugh AC, Janda JM, Strockbine NA, et al. Evolutionary genetics of a new pathogenic *Escherichia* species: *Escherichia albertii* and related *Shigella boydii* strains. *J Bacteriol*. 2005 Jan 1; 187(2): p. 619-28.
75. Beutin L, Strauch E, Fischer I. Isolation of *Shigella sonnei* lysogenic for a bacteriophage encoding gene for production of Shiga toxin. *Lancet*. 1999 May 1; 353(9163): p. 1498.
76. Vargas M, Gascon J, Jimenez De Anta MT, Vil J. Prevalence of *Shigella* enterotoxins 1 and 2 among *Shigella* strains isolated from patients with traveler's diarrhea. *J Clin Microbiol*. 1999 Nov; 37(11): p. 3608-11.
77. Grotiuz G, Sirok A, Gadea P, Varela G, Schelotto F. Shiga toxin 2-producing *Acinetobacter haemolyticus* associated with a case of bloody diarrhea. *J Clin Microbiol*. 2006 Oct; 44(10): p. 3838-41.
78. Probert WS, McQuaid C, Schrader K. Isolation and identification of an *Enterobacter cloacae* strain producing a novel subtype of Shiga toxin type 1. *J Clin Microbiol*. 2014 Jul; 52(7): p. 2346-51.
79. Schmidt H, Montag M, Bockemühl J, Heesemann J, Karch H. Shiga-like toxin II-related cytotoxins in *Citrobacter freundii* strains from humans and beef samples. *Infect Immun*. 1993 Feb; 61(2): p. 534-43.

80. Schroeder CM, Zhao C, DebRoy C, Torcolini J, Zhao S, White DG, et al. Antimicrobial resistance of *Escherichia coli* O157 isolated from humans, cattle, swine, and food. *Appl Environ Microbiol*. 2002 Feb; 68(2): p. 576-81.
81. Chakraborty S, Khan A, Kahali S, Faruque SM, Yamasaki S, Ramamurthy T. Infantile diarrhoea associated with sorbitol-fermenting, non-shiga toxin-producing *Escherichia coli* O157:H-. *Eur J Clin Microbiol Infect Dis*. 2003 May; 22(5): p. 324-6.
82. Lamba K, Nelson JA, Kimura AC, Poe A, Collins J, Kao AS, et al. Shiga Toxin 1-Producing *Shigella sonnei* Infections, California, United States, 2014-2015. *Emerging infectious diseases*. 2016; 22(4): p. 679–686.
83. Brandal LT, Tunsjø HS, Ranheim TE, Løbersli I, Lange H, Wester AL. Shiga toxin 2a in *Escherichia albertii*. *Journal of clinical microbiology*. 2015; 53(4): p. 1454–1455.
84. Tschäpe H, Prager R, Streckel W, Fruth A, Tietze E, Böhme G. Verotoxinogenic *Citrobacter freundii* associated with severe gastroenteritis and cases of haemolytic uraemic syndrome in a nursery school: green butter as the infection source. *Epidemiol Infect*. 1995; 114(3): p. 441-450.
85. Paton AW, Paton JC. *Enterobacter cloacae* producing a Shiga-like toxin II-related cytotoxin associated with a case of hemolytic-uremic syndrome. *J Clin Microbiol*. 1996; 34(2): p. 463-465.
86. Alperi A, Figueras M. Human isolates of *Aeromonas* possess Shiga toxin genes (*stx1* and *stx2*) highly similar to the most virulent gene variants of *Escherichia coli*. *Clinical Microbiology and Infection*. 2010 Oct; 16(10): p. 1563–1567.
87. EFSA BIOHAZ Panel , Koutsoumanis K, Allende A, Alvarez-Ordóñez A, Bover-Cid S, Chemaly M, et al. Scientific Opinion on the pathogenicity assessment of Shiga toxin-producing *Escherichia coli* (STEC) and the public health risk posed by contamination of food with STEC. *EFSA Journal*. 2020 Jan 29; 18(1): p. 5967, 105 pp.

88. Khan SM, Witola WH. Past, current, and potential treatments for cryptosporidiosis in humans and farm animals: A comprehensive review. *Front Cell Infect Microbiol.* 2023 Jan 24; 13: p. 1115522.
89. CDC. Centers for Disease Control and Prevention. [Online].; 2025 [cited 2025 Oct 16]. Available from: <https://www.cdc.gov/amebiasis/about/index.html>

# Document Revision History

Revision	Description
R2, June 2024	Inclusion of QIAstat-Dx Analyzer 2.0
R3, July 2025	The handbook has been revised to address some inaccuracies identified in the previous version.
R4, October 2025	<p>Removal of Analyzer 1.0 in the document</p> <p>Removal of sample type selection for Para-Pak C&amp;S and FecalSwab (sample type will be defined as "Preserved Stool", independently of the sample collection method) and removal of limitation of reporting of Enteropathogenic <i>E. coli</i> (EPEC), STEC <i>stx1/stx2</i> and STEC O157 for FecalSwab</p> <p>Addition of clinical performance data in FecalSwab specimens for EPEC, STEC, and <i>E. coli</i> O157</p> <p>Minor editorial/grammatical changes/corrections/updates throughout the document</p>

## Limited License Agreement for QIAstat-Dx® Gastrointestinal Panel 2

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this Instructions for Use and for use with components contained in the panel only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this panel with any components not included within this panel except as described in the protocols provided with the product, this Instructions for Use, and additional protocols available at [www.qiagen.com](http://www.qiagen.com). Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this panel and/or its use(s) do not infringe the rights of third-parties.
3. This panel and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the panel agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the panel and/or its components.

For updated license terms, see [www.qiagen.com](http://www.qiagen.com).

Trademarks: QIAGEN®, Sample to Insight®, QIAstat-Dx® (QIAGEN Group). Registered names, trademarks, etc., used in this document, even when not specifically marked as such are not to be considered unprotected by law.

10/2025 HB-3018-004 © 2025 QIAGEN, all rights reserved.

