

QuantiTect® Probe PCR Kit

The QuantiTect Probe PCR Kit (cat. nos. 204341, 204343 and 204345) should be stored immediately upon receipt at -30 to -15°C in a constant-temperature freezer and protected from light. 2x QuantiTect Probe PCR Master Mix can also be stored protected from light at 2 – 8°C for up to 6 months, depending on the expiration date.

Further information

- *QuantiTect Probe PCR Handbook*: www.qiagen.com/HB-0233
- Safety Data Sheets: www.qiagen.com/safety
- Technical assistance: support.qiagen.com

Notes before starting

- This protocol is optimized for quantification of gDNA and cDNA targets using dual-labeled probes (including TaqMan® probes, FRET probes and Molecular Beacons) with most real-time cyclers. PCR is carried out in the presence of ROX passive reference dye, which is included in 2x QuantiTect Probe PCR Master Mix and is necessary for all real-time cyclers from Applied Biosystems. The presence of ROX dye does not interfere with real-time PCR on any other instrument.
- For the highest efficiency in real-time PCR using sequence-specific probes, targets should ideally be 100–150 bp in length.
- 2x QuantiTect Probe PCR Master Mix contains dUTP, which allows the use of a uracil-N-glycosylase (UNG) pretreatment of the reaction if contamination with carried-over PCR products is suspected.
- Always start with the cycling conditions specified in this protocol, even if using previously established primer–probe systems.
- The PCR must start with an initial incubation step of 15 min at 95°C to activate HotStarTaq® DNA Polymerase.
- Always readjust the threshold value for analysis of every run.

1. Thaw 2x QuantiTect Probe PCR Master Mix (if stored at -30 to -15°C), template gDNA or cDNA, primer and probe solutions and RNase-free water. Mix the individual solutions.
2. Prepare a reaction mix according to Table 1. Due to the hot start, it is not necessary to keep samples on ice during reaction setup or while programming the realtime cycler.

Note: We strongly recommend starting with an initial Mg^{2+} concentration of 4 mM as provided by 2x QuantiTect Probe PCR Master Mix. For a very limited number of targets, reactions may be improved by using Mg^{2+} concentrations of up to 6 mM.

3. Mix the reaction mix thoroughly, and dispense appropriate volumes into PCR tubes, PCR capillaries or the wells of a PCR plate.

Table 1. Reaction setup

Component	LightCycler® 1.x and 2.0		Other cyclers*	
	Volume	Final conc.	Volume	Final conc.
Reaction mix				
2x QuantiTect Probe PCR Master Mix†	10 μl †	1x	25 μl †	1x
Primer A	Variable	0.5 μM	Variable	0.4 μM
Primer B	Variable	0.5 μM	Variable	0.4 μM
Probe	Variable	0.1–0.2 μM	Variable	0.1–0.2 μM
QuantiTect RT Mix	0.2 μl †		0.5 μl †	
RNase-free water	Variable	–	Variable	–
Optional: Uracil-N-glycosylase, heatHabile	Variable	0.2 units/reaction	Variable	0.5 units/reaction
Template gDNA or cDNA (added at step 4)	Variable	≤ 1 μg /reaction	Variable	≤ 500 ng/reaction
Total reaction volume	20 μl		50 μl	

* Includes Rotor-Gene® cyclers and instruments from Applied Biosystems, Bio-Rad, Cepheid, Eppendorf, Roche and Agilent/Stratagene.

† Provides a final concentration of 4 mM MgCl_2 .

‡ If using a total reaction volume other than indicated, adjust the volume of 2x QuantiTect Probe PCR Master Mix accordingly.

- Add template gDNA or cDNA to the individual PCR tubes, capillaries or wells containing the reaction mix.

Note: For two-step RT-PCR, the volume of the cDNA added (from the undiluted reverse-transcription reaction) should not exceed 10% of the final PCR volume.

- Program the real-time cycler according to Table 2 or, if using FRET probes on the LightCycler 1.x or LightCycler 2.0, Table 3.

Note: Data acquisition should be performed during the combined annealing/extension (2-step cycling) or annealing (3-step cycling) step.

- Place the PCR tubes or plates in the real-time cycler, and start the cycling program.

Note: If using the Applied Biosystems® 7500, it is necessary to adjust the preset threshold value to a lower value. Use a value of 0.01 as a starting point.

Note: If using the LightCycler 1.x or 2.0, we recommend using the “second derivative maximum” method for data analysis. Always readjust the noise band for analysis of every run if using the “fit-point” method for data analysis.

Table 2. Cycling conditions for dual-labeled probes

Step	LightCycler® 1.x and 2.0		Other cyclers	
	Time*	Temperature	Time	Temperature
Optional: UNG pretreatment (if UNG included in reaction mix)	2 min	50°C	2 min	50°C
PCR initial heat activation	15 min	95°C	15 min	95°C
2-step cycling:				
Denaturation	0 s	95°C	15 s†	94°C
Combined annealing/extension	60 s	60°C	60 s	60°C
Number of cycles	35–55‡		35–45‡	

* Ramp rate: 20°C/s.

† SmartCycler® users can reduce denaturation time to 1 s to take advantage of cycling capacities.

‡ The number of cycles depends on the amount of template DNA.

Table 3. Cycling conditions for FRET probes on the LightCycler 1.x and 2.0

Step	Time	Temperature	Ramp rate
Optional: UNG pretreatment (if UNG included in reaction mix)	2 min	50°C	20°C/s
PCR initial heat activation	15 min	95°C	20°C/s
3-step cycling:			
Denaturation	0 s	95°C	20°C/s
Annealing	30 s	50–60°C	20°C/s
Extension	30 s	72°C	2°C/s
Number of cycles	35–55*		

* The number of cycles depends on the amount of template DNA.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, HotStarTaq®, QuantiTect®, Rotor-Gene® (QIAGEN Group); Applied Biosystems® (Applied Biosystems Corporation or its subsidiaries); LightCycler®, TaqMan® (Roche Group); SmartCycler® (Cepheid) 1101228 03/2016 HB-0735-002 © 2016 QIAGEN, all rights reserved.