

# Ni-NTA Spin Column Purification of 6xHis-Tagged Proteins under Native Conditions from *E. coli* Cell Lysates

The Ni-NTA Spin Kit (cat. no. 31314) and Ni-NTA Spin Columns (cat. no. 31014) can be stored at 2–8°C for up to 18 months if not otherwise stated on label.

## Further information

- *Ni-NTA Spin Kit Handbook*: [www.qiagen.com/HB-0883](http://www.qiagen.com/HB-0883)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](mailto:support.qiagen.com)

## Notes before starting

- Prepare bacterial cell pellet (stored at –20°C or –80°C).
  - Buffer compositions are provided in the appendix of the *Ni-NTA Spin Kit Handbook*.
  - Lysis requires Benzonase® nuclease (e.g., Novagen cat. no. 70664-3).
1. Resuspend a pellet derived from 5 ml cell culture volume in 630 µl Lysis Buffer (NPI-10). Add 70 µl Lysozyme Stock Solution (10 mg/ml) and add 3 Units/ml culture volume Benzonase Nuclease (i.e., for cell pellets from 5 ml cultures, add 15 Units Benzonase Nuclease).
  2. Incubate on ice for 15–30 min.
  3. Centrifuge lysate at 12,000 x g for 15–30 min at 4°C. Collect supernatant.  
**Note:** Save 20 µl of the cleared lysate for SDS-PAGE analysis.

4. Equilibrate the Ni-NTA spin column with 600  $\mu$ l Buffer NPI-10. Centrifuge for 2 min at 890  $\times$   $g$  (approx. 2900 rpm).

**Note:** The spin columns should be centrifuged with an open lid to ensure that the centrifugation step is completed under 2 min.

5. Load up to 600  $\mu$ l of the cleared lysate containing the 6xHis-tagged protein onto the pre-equilibrated Ni-NTA spin column. Centrifuge for 5 min at 270  $\times$   $g$  (approx. 1600 rpm), and collect the flow-through.

**Note:** To ensure efficient binding, do not exceed 270  $\times$   $g$  when centrifuging Ni-NTA spin columns. The spin columns can be centrifuged with an open lid.

**Note:** Save the flow-through for analysis by SDS-PAGE.

6. Wash the Ni-NTA spin column twice with 600  $\mu$ l Buffer NPI-20. Centrifuge for 2 min at 890  $\times$   $g$  (approx. 2900 rpm).

**Note:** Save the flow-through for analysis by SDS-PAGE.

7. Elute the protein twice with 300  $\mu$ l Buffer NPI-500. Centrifuge for 2 min at 890  $\times$   $g$  (approx. 2900 rpm), and collect the eluate.

**Note:** Most of the 6xHis-tagged protein (>80%) should elute in the first 300  $\mu$ l eluate. The remainder will elute in the second 300  $\mu$ l.



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup> (QIAGEN Group); Benzonase<sup>®</sup> (Merck KGaA). 1101274 03/2016 HB-0683-002 © 2016 QIAGEN, all rights reserved.