

## Quick-Start Protocol

# AllPrep<sup>®</sup> PowerFecal<sup>®</sup> Pro DNA/RNA Kit

Solution CD2 should be stored at 2–8°C upon arrival. All other reagents and kit components should be stored at room temperature (15–25°C) until the expiry date printed on the box label.

AllPrep PowerFecal Pro DNA/RNA Kit is for the simultaneous isolation of microbial genomic DNA and total RNA from stool samples in two separate eluates.

### Further information

- *AllPrep PowerFecal Pro DNA/RNA Kit Handbook*: [www.qiagen.com/HB-2851](http://www.qiagen.com/HB-2851)
- Safety Data Sheets: [www.qiagen.com/safety](http://www.qiagen.com/safety)
- Technical assistance: [support.qiagen.com](http://support.qiagen.com)

### Notes before starting

- Phenol-chloroform-isoamyl alcohol (25:24:1, pH 6.5-8.0) is required for this protocol.
- Ensure that the PowerBead Pro Tubes rotate freely in the centrifuge without rubbing.
- Perform all centrifugation steps at room temperature.
- Refer to the *AllPrep PowerFecal Pro DNA/RNA Kit Handbook* for optimal homogenization method in step 3.

---

## Procedure

### Sample Pretreatment

1. Spin the PowerBead Pro Tube briefly to ensure that the beads have settled at the bottom.
2. Add 50–100 mg of stool, 650  $\mu$ l of Solution CD1 and 100  $\mu$ l phenol–chloroform–isoamyl alcohol (25:24:1, pH 6.5–8.0) to the PowerBead Pro Tube and vortex briefly to mix.
3. Secure the PowerBead Pro Tube horizontally on a vortex adapter for 1.5–2 ml tubes (cat. no. 13000-V1-24). Orient tube caps to point toward the center of the Vortex Adapter. Vortex at maximum speed for 10 min.

**Note:** If using the vortex adapter for more than 12 preps simultaneously, increase the vortex time by 5–10 min.

For more information about other bead beating methods, see the “Protocol: Detailed” section of *AllPrep PowerFecal Pro DNA/RNA Kit Handbook*.

4. Centrifuge the PowerBead Pro Tube at 15,000  $\times g$  for 1 min. Transfer the supernatant to a clean 2 ml microcentrifuge tube (provided).

**Note:** Expect a supernatant volume of 500–600  $\mu$ l. The supernatant may still contain some stool particles.

5. Add 200  $\mu$ l Solution CD2 and vortex for 5 s. Centrifuge at 15,000  $\times g$  for 1 min at room temperature.
6. Avoiding the pellet, transfer 300  $\mu$ l of supernatant to a clean 2 ml microcentrifuge tube (provided).

**Note:** It is feasible to use higher supernatant volumes, please refer to the HB for detailed information.

7. Add 300  $\mu$ l of Solution CD3. Vortex briefly to mix.

---

## DNA binding

8. Load 600  $\mu\text{l}$  supernatant-CD3 mix into a MB DNA Spin Column (white) and centrifuge at 15,000  $\times g$  for 1 min collecting the flow-through in a 2 ml tube for RNA purification, then place the spin column in a new 2 ml collection tube.

## RNA binding

9. Add 300  $\mu\text{l}$  96–100% ethanol to the flow-through from step 8 and mix by pipetting up and down.

**Note:** If you used a higher volume of the supernatant in step 6, you would have to adjust the binding conditions according to the instructions in the handbook.

10. Transfer up to 700  $\mu\text{l}$  of the mix to a MB RNA Spin Column (pink) placed in a 2 ml collection tube. Centrifuge at 15,000  $\times g$  for 1 min. Discard the flow-through.

**Note:** If the volume of the mixture exceeds 700  $\mu\text{l}$ , centrifuge successive aliquots in the same MB RNA Spin Column. Discard the flow-through after each centrifugation.

## DNA and RNA washing and elution

11. Add 650  $\mu\text{l}$  Solution EA to a MB DNA Spin Column (white) and a MB RNA Spin Column (pink) and centrifuge at 15,000  $\times g$  for 1 min. Discard the flow-through.
12. Add 500  $\mu\text{l}$  Solution C5. Centrifuge at 15,000  $\times g$  for 1 min.
13. Discard flow-through and place the MB RNA and DNA Spin Columns into clean 2 ml collection tubes (provided). Centrifuge at 20,000  $\times g$  (or full speed) for 1 min.
14. Place the MB RNA and DNA Spin Columns into clean 1.5 ml Elution Tubes (provided).
15. Add 100  $\mu\text{l}$  RNase-free water to the center of the white filter membrane.
16. Incubate at room temperature for at least 1 min.
17. Centrifuge at 15,000  $\times g$  for 1 min. Discard the MB RNA and DNA Spin Columns. The RNA and DNA are now ready for any downstream applications.

## Document Revision History

Date	Changes
02/2022	Initial release



Scan QR code for handbook.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN<sup>®</sup>, Sample to Insight<sup>®</sup>, AllPrep<sup>®</sup>, PowerFecal<sup>®</sup> (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1126004 02/2022 HB-2950-001 © 2022 QIAGEN, all rights reserved.