




August 2025

QuantiFERON[®]-TB Gold Plus (QFT[®]-Plus) Instructions for Use

Version 1

 2 x 96 (cat. no. 622130)

 20 x 96 (cat. no. 622832)


The whole blood IFN- γ test measuring responses to ESAT-6 and CFP-10 peptide antigens

Rx ONLY

IVD For In Vitro Diagnostic Use

REF 622130, 622832

 www.qiagen.com

 QIAGEN, 19300 Germantown Road, Germantown, MD 20874, USA
Phone: +1-800-426-8157

R11

Table of Contents

Intended Use	4
Summary and Explanation of the Test	5
Principles of the Assay	8
Components and Storage	10
ELISA components	10
Materials Required but Not Provided	11
Warnings and Precautions	12
Warnings	12
Precautions	13
Storage and Handling	16
ELISA kit reagents	16
Reconstituted and unused reagents	16
Directions for Performing ELISA	17
Time required for performing assay	17
Stage 1: Post-incubation of blood collection tubes and harvesting of plasma	17
Stage 2: IFN-g ELISA	18
Calculations and Test Interpretation	24
Generation of standard curve and sample values	24
Quality control of the test	26
Interpretation of results	27
Limitations	30
Performance Characteristics	32
Clinical studies	32
Specificity	32
Sensitivity	34
Performance in subjects with identified risk factors for a <i>M. tuberculosis</i> infection (mixed-risk individuals)	38
Expected Values	41
Observed response distribution - risk stratified	41

Assay Performance Characteristics	46
Technical Information	49
Indeterminate results	49
Clotted plasma samples	49
Lipemic plasma samples	49
Troubleshooting Guide	50
Appendices	52
Appendix A: US Centers for Disease Control and Prevention (CDC) Guidelines	52
Appendix B: Abbreviated ELISA Test Procedure	52
References	54
Symbols	63
Contact Information	65
Document Revision History	66

Intended Use

QuantiFERON®-TB Gold Plus (QFT-Plus) is an in vitro diagnostic test using a peptide cocktail simulating ESAT-6 and CFP-10 proteins to stimulate cells in heparinized whole blood. Detection of interferon- γ (IFN- γ) by enzyme-linked immunosorbent assay (ELISA) is used to identify in vitro responses to these peptide antigens that are associated with *Mycobacterium tuberculosis* infection.

QFT-Plus is an indirect test for *M. tuberculosis* infection (including disease) and is intended for use in conjunction with risk assessment, radiography, and other medical and diagnostic evaluations.

Summary and Explanation of the Test

Tuberculosis is a communicable disease caused by infection with *M. tuberculosis* complex organisms (*M. tuberculosis*, *M. bovis*, *M. africanum*, *M. microti*, *M. canetti* and *M. caprae*), which typically spreads to new hosts via airborne droplet nuclei from patients with pulmonary tuberculosis disease. A newly infected individual can become ill from tuberculosis within weeks to months, but most infected individuals remain well. Latent tuberculosis infection (LTBI), a non-communicable asymptomatic condition, persists in some, who might develop tuberculosis disease months or years later. The main purpose of diagnosing LTBI is to consider medical treatment for preventing tuberculosis disease. For more than 100 years, the tuberculin skin test (TST) was the only available method for diagnosing LTBI. Cutaneous sensitivity to tuberculin develops from 2 to 10 weeks after infection. However, some infected individuals, including those with a wide range of conditions hindering immune functions, but also others without these conditions, do not respond to tuberculin. Conversely, some individuals who are unlikely to have *M. tuberculosis* infection exhibit sensitivity to tuberculin and have positive TST results after vaccination with Bacille Calmette-Guérin (BCG) or infection with mycobacteria other than *M. tuberculosis* complex, or undetermined other factors.

LTBI must be distinguished from tuberculosis disease, a reportable condition that usually involves the lungs and lower respiratory tract but may also affect other organ systems. Tuberculosis disease is diagnosed from historical, physical, radiological, and mycobacteriological findings.

QuantiFERON-TB Gold Plus (QFT-Plus) test is the fourth generation in QuantiFERON-TB testing technology assessing cell-mediated response through a quantitative measurement of IFN- γ in a whole blood sample. QFT-Plus is a test that measures the cell-mediated immune (CMI) responses to peptide antigens that simulate mycobacterial proteins. These proteins, ESAT-6 and CFP-10 are absent from all BCG strains and from most non-tuberculosis mycobacteria with the exception of *M. kansasii*, *M. szulgai* and *M. marinum* (1). Individuals infected with

M. tuberculosis complex organisms usually have lymphocytes in their blood that recognize these and other mycobacterial antigens. This recognition process involves the generation and secretion of the cytokine, IFN- γ . The detection and subsequent quantification of IFN- γ forms the basis of this test.

Tuberculin skin test and IGRA tests are helpful but insufficient for diagnosing *M. tuberculosis* complex infection in sick patients – a positive result can support the diagnosis of tuberculosis disease; however, infections by other mycobacteria (e.g., *M. kansasii*) could also cause positive results. Other medical and diagnostic evaluations are necessary to confirm or exclude tuberculosis disease.

The antigens used in QFT-Plus are a peptide cocktail simulating the proteins ESAT-6 and CFP 10. Numerous studies have demonstrated that these peptide antigens stimulate IFN- γ responses in T cells from individuals infected with *M. tuberculosis* but generally not from uninfected or BCG-vaccinated persons without disease or risk for LTBI (1–32). However, medical treatments or conditions that impair immune functionality can potentially reduce IFN- γ responses. Patients with certain other mycobacterial infections might also be responsive to ESAT-6 and CFP-10, as the genes encoding these proteins are present in *M. kansasii*, *M. szulgai*, and *M. marinum* (1, 23).

In *M. tuberculosis* infection, CD4⁺ T cells play a critical role in immunological control through their secretion of the cytokine IFN- γ . Evidence now supports a role for CD8⁺ T cells participating in the host defense to *M. tuberculosis* by producing IFN- γ and other soluble factors, which activate macrophages to suppress growth of *M. tuberculosis*, kill infected cells, or directly lyse intracellular *M. tuberculosis* (33–35). IFN- γ producing *M. tuberculosis*-specific CD8⁺ cells have been detected in subjects with LTBI and with active TB (36–39). Moreover, ESAT-6 and CFP-10 specific CD8⁺ T lymphocytes are described as being more frequently detected in subjects with active TB disease versus LTBI and may be associated with a recent *M. tuberculosis* exposure (40–42). In addition, *M. tuberculosis*-specific CD8⁺ T cells producing

IFN- γ have also been detected in active TB subjects with HIV co-infection (43, 44) and in young children with TB disease (45).

QFT-Plus has two distinct TB antigen tubes: TB Antigen Tube 1 (TB1) and TB Antigen Tube 2 (TB2). Both tubes contain peptide antigens from the *M. tuberculosis*-complex-associated antigens, ESAT-6 and CFP-10. Both the TB1 tube and TB2 tubes contain peptides from ESAT-6 and CFP-10 that are designed to elicit CMI responses from CD4⁺ T-helper lymphocytes; the TB2 tube contains an additional set of peptides targeted to the induction of CMI responses from CD8⁺ cytotoxic T lymphocytes.

Diagnostic testing for *M. tuberculosis* using Interferon Gamma Release Assays should follow applicable published guidelines. QFT-Plus has been studied in adults and children identified at higher risk for TB infection, including individuals recently exposed to persons with active tuberculosis, persons living with HIV or otherwise immunocompromised, immigrants from high-burden countries, and residents or employees of high-risk congregate settings. (See "Warnings and Precautions" on page 12 regarding use in co-morbid conditions that may affect immune function). (46–84)

Links to the most recent American Thoracic Society/Infectious Diseases Society of America/Centers for Disease Control and Prevention Clinical Practice Guidelines: Diagnosis of Tuberculosis in Adults and Children, as well as to other information regarding diagnostic testing for tuberculosis, are available at [Tuberculosis \(TB\)](#).

Principles of the Assay

The QFT-Plus assay uses specialized blood collection tubes, containing peptide antigens that simulate *M. tuberculosis* proteins, which are used to collect whole blood. Incubation of the blood occurs in the tubes for 16 to 24 hours, after which, plasma is harvested and tested for the presence of IFN- γ produced in response to the peptide antigens.

The QFT-Plus test is performed in two stages. First, whole blood is collected into each of the QFT-Plus Blood Collection Tubes, which include a Nil tube, TB1 tube, TB2 tube, and a Mitogen tube. Alternatively, blood may be collected in a single blood collection tube that contains lithium or sodium-heparin as the anticoagulant, and then transferred to QFT-Plus Blood Collection Tubes.

The QFT-Plus Blood Collection Tubes are shaken to mix antigen with the blood and should be incubated at $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$ as soon as possible, and within 16 hours of collection. Following a 16 to 24 hour incubation period, the tubes are centrifuged, the plasma is removed, and the amount of IFN- γ (IU/mL) is measured by ELISA. The QFT-Plus ELISA uses a recombinant human IFN- γ standard, which has been assayed against a reference IFN- γ preparation (NIH Ref: Gxg01-902-535). Results for test samples are reported in International Units per mL (IU/mL) relative to a standard curve prepared by testing dilutions of the standard supplied with the kit.

Heterophile (e.g., human anti-mouse) antibodies in serum or plasma of certain individuals are known to cause interference with immunoassays. The effect of heterophile antibodies in the QFT-Plus ELISA is minimized by the addition of normal mouse serum to the Green Diluent and the use of F(ab')₂ monoclonal antibody fragments as the IFN- γ capture antibody coated to the microplate wells.

A QFT-Plus assay is considered positive with an IFN- γ response to either TB antigen tube that is significantly above the Nil IFN- γ IU/mL value. The plasma sample from the Mitogen tube serves as an IFN- γ positive control for each specimen tested. A low response to Mitogen (<0.5

IU/mL) indicates an indeterminate result when a blood sample also has a negative response to the TB antigens. This pattern may occur with insufficient lymphocytes, reduced lymphocyte activity due to improper specimen handling, filling/mixing of the Mitogen tube, or inability of the patient's lymphocytes to generate IFN- γ . Elevated levels of IFN- γ in the Nil sample may occur with the presence of heterophile antibodies, or to intrinsic IFN- γ secretion. The Nil tube adjusts for background (e.g., elevated levels of circulating IFN- γ or presence of heterophile antibodies). The IFN- γ level of the Nil tube is subtracted from the IFN- γ level for the TB antigen tubes and the Mitogen tube.

Components and Storage

ELISA components

ELISA components	2-plate kit	Reference Lab Pack
Catalog no.	622130	622832
Microplate strips (12 x 8 wells) coated with murine anti- human IFN- γ monoclonal antibody	2 sets of 12 x 8 Microplate Strips	20 sets of 12 x 8 Microplate Strips
IFN- γ Standard, lyophilized (contains recombinant human IFN-g, bovine casein, 0.01% w/v Thimerosal)	1 x vial (8 IU/mL when reconstituted)	10 x vials (8 IU/mL when reconstituted)
Green Diluent (contains bovine casein, normal mouse serum, 0.01% w/v Thimerosal)	1 x 30 mL	10 x 30 mL
Conjugate 100x Concentrate, lyophilized (murine anti- human IFN- γ HRP, contains 0.01% Thimerosal)	1 x 0.3 mL (when reconstituted)	10 x 0.3 mL (when reconstituted)
Wash Buffer 20x Concentrate (pH 7.2, contains 0.05% v/v ProClin® 300)	1 x 100 mL	10 x 100 mL
Enzyme Substrate Solution (contains H ₂ O ₂ , 3,3',5,5' Tetramethylbenzidine)	1 x 30 mL	10 x 30 mL
Enzyme Stopping Solution (contains 0.5 M H ₂ SO ₄)*	1 x 15 mL	10 x 15 mL

* See "Precautions" on page 13

Materials Required but Not Provided

QuantiFERON-TB Gold Plus Blood Collection Tubes – Please refer to the *QuantiFERON-TB Gold Plus (QFT-Plus) Blood Collection Tubes Instructions for Use* (HB-3304).

Important: Make sure that all instruments used in this procedure have been checked and calibrated according to the manufacturer’s recommendations.

- 37°C ± 1°C incubator (with or without CO₂)
- Calibrated variable-volume pipets for delivery of 10 µL to 1000 µL with disposable tips
- Calibrated multichannel pipet capable of delivering 50 µL to 100 µL with disposable tips
- Deionized or distilled water, 2 liters
- **Optional:** 1 mL microtubes with caps in 96-well format racks or uncoated microplates with plastic seals for plasma storage (22 patients/rack or plate)
- Centrifuge capable of centrifuging the blood tubes at least to 2000 RCF (*g*)
- Microplate shaker capable of speeds between 500 and 1000 rpm
- Microplate washer (for safety in handling plasma samples, an automated washer is recommended)
- Microplate reader fitted with 450 nm filter and 620 nm to 650 nm reference filter
- Variable speed vortex
- Timer
- Graduated cylinder, 1 liter or 2 liters
- Reagent reservoirs
- Plate lid

Warnings and Precautions

For in vitro diagnostic use only.

For prescription use only. Federal law restricts this device to sale by or on the order of a physician.

Warnings

- A negative QFT-Plus result does not preclude the possibility of *M. tuberculosis* infection or tuberculosis disease: false-negative results can be due to stage of infection (e.g., specimen obtained prior to the development of cellular immune response), co-morbid conditions that affect immune function, incorrect handling of the blood collection tubes following venipuncture, incorrect performance of the assay, or other individual immunological variables. Heterophile antibodies or non-specific IFN- γ production from other inflammatory conditions may mask specific responses to ESAT-6 or CFP-10 peptides.
- A positive QFT-Plus result should not be the sole or definitive basis for determining infection with *M. tuberculosis*. Incorrect performance of the assay may cause false positive QFT-Plus results.
- A positive QFT-Plus result should be followed by further medical evaluation for active tuberculosis disease (e.g., Acid Fast Bacilli smear and culture, chest X-ray).
- While ESAT-6 and CFP-10 are absent from all BCG strains and from most known nontuberculous mycobacteria, it is possible that a positive QFT-Plus result may be due to infection by *M. kansasii*, *M. szulgai*, or *M. marinum*. If such infections are suspected, alternative tests should be performed.
- A false-negative QFT-Plus result can be caused by incorrect blood sample collection or improper handling of the specimen affecting lymphocyte function. Please refer to

QuantiFERON-TB Gold Plus (QFT-Plus) Blood Collection Tubes Instructions for Use (HB-3304), for correct handling of the blood specimens. Delay in incubation may cause false negative or indeterminate results, and other technical parameters may affect ability to detect a significant IFN- γ response.

- The effect of lymphocyte count on reliability of QFT-Plus results is unknown. Lymphocyte counts may vary over time for any individual person, and from person to person. The minimum number of lymphocytes required for a reliable test result has not been established and may also be variable.
- A positive QFT-Plus result can suggest and support the diagnosis of tuberculosis disease. ESAT-6 and CFP-10 are present in *M. tuberculosis*, but infections by other mycobacteria, including *M. kansasii*, *M. szulgai*, and *M. marinum* may also cause positive results. Other diagnostic evaluations (e.g., AFB smear and culture, chest X-ray) besides QFT-Plus are needed to confirm tuberculosis disease.

The predictive value of a negative QFT-Plus result in immunosuppressed persons has not been determined.

Precautions

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN kit and kit component.

CAUTION: Handle human blood as if potentially infectious.



Observe relevant blood handling guidelines. Dispose of samples and materials in contact with blood or blood products in accordance with federal, state, and local regulations.

QuantiFERON Enzyme Stopping Solution



Contains: sulfuric acid. Warning! May be corrosive to metals. Causes skin irritation. Causes serious eye irritation. Wear protective gloves/protective clothing/eye protection/face protection.

QuantiFERON Enzyme Substrate Solution

Warning! Causes mild skin irritation. Wear protective gloves/ protective clothing/ eye protection/ face protection.

QuantiFERON Green Diluent



Contains: tartrazine. Warning! May cause an allergic skin reaction. Wear protective gloves/ protective clothing/ eye protection/ face protection.

QuantiFERON Wash Buffer 20x Concentrate

Harmful to aquatic life with long lasting effects. Avoid release to the environment.

Further information

Safety Data Sheets: www.qiagen.com/safety

- Thimerosal is used as a preservative in some QFT-Plus reagents. It may be toxic upon ingestion, inhalation or skin contact.
- Deviations from the *QuantiFERON-TB Gold Plus (QFT-Plus) Instructions for Use* may yield erroneous results. Please read the instructions carefully before use.
- Do not use kit if any reagent bottle shows signs of damage or leakage prior to use.
- **Important:** Inspect vials prior to use. Do not use Conjugate or IFN- γ Standard vials that show signs of damage or if the rubber seal has been compromised. Do not handle broken

vials. Take the appropriate safety precautions to dispose of vials safely.

Recommendation: Use a vial de-crimper to open the Conjugate or IFN- γ Standard vials to minimize risk of injury from the metal crimp cap.

- Do not mix or use the Microplate strips, IFN- γ Standard, Green Diluent, or Conjugate 100x Concentrate from different QFT-Plus kit batches. Other reagents (Wash Buffer 20x Concentrate, Enzyme Substrate Solution, and Enzyme Stopping Solution) can be interchanged between kits providing the reagents are within their expiration periods and lot details recorded.
- Discard unused reagents and biological samples in accordance with Local, State, and Federal regulations.
- Do not use the QFT-Plus ELISA kit after the expiration date.
- Correct laboratory procedures should be adhered to at all times.
- Make sure that laboratory equipment such as plate washers and readers have been calibrated/validated before use.

Storage and Handling

Refer to the *QuantiFERON-TB Gold Plus (QFT-Plus) Blood Collection Tubes Instructions for Use* (HB-3304) for details on the blood collection workflow for the QFT-Plus Test.

ELISA kit reagents

- Store kit at 2–8°C.
- Always protect Enzyme Substrate Solution from direct sunlight.

Reconstituted and unused reagents

- For instructions on how to reconstitute the reagents, refer to "Stage 2: IFN-g ELISA" on page 18.
- The reconstituted kit standard may be kept for up to 3 months if stored at 2–8°C.

Note the date the kit standard was reconstituted.

- The reconstituted Conjugate 100x Concentrate must be returned to storage at 2–8°C and must be used within 3 months.

Note the date the conjugate was reconstituted.

- Working strength conjugate must be used within 6 hours of preparation.
- Working strength wash buffer may be stored at room temperature ($22 \pm 5^\circ\text{C}$) for up to 2 weeks.

Directions for Performing ELISA

Time required for performing assay

In order to obtain valid results from the QFT-Plus assay, the operator needs to perform specific tasks within set times. Prior to performing the assay, it is recommended to plan each stage of the assay carefully to allow adequate time between steps. The time required is estimated below; the time of testing multiple samples when batched is also indicated.

ELISA plate

- Approximately 3 hours for one ELISA Plate
- <1 hour labor
- Add 10 to 15 minutes for each extra plate

Stage 1: Post-incubation of blood collection tubes and harvesting of plasma

Prior to harvesting plasma, samples in QFT-Plus Blood Collection Tubes must have been incubated at 37°C for 16–24 hours. The incubator does not require CO₂ or humidification.

Procedure

1. After incubation of the blood collection tubes at 37 ± 1°C, tubes may be held between 4–27°C for up to 3 days pre or post centrifugation.
2. After incubation of the tubes at 37 ± 1°C, harvesting of the plasma is facilitated by centrifuging tubes for 5–15 minutes at 2000 to 3000 RCF (*g*). The gel plug will separate the cells from the plasma. If this does not occur, the tubes should be re-centrifuged.

3. It is possible to harvest the plasma without centrifugation, but additional care is required to remove the plasma without disturbing the cells.
4. Plasma samples should only be harvested using a pipet.

Important: After centrifugation, avoid pipetting plasma up and down or mixing plasma by any means prior to harvesting. At all times, take care not to disturb material on the surface of the gel.

Plasma samples can be loaded directly from centrifuged blood collection tubes into the QFT-Plus ELISA plate.

Plasma samples can be stored in centrifuged QFT-Plus Blood Collection Tubes for up to 28 days at 2–8°C. Harvested plasma samples can be stored for up to 28 days at 2–8°C or stored below –20°C (preferably less than –70°C) for extended periods.

For adequate test samples, harvest at least 150 µL of plasma.

Stage 2: IFN-γ ELISA

Refer to "Components and Storage" on page 10 and "Materials Required but Not Provided" on page 11, for materials required to perform ELISA.

Procedure

1. All plasma samples and reagents, except for Conjugate 100x Concentrate, must be brought to room temperature ($22 \pm 5^\circ\text{C}$) before use. Allow at least 60 minutes for equilibration.
2. Remove ELISA plate strips that are not required from the frame, reseal in the foil pouch, and return to the refrigerator for storage until required.
3. Allow at least 1 strip for the QFT-Plus standards and sufficient strips for the number of subjects being tested (refer to Figure 2 for the recommended plate format). After use,

retain the frame and lid for use with the remaining strips.

- a. Reconstitute the IFN- γ Standard with the volume of deionized or distilled water indicated on the label of the vial. Mix gently to minimize frothing and ensure that the entire content of the vial is completely dissolved. Reconstitution of the IFN- γ standard to the correct volume will produce a solution with a concentration of 8.0 IU/mL.
- b. Using the reconstituted standard, prepare a dilution series of 4 IFN- γ concentrations (refer to Figure 1).
- c. A standard curve should be generated with the following IFN- γ concentrations:
 - S1 (Standard 1) contains 4.0 IU/mL
 - S2 (Standard 2) contains 1.0 IU/mL
 - S3 (Standard 3) contains 0.25 IU/mL
 - S4 (Standard 4) contains 0 IU/mL (Green Diluent [GD] alone)
- d. The standards must be assayed at least in duplicate.
- e. Prepare fresh dilutions of the kit standard for each ELISA session.

Table 1. Example of procedure for duplicate standards

A	Label 4 tubes: S1, S2, S3, S4
B	Add 150 μ L of GD to S1, S2, S3, S4
C	Add 150 μ L of the kit standard to S1 and mix thoroughly
D	Transfer 50 μ L from S1 to S2 and mix thoroughly
E	Transfer 50 μ L from S2 to S3 and mix thoroughly
F	GD alone serves as the zero standard (S4)




Figure 1. Preparation of standard curve dilution series.

4. Reconstitute lyophilized Conjugate 100x Concentrate with 0.3 mL of deionized or distilled water. Mix gently to minimize frothing and ensure that the entire content of the vial is completely dissolved.
 - a. Working strength conjugate is prepared by diluting the required amount of reconstituted Conjugate 100x Concentrate in Green Diluent (Table 2).
 - b. Working strength conjugate should be used within 6 hours of preparation.
 - c. Return any unused Conjugate 100x Concentrate to 2–8°C immediately after use.

Table 2. Conjugate preparation (working strength)

Number of strips	Volume of conjugate (100x concentrate) (µL)	Volume of Green Diluent (mL)
2	10	1.0
3	15	1.5
4	20	2.0
5	25	2.5
6	30	3.0
7	35	3.5

Table 2. Conjugate preparation (working strength) (continued)

Number of strips	Volume of conjugate (100x concentrate) (µL)	Volume of Green Diluent (mL)
8	40	4.0
9	45	4.5
10	50	5.0
11	55	5.5
12	60	6.0

- For plasma samples harvested from blood collection tubes and subsequently stored (refrigerated or frozen), thoroughly mix the stored sample before addition to the ELISA well.

Important: If plasma samples are to be added directly from the centrifuged QFT-Plus Blood Collection Tubes, any mixing of the plasma should be avoided. At all times, take care not to disturb material on the surface of the gel.

- Add 50 µL of freshly prepared working strength conjugate to each ELISA plate well.
- Add 50 µL of test plasma sample to appropriate wells (refer to recommended ELISA plate layout in Figure 2).
- Finally, add 50 µL each of the Standards 1 to 4 to the appropriate plate wells (refer to recommended ELISA plate layout in Figure 2). The standards should be assayed in at least duplicate.

	1	2	3	4	5	6	7	8	9	10	11	12
A	1 N	3 N	5 N	7 N	9 N	S1	S1	13 N	15 N	17 N	19 N	21 N
B	1 TB1	3 TB1	5 TB1	7 TB1	9 TB1	S2	S2	13 TB1	15 TB1	17 TB1	19 TB1	21 TB1
C	1 TB2	3 TB2	5 TB2	7 TB2	9 TB2	S3	S3	13 TB2	15 TB2	17 TB2	19 TB2	21 TB2
D	1 M	3 M	5 M	7 M	9 M	S4	S4	13 M	15 M	17 M	19 M	21 M
E	2 N	4 N	6 N	8 N	10 N	11 N	12 N	14 N	16 N	18 N	20 N	22 N
F	2 TB1	4 TB1	6 TB1	8 TB1	10 TB1	11 TB1	12 TB1	14 TB1	16 TB1	18 TB1	20 TB1	22 TB1
G	2 TB2	4 TB2	6 TB2	8 TB2	10 TB2	11 TB2	12 TB2	14 TB2	16 TB2	18 TB2	20 TB2	22 TB2
H	2 M	4 M	6 M	8 M	10 M	11 M	12 M	14 M	16 M	18 M	20 M	22 M

Figure 2. Recommended ELISA plate layout. S1 (Standard 1), S2 (Standard 2), S3 (Standard 3), S4 (Standard 4). 1N (Sample 1. Nil Control plasma), 1 TB1 (Sample 1. TB1 plasma), 1 TB2 (Sample 1. TB2 plasma), 1M (Sample 1. Mitogen plasma).

9. Cover the ELISA plate and mix the conjugate and plasma samples/standards thoroughly using a microplate shaker for 1 minute at 500 to 1000 rpm. Avoid splashing.
10. Cover the ELISA plate and incubate at room temperature ($22 \pm 5^{\circ}\text{C}$) for 120 ± 5 minutes. The ELISA plate should not be exposed to direct sunlight during incubation. Deviation from the specified temperature range can lead to erroneous results.
11. During the ELISA plate incubation prepare working strength wash buffer. Dilute one part Wash Buffer 20x Concentrate with 19 parts deionized or distilled water and mix thoroughly. Sufficient Wash Buffer 20x Concentrate has been provided to prepare 2 liters of working strength wash buffer.
12. When ELISA plate incubation is complete, wash the ELISA plate wells with 400 μL of working strength wash buffer. Perform the wash step at least 6 times. An automated plate washer is recommended for safety reasons when handling plasma samples.

Thorough washing is very important to the performance of the assay. Make sure each well is completely filled with wash buffer to the top of the well for each wash cycle. A soak period of at least 5 seconds between each cycle is recommended.

Standard laboratory disinfectant should be added to the effluent reservoir, and established procedures followed for the decontamination of potentially infectious material.

13. Tap the ELISA plate face down on an absorbent (low lint) towel to remove residual wash buffer. Add 100 μ L of Enzyme Substrate Solution to each plate well, cover the plate, and mix thoroughly for 1 minute at 500 to 1000 rpm using a microplate shaker.
14. Cover the ELISA plate and incubate at room temperature ($22 \pm 5^{\circ}\text{C}$) for 30 minutes. The ELISA plate should not be exposed to direct sunlight during incubation.
15. Following the 30 minute incubation, add 50 μ L of Enzyme Stopping Solution to each plate well in the same order as the substrate was added and mix thoroughly at 500 to 1000 rpm using a microplate shaker.
16. Measure the Optical Density (OD) of ELISA plate wells within 5 minutes of stopping the reaction using a microplate reader fitted with a 450 nm filter and with a 620 nm to 650 nm reference filter. OD values are used to calculate results.

Calculations and Test Interpretation

QFT-Plus Analysis Software can be used to analyze raw data and calculate results. It is available at www.qiagen.com. Please make sure that the most current version of the QFT-Plus Analysis Software is used.

The software performs a Quality Control assessment of the assay, generates a standard curve and provides a test result for each subject, as detailed in "Interpretation of results" on page 27. The software reports all concentrations greater than 10 IU/mL as ">10" as such values fall beyond the validated linear range of the ELISA.

As an alternative to using the QFT-Plus Analysis Software, results can be determined according to the following method.

Generation of standard curve and sample values

If QFT-Plus Analysis Software is not used

Determination of the standard curve and determination of sample IU/mL values require a spreadsheet program, such as Microsoft® Excel®, if the QFT-Plus software is not used.

Using a spreadsheet program:

1. Determine the mean OD values of the kit standard replicates on each plate.
2. Construct a $\log_{(e)} - \log_{(e)}$ standard curve by plotting the $\log_{(e)}$ of the mean OD (y axis) against the $\log_{(e)}$ of the IFN- γ concentration of the standards in IU/mL (x axis), omitting the zero standard from these calculations. Calculate the line of best fit for the standard curve by regression analysis.

- Use the standard curve to determine the IFN- γ concentration (IU/mL) for each of the test plasma samples, using the OD value of each sample.
- These calculations can be performed using software packages available with microplate readers, and standard spreadsheet or statistical software (such as Microsoft Excel). It is recommended that these packages be used to calculate the regression analysis, the coefficient of variation (%CV) for the standards, and the correlation coefficient (r) of the standard curve.

Sample calculation

If the following OD readings were obtained for the standards, the calculations using $-\log_{(e)}$ would follow those in Table 3.

Table 3. Standard curve

Standard	IU/mL	OD values a and b	Mean OD	%CV	Log _(e) IU/mL	Log _(e) Mean (OD)
Standard 1	4	1.089, 1.136	1.113	3.0	1.386	0.107
Standard 2	1	0.357, 0.395	0.376	7.1	0.000	-0.978
Standard 3	0.25	0.114, 0.136	0.125	NA	-1.386	-2.079
Standard 4	0	0.034, 0.037	0.036	NA	NA	NA

The equation of the curve is $y = 0.7885(X) - 0.9837$, where " m " = 0.7885 and " c " = -0.9837. These values are used in the equation $X = (Y-c)/m$ to solve for X . Based on the standard curve, the calculated correlation coefficient (r) = 1.000.

NA: Not applicable.

Using the criteria specified in "Quality control of the test" on the next page, the assay is determined to be valid.

The standard curve (Table 3) is used to convert the Antigen OD responses to International Units per milliliter (IU/mL).

Table 4. Sample calculation

Antigen	OD value	Log _(e) OD value	X	e ^X (IU/mL)	Antigen –Nil (IU/mL)
Nil	0.037	-3.297	-2.934	0.05	–
TB1	1.161	0.149	1.437	4.21	4.15
TB2	1.356	0.305	1.634	5.12	5.07
Mitogen	1.783	0.578	1.981	7.25	7.20

IFN- γ values (in IU/mL) for the TB1, TB2, and Mitogen tubes are corrected for background by subtracting the IU/mL value obtained for the respective Nil control. These corrected values are used for interpretation of the test results.

Quality control of the test

The accuracy of test results is dependent on the generation of an accurate standard curve. Therefore, results derived from the standards must be examined before test sample results can be interpreted.

For the ELISA to be valid:

- The mean OD value for Standard 1 must be ≥ 0.600 .
- The %CV for Standard 1 and Standard 2 replicate values must be $\leq 15\%$.
- Replicate OD values for Standard 3 and Standard 4 must not vary by more than 0.040 optical density units from their mean.
- The correlation coefficient (r) calculated from the mean absorbance values of the standards must be ≥ 0.98 .
- If the above criteria are not met, the run is invalid and must be repeated.
- The mean OD value for the Zero Standard (Green Diluent) should be ≤ 0.150 . If the mean OD value is > 0.150 , the plate washing procedure should be investigated.

The QFT-Plus Analysis Software calculates and reports these quality control parameters.

Each laboratory should determine appropriate types of control materials and frequency of testing in accordance with Local, State, Federal, or other applicable accrediting organizations. External quality assessment and alternative validation procedures should be considered.

Note: Plasmas spiked with recombinant IFN- γ have shown reductions of up to 50% in concentration when stored at either 2–8°C and –20°C. Recombinant IFN- γ is not recommended for establishing control standards.

Interpretation of results

Important: Diagnosing or excluding tuberculosis disease, and assessing the probability of LTBI, requires a combination of epidemiological, historical, medical, and diagnostic findings that should be taken into account when interpreting QFT-Plus results. See general guidance on the diagnosis and treatment of TB disease and LTBI at [Tuberculosis \(TB\)](#) provided under section “US Centers for Disease Control and Prevention (CDC) Guidelines” below.

QFT-Plus results are interpreted using the following criteria (Table 5).

Table 5. Interpretation of QFT-Plus test results

Nil (IU/mL)	TB1 minus Nil (IU/mL)	TB2 minus Nil (IU/mL)	Mitogen minus Nil (IU/mL)*	QFT-Plus Result	Report/interpretation
≤8.0	≥0.35 and ≥25% of Nil	Any	Any	Positive†	<i>M. tuberculosis</i> infection likely
	Any	≥0.35 and ≥25% of Nil			
	<0.35 or ≥0.35 and <25% of Nil	<0.35 or ≥0.35 and <25% of Nil	≥0.50	Negative	<i>M. tuberculosis</i> infection NOT likely
>8.0§	<0.35 or ≥0.35 and <25% of Nil	<0.35 or ≥0.35 and <25% of Nil	<0.50	Indeterminate‡	Likelihood of <i>M. tuberculosis</i> infection cannot be determined
	Any	Any	Any		

* Responses to the Mitogen positive control (and occasionally TB Antigens) can be outside the range of the microplate reader. This has no impact on test results. Values >10 IU/mL are reported by the QFT-Plus software as >10 IU/mL.

† Where *M. tuberculosis* infection is not suspected, initially positive results can be confirmed by retesting the original plasma samples in duplicate in the QFT-Plus ELISA. If repeat testing of one or both replicates is positive, the test result is considered positive.

‡ Refer to "Troubleshooting Guide" on page 50 for possible causes.

§ In clinical studies, less than 0.25% of subjects had IFN-γ levels of >8.0 IU/mL for the Nil value.

The magnitude of the measured IFN-γ level cannot be correlated to stage or degree of infection, level of immune responsiveness, or likelihood for progression to active disease. A positive TB response in persons who are negative to Mitogen is rare but has been seen in patients with TB disease. This indicates the IFN-γ response to TB antigens is greater than that to Mitogen, which is possible as the level of Mitogen does not maximally stimulate IFN-γ production by lymphocytes.




Figure 3. QFT-Plus test interpretation. *For TB1 minus Nil or TB2 minus Nil value to be valid, the $\geq 25\%$ of Nil IU/mL value must be from the same tube as the original ≥ 0.35 IU/mL result.

Limitations

Results from QFT-Plus testing must be used in conjunction with each individual's epidemiological history, current medical status, and other diagnostic evaluations.

Individuals with Nil values greater than 8 IU/mL are classed as "Indeterminate" because a 25% higher response to TB Antigens may be outside the assay measurement range.

- The predictive value of a positive QFT-Plus result in diagnosing *M. tuberculosis* infection depends on the probability of infection, which is assessed by historical, epidemiological, diagnostic, and other findings.
- A diagnosis of LTBI requires that tuberculosis disease must be excluded by medical evaluation including an assessment of current medical and diagnostic tests for disease as indicated.
- A negative result must be considered with the individual's medical and historical data relevant to probability of *M. tuberculosis* infection and potential risk of progression to tuberculosis disease, particularly for individuals with impaired immune function. Negative predictive values are likely to be low for persons suspected to have *M. tuberculosis* disease and should not be relied on to exclude disease.

Unreliable or indeterminate results may occur due to:

- Deviations from the procedure described in the package insert
- Incorrect transport/handling of blood specimen
- Elevated levels of circulating IFN- γ or presence of heterophile antibodies
 - Exceeding validated blood times from blood specimen draw to incubation. Refer to the *QuantiFERON-TB Gold Plus (QFT-Plus) Blood Collection Tubes Instructions for Use* (HB-3304).

- Blood samples collected directly into QFT-Plus Blood Collection Tubes stored longer than 16 hours at room temperature (17–25°C).
- Blood samples collected in lithium or sodium-heparin tube stored longer than 12 hours at room temperature (17–25°C) prior to transfer to QFT-Plus Blood Collection Tubes.
- Blood samples collected in lithium or sodium- heparin tube for refrigeration stored outside temperature and time ranges.

Performance Characteristics

Clinical studies

As there is no definitive standard test for confirming or excluding the diagnosis of LTBI, an estimate of sensitivity and specificity for QFT-Plus cannot be practically evaluated. Specificity of QFT-Plus was approximated by evaluating false-positive rates in persons with low risk (no known risk factors) of tuberculosis infection. Sensitivity was approximated by evaluating groups of study subjects with culture-confirmed active TB disease. In addition, assay performance was evaluated for positive and negative rate in a population of healthy subjects with identified risk factors for tuberculosis infection (a mixed-risk population).

Specificity

A multi-center study evaluating the clinical specificity of QFT-Plus was performed including 733 study subjects who were considered to have either low risk of *M. tuberculosis* infection or no risk factors for exposure to infection or disease. Demographic information and risk factors for TB exposure were determined using a standardized survey at the time of testing. The study was conducted at four independent sites, including one in the United States, two in Japan, and one in Australia. The QFT-Plus test was compared to the QuantiFERON-TB Gold-In-Tube (QFT-Gold) test. A summary of the clinical specificity performance data, stratified by study site and region is provided in Table 6. The performance results are based on the total number of valid tests. There were no indeterminate results.

Table 6. Interpretation of QFT-Plus specificity clinical results

Site	N	Positive		Negative		Indeterminate		Specificity (95% CI)	
		QFT-Gold	QFT-Plus	QFT-Gold	QFT-Plus	QFT-Gold	QFT-Plus	QFT-Gold	QFT-Plus
United States									
(#1) USA-4	212	2	4	210	208	0	0	99.06% (210/212) (96.63–99.74)	98.11% (208/212) (95.25–99.26)
Japan									
(#2) JPN-3	106	1	2	105	104	0	0	99.06% (105/106) (94.85–99.83)	98.11% (104/106) (93.38–99.48)
(#3) JPN-1	216	3	5	213	211	0	0	98.61% (213/216) (96.00–99.53)	97.69% (211/216) (94.70–99.01)
Total Japan	322	4	7	318	315	0	0	98.76% (318/322) (96.85–99.52)	97.83% (315/322) (95.6–98.9)
Australia									
(#4) AU-3	199	8	9	191	190	0	0	95.98% (191/199) (92.27–97.95)	95.48% (190/199) (91.63–97.60)

The overall specificity of QFT-Plus was 97.27% (713/733). The specificity of QFT-Plus in the United States was 98.11%, 97.83% in Japan and 95.48% in Australia. The specificity of QFT-Gold in the United States region was 99.06%, 98.76% in Japan and 95.98% in Australia. The overall specificity of QFT-Gold was 98.09% (719/733).

A breakdown of the results by TB antigen tube type and combinations thereof is shown to provide an example of expected results in a low-risk population (Table 7).

Table 7. QFT-Plus specificity study results by TB antigen tube

Interpretation based on TB Antigen-Nil

IU/mL in	TB1	TB2	QFT-Plus (positive by TB1 and/or TB2)*	Concordant positive TB1 and TB2 (alternate analysis)†
Positive	10	18	20	8
Negative	723	715	713	725
Indeterminate	0	0	0	0
Specificity (95% CI)	–	–	97.3% (713/733) (95.8–98.2)	–
Negativity rate (95% CI)	98.6% (723/733) (97.5–99.3)	97.5% (715/733) (96.2–98.4)	–	98.9% (725/733) (97.9–99.5)

* Interpretation based on a TB antigen – Nil value ≥ 0.35 IU/mL in both (TB1 and TB2) or either TB tube to fit the interpretation criteria for the QFT-Plus (TB1 or TB2) to be determined positive.

† Alternate analysis provided for information only.

In subjects with low risk for TB infection, a total of 20/733 subjects returned a positive result. Of these, only 8 subjects returned a value of >0.35 IU/mL in both TB1 and TB2 tubes. A comparison of the QFT-Gold versus QFT-Plus assays was performed in the low-risk study cohort and showed an overall percent agreement of 97.5% (715/733), and a negative percent agreement of 98.3% (707/719).

Sensitivity

While there is no definitive standard test for LTBI, a suitable surrogate is microbiological culture of *M. tuberculosis* because infection with TB is a necessary precursor to disease.

A multi-center study evaluating the clinical sensitivity of QFT-Plus was performed including 434 study subjects who presented with signs and symptoms of active *M. tuberculosis* disease confirmed by culture and/or PCR and were not receiving TB treatment or with ≤ 14 days of treatment prior to blood collection. The study was performed at 7 independent sites including three sites in the United States, three sites in Japan, and one site in Australia. The QFT-Plus test was compared to the QFT-Gold test. A summary of the clinical sensitivity performance data, stratified by study site and region is provided in Table 8. The performance results are based on the total number of valid tests. The frequency of indeterminate results for QFT-Gold and QFT-Plus was 2.3% (10/434) and 2.5% (11/434), respectively.

Table 8. Clinical sensitivity study performance summary stratified by site

Site	N	Positive		Negative		Indeterminate		Sensitivity (95% CI)	
		QFT-Gold	QFT-Plus	QFT-Gold	QFT-Plus	QFT-Gold	QFT-Plus	QFT-Gold	QFT-Plus
United States									
(#1) USA-5	15	13	13	2	2	0	0	86.67% (13/15) (62.12–96.26)	86.67% (13/15) (62.12–96.26)
(#2) USA-1	33	29	29	4	4	0	0	87.88% (29/33) (72.67–95.18)	87.88% (29/33) (72.67–95.18)
(#3) USA-4	5	5	5	0	0	0	0	100.0% (5/5) (56.55–100.0)	100.0% (5/5) (56.55–100.0)
Total United States	53	47	47	6	6	0	0	88.7% (47/53) (77.4–94.7)	88.7% (47/53) (77.4–94.7)
Japan									
(#4) JPN-2	76	72	67	1	3	3	6	98.63% (72/73) (92.64–99.76)	95.71% (67/70) (88.14–98.53)
(#5) JPN-3	99	97	98	2	1	0	0	97.98% (97/99) (92.93–99.44)	98.99% (98/99) (94.50–99.82)
(#6) JPN-1	177	159	157	12	15	6	5	92.98% (159/171) (88.14–95.94)	91.28% (157/172) (86.11–94.64)
Total Japan	352	328	322	15	19	9	11	95.63% (328/343) (92.91–97.33)	94.43% (322/341) (91.5–96.4)
Australia									

Table 8. Clinical sensitivity study performance summary stratified by site (continued)

Site	N	Positive		Negative		Indeterminate		Sensitivity (95% CI)	
		QFT-Gold	QFT-Plus	QFT-Gold	QFT-Plus	QFT-Gold	QFT-Plus	QFT-Gold	QFT-Plus
(#7) AU-2	29	27	29	1	0	1	0	96.43% (27/28) (82.29–99.37)	100.0% (29/29) (88.30–100.0)

The analysis above in the table does not include indeterminate results.

The overall sensitivity of QFT-Plus was 94.09% (398/423). The sensitivity of QFT-Plus in the United States was 88.7%, 94.43% in Japan, and 100.0% in Australia. The sensitivity of QFT-Gold in the United States was 88.7%, 95.63% in Japan, and 96.43% in Australia. The overall sensitivity of QFT-Gold was 94.81% (402/424).

A breakdown of the results by TB antigen tube type and combinations of tubes is shown to provide an example of expected results in a confirmed TB infected population (Table 9).

Table 9. QFT-Plus sensitivity study results by TB antigen tube

Interpretation based on TB Antigen-Nil

IU/mL in	TB1	TB2	QFT-Plus (positive by TB1 and/or TB2)
Positive	388	397	398
Negative	32	26	25
Indeterminate	14	11	11
Sensitivity* (95% CI)	–	–	94% (398/423) (91.4–96.0)
Positivity rate* (95% CI)	92.4% (388/420) (89.4–94.6)	93.9% (397/423) (91.1–95.8)	–

* Excluding indeterminate values.

A comparison of the QFT-Gold and QFT-Plus assays was assessed in a culture confirmed active TB cohort (sensitivity study cohorts) and showed an overall percent agreement of 95.9% and a positive percent agreement of 97.3% (391/402).

Performance in subjects with identified risk factors for a *M. tuberculosis* infection (mixed-risk individuals)

A cohort of 601 individuals with mixed risk factors for TB infection (e.g., HIV positivity, history of treatment for active or latent TB, exposure to active TB case, HCW status, etc.) was assessed with both QFT-Gold and QFT-Plus tests. Risk factors were identified using a standardized survey and individuals displayed no symptoms associated with active TB at the time of recruitment. Demographics and risk factors are reported in Table 10. In this population, 68/601 (11.3%) subjects returned a positive QFT-Plus result, with a positive percent agreement of 98.44% and negative percent agreement of 99.07% (Table 11). In this cohort of 68 QFT-Plus positive subjects, a total of 62 subjects were positive by both TB1 and TB2 tubes, 2 subjects were positive by TB1 only, and 4 subjects were positive by TB2 only. No indeterminate results (0/601) were observed.

Table 10. Demographics and factors associated with risk of TB infection in a mixed cohort

Total subjects (601)		Number	Percentage (%)
Gender	Male	539	89.7
	Female	62	10.3
Age (years)	Range	18–70	–
	Mean	46.7	
BCG vaccinated	Yes	15	2.5
	No	586	97.5
HIV positive or tested positive for HTLV viruses	Yes	12	2.0
	No	589	98

Table 10. Demographics and factors associated with risk of TB infection in a mixed cohort (continued)

Total subjects (601)		Number	Percentage (%)
Previously diagnosed with active TB	Yes	11	1.8
	No	590	98.2
Had a positive Tuberculin Skin Test (TST)/Mantoux test for TB	Yes	47	7.8
	No	554	92.2
Ever been treated for active or latent TB	Yes	35	5.8
	No	566	94.2
Lived, worked or volunteered (>1 month) in a jail or prison	Yes	373	62.1
	No	228	37.9
Lived, worked or volunteered (>1 month) in a homeless shelter	Yes	525	87.4
	No	76	12.6
Healthcare worker	Yes	8	1.3
	No	593	98.7
Close contact of someone with or suspected of having active TB disease	Yes	9	1.5
	No	592	98.5

Table 11. Summary performance of QFT-Plus versus QFT-Gold in subjects with known risk factors for latent TB infection

		QFT-Gold		
		Positive (+)	Negative (-)	Total
QFT-Plus	Positive (+)	63	5*	68
	Negative (-)	1*	532	533
	Total	64	537	601

*All 6 discordant samples had IFN- γ levels of the TB Antigen tubes that were close to the assay cut-off.

The positive percent agreement (PPA) and negative percent agreement (NPA) between the results of QFT-Gold and the QFT-Plus were as follows:

- PPA: 98.44% (63/64), 95%CI (91.67, 99.72)
- NPA: 99.07% (532/537), 95% CI (97.84, 99.60)

Table 12 below illustrates the performance of QFT-Plus as compared to the QFT-Gold test in BCG vaccinated study subjects.

Table 12. Performance of QFT-Plus as compared to the QFT-Gold Test in BCG-vaccinated study subjects (combined data from sensitivity, specificity, and LTBI study subjects)

		QFT-Gold		
		Positive (+)	Negative (-)	Total
QFT-Plus	Positive (+)	69	4	73
	Negative (-)	3	267	270
	Total	72	271	343*

*Two Sensitivity Study Subjects were excluded from the analysis due to indeterminate results.

- PPA = 95.8% (69/72), 95%CI (88.5, 98.6)
- NPA = 98.5% (267/271), 95%CI (96.3, 99.4)

Expected Values

Observed response distribution - risk stratified

A range of IFN- γ responses to TB1, TB2, and control tubes (Nil and Mitogen) were observed in clinical trials and stratified by risk of *M. tuberculosis* infection (Figure 4 through Figure 7). The mixed risk group consists of subjects representative of a general testing population, including subjects with and without risk factors for TB exposure, and where active TB is unlikely (i.e., LTBI).




Figure 4. Distribution of Nil. A) Distribution of Nil values in a low-risk population (n=733). B) Distribution of Nil values in a mixed-risk population (n=601). C) Distribution of Nil values in a population with culture-confirmed *M. tuberculosis* infection (n=434).




Figure 5. Distribution of TB1 and TB2 (Nil subtracted). **A)** Distribution of TB1 and TB2 (Nil subtracted) values in a low-risk population (n=733). **B)** Distribution of TB1 and TB2 (Nil subtracted) values in a mixed-risk population (n=601). **C)** Distribution of TB1 and TB2 (Nil subtracted) values in a population with culture-confirmed *M. tuberculosis* infection (n=434).




Figure 6. Distribution of Mitogen (Nil subtracted). **A)** Distribution of Mitogen (Nil subtracted) values in a low-risk population (n=733). **B)** Distribution of Mitogen (Nil subtracted) values in a mixed-risk population (n=601). **C)** Distribution of Mitogen (Nil subtracted) values in a population with culture-confirmed *M. tuberculosis* infection (n=434).




Figure 7. Observed difference between TB1 and TB2 values (Nil subtracted), stratified by risk. Includes data from the Mixed-risk cohort study to show differences between low risk, mixed-risk, and active TB cohorts. This data analysis included a mixed-risk cohort with known risk factors. From the low risk cohort n=733, from the mixed-risk cohort n=588 and from the active TB cohort n=357. The quantitative difference in IU/mL for each subject was obtained by subtracting the TB1 value from the TB2 value.

Assay Performance Characteristics

A study was conducted to assess the linearity of the QFT-Plus ELISA using a 4-point standard curve. Two panels each of 11 contrived plasma samples were tested. The contrived samples were prepared by spiking negative plasma samples with either a 10 IU/mL IFN- γ concentrated plasma sample (High Pool), or a 1.5 IU/mL IFN- γ concentrated plasma sample (Low Pool) in a dilution series.

Weighted linear regression analysis of the calculated mean across replicates versus expected values (based on the dilution factor) was performed for both High and Low Pools.




Figure 8. Illustration of Linearity Study Regression Analysis – High Pool Mean = $-0.24 + 0.9964 \cdot \text{Expected}$.

Analysis of samples of the High Pool demonstrated linearity across the range of 0.79 IU/mL to 10 IU/mL with deviation from linearity $\leq 6.1\%$.




Figure 9. Illustration of Linearity Study Regression Analysis – Low Pool Mean = 0.016 + 1.003•Expected.

Analysis of samples of the Low Pool demonstrated linearity across the range of 0.17 IU/mL to 1.52 IU/mL with deviation from linearity $\leq 6.2\%$.

The combined study data demonstrated linearity across the range of 0.17 IU/mL to 10 IU/mL with deviation $\leq 6.2\%$. Deviation from linearity using the 4-point standard curve was observed to be less than 7%.

A multi-center study reproducibility study was conducted to evaluate performance of QFT-Plus across study sites with multiple operators. This was a prospective study conducted at three external testing sites and one collection site. A total of 32 positive and 34 negative (determined by the QFT-Gold test) study subjects were enrolled. The study subjects were comprised of healthcare workers in the United States. The study subjects represented group

with mixed risk for TB exposure due to their occupation or as foreign born healthcare workers originating from a location with at TB rate exceeding 50/100,000.

Three lithium heparin (Li Hep) blood collection tubes were obtained from each study subject at the collection site. The Li Hep blood collection tubes were then transferred to each of three testing sites where they were aliquoted into two sets of QFT-Plus Blood Collection Tubes (QFT-Plus TB1, TB2, Mitogen, and Nil) then tested in accordance with the QFT-Plus assay procedure. At each site at least two operators ran the two tests per study subject independently. Each operator was blinded to the results obtained by the other operator and blinded to the QFT-Gold test result of the study subject.

There were six results generated across all three testing sites per each of the 66 study subjects, resulting in a total of 396 data points. A summary of the reproducibility summary results are provided in Table 13.

Table 13. Multi-site reproducibility study results summary samples

Site 1 – 2 Operators	Site 2 – 2 Operators	Site 3 – 3 Operators
64/66 = 96.97%	64/66 = 96.97%	59/66 = 89.39%
Agreement of Qualitative Results of Tube Set 1 and Tube Set 2	Agreement of Qualitative Results of Tube Set 1 and Tube Set 2	Agreement of Qualitative Results of Tube Set 1 and Tube Set 2

The qualitative overall percent agreement across the three study sites was 94.7% (375/396). In this calculation, the total number of test results in agreement (375) included those instances where there was agreement of all 6 results, agreement of 5 out of 6 results, agreement of 4 out of 6 results, and agreement of 3 out of 6 results, combined.

Technical Information

Indeterminate results

Indeterminate results are uncommon and may relate to the immune status of the individual being tested (5) but may also be related to a number of technical factors (e.g., inappropriate handling/storage of blood collection tubes, incomplete ELISA plate washing) if the above instructions for use are not followed.

If technical issues are suspected with the reagent storage, blood collection, or handling of the blood samples, repeat the entire QFT-Plus test with new blood specimens. Repeating the ELISA testing using the stimulated plasmas can be performed if inadequate washing or other procedural deviation with the ELISA test is suspected. Physicians may choose to redraw a specimen or perform other procedures as appropriate.

Clotted plasma samples

Should fibrin clots occur with long-term storage of plasma samples, centrifuge samples to sediment clotted material and facilitate pipetting of plasma.

Lipemic plasma samples

Care should be exercised when pipetting lipemic samples as fatty deposits can block pipet tips.

Troubleshooting Guide

This troubleshooting guide may be helpful in solving any problems that may arise. For more information, see also the Frequently Asked Questions page at our Technical Support Center: www.qiagen.com/FAQ/FAQList.aspx (for contact information, visit www.qiagen.com).

Comments and suggestions

ELISA troubleshooting

Non-specific color development

- | | |
|--|---|
| a. Incomplete washing of the plate | Wash the plate at least 6 times with 400 μ L/well of wash buffer. More than 6 washing cycles may be required depending on the washer being used. A soak time of at least 5 seconds between cycles should be used. |
| b. Cross-contamination of ELISA wells | Take care while pipetting and mixing sample to minimize risk. |
| c. Kit/components have expired | Ensure kit is used before the expiry date. Ensure reconstituted standard and Conjugate 100x Concentrate are used within three months of the reconstitution date. |
| d. Enzyme Substrate Solution is contaminated | Discard substrate if blue coloration exists. Ensure clean reagent reservoirs are used. |
| e. Mixing of plasma in QFT-Plus Blood Collection Tubes before harvesting | After centrifugation, avoid pipetting up and down or mixing plasma by any means prior to harvesting. At all times, take care not to disturb material on the surface of the gel. |

Low optical density readings for standards

- | | |
|-----------------------------------|---|
| a. Standard dilution error | Ensure dilutions of the Kit Standard are prepared correctly as per this Package Insert. |
| b. Pipetting error | Ensure pipets are calibrated and used according to manufacturer's instructions. |
| c. Incubation temperature too low | Incubation of the ELISA should be performed at room temperature ($22 \pm 5^{\circ}\text{C}$). |

Comments and suggestions

- | | |
|---------------------------------------|---|
| d. Incubation time too short | Incubation of the plate with the conjugate, standards and samples should be for 120 ± 5 minutes. The Enzyme Substrate Solution should be incubated on the plate for 30 minutes. |
| e. Incorrect plate reader filter used | Plate should be read at 450 nm with a reference filter of between 620 and 650 nm. |
| f. Reagents are too cold | All reagents, with the exception of the Conjugate 100x Concentrate, must be brought to room temperature prior to commencing the assay. This takes approximately 1 hour. |
| g. Kit/components have expired | Ensure that the kit is used before the expiry date. Ensure reconstituted Standard and Conjugate 100x Concentrate are used within 3 months of the reconstitution date. |

High background

- | | |
|--|--|
| a. Incomplete washing of the plate | Wash the plate at least 6 times with 400 μ L/well of wash buffer. More than 6 washing cycles may be required. A soak time of at least 5 seconds between cycles should be used. |
| b. Incubation temperature too high | Incubation of the ELISA should be performed at room temperature ($22 \pm 5^\circ\text{C}$). |
| c. Kit/components have expired | Ensure that the kit is used within the expiry date. Ensure reconstituted standard and Conjugate 100x Concentrate are used within three months of the reconstitution date. |
| d. Enzyme Substrate Solution is contaminated | Discard substrate if blue coloration exists. Ensure clean reagent reservoirs are used. |

Nonlinear standard curve and duplicate variability

- | | |
|--|--|
| a. Incomplete washing of the plate | Wash the plate at least 6 times with 400 μ L/well of wash buffer. More than 6 washing cycles may be required. A soak time of at least 5 seconds between cycles should be used. |
| b. Standard dilution error | Ensure dilutions of the standard are prepared correctly as per this Instructions for Use. |
| c. Poor mixing | Mix reagents thoroughly by inversion or gentle vortexing prior to their addition to the plate. |
| d. Inconsistent pipetting technique or interruption during assay setup | Sample and standard addition should be performed in a continuous manner. All reagents should be prepared prior to commencing the assay. |

Appendices

Appendix A: US Centers for Disease Control and Prevention (CDC) Guidelines

Diagnostic testing for *Mycobacterium tuberculosis* using Interferon Gamma Release Assays should follow applicable published guidelines, including when testing in populations such as children, pregnant women, and HIV-infected or otherwise immunocompromised individuals. (Also see "Warnings and Precautions" on page 12 regarding use in co-morbid conditions which affect immune function). Links to the most recent American Thoracic Society/Infectious Diseases Society of America/Centers for Disease Control and Prevention Clinical Practice Guidelines: Diagnosis of Tuberculosis in Adults and Children, as well as to other information regarding diagnostic testing for tuberculosis, are available at www.cdc.gov/tb/testing/index.html.

Appendix B: Abbreviated ELISA Test Procedure

1. Equilibrate ELISA components, with the exception of the Conjugate 100x Concentrate, to room temperature ($22 \pm 5^\circ\text{C}$) for at least 60 minutes.



2. Reconstitute the Kit Standard to 8.0 IU/mL with distilled or deionized water. Prepare four (4) standard dilutions.



3. Reconstitute freeze-dried Conjugate 100x Concentrate with distilled or deionized water.



4. Prepare working strength conjugate in Green Diluent and add 50 μL to all wells.


5. Add 50 μL of test plasma samples and 50 μL standards to appropriate wells. Mix using shaker.



6. Incubate for 120 minutes at room temperature ($22 \pm 5^\circ\text{C}$).



7. Wash wells at least 6 times with 400 μL /well of wash buffer.



8. Add 100 μL Enzyme Substrate Solution to wells. Mix using shaker.




9. Incubate for 30 minutes at room temperature ($22 \pm 5^\circ\text{C}$).



10. Add 50 μL Enzyme Stopping Solution to all wells. Mix using shaker.



11. Read results at 450 nm with a 620 to 650 nm reference filter.



12. Analyze results.



References

1. Andersen, P. et al. (2000) Specific immune-based diagnosis of tuberculosis. *Lancet* 356, 1099.
2. Balcells, M.E. et al. (2008) A comparative study of two different methods for the detection of latent tuberculosis in HIV-positive individuals in Chile. *Int. J. Infect. Dis.* 12, 645.
3. Bartalesi, F. et al. (2009) QuantiFERON-TB Gold and TST are both useful for latent TB screening in autoimmune diseases. *Eur. Respir. J.* 33, 586.
4. Bocchino, M. et al. (2008) Performance of two commercial blood IFN-gamma release assays for the detection of *Mycobacterium tuberculosis* infection in patient candidates for anti-TNF-alpha treatment. *Eur. J. Clin. Microbiol. Infect. Dis.* 27, 907.
5. Brock, I. et al. (2006) Latent tuberculosis in HIV positive, diagnosed by the *M. tuberculosis* specific interferon-gamma test. *Respir. Res.* 7, 56.
6. Chun, J.K. et al. (2008) The role of a whole blood interferon gamma assay for the detection of latent tuberculosis infection in Bacille Calmette-Guerin vaccinated children. *Diagn. Microbiol. Infect. Dis.* 62, 389.
7. Connell, T.G. et al. (2008) A three-way comparison of tuberculin skin testing, QuantiFERON-TB Gold and T-SPOT.TB in children. *PLoS ONE* 3, e2624.
8. Detjen, A.K. et al. (2007) Interferon-gamma release assays improve the diagnosis of tuberculosis and nontuberculous mycobacterial disease in children in a country with a low incidence of tuberculosis. *Clin. Infect. Dis.* 45, 322.
9. Diel, R. et al. (2009) Comparative performance of tuberculin skin test, QuantiFERON-TB Gold In Tube assay, and T-Spot.TB test in contact investigations for tuberculosis. *Chest* 135, 1010.

10. Diel, R. et al. (2008) Predictive value of a whole-blood IFN- γ assay for the development of active TB disease. *Am. J. Respir. Crit. Care Med.* 177, 1164.
11. Diel, R. et al. (2006) Tuberculosis contact investigation with a new, specific blood test in a low-incidence population containing a high proportion of BCG-vaccinated persons. *Respir. Res.* 7, 77.
12. Dogra, S. et al. (2007) Comparison of a whole blood interferon-gamma assay with tuberculin skin testing for the detection of tuberculosis infection in hospitalized children in rural India. *J. Infect.* 54, 267.
13. Drobniewski, F. et al. (2007) Rates of latent tuberculosis in health care staff in Russia. *PLoS Med.* 4, e55.
14. Gerogianni, I. et al. (2008) Whole-blood interferon-gamma assay for the diagnosis of tuberculosis infection in an unselected Greek population. *Respirology* 13, 270.
15. Harada, N. et al. (2008) Comparison of the sensitivity and specificity of two whole blood interferon-gamma assays for *M. tuberculosis* infection. *J. Infect.* 56, 348.
16. Higuchi, K. et al. (2009) Comparison of performance in two diagnostic methods for tuberculosis infection. *Med. Microbiol. Immunol.* 198, 33.
17. Kang, Y.A. et al. (2005) Discrepancy between the tuberculin skin test and the whole-blood interferon gamma assay for the diagnosis of latent tuberculosis infection in an intermediate tuberculosis-burden country. *JAMA* 293, 2756.
18. Katiyar, S. K. et al. (2008) Use of the QuantiFERON-TB Gold In-Tube test to monitor treatment efficacy in active pulmonary tuberculosis. *Int. J. Tuberc. Lung Dis.* 12, 1146.
19. Kipfer, B. et al. (2008) Tuberculosis in a Swiss army training camp: contact investigation using an Interferon gamma release assay. *Swiss Med. Wkly.* 138, 267.
20. Luetkemeyer, A. et al. (2007) Comparison of an interferon-gamma release assay with tuberculin skin testing in HIV-infected individuals. *Am. J. Respir. Crit. Care Med.* 175,

737.

21. Mackensen, F. et al. (2008) QuantiFERON-TB Gold – a new test strengthening long-suspected tuberculous involvement in serpiginous-like choroiditis. *Am. J. Ophthalmol.* 146, 761.
22. Manuel, O. et al. (2007) Comparison of QuantiFERON-TB Gold with tuberculin skin test for detecting latent tuberculosis infection prior to liver transplantation. *Am. J. Transplant.* 7, 2797.
23. Matulis, G. et al. (2007) Detection of latent tuberculosis in immunosuppressed patients with autoimmune diseases performance of a Mycobacterium tuberculosis antigen specific IFN-gamma assay. *Ann. Rheum. Dis.* 67, 84.
24. Mirtskhulava, V. et al. (2008) Prevalence and risk factors for latent tuberculosis infection among health care workers in Georgia. *Int. J. Tuberc. Lung Dis.* 12, 513.
25. Nakaoka, H. et al. (2006) Risk for tuberculosis among children. *Emerg. Infect. Dis.* 12, 1383.
26. Pai, M. et al. (2005) Mycobacterium tuberculosis infection in health care workers in rural India: comparison of a whole-blood, interferon-gamma assay with tuberculin skin testing. *JAMA* 293, 2746.
27. Ponce, de Leon, D. et al. (2008) Comparison of an interferon-gamma assay with tuberculin skin testing for detection of tuberculosis (TB) infection in patients with rheumatoid arthritis in a TB-endemic population. *J. Rheumatol.* 35, 776.
28. Richeldi, L. et al. (2008) Prior tuberculin skin testing does not boost QuantiFERON-TB results in paediatric contacts. *Eur. Respir. J.* 32, 524.
29. Rothel, J.S., Andersen, P. (2005) Diagnosis of latent Mycobacterium tuberculosis infection: Is the demise of the Mantoux test imminent? *Expert Rev. Anti Infect. Ther.* 3, 981.

30. Schoepfer, A.M. et al. (2008) Comparison of interferon-gamma release assay versus tuberculin skin test for tuberculosis screening in inflammatory bowel disease. *Am. J. Gastroenterol.* 103, 2799.
31. Silverman, M.S. et al. (2007) Use of an interferon-gamma based assay to assess bladder cancer patients treated with intravesical BCG and exposed to tuberculosis. *Clin. Biochem.* 40, 913.
32. Stebler, A. et al. (2008) Whole-blood interferon-gamma release assay for baseline tuberculosis screening of healthcare workers at a Swiss university hospital. *Infect. Control Hosp. Epidemiol.* 29, 681.
33. Turner, J. et al. (1996) Stimulation of human peripheral blood mononuclear cells with live *Mycobacterium bovis* BCG activates cytolytic CD8⁺ T cells in vitro. *Immunology* 87, 339.
34. Brookes, R.H. et al. (2003) CD8⁺ T cell-mediated suppression of intracellular *Mycobacterium tuberculosis* growth in activated human macrophages. *Eur. J. Immunol.* 33, 3293.
35. Stenger, S. et al. (1998) An antimicrobial activity of cytolytic T cells mediated by granulysin. *Science* 282, 121.
36. Lalvani, A. et al. (1998) Human cytolytic and interferon gamma-secreting CD8⁺ T lymphocytes specific for *Mycobacterium tuberculosis*. *Proc. Natl. Acad. Sci. U.S.A.* 95, 270.
37. Lewinsohn, D.M. et al. (2001) Classically restricted human CD8⁺ T lymphocytes derived from *Mycobacterium tuberculosis*-infected cells: definition of antigenic specificity. *J. Immunol.* 166, 439.
38. Lewinsohn, D.A. et al. (2007) Immunodominant tuberculosis CD8 antigens preferentially restricted by HLA-B. *PLoS Pathol.* 3, 1240.

39. Barcellini, L. et al. (2016) First independent evaluation of QuantiFERON-TB Plus performance. *Eur. Respir. J.* 47, 1587.
40. Day, C.L. et al. (2011) Functional capacity of *Mycobacterium tuberculosis*-specific T cell responses in humans is associated with mycobacterial load. *J. Immunol.* 187, 2222.
41. Rozot, V. et al. (2013) *Mycobacterium tuberculosis*-specific CD8+ T cells are functionally and phenotypically different between latent infection and active disease. *Eur. J. Immunol.* 43, 1568.
42. Nikolava, M. et al. (2013) Antigen-specific CD4- and CD8-positive signatures in different phases of *Mycobacterium tuberculosis* infection. *Diagn. Microbiol. Infect. Dis.* 75, 277.
43. Chicchio, T. et al. (2014) Polyfunctional T-cells and effector memory phenotype are associated with active TB in HIV-infected patients. *J. Infect.* 69, 533.
44. Ongaya, A. et al. (2013) *Mycobacterium tuberculosis*-specific CD8+ T cell recall in convalescing TB subjects with HIV co-infection. *Tuberculosis* 93, S60.
45. Lanicioni, C. et al. (2012) CD8+ T cells provide an immunologic signature of tuberculosis in young children. *Am. J. Respir. Crit. Care Med.* 185, 206.
46. Mazurek, G.H., et al.; IGRA Expert Committee; Centers for Disease Control and Prevention (CDC). (2010) Updated guidelines for using Interferon Gamma Release Assays to detect *Mycobacterium tuberculosis* infection – United States, 2010. *MMWR Recomm. Rep.* 59, 1.
47. Chiappini, E. et al. (2014). Interferon-gamma release assay sensitivity in children younger than 5 years is insufficient to replace the use of tuberculin skin test in western countries. *Pediatric Infectious Disease Journal*, 33(12), 1291–1293.
48. Garazzino, S. et al. (2014). Performance of interferon- γ release assay for the diagnosis of active or latent tuberculosis in children in the first 2 years of age: a multicenter study of the Italian Society of Pediatric Infectious Diseases. *Pediatric Infectious Disease Journal*, 33(9).

49. Rose, W. et al. (2014). Quantiferon gold-in-tube assay for TB screening in HIV infected children: Influence of Quantitative Values. *BMC Infectious Diseases*, 14:516.
50. Pavić, I. et al. (2015). Discordance between tuberculin skin test and interferon- γ release assay in children younger than 5 years who have been vaccinated with bacillus calmette-guérin. *Laboratory Medicine*, 46(3), 200–206.
51. Howley, M. M. et al. (2015). Evaluation of Quantiferon-TB gold in-tube and tuberculin skin tests among immigrant children being screened for latent tuberculosis infection. *Pediatric Infectious Disease Journal*, 34(1), 35–39.
52. Rose, W. et al. (2014). Relating tuberculosis (TB) contact characteristics to quantiferon-tb-gold and tuberculin skin test results in the Toronto Pediatric TB Clinic. *Journal of the Pediatric Infectious Diseases Society*, 4(2), 96–103.
53. Grinsdale, J. A. et al. (2014). Interferon-gamma release assays and pediatric public health tuberculosis screening: The San Francisco Program Experience 2005 to 2008. *Journal of the Pediatric Infectious Diseases Society*, 5(2), 122–130.
54. Gudjónsdóttir, M. J. et al. (2016). Relation between BCG vaccine scar and an interferon-gamma release assay in immigrant children with “positive” tuberculin skin test (≥ 10 mm). *BMC Infectious Diseases*, 16(1):540.
55. Nejat, S. & Bennet, R. (2016). Interferon-gamma release assays can effectively screen migrants for the tuberculosis infection, but urgent, active cases need clinical recognition. *Acta Paediatrica*, 105(6), 671–675.
56. Dehority, W. et al. (2017). Comparison of the Quantiferon TB gold in-tube assay with tuberculin skin test for the diagnosis of latent tuberculosis infection among HIV-infected and uninfected children. *Pediatric Infectious Disease Journal*, 36(12): e317-e321.
57. Marino, A. et al. (2017). Prebiologic therapy tuberculosis screening experience in a pediatric rheumatology center: TST and IGRA are both necessary. *Pediatric Infectious Disease Journal*, 36(4), 440–441.










58. Petrucci, R. et al. (2017). QuantiFERON-TB gold in-tube improves tuberculosis diagnosis in children. *Pediatric Infectious Disease Journal*, 36(1), 44–49.
59. Hormi, M. et al. (2018). Performance of the quantiFERON-TB gold assay among HIV-infected children with active tuberculosis in France. *Pediatric Infectious Disease Journal*, 37(4), 339–344.
60. Kay, A. W. et al. (2018). Interferon- γ release assay performance for tuberculosis in childhood. *Pediatrics*, 141(6):e20173918.
61. Stout, J.E. et al. (2018). Tuberculosis Epidemiologic Studies Consortium. Evaluating latent tuberculosis infection diagnostics using latent class analysis. *Thorax*, 73(11):1062-1070.
62. Chiappini, E. et al. (2019). QuantiFERON-TB Gold In-Tube test performance in a large pediatric population investigated for suspected tuberculosis infection. *Paediatr Respir Rev*. 32:36-47.
63. Debulpaep, S. et al. (2019). Contribution of QuantiFERON-TB Gold-in-Tube to the Diagnosis of *Mycobacterium tuberculosis* Infection in Young Children in a Low TB Prevalence Country. *Front Pediatr*. 18;7:291.
64. Lombardi, G. et al. (2019). QuantiFERON-TB Performs Better in Children, Including Infants, than in Adults with Active Tuberculosis: a Multicenter Study. *J Clin Microbiol*. 24;57(10):e01048-19.
65. Mastrolia, M.V. et al. (2019). Utility of tuberculin skin test and IGRA for tuberculosis screening in internationally adopted children: Retrospective analysis from a single center in Florence, Italy. *Travel Med Infect Dis*. 28:64-67.
66. Venkatappa, T.K. et al. (2019). Comparing QuantiFERON-TB Gold Plus with Other Tests To Diagnose *Mycobacterium tuberculosis* Infection. *J Clin Microbiol*. 23;57(11):e00985-19.

67. Ahmed, A. et al. (2020). Tuberculosis Epidemiologic Studies Consortium. Interferon- γ Release Assays in Children <15 Years of Age. *Pediatrics*. 145(1):e20191930.
68. Velasco-Arnaiz, E. et al. (2020). Impact of Baseline Tuberculin Skin Test and Isoniazid Chemoprophylaxis on Subsequent Quantiferon-TB Gold In-Tube Performance in Young Children Assessed After Tuberculosis Contact in Catalonia. *Pediatr Infect Dis J*. 39(2):e22-e25.
69. Buonseso, D. et al. (2020). Accuracy of QuantiFERON-TB Gold Plus Test for Diagnosis of Mycobacterium tuberculosis Infection in Children. *J Clin Microbiol*. 26;58(6):e00272-20.
70. Quintana-Ortega, C. et al. (2020). False-positive Results of Quantiferon-Tb-Gold Assay in Children. *Pediatr Infect Dis J*. 39(7):620-623.
71. Worjloh, A. et al. (2011). Interferon gamma release assay compared with the tuberculin skin test for latent tuberculosis detection in pregnancy. *Obstet Gynecol*. 118(6):1363-70.
72. Lighter-Fisher, J. & Surette, A.M. (2012). Performance of an interferon-gamma release assay to diagnose latent tuberculosis infection during pregnancy. *Obstet Gynecol*. 119(6):1088-95.
73. Aichelburg, MC. et al. (2014). Reversion and conversion of interferon- γ release assay results in HIV-1-infected individuals. *J Infect Dis*. 209(5):729-33.
74. Aichelburg, MC. et al. (2014). The association of smoking with IGRA and TST results in HIV-1-infected subjects. *Int J Tuberc Lung Dis*. 18(6):709-16.
75. Pullar, ND. et al. (2014). Low prevalence of positive interferon-gamma tests in HIV-positive long-term immigrants in Norway. *Int J Tuberc Lung Dis*. 18(2):180-7.
76. Pullar, ND. et al. (2014). HIV patients with latent tuberculosis living in a low-endemic country do not develop active disease during a 2 year follow-up; a Norwegian prospective multicenter study. *BMC Infect Dis*. 14:667.

77. Sauzullo, I. et al. (2014). Interferon- γ release assay in HIV-infected patients with active tuberculosis: impact of antituberculous drugs on host immune response. *New Microbiol.* 37(2):153-61.
78. Sester, M. et al. (2014). Risk assessment of tuberculosis in immunocompromised patients. A TBNET study. *Am J Respir Crit Care Med.* 190(10):1168-76.
79. Bourgarit, A. et al. (2015). IGRAVIH Study Group. Latent Tuberculosis Infection Screening and 2-Year Outcome in Antiretroviral-Naive HIV-Infected Patients in a Low-Prevalence Country. *Ann Am Thorac Soc.* 12(8):1138-45.
80. Parrella, R. et al. (2015). Interferon gamma release assays and tuberculin skin test performance in different settings of HIV immunodeficiency. *In Vivo.* 29(1):137-40.
81. Kurti, Z. et al. (2015). Tuberculin Skin Test and Quantiferon in BCG Vaccinated, Immunosuppressed Patients with Moderate-to-Severe Inflammatory Bowel Disease. *J Gastrointest Liver Dis.* 24(4):467-72.
82. Perifanou, D. et al. (2018). Screening for Latent Tuberculosis Infection in Patients with Autoimmune Diseases Before Initiating TNF- α Inhibitors Therapy. *Mater Sociomed.* 30(1):32-37.
83. Chiacchio, T. et al. (2019). Characterization of QuantiFERON-TB-Plus results in latent tuberculosis infected patients with or without immune-mediated inflammatory diseases. *J Infect.* 79(1):15-23.
84. Diamantopoulos, P.T. et al. (2017). Evaluation of the performance of tuberculin skin test and QuantiFERON-TB gold in tube test in patients with hematologic malignancies. *Infect Dis (Lond).* 49(7):545-548.

Symbols

The following symbols appear in the instructions for use or on the packaging and labeling:

Required Symbol	Symbol title / Number	Description
	Rx Only / N/A	Federal law restricts this device to sale by or on the order of a physician*
	Manufacturer / 5.1.1	Indicates the medical device manufacturer**
	Use-by date / 5.1.4	Indicates the date after which the medical device is not to be used**
	Batch code / 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified**
	Catalog number / 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified**
	Keep away from sunlight / 5.3.2	Indicates a medical device that needs protection from light sources**
	Temperature limit / 5.3.7	Indicates the temperature limits to which the medical device can be safely exposed**
		Global Trade Item Number
	Consult instructions for use or consult electronic instructions for use / 5.4.3	Indicates the need for the user to consult the instruction for use**

**Required
Symbol**

Symbol title / Number Description



Caution / 5.4.4

Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences**



In vitro diagnostic medical device / 5.5.1

Indicates a medical device that is intended to be used as an in vitro diagnostic medical device**



Contains sufficient for <n> tests / 5.5.5

Indicates the total number of tests that can be performed with the medical device**



Unique device Identifier / 5.7.10

Indicates a carrier that contains unique device identifier information**

*Regulation: 21 CFR 809.10 (a)(4)

**Regulation: ISO 15223-1: Medical devices – Symbols to be used with information to be supplied by the manufacturer

Contact Information

For technical assistance and more information, please see our Technical Support Center at www.qiagen.com/Support, call 1-800-362-7737, or contact one of the QIAGEN Technical Service Departments or local distributors (see back cover or visit www.qiagen.com).

Document Revision History

Significant changes in this edition of the QFT-Plus Package Insert are summarized in the table below:

Revision	Description
R8, June 2022	Updated Components and Storage section to remove references to package insert (will not be online)
R9, January 2023	Added sodium-heparin tube as an alternative blood collection tube in Specimen Collection and Handling section
R10, January 2025	Clarification in the Directions for Performing ELISA section regarding pre or post centrifugation stability up to 3 days Updated centrifugation time to 5–15 minutes
R11, August 2025	Removal of blood collection tube workflow, general formatting

Limited License Agreement for QuantiFERON®-TB Gold Plus (QFT-Plus)

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this Instructions for Use and for use with components contained in the panel only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this panel with any components not included within this panel except as described in the protocols provided with the product, this Instructions for Use, and additional protocols available at www.qiagen.com. Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this panel and/or its use(s) do not infringe the rights of third-parties.
3. This panel and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the panel agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the panel and/or its components.

For updated license terms, see www.qiagen.com.

Trademarks: QIAGEN®, Sample to Insight®, QuantiFERON® (QIAGEN Group); Microsoft®, Excel® (Microsoft), ProClin® (Rohm and Haas Co.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

08/2025 HB-3303-011 © 2025 QIAGEN, all rights reserved.

