



Quick-Start Protocol

April 2025

QIAseq[®] xHYB Custom Viral Analysis

Further information


- QIAseq xHYB Microbial Panels: qiagen.com/ngs/microbial
- GeneGlobe: geneglobe.qiagen.com
- Technical assistance: support.qiagen.com

Important points before starting

- Download the custom viral panel Taxonomy List Excel[®] file from GeneGlobe[®] Custom Designer.
- Ensure that the latest version of CLC Genomics Workbench Microbial Genomics Module is installed.
- Ensure that the "Analyze QIAseq xHYB Viral Panel Data (Human Host)" workflow is installed.

Create custom taxonomic profiling index

1. Open “Download Custom Microbial Reference Database” tool through the following path:
Tool → Microbial Genomics Module → Databases → Taxonomic Analysis → Download Custom Microbial Reference Database.
2. Select the following parameters.




The screenshot shows a web interface for configuring the 'Download Custom Microbial Reference Database' tool. It contains two sections: 'Plasmid handling' and 'Skip manual selection'. In the 'Plasmid handling' section, there are three radio buttons: 'Include all' (unselected), 'Include only plasmids' (unselected), and 'Exclude all plasmids' (selected). In the 'Skip manual selection' section, there is a checked checkbox for 'Skip Database Builder', an unchecked checkbox for 'Include all annotation tracks (CDS, genes, etc.)', and a text input field for 'Minimum contig length' with the value '100'.

3. Click **Next** and then select **Build database from accessions or TaxIDs.**
4. From the “Taxonomy List” Excel file downloaded from GeneGlobe Custom Builder, copy values from “Species Taxonomy” in column B. These are species level taxonomy IDs for the targets from the custom panel.
5. Paste the “Species Taxonomy” values in the “ID matching” field.
6. Click **Next** and then select **Save.**
7. After downloading sequences, select “Create Taxonomic Profiling Index” tool using the following pathway: **Tool → Microbial Genomics Module → Databases → Taxonomic Analysis → Create Taxonomic Profiling Index.**
8. Click **Next.** Select **Save**, then **Finish.**

Create “Find Best Reference Sequence” list

1. Open the “Search for Sequences at NCBI” tool.
2. Select **Nucleotide** as the database and **Accession** as the search parameter.
3. From the “Taxonomy List” Excel file downloaded from GeneGlobe Custom Builder, copy an accession number from “Reference” in column E.

4. Paste the accession number in the search field and start search.
5. Save sequence (suggest to make a new folder to store sequences) and then continue until all sequences from the Reference column are downloaded.
6. Select all downloaded sequences, and then click **New** → **Sequence List tool**.




7. Select **Next** and then **Save**.

Viral analysis

1. Select Analyze QIAseq xHYB Viral Panel Data (Human host) workflow and select samples to be analyzed. Click **Next**.
2. Under “Specify reference data handling”, select **Use specified data elements**.



3. The analysis is the same as cataloged panels but with the following 2 changes:
 - a. In the Viral Taxonomic Profiling Index step, replace with the custom Taxonomic Profiling Index that was created under “Create Custom Taxonomic Profiling Index”.



- b. In the Viral Reference Database step, replace with the Reference Sequence List that

was created under “Create Find Best Reference Sequence List”.

Viral Reference Database

Workflow Input



Document Revision History

Date	Description
04/2025	Initial release

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

Trademarks: QIAGEN®, Sample to Insight®, QIAseq®, GeneGlobe® (QIAGEN GmbH); Excel® (Microsoft Corporation). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

04/2025 HB-3694-001 © 2025 QIAGEN, all rights reserved.