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Albumin/IgG Depletion Facilitates Analysis of Low-Abundance Proteins


Figure 1 Coomassie® stained 2-D PAGE gels showing non-depleted (upper panel) and depleted (lower panel) plasma samples.

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Albumin Depletion Facilitates Analysis of Low-Abundance Proteins



Figure 2 Coomassie stained 2D PAGE gels showing non-depleted(upper panel) and depleted (lower panel) plasma samples

Cell Compartment Kit Fractionation Procedure

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Reproducible, Efficient Separation of Marker Proteins


Figure 3 Western blots of fractionated NIH 3T3 cells. Protein (20 μ g) from the cytosolic (C), membrane (M), and nuclear (N) fractions was separated by SDS-PAGE. After Western blotting, proteins specific to each fraction were detected using **A** annexin, **B** TIM23, and **C** lamin antibodies, and an HRP-conjugated secondary antibody with chemiluminescent detection.

Mitochondrial Fractions Virtually Free of Cross-Contamination


Figure 4 Western blots of mitochondrial preparations probed with antibodies specific for organelle-/cell-compartment specific marker proteins. PM: plasma membrane; ER: endoplasmic reticulum.

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Qproteome Plasma Membrane Protein Kit

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Plasma Membrane Fractions free from Cross-Contamination



Figure 5 Plasma membrane proteins were purified from either HeLa or NIH3T3 cell cultures using the Qproteome Plasma Membrane Protein Kit. Cell lysates (CL) and elution fractions (E) were separated by SDS-PAGE and transferred to a nitrocellulose membrane by western blotting. Proteins regarded as markers for different cell compartments were detected using protein-specific antibodies and an HRP-conjugated secondary antibody with chemiluminescent detection.

Qproteome Nuclear Protein Procedure

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PhosphoProtein Purification Kit PhosphoProtein Purification Cartridge

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Highly Specific Separation of Phosphorylated Proteins


Figure 6 Protein-specific immunodetection of **A** unphosphorylated HSP-60 protein, and **B** phosphorylated p44 and p42 mitogen-activated protein kinase (MAPK) proteins. F: flow-through; E: eluate fractions. The antibody used to detect MAPK recognizes an epitope containing phosphorylated residues at Thr202 and Tyr204 in the p44 (upper band) and p42 (lower band) MAPK amino acid sequences. The absence of unphosphorylated HSP-60 in the eluate fraction and the absence of phosphorylated MAPK in the flow-through fraction demonstrate the complete separation of phosphorylated proteins using the PhosphoProtein Purification Kit.

Highly Specific Glycoprotein Fractionation Using Lectin Spin Columns


Figure 7 Glycoproteins were fractionated from serum using the different lectin spin columns in glycoprotein fractionation kits and analyzed by SDS-PAGE followed by silver staining. **A** Elution steps 1–3 and 4–6 from Total Lectin Spin Columns in the Total Glycoprotein Kit. **B** Eluted glycoproteins from ConA, GNA, and LCH Spin Columns in the Mannose Glycoprotein Kit. **C** Eluted glycoproteins from WGA, SNA, and MAL Spin Columns in the Sialic Glycoprotein Kit. **D** Eluted glycoproteins from AIL and PNA Spin Columns in the O-Glycan Glycoprotein Kit.

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Qproteome FFPE Tissue Kit Qproteome FFPE Tissue 2D PAGE Kit ^{NEW}

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Qproteome Glycoprotein Kit and Lectin Cartridges

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
Figure 8 Protein yields from FFPE tissues using the Qproteome kit are comparable to these from frozen material. Protein yields from frozen and FFPE tissues were compared by processing material in parallel, separation by SDS-PAGE, and western blotting. Controls were lysates from SW480 colorectal adenocarcinoma cells (colon), Huh7 human hepatoma cells (liver), or SH-SY5Y neuroblastoma cells (rat brain).

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