

# QIAGEN Supplementary Protocol

---


## Rotor-Gene<sup>®</sup> software setup for the Rotor-Gene SYBR<sup>®</sup> Green RT-PCR Kit

This protocol shows the necessary parameters that need to be entered into the Rotor-Gene software (version 6.0) when using the Rotor-Gene SYBR Green RT-PCR Kit.

**IMPORTANT:** Please read the *Rotor-Gene SYBR Green Handbook*, paying careful attention to the “Safety Information” section, before beginning this procedure.

### Procedure

1. Launch the software and set up the program as described in the next few steps.
2. Select “Hold” and set up the parameters for the reverse-transcription step as shown below.



3. Select "Hold 2" and set up the parameters for the initial activation step as shown below.

The run will take approximately 71 minute(s) to complete. The graph below represents the run to be performed :



Click on a cycle below to modify it :

|               |                  |
|---------------|------------------|
| Hold          | Insert after...  |
| <b>Hold 2</b> | Insert before... |
| Cycling       | Remove           |
| Melt          |                  |

Hold Temperature :  deg.

Hold Time :  mins  secs

< Back   Save Template   **Start Run**   Cancel

4. Select "Cycling" and set up the parameters for PCR cycling as shown below.

The run will take approximately 71 minute(s) to complete. The graph below represents the run to be performed :




Click on a cycle below to modify it :

|                |                  |
|----------------|------------------|
| Hold           | Insert after...  |
| Hold 2         | Insert before... |
| <b>Cycling</b> | Remove           |
| Melt           |                  |

This cycle repeats  time(s).


Click on one of the steps below to modify it, or press + or - to add and remove steps for this cycle.

|                                     |  |  |  |  |
|-------------------------------------|--|--|--|--|
| Timed Step                          |  |  |  |  |
| 60 deg.                             |  |  |  |  |
| 10 seconds                          |  |  |  |  |
| Acquiring to Cycling A              |  |  |  |  |
| on Green                            |  |  |  |  |
| <input type="checkbox"/> Long Range |  |  |  |  |
| <input type="checkbox"/> Touchdown  |  |  |  |  |



< Back   Save Template   **Start Run**   Cancel

**5. Select "Melt" and set up the parameters for melting curve analysis as shown below.**



**6. Load your PCR tubes and start the program.**

QIAGEN handbooks can be requested from QIAGEN Technical Service or your local QIAGEN distributor. Selected handbooks can be downloaded from [www.qiagen.com/literature](http://www.qiagen.com/literature).

Material safety data sheets (MSDS) for any QIAGEN product can be downloaded from [www.qiagen.com/Support/MSDS.aspx](http://www.qiagen.com/Support/MSDS.aspx).

The Rotor-Gene SYBR Green RT-PCR Kit is intended for research use. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

Trademarks: QIAGEN® (QIAGEN Group); Rotor-Gene® (Corbett Research); SYBR® (Molecular Probes, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Purchase of this product (Rotor-Gene SYBR Green RT-PCR Kit) is accompanied by a limited, non-transferable immunity from suit to use it with detection by a dsDNA-binding dye as described in U.S. Patents Nos. 5,994,056 and 6,171,785 and corresponding patent claims outside the United States for the purchaser's own internal research. No real-time apparatus or system patent rights or any other patent rights, and no right to use this product for any other purpose are conveyed expressly, by implication or by estoppel.

The purchase of this product (Rotor-Gene SYBR Green RT-PCR Kit) includes a limited, non-transferable license under U.S. Patent No. 5,871,908 and all continuations and divisionals, and corresponding claims in patents and patent applications outside the United States, owned by Roche Diagnostics GmbH, for internal research use or for non-*in vitro* diagnostics applications with authorized reagents with regard to Melting Curve Analysis. No right is conveyed, expressly, by implication or estoppel, under any other patent or patent claims owned by Roche Diagnostics GmbH, or by any other Party.

Purchase of this product (Rotor-Gene SYBR Green RT-PCR Kit) is accompanied by a limited, non-transferable license under RT and Reverse Transcription-PCR patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche Ltd to use it for the purchaser's own internal research. No real-time patent rights of any kind, no right under any other patent claims (such as apparatus or system claims), and no right to use this product for any other purpose is hereby granted expressly, by implication or by estoppel.

The QIAGEN silver logo is exclusively licensed to Corbett Research.

#### NOTICE TO PURCHASER:

The purchase of this product (Rotor-Gene Q, Rotor-Disc) includes a limited, non-transferable license to certain patents (see details below) surrounding rapid polymerase chain reaction (PCR) methods and instrumentation, the use of SYBR Green I in PCR reactions, melting curve analysis, analysis methods of DNA melting data, specifically high resolution melting (HRM) and others.

The purchase of this product (Rotor-Gene Q, Rotor-Disc) includes a limited, non-transferable license to one or more of US Patents Nos 6,787,338; 7,238,321; 7,081,226; 6,174,670; 6,245,514; 6,569,627; 6,303,305; 6,503,720; 5,871,908; 6,691,041; 7,387,887; and U.S. Patent Applications Nos. 2003-0224434 and 2006-0019253 and all continuations and divisionals, and corresponding claims in patents and patent applications outside the United States, owned by the University of Utah Research Foundation, Idaho Technology, Inc., and/or Roche Diagnostics GmbH, for internal research use or for non-*in vitro* diagnostics applications. No right is conveyed, expressly, by implication or estoppel, for any reagent or kit, or under any other patent or patent claims owned by the University of Utah Research Foundation, Idaho Technology, Inc., and/or Roche Diagnostics GmbH, or by any other Party. For information on purchasing licenses for in-vitro diagnostics applications or reagents, contact Roche Molecular Systems, 4300 Hacienda Drive, Pleasanton, CA 94588, USA.

#### ADDITIONAL NOTICE TO PURCHASER IN USA, EUROPE, AND AUSTRALIA ONLY:

Practice of the polymerase chain reaction (PCR) process may require a license depending on the reagents and methods applied. Purchase of this instrument (Rotor-Gene Q) does not itself convey to the purchaser a complete license or right to perform the PCR process. The use of this instrument (Rotor-Gene Q) with Authorized Reagents such as Authorized QIAGEN kits and assays provides a limited PCR license in accordance with the label rights accompanying such reagents.

Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA and/or Roche Molecular Systems, 4300 Hacienda Drive, Pleasanton, CA 94588, USA. No right is conveyed expressly, by implication or by estoppel under any other patent claim, such as claims to apparatus, reagents, kits, or methods such as 5' nuclease methods.

PCR107 Jan-09 © 2009 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)

Canada = 800-572-9613

Ireland = 1800 555 049

Norway = 800-18859

China = 0086-21-3865-3865

Italy = 800-787980

Singapore = 65-67775366

Denmark = 80-885945

Japan = 03-6890-7300

Spain = 91-630-7050

Australia = 1-800-243-800

Finland = 0800-914416

Korea (South) = 1544 7145

Sweden = 020-790282

Austria = 0800/281010

France = 01-60-920-930

Luxembourg = 8002 2076

Switzerland = 055-254-22-11

Belgium = 0800-79612

Germany = 02103-29-12000

Mexico = 01-800-7742-639

UK = 01293-422-911

Brazil = 0800-557779

Hong Kong = 800 933 965

The Netherlands = 0800 0229592

USA = 800-426-8157

