

December 2017

QIAGEN GeneReader® User Manual

For use with Advanced Process Flow (APF/HP) Instrument
Configuration and GeneReader™ Software version 1.6

For Research Use Only. Not for use in diagnostic procedures.

Contents

Instrument Configuration History	6
1 Introduction	8
1.1 About this user manual	8
1.2 General information	9
1.2.1 Technical assistance	9
1.2.2 Policy statement	9
1.3 Intended use of the GeneReader	9
1.4 Requirements for GeneReader users	10
2 Safety Information	11
2.1 Proper use	12
2.2 Electrical safety	14
2.3 Environment	15
2.3.1 Operating conditions	15
2.4 Chemicals	16
2.5 Waste disposal	16
2.6 Mechanical hazards	17
2.7 Maintenance safety	17
2.8 Symbols on the GeneReader	18
3 General Description	20
3.1 QIAGEN GeneReader Sample to Insight® NGS workflow	20
3.2 GeneReader principle	21
3.3 External features of the GeneReader	22
3.3.1 Status lights	22
3.3.2 Flow Cell door	22
3.3.3 Hood	23
3.3.4 Fluidic drawer	24
3.3.5 USB ports	24
3.3.6 Power switch	24
3.3.7 Power cord socket	24

3.3.8	Cooling air outlet.....	24
3.3.9	Workstation equipment	25
3.4	Internal features of the GeneReader	25
3.4.1	GeneReader Flow Cell	25
4	Installation Procedures	27
4.1	System delivery and installation.....	27
4.2	Site requirements	27
4.3	AC Power and USB cable connections	28
4.3.1	Power requirements	28
4.3.2	Grounding requirements.....	29
4.3.3	Installation of AC power cords, USB cables and workstation	29
4.3.4	Installation of AC power cords, USB cables and workstation using an optional UPS	29
4.4	Workstation requirements.....	30
4.5	Getting started	31
4.5.1	Powering ON the GeneReader and workstation.....	31
4.5.2	Software upgrade.....	32
5	Operating Procedures.....	34
6	GeneReader Software	35
6.1	Use of the GeneReader software	35
6.1.1	Start-up	35
6.1.2	Software workflow.....	35
6.1.3	User interface.....	36
6.1.4	Status of Flow Cells view	37
6.1.5	Status of Fluids view.....	38
6.1.6	Status of Configuration view	39
6.1.7	File handling.....	41
6.2	Workflow procedures	41
6.2.1	Flow Cell setup.....	41
6.2.2	Loading and running the GeneReader.....	43
6.2.3	Staggered loading of Flow Cells	43

6.2.4	Unloading reagents and Flow Cells	43
6.2.5	Run finished	43
7	Maintenance	48
7.1	Equipment and reagents to be supplied by user	48
7.2	Maintenance wash	48
7.2.1	Routine maintenance (DI water wash).....	48
7.2.2	Weekly maintenance	51
7.2.3	Preparing the GeneReader fluidics lines for long-term storage	58
7.2.4	Monthly cleaning procedure	59
7.3	General cleaning procedures.....	59
7.3.1	Cleaning agents.....	59
7.3.2	General instructions	60
7.3.3	Servicing	61
7.4	Cleaning the workstation hard disk.....	61
7.4.1	Procedure	61
7.5	Creating a service archive folder for troubleshooting	62
7.5.1	Procedure	62
8	Troubleshooting	64
8.1	Hardware and software errors	64
8.1.1	Application module	64
8.1.2	Manage Files module.....	64
8.1.3	Flow Cells.....	65
8.1.4	Run module.....	65
8.1.5	Start Run module	65
9	Technical Data.....	66
9.1	Environmental conditions – operating conditions	66
9.2	Mechanical data and hardware features	66
9.3	Workstation specifications (hardware and software)	66
9.3.1	Workstation	66
9.3.2	Software.....	67
	Appendix A.....	68

Declaration of Conformity	68
License Terms.....	69
Waste Electrical and Electronic Equipment (WEEE)	70
FCC Declaration	71
Liability Clause	72
Appendix B	73
Ordering Information.....	73
Appendix C	74
Safety Information Translated into German and French	74
Sicherheitsinformationen	74
Informations de Sécurité.....	74
Index	86

Instrument Configuration History

Document title	Document number	Date	Description of changes and specifications	Compatible software and sequencing kits
QIAGEN GeneReader User Manual for use with APF/HP Instrument Configuration	HB-2503- 001	December 2017	<p>Supports Advanced Process Flow (APF) workflow and UMI Advanced (HP) chemistry compatible with:</p> <ul style="list-style-type: none"> • GeneRead™ QIAact Actionable Insights Tumor Panel, cat. no. 181910 • GeneRead QIAact BRCA 1/2 Panel, cat. no. 181920 • GeneRead QIAact Lung DNA UMI Panel, cat. no. 181931 • GeneRead QIAact Lung Fusion UMI Panel, cat. no. 181936 • Any future GeneRead QIAact panel with UMI technology <p>Updated the following specifications:</p> <ul style="list-style-type: none"> • Higher data output, >7% • Enhanced maintenance wash procedure <p>GeneReader instruments released prior to the Advanced Process Flow instrument configuration can be upgraded to the APF/HP instrument configuration.</p>	<p>GeneReader Software version 1.6</p> <p>QIAGEN GeneRead® UMI Advanced Sequencing Q Kit (3), cat. no. 185251, which consists of:</p> <ul style="list-style-type: none"> • GeneRead UMI Advanced Sequencing Q Add-Ons (3), cat. no. 1108488 • GeneRead UMI Advanced Sequencing Q Flow Cell (3), cat. no. 1108489 • GeneRead UMI Advanced Sequencing Q Buffers (3), cat. no. 1108487 <p>GeneRead UMI Advanced Sequencing Q Wash Buffers (9), cat. no. 185905, purchased separately</p>

Document title	Document number	Date	Description of changes and specifications	Compatible software and sequencing kits
QIAGEN GeneReader User Manual for use with APF Instrument Configuration	HB-2325-001	December 2016	<p>Supports Advanced Process Flow (APF) workflow compatible with:</p> <ul style="list-style-type: none"> GeneRead™ QIAact Actionable Insights Tumor Panel, cat. no. 181910 GeneRead QIAact Lung DNA Panel, cat. no. 181930 GeneRead QIAact BRCA 1/2 DNA Panel <p>Updated the following specifications:</p> <ul style="list-style-type: none"> Higher read length, >30% Higher data output, >50% Higher sample multiplexing capacity per flow cell for the GeneRead QIAact Actionable Insights Tumor Panel (ATP) Approximately 20% less reagent consumption per cycle Equal or less sequencing time Parallel or staggered processing of up to 3 flow cells <p>GeneReader instruments released prior to the Advanced Process Flow instrument configuration can be upgraded to the APF instrument configuration</p>	<p>GeneReader Software version 1.4.0</p> <p>QIAGEN GeneRead Advanced Sequencing Q Kit (3), cat. no. 185231</p> <p>GeneRead Sequencing Buffer Q Kit (16), cat. no. 185901</p> <p>IMPORTANT: When used in combination with the QIAGEN GeneRead Advanced Sequencing Q Kit, the GeneRead Sequencing Buffer Q Kit (16) will allow the user to process 6 Flow Cells.</p>
GeneReader User Manual	HB-2023-001	December 2015	<p>Specifications:</p> <ul style="list-style-type: none"> Compatible with GeneRead QIAact Panels, Powered by QCI™, cat. no. 181910 Approximately 45 hours sequencing time Parallel or staggered processing of up to 4 Flow Cells 	<p>GeneReader Software up to version 1.1.3</p> <p>GeneRead Sequencing Q Kit (1), cat. no. 185200</p> <p>GeneRead Sequencing Q Kit (4), cat. no. 185201</p> <p>GeneRead Sequencing Buffer Q Kit, (16) cat. no. 185901</p>

1 Introduction

Thank you for choosing the GeneReader. We are confident it will become an integral part of your laboratory.

Before using the GeneReader, it is essential that you read this user manual carefully and pay particular attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition.

1.1 About this user manual

This user manual provides information about the GeneReader in the following sections:

- Instrument Configuration History
- Introduction
- Safety Information
- General Description
- Installation Procedures
- Operating Procedures
- Maintenance
- Troubleshooting
- Technical Data
- Appendices
- Index

The appendices contain the following information:

- Declaration of Conformity
- License Terms
- Waste Electrical and Electronic Equipment (WEEE)
- FCC Declaration
- Liability Clause
- GeneReader Accessories
- Safety Information Translated into German and French

1.2 General information

1.2.1 Technical assistance

At QIAGEN, we pride ourselves on the quality and availability of our technical support. Our Technical Services Departments are staffed by experienced scientists with extensive practical and theoretical expertise in molecular biology and the use of QIAGEN® products. If you have any questions or experience any difficulties regarding the GeneReader or QIAGEN products in general, do not hesitate to contact us.

QIAGEN customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at QIAGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance, contact QIAGEN Technical Services (see back cover).

1.2.2 Policy statement

It is the policy of QIAGEN to improve products as new techniques and components become available. QIAGEN reserves the right to change specifications at any time. In an effort to produce useful and appropriate documentation, we appreciate your comments on this user manual. Please contact QIAGEN Technical Services (see back cover).

1.3 Intended use of the GeneReader

The GeneReader is designed to perform next-generation sequencing (NGS) applications by integrating highly parallel fluorescence-based sequencing chemistry with detection of the corresponding fluorescent signals from templates that have been clonally amplified using the GeneRead QIAcube®.

GeneReader software provides a wizard for setting up the sequencing, data storage management, and the functionality for base calling and generation of FASTQ files.

The GeneReader is intended to be used only in combination with QIAGEN kits indicated for use with the GeneReader for applications described in the respective QIAGEN kit product sheets or handbooks.

The GeneReader is intended for Research Use Only. Not for use in diagnostic procedures.

The GeneReader is intended for use by professional users trained in molecular biological techniques and in the operation of the GeneReader.

1.4 Requirements for GeneReader users

Table 1 covers the general level of competence and training necessary for transportation, installation, use, maintenance and servicing of the GeneReader.

Table 1. Recommended training and skill proficiency to use, maintain and service GeneReader.

Task	Personnel	Training and experience
Transportation	No special requirements	No special requirements
Installation	QIAGEN Field Service Specialists only	Special training required
System relocation	QIAGEN Field Service Specialists only	Special training required
Routine use (running protocols)	Laboratory technicians or equivalent	Appropriately trained and experienced personnel familiar with use of computers and automation in general
Regular and weekly maintenance	Laboratory technicians or equivalent	Appropriately trained and experienced personnel familiar with use of computers and automation in general
Annual preventative maintenance and servicing	QIAGEN Field Service Specialists only	Special training required

2 Safety Information

Before using the GeneReader, it is essential that you read this user manual carefully and pay particular attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of the instrument and to maintain the condition of the instrument to ensure its continued safe operation.

IMPORTANT: safety information in German (*Sicherheitsinformationen*) and French (*Informations de Sécurité*) is provided in Appendix C, page 74.

The following types of safety information appear throughout the *QIAGEN GeneReader User Manual for Advanced Process Flow (APF/HP) Instrument Configuration*.

WARNING 	The term WARNING is used to inform you about situations that could result in personal injury to you or others. Details about these circumstances are given in a box like this one.
---	--

CAUTION 	The term CAUTION is used to inform you about situations that could result in damage to an instrument or other equipment. Details about these circumstances are given in a box like this one.
---	--

The guidance provided in this manual is intended to supplement, not supersede, the normal safety requirements prevailing in the user's country.

2.1 Proper use

WARNING 	Risk of personal injury and material damage (W1) Improper use of the GeneReader may cause personal injuries or damage to the instrument. The GeneReader must only be operated by qualified personnel who have been appropriately trained. Servicing of the GeneReader instrument must only be performed by a QIAGEN Field Service Specialist.
---	---

Perform the maintenance as described in Section 7. QIAGEN charges for repairs that are required due to incorrect maintenance.

WARNING 	Risk of personal injury and material damage (W2) The GeneReader is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone. Contact QIAGEN Technical Services to relocate the instrument.
--	---

WARNING 	Risk of personal injury and material damage (W3) Do not attempt to move the GeneReader during operation.
---	--

WARNING 	Risk of personal injury and material damage (W4) Load flow cell only in accordance with step-by-step instructions provided by GeneReader software. Beware of moveable parts.
---	--

WARNING 	Risk of personal injury and material damage (W5) Do not stare into the beam of the flow cell bar code reader.
---	---

CAUTION 	Risk of material damage (C1) Avoid moving the workbench and causing vibrations to the GeneReader during operation to prevent disturbing sensitive optical measurements.
---	---

CAUTION 	Damage to the instrument (C2) Avoid spilling water or chemicals onto the GeneReader. Damage caused by water or chemical spillage will void your warranty.
---	---

CAUTION 	Risk of material damage (C3) Do not place any items on top of the instrument.
---	---

In case of emergency, power OFF the GeneReader using the power switch at the right, rear panel of the instrument and unplug the power cord from the power outlet.

CAUTION 	Damage to the instrument (C4) Only use QIAGEN consumables with the GeneReader. Damage caused by use of other types of consumables will void your warranty.
---	--

CAUTION 	Damage to the instrument (C5) Make sure that the flow cell is inserted in the correct position. Incorrect insertion of the flow cell can damage the instrument.
---	---

WARNING 	Fire hazard (W6) Empty the liquid waste bottle before each run and make sure to place it in the correct orientation back in the GeneReader instrument. Spilling of liquid-waste may cause an electrical short-circuit and fire.
---	---

2.2 Electrical safety

Disconnect the line power cord from the power outlet before servicing.

<p>WARNING</p> 	<p>Electrical hazard (W7)</p> <p>Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.</p> <p>Intentional interruption is prohibited.</p> <p>Lethal voltages inside the instrument</p> <p>When the instrument is connected to line power, terminals may be live and opening covers or removing parts is strongly discouraged as doing so may result in electric shock or death.</p>
---	---

<p>WARNING</p> 	<p>Damage to electronics (W8)</p> <p>Before powering ON the instrument make sure that the correct supply voltage is used.</p> <p>Incorrect use of supply voltage may cause damage to electronics.</p> <p>See specifications indicated on the type plate of the instrument.</p>
---	---

<p>WARNING</p> 	<p>Risk of electric shock (W9)</p> <p>Do not open any panels on the GeneReader.</p> <p>Risk of personal injury and material damage</p> <p>Only perform maintenance that is specifically described in this user manual.</p>
---	--

To ensure satisfactory and safe operation of the GeneReader, follow the advice below:

- The line power cord must be connected to a line power outlet that has a protective conductor (earth/ground).
- Place instrument in a location so that the power cord is accessible and can be connected/disconnected.
- Use only the power cord delivered by QIAGEN.
- Do not adjust or replace internal parts of the instrument.
- Do not operate the instrument with any covers or parts removed.
- If liquid has spilled inside the instrument, switch off the instrument, disconnect it from the power outlet and contact QIAGEN Technical Services.

If the instrument becomes electrically unsafe, prevent other personnel from operating it and contact QIAGEN Technical Services (see back cover).

The instrument may be electrically unsafe when:

- It or the line power cord appears to be damaged.
- It has been stored under unfavorable conditions for a prolonged period.
- It has been subjected to severe transport stresses.

2.3 Environment

2.3.1 Operating conditions

WARNING 	Explosive atmosphere (W10) The GeneReader is not designed for use in an explosive atmosphere.
---	---

WARNING 	Risk of explosion (W11) The GeneReader is intended for use with reagents and substances supplied with QIAGEN kits. Use of other reagents and substances may lead to fire or explosion.
---	--

CAUTION 	<p>Damage to the instrument (C6)</p> <p>Direct sunlight may bleach parts of the instrument and cause damage to plastic parts.</p> <p>The GeneReader must be located out of direct sunlight.</p>
---	--

2.4 Chemicals

WARNING 	<p>Hazardous chemicals (W12)</p> <p>Some chemicals used with this instrument may be hazardous or may become hazardous after completion of the protocol run. Always wear safety glasses, gloves, and a lab coat. The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe and that the instrument operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable Safety Data Sheets (SDSs) or OSHA,* ACGIH[†] or COSHH[‡] documents.</p> <p>Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
---	---

* OSHA: Occupational Safety and Health Administration (United States of America).

[†] ACGIH: American Conference of Government Industrial Hygienists (United States of America).

[‡] COSHH: Control of Substances Hazardous to Health (United Kingdom).

2.5 Waste disposal

Used labware may contain hazardous chemicals. Waste must be collected and disposed of properly according to local safety regulations.

For more information about how to dispose of the GeneReader instrument, see “Waste Electrical and Electronic Equipment (WEEE)” in Appendix A, page 70.

<p>WARNING</p> 	<p>Hazardous chemicals and infectious agents (W13)</p> <p>The waste contains samples and reagents. This waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.</p>
---	--

2.6 Mechanical hazards

The door of the GeneReader must remain closed during operation of the instrument. Only handle the flow cell loading station when the flow cell door has been released by the software.

Note: Only power OFF the instrument if the process has been properly terminated by the software and the flow cell door is closed.

<p>WARNING</p> 	<p>Moving parts (W14)</p> <p>To avoid contact with moving parts during operation of the GeneReader, the instrument must be operated with the flow cell door closed.</p> <p>If the door sensor is not functioning correctly, contact QIAGEN Technical Services.</p>
---	---

<p>WARNING</p> 	<p>Risk of overheating (W15)</p> <p>To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the GeneReader.</p> <p>Slits and openings that ensure the ventilation of the GeneReader must not be covered.</p>
---	--

2.7 Maintenance safety

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W16)</p> <p>Only perform maintenance that is specifically described in this user manual.</p>
---	---

CAUTION 	Damage to the instrument (C7) Do not use bleach, solvents, or reagents containing acids, alkalis, or abrasives to clean the GeneReader.
---	--

2.8 Symbols on the GeneReader

Symbol	Location	Language	Description
	Inside instrument	–	Heat hazard – do not perform maintenance before the system has cooled down.
	On the instrument	–	Mechanical hazard – avoid contact with moving parts.
	On front of the instrument, open door	–	Mechanical hazard – avoid contact with moving parts.
	Inside instrument	–	Electric shock hazard
	On the instrument, right side panel	EN	This product contains a class 2 laser. Do not stare into the beam.
	On the instrument, right side panel	FR	This product contains a class 2 laser. Do not stare into the beam.
	On front of the instrument, open door	–	This product contains a class 2 laser. Do not stare into the beam.
	Type plate on the right side panel	–	WEEE about the disposal of waste electrical and electronic equipment for Europe and rest of the world.
	Type plate on the right side panel	–	Legal manufacturer.

Symbol	Location	Language	Description
	On the instrument, right side panel	–	Consult instructions for use.
	On the instrument, right side panel	EN	Disconnect power supply before servicing.
	Inside instrument	–	Earth (Ground)
	Type plate on the back of the instrument	–	CE mark for Europe
	Type plate on the back of the instrument	–	FCC mark of the United States Federal Communications Commission
	Type plate on the back of the instrument	–	RCM (former C-Tick) for Australia (supplier identification N17965)
	Type plate on the back of the instrument	–	RoHS mark for China (the restriction of the use of certain hazardous substances in electrical and electronic equipment)
	Type plate on the back of the instrument	–	Instrument serial number
	Type plate on the back of the instrument	–	Certification mark Indicates that the product was tested and has met the certification requirements for electrical, plumbing and/or mechanical products
	Sticker on the instrument, right side panel	–	Sticker indicating GeneReader with Advanced Process Flow (APF) and UMI Advanced (HP) instrument configuration

3 General Description

The QIAGEN GeneReader performs fully automated next-generation sequencing (NGS) by integrating highly parallel fluorescence-based sequencing chemistry with detection of the corresponding fluorescent signals on templates that have been clonally amplified using the GeneRead QIAcube.

The GeneReader sequencer consists of the GeneReader, the workstation, the GeneReader software and a handheld bar code scanner that connects to the workstation for scanning bar codes of kits and buffers, which are then automatically entered into the GeneReader software. The initial software installation is performed by a QIAGEN Field Service Specialist. There are two USB connections between the GeneReader instrument and workstation. The GeneReader sequencer includes several additional components, which are listed in Section 4.1.

The GeneReader software provides a FASTQ file of sequence information for each analyzed sample that is ready for QCI Analyze or GeneRead Link. QCI Analyze automatically runs an optimized workflow for GeneReader panels and generates a VCF result file that is ready for upload to QCI Interpret.

3.1 QIAGEN GeneReader Sample to Insight[®] NGS workflow

The QIAGEN GeneReader Sample to Insight workflow provides a streamlined and standardized approach to next-generation sequencing (NGS), from sample preparation to the biological interpretation of sequencing data.

The majority of the workflow is automated, ensuring greater standardization and more accurate results.

QIAGEN Clinical Insight[™] combines analytical tools and integrated human disease content, providing access to current and advanced interpretations of genomic data. QCI Analyze or GeneRead Link automatically sends samples through predefined workflows and provides a web-based result viewer. Afterwards, QCI Interpret (another web-based viewer) provides a biological interpretation.

3.2 GeneReader principle

The workflow includes the following six (6) processes:

- sequencing primer hybridization
- Flow Cell preparation
- reagent preparation
- experiment set-up
- Flow Cell loading and start of run
- post-run maintenance wash

The principle behind a sequencing run on the GeneReader is described below.

The first step of the GeneReader sequencing-by-synthesis technology consists of the incorporation of unique deoxyribonucleotide triphosphates (dNTPs) which are both fluorescently labelled and reversibly terminated. These are known as “labeled nucleotides”. This is followed by the addition of unlabeled reversibly terminated dNTPs or dark nucleotides”.

The GeneReader sequencing chemistry uses four dye colors for labeling each dNTP to differentiate each base (A, C, G or T) that is incorporated onto the DNA fragment. Furthermore, the reversible terminators facilitate the addition of only one engineered nucleotide at a time to the growing strand of all DNA templates. Upon signal detection from each bead, the fluorescent labels as well as the terminators are removed which then allows for a new cycle of incorporation, hence ensuring highly accurate and cost-effective next-generation sequencing.

DNA libraries are clonally amplified on beads using the GeneRead QIAcube to serve as a sequencing template. After hybridization of a sequencing primer, the primer-template carrying beads are immobilized via direct bead-glass interaction to produce a high-density array on a GeneReader Flow Cell. To read out the content of templates on each bead, the array of fragments is first subjected to reagents containing uniquely engineered dNTPs, as described above. These bases are incorporated by the addition of a modified DNA polymerase to the end of the growing strand of DNA in accordance with the base on the complementary strand. The array is subsequently scanned by a high-resolution digital camera and the fluorescent output of each of the four dye colors at each array position is measured and recorded. Finally, the array is exposed to cleavage chemistry to break off the fluorescent dyes and reversible terminators that will then allow additional bases to be incorporated. This cycle is then repeated on the GeneReader.

3.3 External features of the GeneReader




Figure 1. External features and functions of GeneReader.

3.3.1 Status lights

The status lights illuminate in the following patterns:

- When the instrument is not running, the status lights are off.
- When a protocol is running normally, the green light is on.
- When a protocol is running, but a pause has been requested, or when a protocol is paused, the green light blinks.
- When all protocols have been stopped, canceled or an error has occurred, the red light is on.
- When the protocol has finished successfully, the green light blinks.

3.3.2 Flow Cell door

Opening the Flow Cell door

Flow Cells are inserted through the flow cell door. Opening the Flow Cell door is controlled by software.

Note: The Flow Cell door cannot be opened manually.

Closing the Flow Cell door

Push the Flow Cell door manually until the fastener snaps in place to lock the door.

Note: The GeneReader will not work if the Flow Cell door is not locked.

3.3.3 Hood

Opening the hood

Opening the hood is controlled by software. The hood must be opened prior to opening the fluidic drawer.

Note: The hood can be manually released if the GeneReader loses power. Insert your hand, palm facing upward, and feel for the hole located approximately where the Main Hood meets the fluidic drawer. When your fingers are in the hole you will feel a lever on the right and by pulling slightly towards the front the hood will open.



Figure 2. Diagram illustrating how to manually open the hood of GeneReader.

Closing the hood

The fluidic drawer must be closed (pushed all the way in until a clicking sound is heard) prior to manually pulling down the hood.

Note: The Main Hood on the left side of the instrument must only be opened by QIAGEN Field Service Specialists.

3.3.4 Fluidic drawer

Open the fluidic drawer to:

- Load or unload 1 liter wash bottles.
- Load or unload 50 ml tubes.
- Insert or remove the liquid waste bottle.
- Clean the cooling block, remove excessive condensation.
- Clean the drawer.
- Clean the dip sticks.
- Clean the waste level sensor surface.

The fluidic drawer remains locked during a run.

Note: If the hood is closed, the fluidic drawer cannot be opened.

3.3.5 USB ports

The two (2) USB ports are used to connect the GeneReader with the workstation. They are located on the right side panel of the instrument.

3.3.6 Power switch

The power switch is located on the right side panel of the GeneReader.

3.3.7 Power cord socket

The power cord socket is located on the right side panel of the GeneReader. It allows connection of the GeneReader to a power outlet via the supplied power cord.

3.3.8 Cooling air outlet

Cooling air outlets are on the right side and back of the GeneReader. They allow cooling of the internal components. The instrument should not be positioned in close proximity to a wall that could block airflow.

3.3.9 Workstation equipment

The GeneReader system is operated with a workstation.

The workstation specifications are listed in Section 9.3.

3.4 Internal features of the GeneReader



- 1 50 ml conical tubes in cooling compartment
- 2 One liter bottles for sequencing wash buffer or Maintenance Wash Buffer
- 3 Liquid waste bottle

Figure 3. Internal view of the fluidic drawer of GeneReader.

3.4.1 GeneReader Flow Cell

Sequencing beads carrying single strand DNA templates with sequencing primers annealed are deposited into the GeneReader Flow Cell according to the current *QIAGEN GeneRead UMI Advanced Sequencing Q Handbook*. The Flow Cell should be loaded into the GeneReader through the Flow Cell door with the bar code towards the left side of the GeneReader.



Figure 4. Loading the Flow Cell of GeneReader.



Figure 5. View inside the Flow Cell door of GeneReader.

Flow Cell bar code reader

Upon loading of a Flow Cell, the Flow Cell bar code reader scans the bar code on the Flow Cell and provides the information to the GeneReader software. If the Flow Cell has been inserted in the wrong orientation an error message is displayed.

4 Installation Procedures

4.1 System delivery and installation

The unpacking and installation of the GeneReader is performed by a certified QIAGEN Field Service Specialist. A person who is familiar with your laboratory and computer equipment should be present during the installation.

The following items are delivered:

- GeneReader instrument
- *GeneReader User Manual for Advanced Process Flow (APF/HP) Instrument Configuration*
- Workstation
- GeneReader software (will be installed by QIAGEN Field Service during initial set up)
- Additional components: 1 international power cable set, 2 USB cables, 1 handheld bar code scanner, 1 waste container and 4 one liter bottles

Note: The computer recommended for data analysis is the QIAGEN Powerstation, however, any hardware that meets the specifications as defined in the current version of the *QIAGEN Clinical Insight Analyze User Manual* may be used.

4.2 Site requirements

The GeneReader must be located out of direct sunlight, away from heat sources, and away from sources of vibration and electrical interference. Refer to Section 9.1 for the operating conditions (temperature and humidity). The site of installation should be free of excessive drafts, excessive moisture and dust, and not be subject to large temperature fluctuations.

Refer to Section 9.2 for the weight and dimensions of the instrument.

Use an appropriate workbench to accommodate the GeneReader. Ensure that the workbench is dry and clean, and has additional space for accessories. To accommodate the GeneReader instrument with the hood open, 125 cm (49.3 in.) minimum clearance above the workbench is required. Allow at least 10 cm (4 in.) of free space behind the instrument for cables and cooling of the instrument.

The GeneReader must be placed within approximately 1.5 m (59 in.) of a properly grounded (earthed) AC power outlet. The power line to the GeneReader should be voltage regulated and surge protected, and an uninterruptable power supply (UPS) is optional.

Note: Do not place the GeneReader instrument on a vibrating surface or near vibrating objects.

WARNING 	Risk of overheating (W15) To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the GeneReader. Slits and openings that ensure the ventilation of the GeneReader must not be covered.
---	---

WARNING 	Risk of personal injury and material damage (W2) The GeneReader is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.
---	--

4.3 AC Power and USB cable connections

4.3.1 Power requirements

The GeneReader operates at:

- 100–240 V AC, 50/60 Hz, max. 600 VA

The workstation operates at:

- 100–240 V AC, 50/60 Hz, max. 400 VA

Make sure that the voltage rating of the GeneReader and workstation are compatible with the AC voltage available at the installation site.

WARNING 	Damage to electronics (W8) Incorrect use of supply voltage may cause damage to electronics. See specifications indicated on the type plate of the instrument.
---	--

4.3.2 Grounding requirements

To protect operating personnel, the National Electrical Manufacturers' Association (NEMA) recommends that the GeneReader be correctly grounded (earthed). The instrument is equipped with a 3-conductor AC power cord that, when connected to an appropriate AC power outlet, grounds (earths) the instrument. To preserve this protection feature, do not operate the instrument from an AC power outlet that does not have a ground (earth) connection.

4.3.3 Installation of AC power cords, USB cables and workstation

IMPORTANT: Make sure that the GeneReader and workstation power switches are set to the OFF position.

Connect the power cords and USB cables as shown in the diagram below.

Note: It is recommended to plug the instrument and workstation directly into a shared line power outlet because the instrument and workstation are connected with 2 USB cables. This configuration results in less interference caused by ground loops.




Figure 6. Diagram of cable connections of GeneReader.

4.3.4 Installation of AC power cords, USB cables and workstation using an optional UPS

IMPORTANT: Make sure that the GeneReader and workstation power switches are set to the OFF position.

Connect the GeneReader and workstation to the UPS as shown in the diagram below. The UPS should be configured in such a way that a Windows® shutdown command is issued shortly before backup power is lost. The UPS should be able to supply the GeneReader and workstation with a maximum 1000 VA for the desired time.




Figure 7. Diagram of cable connections using an optional UPS.

4.4 Workstation requirements

The workstation, supplied with the GeneReader instrument, fulfills the requirements of the GeneReader software as detailed in Table 2.

Table 2. Workstation system requirements

Description	Minimum requirement
Operating system	Microsoft® Windows 7 Professional 64-bit (D/F/E/I)
Processor	Intel® Xeon® E5-1650 v4 (3.60 GHz 6 Core)
Main memory	32 GB (2x16 GB) DDR4-2400
Hard disk space	1 x 256 GB SATA SSD 3 x 4.0 TB SATA 7200 rpm HDD
Graphic card	NVIDIA NVS 310 1GB 1st GFX
Monitor	Resolution at least 1920 x 1080 pixels
Ports	Front: 4x USB 3.0 Rear: 2x USB 2.0, 4x USB 3.0 Integrated Intel I-218 Gbit LAN

Adobe® Reader® software is preinstalled on the computer to view reports generated in PDF format.

IMPORTANT: Microsoft Windows 7 Professional comes pre-installed on the GeneReader workstation. By default, the following Microsoft background activities are disabled and any change to the Microsoft Windows configuration by an operator may affect sequencing performance:

- Automatic Windows update
- Background indexing
- Automatic background defragmentation
- Automatic backup
- Energy saver mode

4.5 Getting started

4.5.1 Powering ON the GeneReader and workstation

1. Power ON the GeneReader using the power switch located on the right side panel of the GeneReader.
2. Power ON the workstation and monitor.
3. Log in to the system using your password.
4. Start the GeneReader software.

When the GeneReader software has been started, the following window will appear:





IMPORTANT: If virus scan software is installed on the GeneReader workstation, make sure that active scanning or anti-virus software updates are not performed during a sequencing run.

Note: After the instrument is powered ON and the software is running, the temperature regulation of the GeneReader starts. The specified temperature will be reached in 30 minutes.

4.5.2 Software upgrade

The GeneReader software is not preinstalled. QIAGEN Field Service will perform the installation. To upgrade software, follow the steps below:

1. Confirm that the following requirements are met before you begin the upgrade process:
 - a. A run is not being performed.
 - b. There are no loaded Flow Cells on the instrument.
 - c. The GeneReader software is closed.
2. Create a folder named "Upgrade" on the desktop. Copy the **GeneReaderInstaller** zip file that is provided by QIAGEN to the newly created "Upgrade" folder.
3. Extract the **GeneReaderInstaller** zip file into the "Upgrade" folder by right-clicking the file and selecting **Extract all...** from the content menu.
4. Go to the "Upgrade" folder.
5. Go to the "GeneReaderInstaller" subfolder and run the **setup.exe** file.
6. Follow the on-screen instructions to confirm all license agreements and install all parts of the software. The sequence of this process may differ depending on the existing and new software version.



7. The computer may require a restart after the installation is complete. Follow the on-screen instructions and restart the computer.

-
8. The installer places shortcut icons for starting the GeneReader software on the desktop and in the start menu. In the start menu additional shortcuts are created to start the software in service mode, but they are intended for service purpose only.
 9. Start the GeneReader software by clicking the **GeneReader** shortcut icon. If necessary, the GeneReader firmware updater appears. If this occurs click **OK** to perform a firmware update. This may take several minutes to finish.

Note: Make sure that a required firmware update is performed successfully.

10. After the firmware update is complete, wait until the software has finished the initialization and homing processes. The GeneReader software will re-start automatically.
11. Delete the "Upgrade" folder that was created on the desktop in step 2. The software is now ready for use.

Note: To uninstall, run the **setup.exe** file in the "GeneReaderInstaller" subfolder.

5 Operating Procedures

Before proceeding, we recommend that you familiarize yourself with the features of the instrument by referring to Section 3.

CAUTION 	Damage to the instrument (C4) Only use QIAGEN flow cells and consumables with the GeneReader. Damage caused by use of other types of flow cell or consumable will void your warranty.
---	---

CAUTION 	Risk of material damage (C1) Avoid moving the workbench and causing vibrations to the GeneReader during operation to prevent disturbing sensitive optical measurements.
---	---

6 GeneReader Software

6.1 Use of the GeneReader software

6.1.1 Start-up

During start-up, the software will perform a self-test and take an inventory of its own configuration files and any loaded Flow Cells and reagents. If Flow Cells are still loaded, you will be asked to unload them. Should a configuration file be missing, or any failure occur during the start-up procedure, contact QIAGEN Technical Services (see back cover).


6.1.2 Software workflow

The GeneReader software provides a sample-to-result workflow, which includes the following stages:

- Definition of flow cell parameters, including multiplexing information.
- Step-by-step loading instructions.
- Sequence detection.
- Generation of FASTQ data file.

The software also provides step-by-step instructions for the maintenance wash that is part of routine GeneReader maintenance (see Section 7).

6.1.3 User interface




The GeneReader software is divided into four distinct areas (see image above). At the top, the menu offers access to general functionalities. Below the menu bar, the navigation area contains quick links to the commonly used views and workflows. The main portion of the screen is dedicated to the work area, which displays the current setup of flow cells and status of a protocol run, and on the right side the message area displays messages relating to protocol function, setup and errors.

Menu and navigation



Select **File** and **Exit** in the menu to close the GeneReader software. This cannot be done while a protocol is running or if any flow cells are still loaded.




Select **Maintenance** and the **Manage Files** command to start the Manage Files wizard, which can be used to clear space on the hard drive prior to starting a run. Any selected data for clearing will be permanently deleted. If data is still needed, make sure to create a proper backup on an external drive. Select **Maintenance** and the **Service Archive** command to create a service archive if requested by QIAGEN Technical Services.



Use the **Help** menu to access version and copyright information in the **About** dialog.



The **Flow Cells** button displays the status of all flow cells inside the GeneReader. Refer to Section 6.1.4 for more information




The **Fluids** button displays the status of all fluids in the main work area. Refer to Section 6.1.5 for more information.



The **Configuration** button opens the view for editing several configuration values. Refer to Section 6.1.6 for more information




Use the **Run Setup** command to open the workflow and to setup and start a run. Refer to Sections 6.2.1 and 6.2.2 for more information.



Use the **Add** command to load flow cells during a run in staggered mode once the first cycle of sequencing is complete, and prior to 24 hours after the run has started. Refer to Section 6.2.3 for more information.



Use the **Unload** command to unload flow cells from the device. Refer to Section 6.2.4 for more information.



Use the **Wash** command to open the maintenance workflow. Refer to Section 7 for more information.

6.1.4 Status of Flow Cells view

Status of Flow Cells is displayed in the main work area and displays an overview of all Flow Cells loaded inside the instrument.

No.	Barcode	Flow Cell Name	Station Name	Current Step	Cycle	Total Cycles	Progress	Estimated Remaining Time	Start Time	Estimated End Time	Action
1	8130831245678912150...	FlowCell1	19_Dwell	Dwelling	Sequencing 1 of 127	127	1 day 15hr	07-12 01:25 PM	07-14 05 AM		

Signal Information

Data not available

File Transfer

The current status of all loaded flow cells are shown, as well as estimations for remaining run times. In the **Action** column on the right-hand side of the screen, there are icons for performing specific actions regarding the particular flow cell. The icons may not be enabled if the action is not available for the current status of the flow cell.



The run can be interrupted by pressing the **Stop** icon. Be aware that the protocol does not stop immediately. The protocol will stop after the currently running process step is complete. This button is available if the Flow Cell is running, after first cycle of sequencing is complete.



The **Open Flow Cell Data Folder** icon opens the run folder of the corresponding Flow Cell. This button is available if the Flow Cell is either finished or stopped, and the analysis is not currently running.



The **Flow Cells Results** icon opens the folder button that contains the analysis result data of the corresponding Flow Cell. This button is available if the Flow Cell is finished and the analysis is not currently running, but not if the Flow Cell was stopped.

Note: The number of cycles for a loaded Flow Cell will be slightly higher than the amount of cycles selected for sequencing. This is due to extra cycles being required for lead/lag correction and sample index sequencing.

To view the multiplexing details of a flow cell, press the **(+)** icon at the beginning of each row.

6.1.5 Status of Fluids view

The Status of Fluids view is displayed in the main work area and provides an overview of all fluids currently loaded on the GeneReader.

Status Of Fluids						
Reagents						
Kit Name	Lot Number	Installation Date	Expiration Date	Temperature Zone	Current Temperature	Expected Temperature
GeneRead UMI Sequencing Q Buffers (3)	1234567890	12.07.2017 13:24	31.12.2040 23:59	Ambient	5.0 °C	5.0 °C
GeneRead UMI Sequencing Q Add-Ons (3)	1234567890	12.07.2017 13:24	31.12.2040 23:59	Cooled	5.0 °C	5.0 °C

Available cycles for reagents	347
-------------------------------	-----

Wash Buffer					
Name	Id	Lot Number	Installation Date	Expiration Date	Temperature Zone
Wash Buffer 9	81700891234567890124000001	1234567890	12.07.2017 13:24	31.12.2040 23:59	Ambient
Wash Buffer 11	81700901234567890124000003	1234567890	12.07.2017 13:24	31.12.2040 23:59	Ambient

Waste Container	
Fill Level	Waste Bottle State
 53 %	Present

The screen is partitioned into three areas and provides information on the levels of reagents (top), wash bottles (middle) and the waste bottle (bottom). This view provides information only and does not allow the user to perform any actions.

6.1.6 Status of Configuration view

Configuration settings can be changed in The Status of Configuration view. This view is divided into two areas (**Connection Settings and Backup/Restore System Settings**) which can be selected by tabs.

Connection Settings

The **Connection Settings** tab allows the setup of connections to external systems for experiment planning and analyzes (i.e., QCI Analyze, GeneRead Link) . Depending on the local setup, there are three options to select:

- **Not connected to external system:** Select this if there is no external experiment planning system available. A checkbox is displayed to enable the automatic FASTQ file transfer. The destination folder for the FASTQ files can be specified using the “Destination Directory” field.
- **Connected to QCI Analyze:** Select this if QCI Analyze is used as the experiment planning system. The user must provide the network share directory that the QCI Analyze experiment planner files will be stored, and the directory that the result data files will be saved to. For more information refer to the current version of the *QIAGEN Clinical Insight Analyze User Manual*.
- **Connected to GeneRead Link:** Select this if GeneRead Link is used as the experiment planning system. The user must provide several configuration values. First, enter the network share directory where the result data files will be saved. Second, enter the URL, user name and password of the GeneRead Link server. To test if the GeneRead Link is working, click the **Go**

button. Third, the user must provide the time (in minutes) to wait for a retry if the data upload failed. For more information refer to the *GeneRead Link User Manual*.

Note: Previous settings for the connection of the GeneReader software to an external experiment planner system are saved after connection settings are changed.

The **Connection Settings** tab also allows the setup of a connection to GeneRead Databank as an external data management tool. If GeneRead Databank is used to manage files, select the **Connected to Customer Data Management** check box.

- **Connected to Customer Data Management:** If GeneRead Databank is used as the external system, the **Manage Files** command in the **Maintenance** menu is deactivated. For more information about GeneRead Databank (Customer Data Management) refer to the *GeneRead Databank User Manual*.

Status Of Configuration

Connection Settings | Backup / Restore System Settings

Settings for connection to external systems

Not connected to external systems

Transfer FastQ files automatically

Destination Directory
C:\Temp

Connected to QCI Analyze

Destination Directory
Z:\runs

GeneReader Planner Directory
Z:\flowcell

Connected to GeneRead Link

Destination Directory

GeneRead Link URL

User name for GeneRead Link

Password

Test Connection

Retry interval in minutes for failed transfers and uploads
10


Settings for connection to GeneRead Databank

Connected to GeneRead Databank

<Backup/Restore System Settings

The **Backup/Restore Systems Settings** tab displays two buttons: one to save all GeneReader instrument settings to a **Backup** file and one to **Restore** the settings by loading a backup file that

was previously saved. **Backup** and **Restore** may be disabled depending on the current state of the GeneReader instrument and software.



6.1.7 File handling

Required files

The GeneReader software performs an inventory of its own files and configuration upon start-up. If a message is shown about missing or invalid files, please contact QIAGEN Technical Services (see back cover).

6.2 Workflow procedures

6.2.1 Flow Cell setup

There are 3 options for generating a Flow Cell setup in the GeneReader NGS System:

- QCI Analyze for importing Flow Cells.
- GeneRead Link for importing Flow Cells.
- GeneReader software for adding Flow Cells manually.

QCI Analyze (web interface)

QCI Analyze is a browser-based system for analysing NGS data. Underneath QCI Analyze, a CLC Genomics Server is running the analyzes, storing data and handling various processes such as queuing. QCI Analyze can import NGS data, analyse it, display results for inspection and export result data in VCF format. In addition, it can connect directly to QCI Interpret for interpretation and reporting. For more information on how to generate a Flow Cell setup using QCI Analyze refer to the *QIAGEN Clinical Insight Analyse User Manual*.

See “Protocol: Loading and Running the GeneReader” in the current *QIAGEN GeneRead UMI Advanced Sequencing Q Handbook* for information on how to import a Flow Cell setup from QCI Analyze.

GeneRead Link (web interface)

GeneRead Link is QIAGEN middleware that directs planning and execution of the GeneReader NGS System workflow for samples. It provides bidirectional connectivity with a Laboratory Information Management System (LIMS), the work instructions for laboratory users and automatically interfaces to specific QIAGEN instruments. For more information on how to generate a Flow Cell setup refer to the *GeneRead Link User Manual*.

See “Protocol: Loading and Running the GeneReader” in the current *QIAGEN GeneRead UMI Advanced Sequencing Q Handbook* for information on how to import a Flow Cell setup from GeneRead Link.

GeneReader software

GeneReader software can be directly used for a Flow Cell setup. See “Protocol: Loading and Running the GeneReader” in the current *QIAGEN GeneRead UMI Advanced Sequencing Q Handbook* for a step-by-step guide on how to add a Flow Cell manually using the GeneReader software user interface.

IMPORTANT: If GeneReader software is used for Flow Cell setup there will be no automatic data analysis of FASTQ files generated by the sequencer through QCI Analyze. In *QCI Analyze* it is now possible to manually import FASTQ files for secondary analysis: for more information on how to manually upload FASTQ files using QCI Analyze refer to the *QIAGEN Clinical Insight Analyse User Manual*. In order to integrate data analysis in QCI Analyze use the QCI Analyze or the GeneRead Link web interfaces for Flow Cell setup.

Note: If GeneReader software is used for Flow Cell setup with staggered loading, ensure that sufficient reagents are loaded at the beginning of the run, for all of the planned Flow Cells

6.2.2 Loading and running the GeneReader

See “Protocol: Loading and Running the GeneReader” in the current *QIAGEN GeneRead UMI Advanced Sequencing Q Handbook* for a step-by-step guide on how to load one or more Flow Cells using the GeneReader software user interface.

6.2.3 Staggered loading of Flow Cells

See “Protocol: Staggered Loading of Flow Cells” in the current *QIAGEN GeneRead UMI Advanced Sequencing Q Handbook* for a step-by-step guide on how to load one or more Flow Cells in staggered mode after a sequencing run has been started in the GeneReader software.

6.2.4 Unloading reagents and Flow Cells

See “Protocol: Unloading Reagents and Flow Cells” in the current *QIAGEN GeneRead UMI Advanced Sequencing Q Handbook* for a step-by-step guide on how to unload reagents and Flow Cells in the GeneReader software.

6.2.5 Run finished

When the sequencing and analysis of all Flow Cells is completed, the Run Finished wizard will allow:

- Transfer of FASTQ and report files to another local or network directory by clicking **Transfer**.
- Viewing of result directory.
- Unloading of flow cells and fluids.
- A deionized (DI) water or maintenance wash to be performed.
- Running the “Manage Files” wizard.

A run report is generated automatically after the analysis of the Flow Cells. A separate report is created for every Flow Cell in the run and it is saved to the Flow Cell dDirectory within the Run Directory. The report contains sections on the run setup, the quality score distribution, data on bead quality and the sequence quality score, as well as multiplex data. Additionally, run reports contain any errors or warnings that may have been logged/reported during the course of the sequencing and analysis processes.

Log files


- The GeneReader software writes an **EventLog.txt** file to the folder:
C:\ProgramData\QIAGEN\GeneReader\GeneReader. This log file contains warnings and errors and is intended to provide information in case of an issue. To prevent this log file from getting too large, a new log file is created once a month.
- During a run, a run-specific event log file is created within the run folder. Refer to this file for any issues that occurred during the execution of a sequencing protocol.

Final report

- When the analysis software is completed, a Flow Cell report is generated as a PDF file and written to the path:
\An example file is as follow:
D:\QIAGEN\GRSData\2016_10_06_133222_SampleRun\FC_FC1_2c0sbgi977zjkqkbtbg1vuj71\Analyze\1610061332_FC1_FinalReport.pdf

The **FinalReport.pdf** file contains important run information such as:

- General information about the run



The screenshot shows the header of a PDF report titled "GeneReader Flow Cell Report" with the QIAGEN logo. Below the header is a table with the following data:

Run Name	170603_PQ_fullP_GR821	Instrument	1508006
Flow Cell Name	24040		
Lot	0000000000	Id	24040
Number of Cycles	127		
Mode	Multi Flow Cell Run / Standard		
Read Design	Default(V2)		
Note			
Start Time	06/3/2017 4:56 PM	End Time	06/5/2017 9:19 AM
Analysis End Time	06/5/2017 1:39 PM		

- Sample information

Sample Information

Sample Id	Adapter Q	Adapter Q Sequence	GenePanel
1	BC1	ATCACG	
10	BC10	TAGCTT	
2	BC2	CGATGT	
3	BC3	TTAGGC	
4	BC4	TGACCA	
5	BC5	ACAGTG	
6	BC6	GCCAAT	
7	BC7	CAGATC	
8	BC8	ACTTGA	
9	BC9	GATCAG	
unindexed	unindexed	NNNNNN	

- Reagent kits

Reagent Kits

Sequencing Kit	Lot	Material	Exp date
GeneRead UMI Sequencing Q Buffers (3)	0154015256	9904201	12/31/2060
GeneRead UMI Sequencing Q Add-Ons (3)	0154015256	9904301	12/31/2020

- Wash kits

Wash Kits


Wash Kit	Lot	Material	Exp date
Wash Buffer 9	1234567890	8170089	12/31/2040
Wash Buffer 11	1234567890	8170090	12/31/2040

- Software version

Software Version	GeneReader Software 1.6.0
-------------------------	---------------------------

- Quality score distribution

Quality Score Distribution



- Run metrics

Reads Past Filtering	52.04956%	20,052,002
Reads Past Demultiplexing	19,368,738	
Yield [bases]	2,649,561,215	
Yield Past Demultiplexing	2,561,004,714	


Reads Past Filtering reports the percentage and the count of reads that passed Q-score-based quality filtering before demultiplexing (i.e., the count of reads as written in the FASTQ file).

Reads Past Demultiplexing reports the count of reads for which a sample bar code could be identified.

Yield [bases] reports the total number of base calls before demultiplexing.

Yield Past Demultiplexing reports the total number of base calls for which a sample bar code could be identified.

- Sequence quality score



- Reads information

Reads

Sample Id	Sample Validity	Reads Past Filtering	Yield [bases]	FastQ file
1	Valid	606,885	79,290,637	20170603145224_8130063000000000012152404_0_1_BC1_24040.fastq
10	Valid	1,307,705	170,616,986	20170603145224_8130063000000000012152404_0_10_BC10_24040.fastq
2	Valid	1,003,902	131,250,479	20170603145224_8130063000000000012152404_0_2_BC2_24040.fastq
3	Valid	701,861	91,831,754	20170603145224_8130063000000000012152404_0_3_BC3_24040.fastq
4	Valid	699,644	91,434,528	20170603145224_8130063000000000012152404_0_4_BC4_24040.fastq
5	Valid	10,716,623	1,426,376,672	20170603145224_8130063000000000012152404_0_5_BC5_24040.fastq
6	Valid	1,747,901	232,219,492	20170603145224_8130063000000000012152404_0_6_BC6_24040.fastq
7	Valid	888,434	116,272,025	20170603145224_8130063000000000012152404_0_7_BC7_24040.fastq
8	Valid	670,822	87,567,851	20170603145224_8130063000000000012152404_0_8_BC8_24040.fastq
9	Valid	1,024,961	134,144,290	20170603145224_8130063000000000012152404_0_9_BC9_24040.fastq
unindexed	Valid	683,264	88,556,501	20170603145224_8130063000000000012152404_0_unindexed_24040.fastq

Validity indicates if samples analyzed in a flow cell are “Valid” or “Invalid”. The GeneReader software calculates in the background the presence of unexpected bar codes in the pool of samples analyzed in each flow cell. If the number of reads belonging to unexpected bar codes exceeds 5% of the total number of all indexed reads, all samples of this Flow Cell are labeled as “Invalid” due to suspected cross-contamination.

7 Maintenance

This section describes routine and weekly maintenance procedures.

7.1 Equipment and reagents to be supplied by user

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, consult the appropriate safety data sheets (SDSs) available from the product supplier.

The following material is required to perform maintenance of the GeneReader instrument:

- Deionized water or laboratory-grade pure water (referred to as DI water)

7.2 Maintenance wash

This manual describes 2 required maintenance procedures:


- **Routine maintenance (DI water wash):** must be performed when the GeneReader will be idle and/or when a run has completed and no additional runs are planned to follow immediately
- **Weekly maintenance (using the Maintenance Wash Buffer):** the maintenance wash must be performed every 4 runs or once a week (whichever occurs first)

Performing these procedures ensures that the GeneReader is kept in good operating condition and helps prevent clogging of the fluidics system.


7.2.1 Routine maintenance (DI water wash)

After a sequencing run has finished and the GeneReader will be idle or no other sequencing run is planned, a DI water wash will prevent the formation of air pockets and the drying out of the fluidic lines inside the instrument.


1. Open the GeneReader software. Click **Wash** in the navigation panel.
2. Select **DI Water Wash** and click **Next**.
3. The hood will open. Load 2 wash bottles (QIAGEN standard 1 liter bottle) containing 500 ml DI water into positions 9 and 11.
4. Load reagent tubes (QIAGEN standard 50 ml conical tubes) containing 30–35 ml DI water into positions 1, 2, 4, 5, 6 and 8 of the cooling compartment. Click **Done** twice and then click **Next**.




5. Empty the waste bottle and dispose of the contents according to local safety regulations. Reload the empty waste bottle into the instrument. Close the hood and click **Next**.
IMPORTANT: Make sure to reload the waste bottle into the instrument in the correct orientation.




6. The Flow Cell door will open. Load an empty or used Flow Cell (priming Flow Cell) to allow liquid flow. Close the Flow Cell door and click **Next**.



- Instrument will now prime with DI water. The Routine Maintenance wash will take about 15 minutes. Once complete, click **Next**.




- The Flow Cell door will open. Unload the priming Flow Cell. Close the Flow Cell door and click **Next**.



- The reagent hood will open. Unload the reagent and wash kit. Empty the waste bottle. Close the hood and click **Done**.

Note: If the GeneReader is not used immediately for another run, leave the reagent and wash bottles containing DI water in the system. This is important to ensure that the fluidic lines of the system do not dry out.

Note: If no run is planned for **longer than 1 month**, perform the long-term storage procedure to prepare the GeneReader instrument for dry storage conditions (see Section 7.2.3).



- Click **Finish** to close the wizard.

7.2.2 Weekly maintenance

To make sure reagent lines remain clear and to prevent buildup of material that could potentially clog the lines, a maintenance wash must be performed every 4 runs or once per week (whichever occurs first). The maintenance wash uses the Maintenance Wash Buffer provided in the GeneRead

UMI Sequencing Q Wash Buffer Kit (cat. no. 185905). Maintenance buffer preparation and the maintenance wash procedure are described below.

Preparation of Maintenance Wash Buffer (sufficient volume for 9 maintenance washes)

1. Take out the Maintenance Wash Buffer bottle (which contains 900 ml Maintenance Wash Buffer) from the GeneRead UMI Sequencing Q Wash Buffer Kit.
2. Add 100 ml DI water to the Maintenance Wash Buffer for a final volume of 1000 ml.
3. Label a second empty 1000 ml bottle (QIAGEN GeneReader 1 liter bottle) with a Maintenance Wash Buffer sticker.
4. Transfer 500 ml Maintenance Wash Buffer to the second bottle, so that each bottle contains 500 ml Maintenance Wash Buffer.
5. Label 6 reagent tubes (QIAGEN standard 50 ml conical tubes) with Maintenance Wash Buffer stickers.
6. Aliquot approximately 30–35 ml Maintenance Wash Buffer from one bottle to each of the 6 reagent tubes. This initial volume will be sufficient for several maintenance wash runs. **A re-fill of the maintenance wash reagent tubes will be required when the liquid level in the maintenance wash reagent tubes falls below the conical mark.** Refer to Section 7.2.2 for re-fill instructions.
7. Load the Maintenance Wash Buffer bottle that was used for aliquoting in the previous step (containing approximately 300 ml remaining volume) onto position 11 in the GeneReader. Load the bottle containing 500 ml Maintenance Wash Buffer onto position 9 (Figure 8).




Figure 8. Loading 1 liter bottles containing Maintenance Wash Buffer.

8. Load the 6 reagent tubes containing the aliquoted Maintenance Wash Buffer (approximately 30–35 ml fill volume each) into the reagent cooler according to Figure 9 (positions 1, 2, 4, 5, 6 and 8), leaving positions 3 and 7 empty.





Figure 9. Loading 50 ml tubes containing Maintenance Wash Buffer.

Weekly maintenance procedure

1. If the sequencing run has just completed, make sure that all Flow Cells are removed, all reagents have been removed from the cooling compartment and discarded, and all wash buffers are removed and stored at room temperature (15–25°C).
2. If GeneReader software is not running, double-click the GeneReader icon to start the software. The GeneReader instrument will begin to scan for Flow Cells. If prompted by the GeneReader, remove any Flow Cells, and close the Flow Cell door after removal (depending on the previous run, up to 3 Flow Cells may prompt for their removal from the GeneReader instrument).
3. Click **Wash** in the navigation panel.
4. Select **Maintenance Wash** and click **Next**. The reagent hood will open.
5. Load a 1 liter bottle (QIAGEN GeneReader 1 liter bottle) containing approximately **300 ml Maintenance Wash Buffer onto position 11** in the fluidic drawer.
6. Load a 1 liter bottle (QIAGEN GeneReader 1 liter bottle) containing **500 ml Maintenance Wash Buffer onto position 9** in the fluidic drawer.


7. Load 50 ml conical tubes (QIAGEN standard 50 ml conical tubes) **containing approximately 30–35 ml Maintenance Wash Buffer onto positions 1, 2, 4, 5, 6 and 8** of the reagent cooler, leaving positions 3 and 7 empty, according to Figure 9.

Note: A re-fill of the maintenance wash reagent tubes will be required when the liquid level in the maintenance wash reagent tubes falls below the conical mark. Use Maintenance Wash Buffer bottle on position 11 for this re-fill.




8. Click **Done** twice when finished.
9. Click **Next**.
10. Empty the waste bottle and dispose the contents according to local safety regulations. Reload the empty waste bottle into the instrument. Close the reagent hood and click **Next**.


IMPORTANT: Make sure to reload the waste bottle into the instrument in the correct orientation. An incorrect orientation may damage the instrument. Do not use force to close the hood. If blocked, check the waste bottle orientation and correct if necessary before retrying to close the reagent hood.



11. Load an empty or used Flow Cell (priming Flow Cell) to allow liquid flow, close the Flow Cell door and click **Next**.



12. The instrument will now start the wash. The maintenance wash will take about 40 minutes. The software will notify when the maintenance wash is complete, click **Next**.




13. The reagent hood will open.
14. Remove the reagent tubes containing the Maintenance Wash Buffer, re-cap the reagent tubes, and store them at room temperature for subsequent maintenance wash procedure uses.
Note: The initial volume of 30–35 ml will be sufficient for 6–7 maintenance washes, and a re-fill will be required when the liquid level in the maintenance wash reagent tubes falls below the conical mark. Use Maintenance Wash Buffer bottle on position 11 for this re-fill: prepare new when the liquid level in the maintenance wash reagent tubes falls below the conical mark.

15. Remove both 1 liter Maintenance Wash Buffer bottles, re-cap the bottles, and store them at room temperature for subsequent maintenance wash procedure uses.


Note: The initial volumes of 300 ml (wash buffer position 11) and 500 ml (wash buffer position 9), respectively, will be sufficient for 9 maintenance washes: prepare new when the Maintenance Wash Buffer bottle volume is less than 100 ml.

16. Load 50 ml conical tubes containing 30–35 ml DI water onto positions 1, 2, 4, 5, 6 and 8 of the cooling compartment.




17. Click **Done** twice when finished.
18. Click **Next**.
19. Empty the waste bottle and dispose the contents according to local safety regulations. Reload the empty waste bottle into the instrument. Close the reagent hood and click **Next**.


IMPORTANT: Make sure to reload the waste bottle into the instrument in the correct orientation. An incorrect orientation may damage the instrument. Do not use force to close the reagent hood. If blocked, check the waste bottle orientation and correct if necessary before trying to close the hood again.



20. The instrument will now start the water wash. The DI water wash will take about 40 minutes. After the final wash is complete, click **Next**.




21. The Flow Cell door will open. Unload the Flow Cell, close the Flow Cell door, and click **Next**.



22. The reagent hood will open. Unload the reagent and wash kit. Empty the waste bottle. Close the reagent hood and click **Done**.

Note: If the GeneReader is not used immediately for another run, leave the reagent and wash bottles containing DI water in the system. This is important to ensure that the fluidic lines of the system do not dry out.

Note: If no run is planned for longer than 1 month, perform the long-term storage procedure, to prepare the GeneReader instrument for dry storage conditions (see Section 7.2.3).



23. Click **Finish** to close the wizard.

7.2.3 Preparing the GeneReader fluidics lines for long-term storage

If the GeneReader instrument will be stored for a prolonged period of time (i.e., no runs planned for longer than 1 month), the GeneReader must be prepared for storage by executing a DI water wash and the emptying of the fluidic lines. The following procedure describes this process:

1. Fill all reagent containers and buffer bottles with DI water.
2. Perform a routine DI water wash and wait until the run has finished.
3. Remove all liquid from all buffer bottles and reagent containers, empty the waste bottle, and place the empty bottles, the empty reagent containers, and the empty waste bottle back onto the instrument.
4. Run the routine DI water wash protocol again, but with empty bottles, empty reagent containers, and empty waste bottle, and wait until the run has finished.

The instrument will now be free of liquids and can be stored in this state.

7.2.4 Monthly cleaning procedure

Thoroughly wipe the inside and outside of the GeneReader using the cleaning agents described in Section 7.3.1.

IMPORTANT: Do not use alcohol or alcohol-based disinfectants or detergents to decontaminate the GeneReader hood.

7.3 General cleaning procedures

7.3.1 Cleaning agents

The following disinfectants and detergents are recommended for cleaning the GeneReader instrument:

Mikrozid® Wipes (Schülke & Mayr GmbH; www.schuelke-mayr.com) – moistened with ethanol-based disinfectant for wiping alcohol-insensitive surfaces of the GeneReader instrument.

Mikrozid Sensitive Liquid (Schülke & Mayr GmbH; www.schuelke-mayr.com) – quaternary ammonium salt-based disinfectant for alcohol-sensitive surfaces (consists of 0.26 g quaternary ammonium compounds, benzyl-C12-C16-alkyldimethyl, chlorides; 0.26 g Didecyldimethyl-ammonium chloride and 0.26 g quaternary ammonium compounds, benzyl-C12-C14-alkyl[(ethylphenyl)methyl]dimethyl, chlorides per 100 g Mikrozid Sensitive Liquid).

IMPORTANT: Do not use alcohol or alcohol-based disinfectants to clean the GeneReader instrument hood. Exposure of the GeneReader hood to alcohol or alcohol-based disinfectants will cause surface cracking. Clean the GeneReader hood with distilled water or Mikrozid Sensitive Liquid only.

Note: If you would like to use different disinfectants from those recommended, ensure that their compositions are similar to those described above. A suitable alternative to Mikrozid Liquid is Incidin® Liquid (Ecolab; www.ecolab.com).

7.3.2 General instructions

- Do not use spray bottles to spray cleaning liquids onto surfaces of the GeneReader instrument.
- If solvents or saline, acidic or alkaline solutions are spilled on the GeneReader, or if QIAGEN buffers splash the instrument door, wipe the spilled liquid away immediately.
- Follow manufacturer's safety instruction for handling cleaning agents.
- Follow manufacturer's instruction for concentration of the cleaning agents.
- Do not use alcohol or alcohol-based disinfectants to clean the GeneReader hood. Exposing the GeneReader hood, Main Hood, and flow cell door to alcohol or alcohol-based disinfectants will cause surface cracking. Clean the GeneReader hood with distilled water or Mikrozid Sensitive Liquid only.

<p>WARNING</p> 	<p>Risk of electric shock (W9)</p> <p>Do not open any panels on the GeneReader.</p> <p>Risk of personal injury and material damage</p> <p>Only perform maintenance that is specifically described in this user manual.</p>
---	--

<p>WARNING</p> 	<p>Hazardous chemicals and infectious agents (W13)</p> <p>The waste contains samples and reagents. This waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.</p>
---	--

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W1)</p> <p>Improper use of the GeneReader may cause personal injuries or damage to the instrument.</p> <p>The GeneReader must only be operated by qualified personnel who have been appropriately trained.</p> <p>Servicing of the GeneReader instrument must only be performed by a QIAGEN Field Service Specialist.</p>
---	--

7.3.3 Servicing


Contact QIAGEN Technical Services or your local distributor for more information about flexible Service Support Agreements from QIAGEN (see back cover).

7.4 Cleaning the workstation hard disk

A sequencing run requires a large amount of hard disk space. Therefore, the GeneReader software provides a “Manage Files” Wizard to assist with backup tasks in case GeneRead Databank is not used as an external data management tool. To perform the “manage files” process, additional free disk space for the backup is required, such as an external hard drive or network drive.


7.4.1 Procedure

1. Click **Manage Files** in the **Maintenance** menu to open the wizard. The wizard will display all previously run directories.




2. Select the folders you wish to backup by setting the respective action in the **Copy ...** column. If any folder is chosen for backup, the destination drive needs to be set. Select the folders to be cleaned up by setting the respective action in the **Delete ...** column.

IMPORTANT: It is possible to select folders for clean up without adding them to the list of files which will be backed up. Loss of data may occur if a manual backup has not been performed first!



- The wizard will then display a summary of all copy and delete tasks as shown below. Check the estimated free disk space. You may return to the previous screen to change your settings by clicking **Back**.




Note: “Manage Files” is deactivated in the **Maintenance** menu if GeneRead Databank is used as an external data management tool.


7.5 Creating a service archive folder for troubleshooting

7.5.1 Procedure

- Click **Service Archive** in the **Maintenance** menu to open the wizard for creating a service archive. The wizard will display all flow cells for which a service archive can be created.



- Select **Select For Archiving** for each flow cell for which a service archive will be created. All mandatory files for troubleshooting, except images, will be archived. To include images in the service archive, select either **All Images** to archive all images or **Special Images only** to archive only specific images of interest for troubleshooting.



3. Select a destination folder for the service archive by clicking the **Browse...** button.
4. Start the service archive creation by clicking the **Start Archiving** button.
5. Click **Close** to close the wizard.

8 Troubleshooting

This section provides information about what to do if an error occurs when using the GeneReader system.

If you need to contact QIAGEN Technical Services about an error, note the steps leading to the error and the information from any dialog boxes. This will help the QIAGEN Technical Service Specialist to resolve the error.

QIAGEN Technical Services may request a service archive folder. The service archive wizard can be started by selecting **Maintenance** and **Service Archive** (see Section 7.5).

8.1 Hardware and software errors

8.1.1 Application module

Error message	Comments and suggestions
The application cannot be closed because a run is ongoing.	Stop the run and unload all Flow Cells.
The application cannot be closed, because there are still Flow Cells in the instrument.	Unload all Flow Cells.

8.1.2 Manage Files module

Error message	Comments and suggestions
An error occurred while performing the operation.	Verify that the destination path is still available and that you are running the software with the required permissions to write data there.
Not all files that are selected for deletion are selected for backing up.	Verify that you wish to delete files without creating a backup.
The destination drive is too small to copy the data.	Choose a destination drive with sufficient available space or copy less data.
The disk drive for the source and the destination is the same.	If the source and destination drives are identical, this wizard will not be able to release disk space.
The selected destination path does not exist.	Create the specified directories.
The selected destination path is too long to copy the selected files.	Choose a less nested destination path or rename the folders. The path may not exceed 255 characters in length.
The selected destination path is write protected.	Verify that you are running the software with the required permissions to write data to the selected path.

8.1.3 Flow Cells

Error message	Comments and suggestions
Cannot start a run because a maintenance wash has to be performed first.	Perform a maintenance wash.
Cannot start a run because there is not enough space on the hard drive to store run data.	Run the Manage Files Wizard or use GeneRead Databank to free up disk space.
The wash cannot be started during a run.	Complete the run and unload process.
There are still flow cells on the carousel. Please unload these before starting a wash.	Unload all Flow Cells.
Validation of Flow Cell failed. Load the Flow Cell and ensure that the Flow Cell is placed correctly. Close the door to continue.	Verify that the Flow Cell orientation is correct and that the bar code is clean and undamaged. The bar code should be face up and to the left. Refer to the relevant QIAGEN kit handbook for more information.

8.1.4 Run module

Error message	Comments and suggestions
There must be at least one Flow Cell loaded to run a protocol.	Use the Start Run Wizard to load flow cells and fluids before starting a run.
Flow Cell door is open.	Close Flow Cell door.
No Flow Cell detected at the loader but Flow Cell <bar:code> was expected to be there.	Select Yes if this Flow Cell is at the loader. Select No to remove this Flow Cell from the inventory. Contact QIAGEN Technical Services if this message persists.

8.1.5 Start Run module

Error message	Comments and suggestions
This is not a valid bar code. Please re-scan the proper bar code from the reagent kit.	Verify that you are using a GeneRead UMI Sequencing Q Kit and contact QIAGEN Technical Services.
This is not a valid bar code. Please re-scan the proper bottle from the wash kit.	Verify that you are using GeneRead UMI Wash Buffer and contact QIAGEN Technical Services.

9 Technical Data

9.1 Environmental conditions – operating conditions

Power	100–240 V AC ($\pm 10\%$), 50/60 Hz Max. 600 VA (GeneReader) Max. 400 VA (workstation)
Overvoltage category	II
Air temperature	18 to 25°C (64 to 77°F)
Relative humidity	10–75% (noncondensing)
Altitude	Up to 2000 m (6500 ft.)
Place of operation	For indoor use only
Pollution level	2
Environmental class	3K2 (IEC 60721-3-3)

9.2 Mechanical data and hardware features

Dimensions (hoods closed)	Width: 93 cm (36.6 in.) Height: 61 cm (24.0 in.) Depth: 59 cm (23.2 in.)
Dimensions (hood open)	Width: 93 cm (36.6 in.) Height: 112 cm (44.1 in.) Depth: 59 cm (23.2 in.)
Dimensions (Main Hood open)	Width: 93 cm (36.6 in.) Height: 125 cm (49.3 in.) Depth: 59 cm (23.2 in.)
Weight	102 kg (224.9 lb.)
Capacity	Up to 3 Flow Cells simultaneously

9.3 Workstation specifications (hardware and software)

9.3.1 Workstation

- Workstation, type HP z420 or HP z440
- Display resolution of 1920 x 1080 pixels
- 4 x USB 3.0 (front), 2 x USB 2.0 (rear), and 4 x USB 3.0 (rear)
- DVD drive
- Keyboard and mouse
- Monitor
- Handheld bar code scanner

9.3.2 Software

- Microsoft Windows 7 Professional 64-bit (D/F/E/I)
- GeneReader software v1.6

Adobe Reader software is preinstalled on the computer to view reports generated in PDF format.

Appendix A

Declaration of Conformity

Name and address of the legal manufacturer:

QIAGEN GmbH
QIAGEN Strasse 1
40724 Hilden
Germany

An up-to-date Declaration of Conformity can be requested from QIAGEN Technical Services.

License Terms

The license terms for all software used with GeneReader, including QIAGEN software components, commercial software components and open source software components, are provided in the files **licenses.rtf** and **Prerequisite.LicenseAgreements.rtf** located on the GeneReader workstation under the following paths:

C:\ProgramData\QIAGEN\GeneReader\licenses.rtf

C:\ProgramData\QIAGEN\GeneReader\Prerequisite.LicenseAgreements.rtf.

Waste Electrical and Electronic Equipment (WEEE)

This section provides information about disposal of waste electrical and electronic equipment by users.

The crossed-out wheeled bin symbol (see below) indicates that this product must not be disposed of with other waste; it must be taken to an approved treatment facility or to a designated collection point for recycling, according to local laws and regulations.



The separate collection and recycling of waste electronic equipment at the time of disposal helps to conserve natural resources and ensures that the product is recycled in a manner that protects human health and the environment.

Recycling can be provided by QIAGEN upon request at additional cost. In the European Union, in accordance with the specific WEEE recycling requirements and where a replacement product is being supplied by QIAGEN, free recycling of its WEEE-marked electronic equipment is provided.

To recycle electronic equipment, contact your local QIAGEN sales office for the required return form. Once the form is submitted, you will be contacted by QIAGEN either to request follow-up information for scheduling collection of the electronic waste or to provide you with an individual quote.

If you are unclear about how to dispose of your QIAGEN equipment, contact QIAGEN Technical Services (see back cover).

FCC Declaration

The “United States Federal Communications Commission” (USFCC) (in 47 CFR 15. 105) declared that the users of this product must be informed of the following facts and circumstances.

‘This device complies with part 15 of the FCC as stated below:

“Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.”

“This Class B digital apparatus complies with Canadian ICES-0003.”

The following statement applies to the products covered in this manual, unless otherwise specified herein. The statement for other products will appear in the accompanying documentation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and meets all requirements of the Canadian Interference-Causing Equipment Standard ICES-003 for digital apparatus. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that the interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for help.

QIAGEN GmbH Germany is not responsible for any radio television interference caused by unauthorized modifications of this equipment or the substitution or attachment of connection cables and equipment other than those specified by QIAGEN GmbH, Germany. The correction of interference caused by such unauthorized modification, substitution or attachment will be the responsibility of the user.

Liability Clause

QIAGEN shall be released from all obligations under its warranty in the event repairs or modifications are made by persons other than its own personnel, except in cases where the Company has given its written consent to perform such repairs or modifications.

All materials replaced under this warranty will be warranted only for the duration of the original warranty period, and in no case beyond the original expiration date of original warranty unless authorized in writing by an officer of the Company. Read-out devices, interfacing devices, and associated software will be warranted only for the period offered by the original manufacturer of these products. Representations and warranties made by any person, including representatives of QIAGEN, which are inconsistent or in conflict with the conditions in this warranty shall not be binding upon the Company unless produced in writing and approved by an officer of QIAGEN.

Appendix B

Ordering Information

Note: Only use GeneReader products supplied by QIAGEN.

Product	Contents	Cat. no.
GeneReader Platform	Next-generation sequencing instrument: includes installation and training, 1 year warranty on parts and labor	9002312
Related Products		
GeneRead UMI Sequencing Q Kit (3)	Includes reagents, add-ons and 3 Flow Cells for as many as 3 Flow Cell runs, up to 120 cycles each, on the GeneReader	185251
Additional products		
Bar Code Scanner NGS System	Bar code scanner for the NGS system	9022959
Waste Bottle GeneReader	Wash bottle for the NGS system	9243790
Workstation NGS System, V2 W7, 64bit	Workstation for the NGS system	9024359

Appendix C

Safety Information Translated into German and French

This subsection contains translations of the warnings and cautions referenced in this user manual. Each warning or caution has a reference number in square brackets at the top right of its box.

Sicherheitsinformationen

In diesem Unterabschnitt befindet sich die Übersetzung der Warn- und Vorsichtshinweise, die in diesem Handbuch vorkommen. Jeder Warn- oder Vorsichtshinweis ist in der rechten oberen Ecke des Kastens mit einer Referenznummer in eckigen Klammern gekennzeichnet.

Informations de Sécurité

Ce paragraphe contient les traductions des avertissements et mises en garde utilisés dans ce manuel de l'utilisateur. Chaque avertissement ou mise en garde s'accompagne d'un numéro de référence entre crochets, en haut à droite de son encadré.

<p>WARNING</p> 	<p>The term WARNING is used to inform you about situations that could result in personal injury to you or other persons.</p> <p>Details about these circumstances are given in a box like this one.</p>
<p>DE</p>	<p>WARNING (WARNUNG)</p> <p>Warnung weist auf Situationen und Umstände hin, die zu einer Verletzung des Benutzers oder anderer Personen führen können.</p> <p>Nähere Angaben zu der Art der Gefährdung und der Vermeidung solcher Situationen werden in einem Textfeld wie diesem neben der Warnung gemacht.</p>
<p>FR</p>	<p>WARNING (DANGER)</p> <p>La formule WARNING (DANGER) est utilisée pour avertir des situations pouvant occasionner des dommages corporels à l'utilisateur ou à d'autres personnes. Les détails sur ces circonstances sont données dans un encadré semblable à celui-ci.</p>

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W1)</p> <p>Improper use of the GeneReader may cause personal injuries or damage to the instrument.</p> <p>The GeneReader should only be operated by qualified personnel who have been appropriately trained.</p> <p>Servicing of the GeneReader should only be performed by QIAGEN Instrument Service Specialists.</p>
<p>DE</p>	<p>Verletzungsgefahr und Gefahr der Beschädigung des Geräts</p> <p>Die unsachgemäße Bedienung des GeneReader kann zu einer Verletzung des Benutzers oder zur Beschädigung des Geräts führen.</p> <p>Die Bedienung des GeneReader darf nur durch qualifiziertes, entsprechend geschultes Personal erfolgen.</p> <p>Wartungsarbeiten am GeneReader sollten nur durch Mitarbeiter des QIAGEN Kundendienstes durchgeführt werden.</p>
<p>FR</p>	<p>Risque d'accident corporel et de détérioration du matériel</p> <p>L'utilisation inappropriée du GeneReader peut provoquer des accidents corporels ou une détérioration de l'appareil.</p> <p>Le GeneReader ne doit être utilisé que par du personnel qualifié ayant été convenablement formé.</p> <p>Seul un ingénieur du service après-vente QIAGEN est autorisé à effectuer des travaux d'entretien sur le GeneReader.</p>

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W2)</p> <p>The GeneReader is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.</p> <p>Contact QIAGEN Technical Services to relocate the instrument.</p>
<p>DE</p>	<p>Gefahr von Verletzungen und Sachbeschädigungen</p> <p>Der GeneReader ist sehr schwer und darf nicht von einer Person angehoben werden. Heben Sie das Gerät nicht allein an, um eine Verletzung und/oder Beschädigung des Geräts zu vermeiden.</p> <p>Wenn Sie das Gerät umstellen wollen, wenden Sie sich bitte an den technischen Service von QIAGEN.</p>
<p>FR</p>	<p>Risque d'accident corporel et de dommages matériels</p> <p>L'instrument GeneReader est trop lourd pour être soulevé par une seule personne. Afin d'éviter tout accident corporel et toute détérioration du matériel, ne pas soulever l'appareil seul.</p> <p>Contactez les services techniques de QIAGEN pour déplacer l'instrument.</p>

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W3)</p> <p>Do not attempt to move the GeneReader during operation.</p>
<p>DE</p>	<p>Verletzungsgefahr und Beschädigung des Gerätes</p> <p>Den GeneReader während eines Laufes nicht bewegen.</p>
<p>FR</p>	<p>Risque de dommages corporels et matériels</p> <p>Ne pas essayer de bouger le GeneReader pendant son fonctionnement.</p>

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W4)</p> <p>Load flow cell only in accordance with step-by-step instructions provided by GeneReader software. Beware of moveable parts.</p>
<p>DE</p>	<p>Gefahr von Verletzungen und Sachbeschädigungen</p> <p>Beladen Sie die Durchflusszelle nur gemäß der schrittweisen Anweisungen der GeneReader Software. Achten Sie auf bewegliche Teile.</p>
<p>FR</p>	<p>Risque d'accident corporel et de dommages matériels</p> <p>Chargez la cellule de débit en suivant exclusivement les instructions fournies par le logiciel GeneReader. Attention aux pièces mobiles.</p>

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W5)</p> <p>Do not stare into the beam of the flow cell bar code reader.</p>
<p>DE</p>	<p>Gefahr von Verletzungen und Sachbeschädigungen</p> <p>Blicken Sie nicht absichtlich den Strahl des Barcode-Lesers in der Durchflusszelle.</p>
<p>FR</p>	<p>Risque d'accident corporel et de dommages matériels</p> <p>Ne regardez pas directement dans le faisceau du lecteur de code-barres de la cellule de débit.</p>

<p>WARNING</p> 	<p>Fire hazard (W6)</p> <p>Empty the liquid waste bottle before each run and make sure to place it in the correct orientation back in the GeneReader instrument. Spilling of liquid-waste may cause an electrical short-circuit and fire.</p>
--	--

DE	<p>Brandgefahr</p> <p>Leeren Sie die Flüssigabfallflasche vor jedem Lauf und achten Sie darauf, die Flasche korrekt ausgerichtet wieder in das GeneReader Gerät zu setzen. Ein Verschütten von Flüssigabfall kann zu elektrischen Kurzschlüssen und Bränden führen.</p>
FR	<p>Risque d'incendie</p> <p>Videz la bouteille de déchets liquides avant chaque utilisation et assurez-vous de la replacer dans la bonne position dans l'instrument GeneReader. Évitez de renverser les déchets liquides, cela pourrait entraîner un court-circuit électrique et un incendie.</p>

<p>WARNING</p> 	<p>Electrical hazard (W7)</p> <p>Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.</p> <p>Intentional interruption is prohibited.</p> <p>Lethal voltages inside the instrument</p> <p>When the instrument is connected to line power, terminals may be live and opening covers or removing parts is strongly discouraged as doing so may result in electric shock or death.</p>
DE	<p>Stromschlaggefahr</p> <p>Jede Unterbrechung des Schutzleiters (Erdungs- bzw. Masseleiter) im Gerät oder außerhalb des Geräts und jede Abtrennung vom Schutzleiteranschluss erhöht die Gefahr eines Stromschlags.</p> <p>Eine absichtliche Unterbrechung der Schutzleiter ist verboten.</p> <p>Gefährliche Spannung im Gerät</p> <p>Wenn das Gerät an die Stromversorgung angeschlossen ist, sind die Anschlussstellen unter Umständen spannungsführend. Es wird dringend vom Öffnen der Abdeckungen und dem Entfernen von Gehäuseteilen abgeraten, da Stromschlag- und Todesgefahr besteht.</p>

FR	<p>Risque électrique</p> <p>Toute interruption du conducteur de protection (conducteur de terre/de masse) à l'intérieur ou à l'extérieur de l'instrument ou toute déconnexion de la borne du conducteur de protection est susceptible de rendre l'instrument dangereux.</p> <p>Toute interruption intentionnelle est interdite.</p> <p>Tensions mortelles à l'intérieur de l'instrument</p> <p>Lorsque l'instrument est relié à l'alimentation, les bornes peuvent être sous tension et l'ouverture de capots ou le retrait d'éléments est fortement déconseillé car cela risque d'entraîner un choc électrique ou la mort.</p>
----	---

<p>WARNING</p> 	<p>Damage to electronics (W8)</p> <p>Before powering ON the instrument make sure that the correct supply voltage is used.</p> <p>Incorrect use of supply voltage may cause damage to electronics.</p> <p>See specifications indicated on the type plate of the instrument.</p>
--	--

DE	<p>Beschädigung von elektronischen Bauteilen</p> <p>Stellen Sie sicher, dass Sie den für Ihr Land bestimmten Netzteiladapter verwenden, bevor Sie das Netzteil mit dem Controller verbinden.</p> <p>Der Betrieb mit falscher Netzspannung/-Frequenz kann zu Schäden in der Elektronik führen.</p> <p>Bitte beachten Sie auch die entsprechenden Spezifikationen auf dem Typenschild des Geräts</p>
----	---

FR	<p>Endommagement des composants électroniques</p> <p>Veiller à utiliser le bon adaptateur d'alimentation électrique avant de brancher le câble d'alimentation électrique à l'instrument.</p> <p>Une tension d'alimentation incorrecte peut endommager les composants électroniques.</p> <p>Voir les spécifications indiquées sur la plaque signalétique de l'instrument.</p>
----	---

<p>WARNING</p> 	<p>Risk of electric shock (W9)</p> <p>Do not open any panels on the GeneReader.</p> <p>Risk of personal injury and material damage</p> <p>Only perform maintenance that is specifically described in this user manual.</p>
--	--

DE	<p>Gefährdung durch Elektrizität</p> <p>Unter keinen Umständen darf das Gehäuse des GeneReader geöffnet werden.</p> <p>Verletzungsgefahr und Beschädigung des Gerätes</p> <p>Keine Pflege- und Wartungsarbeiten durchführen, die nicht in diesem Handbuch beschrieben sind.</p>
FR	<p>Risque d'électrocution</p> <p>Ne pas ouvrir les panneaux du GeneReader.</p> <p>Risque de dommages personnels et matériels</p> <p>Réaliser uniquement la maintenance décrite spécifiquement dans ce manuel.</p>

<p>WARNING</p> 	<p>Explosive atmosphere (W10)</p> <p>The GeneReader is not designed for use in an explosive atmosphere.</p>
DE	<p>Explosionsfähige Atmosphäre</p> <p>Der GeneReader ist nicht für den Betrieb in einer explosionsfähigen Atmosphäre ausgelegt.</p>
FR	<p>Atmosphère explosive</p> <p>Le GeneReader n'est pas conçu pour être utilisé dans une atmosphère explosive.</p>

<p>WARNING</p> 	<p>Risk of explosion (W11)</p> <p>The GeneReader is intended for use with reagents and substances supplied with QIAGEN kits. Use of other reagents and substances may lead to fire or explosion.</p>
DE	<p>Explosionsgefahr</p> <p>Der GeneReader ist ausschließlich mit Reagenzien und Substanzen aus den QIAGEN Kits zu benutzen. Die Benutzung von anderen Reagenzien oder Substanzen kann Feuer oder eine Explosion auslösen.</p>

FR	<p>Risque d'explosion</p> <p>Le GeneReader a été conçu pour l'utilisation des réactifs et substances fournis par les kits QIAGEN. L'utilisation de réactifs et de substances autres que celles indiquées peut entraîner un risque d'incendie ou d'explosion.</p>
----	---

<p>WARNING</p> 	<p>(W12)</p>	<p>Hazardous chemicals</p> <p>Some chemicals used with this instrument may be hazardous or may become hazardous after completion of the protocol run.</p> <p>Always wear safety glasses, gloves, and a lab coat.</p> <p>The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe and that the instrument operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable Material Safety Data Sheets (MSDSs) or OSHA, ACGIH, or COSHH documents.</p> <p>Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.</p>
DE	<p>Gefährliche Chemikalien</p> <p>Einige der in Verbindung mit diesem Gerät verwendeten Chemikalien sind gesundheitsgefährdend oder können nach Beendigung eines Protokoll-Durchlaufes gesundheitsgefährdend werden.</p> <p>Es sollten immer Sicherheitsbrille, Handschuhe und ein Laborkittel getragen werden.</p> <p>Der Betreiber der Anlage ist für die Gewährleistung der Sicherheit am Arbeitsplatz verantwortlich. Er hat sicherzustellen, dass die Bediener des Gerätes ausreichend geschult sind und nicht gesundheitsgefährdenden Konzentrationen toxischer Substanzen (chemischer oder biologischer) ausgesetzt sind, so wie dies in den Sicherheitsdatenblättern oder in anderen zu beachtenden Dokumenten festgelegt ist.</p> <p>Bei der Behandlung von Abluft und bei der Abfallbeseitigung sind alle gesetzlichen Regelungen zur Gesundheit und Sicherheit auf nationaler, regionaler und lokaler Ebene zu berücksichtigen.</p>	

FR	<p>Substances chimiques dangereuses</p> <p>Certaines substances chimiques utilisées avec cet instrument peuvent être dangereuses ou peuvent le devenir après que le protocole ait été effectué.</p> <p>Toujours porter des lunettes de protection, paire de gants et une blouse de laboratoire.</p> <p>La personne responsable (par exemple le Chef du laboratoire) doit prendre les précautions nécessaires pour assurer la sécurité de l'environnement du poste de travail et pour être sûr que les opérateurs de l'instrument sont suffisamment formés et non exposés à des quantités dangereuses de substances toxiques (chimique ou biologique) comme défini dans "Material Safety Data Sheets (MSDS)" ou des documents "OSHA, ACGIH ou COSHH".</p> <p>L'évacuation des vapeurs et déchets doit être conforme à tous règlements et dispositions légales - au plan national, départemental et local - concernant la santé et la sécurité.</p>
----	--

<p>WARNING</p> 	<p>Hazardous chemicals and infectious agents (W13)</p> <p>The waste contains samples and reagents. This waste may contain toxic or infectious material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.</p>
DE	<p>Gefährliche Chemikalien und infektiöse Agenzien</p> <p>Der Abfall kann gesundheitsgefährdende Reagenzien oder infektiöses Probenmaterial enthalten und muss gemäß den lokalen Sicherheitsvorschriften entsorgt werden.</p>
FR	<p>Agents infectieux et substances chimiques dangereuses</p> <p>Les déchets contiennent des échantillons et des réactifs. Ils peuvent contenir des agents infectieux et doivent être éliminés selon les règles de sécurité de votre laboratoire.</p>

<p>WARNING</p> 	<p>Moving parts (W14)</p> <p>To avoid contact with moving parts during operation of the GeneReader, the instrument must be operated with the flow cell door closed.</p> <p>If the door sensor is not functioning correctly, contact QIAGEN Technical Services.</p>
--	---

DE	<p>Sich bewegende Geräteteile</p> <p>Um beim Betrieb des GeneReader einen Kontakt mit sich bewegenden Teilen zu vermeiden, darf das Gerät nur mit geschlossener Durchflusszellentür betrieben werden.</p> <p>Falls die Türsensoren nicht korrekt funktionieren, kontaktieren Sie bitte den technischen Service von QIAGEN.</p>
FR	<p>Pièces mobiles</p> <p>Pour éviter tout contact avec des pièces en mouvement pendant le fonctionnement du GeneReader, l'instrument doit être utilisé avec le capot de la cellule de débit fermé.</p> <p>Si le capteur du capot ne fonctionne pas correctement, contactez les services techniques de QIAGEN.</p>

<p>WARNING</p> 	<p>Risk of overheating (W15)</p> <p>To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the GeneReader.</p> <p>Slits and openings that ensure the ventilation of the GeneReader must not be covered.</p>
DE	<p>Überhitzung des Gerätes</p> <p>Zur Sicherstellung einer ausreichenden Belüftung des GeneReader muss ein Mindestabstand von 10 cm an den Seiten und an der Rückseite des Gerätes eingehalten werden.</p> <p>Lüftungsschlitze und -öffnungen des Gerätes nicht abdecken.</p>
FR	<p>Risque de surchauffe</p> <p>Laisser un espace d'au moins 10 cm sur les côtés et à l'arrière du GeneReader pour assurer une ventilation efficace.</p> <p>Les grilles et prises d'air assurant la ventilation du GeneReader ne doivent pas être couvertes.</p>

<p>WARNING</p> 	<p>Risk of personal injury and material damage (W16)</p> <p>Only perform maintenance that is specifically described in this user manual.</p>
DE	<p>Verletzungsgefahr und Beschädigung des Gerätes</p> <p>Keine Pflege- und Wartungsarbeiten durchführen, die nicht in diesem Handbuch beschrieben sind.</p>

FR	<p>Risque de dommages personnels et matériels</p> <p>Réaliser uniquement la maintenance décrite spécifiquement dans ce manuel.</p>
----	---

<p>CAUTION</p> 	<p>The term CAUTION is used to inform you about situations that could result in damage to the TissueRuptor II or other equipment.</p> <p>Details about these circumstances are given in a box like this one.</p>
DE	<p>CAUTION (ACHTUNG)</p> <p>ACHTUNG weist auf Situationen und Umstände hin, die zu einer Beschädigung des Gerätes führen können. Um einen Geräteschaden zu vermeiden, muss die genannte Anleitung unbedingt befolgt werden.</p> <p>Nähere Angaben zu der Art der Gefährdung und der Vermeidung solcher Situationen werden in einem Textfeld wie diesem gemacht.</p>
FR	<p>CAUTION (ATTENTION)</p> <p>Le terme CAUTION (Attention) est utilisé pour signaler les situations susceptibles de provoquer des détériorations de l'instrument ou d'autre matériel.</p> <p>Les détails sur ces circonstances figurent dans un encadré semblable à celui-ci.</p>

<p>CAUTION</p> 	<p>Risk of material damage (C1)</p> <p>Avoid moving the workbench and causing vibrations to the GeneReader during operation to prevent disturbing sensitive optical measurements.</p>
DE	<p>Gefahr von Materialbeschädigungen</p> <p>Während des Betriebs des GeneReader sollten Sie eine Bewegung der Arbeitsbank und Vibrationen vermeiden, damit die empfindlichen optischen Messungen nicht gestört werden.</p>
FR	<p>Risque de dommages matériels</p> <p>Évitez de déplacer le plan de travail et de créer des vibrations dans le GeneReader en fonctionnement, afin de ne pas fausser les mesures optiques sensibles.</p>

<p>CAUTION</p> 	<p>Damage to the instrument (C2)</p> <p>Avoid spilling water or chemicals onto the GeneReader. Damage caused by water or chemical spillage will void your warranty.</p>
--	--

DE	<p>Beschädigung des Gerätes</p> <p>Vermeiden Sie es, Wasser oder Chemikalien auf dem GeneReader zu verschütten. Durch verschüttetes Wasser oder verschüttete Chemikalien verursachte Geräteschäden sind nicht durch die Garantie abgedeckt.</p>
FR	<p>Détérioration de l'appareil</p> <p>Eviter de renverser de l'eau ou des substances chimiques sur le GeneReader. Tout dommage causé par de l'eau ou des produits chimiques mettra fin à la garantie.</p>

<p>CAUTION</p> 	<p>Risk of material damage (C3)</p> <p>Do not place any items on top of the instrument.</p>
DE	<p>Gefahr von Materialbeschädigungen</p> <p>Stellen Sie keine Gegenstände auf das Gerät.</p>
FR	<p>Risque de dommages matériels</p> <p>Ne pas placer ou fixer d'objets sur l'instrument.</p>

<p>CAUTION</p> 	<p>Damage to the instrument (C4)</p> <p>Only use QIAGEN consumables with the GeneReader. Damage caused by use of other types of consumables will void your warranty.</p>
DE	<p>Beschädigung des Geräts</p> <p>Verwenden Sie nur QIAGEN Verbrauchsartikel für den GeneReader. Kommt es bei der Verwendung von Verbrauchsartikeln anderer Hersteller zu Geräteschäden, erlischt die Garantie.</p>
FR	<p>Endommagement de l'instrument</p> <p>Utiliser seulement des consommables QIAGEN avec l'instrument GeneReader. Les détériorations causées par l'utilisation d'autres types de consommables annulent la garantie.</p>

<p>CAUTION</p> 	<p>Damage to the instrument (C5)</p> <p>Make sure that the flow cell is inserted in the correct position. Incorrect insertion of the flow cell can damage the instrument.</p>
--	--

DE	<p>Beschädigung des Geräts</p> <p>Achten Sie darauf, dass die Durchflusszelle in die korrekte Position eingesetzt wurde. Ein falsches Einsetzen der Durchflusszelle kann eine Beschädigung des Geräts zur Folge haben.</p>
FR	<p>Détérioration de l'instrument</p> <p>Assurez-vous que la cellule de débit a été insérée dans la bonne position. L'insertion incorrecte de la cellule de débit peut entraîner la détérioration de l'instrument.</p>

<p>CAUTION</p> 	<p>Damage to the instrument (C6)</p> <p>Direct sunlight may bleach parts of the instrument and cause damage to plastic parts.</p> <p>The GeneReader must be located out of direct sunlight.</p>
DE	<p>Beschädigung des Gerätes</p> <p>Direktes Sonnenlicht kann Teile des Gerätes bleichen und Plastikteile schädigen. Der GeneReader darf nicht ins direkte Sonnenlicht gestellt werden.</p>
FR	<p>Détérioration de l'instrument</p> <p>La lumière directe du soleil peut décolorer des parties de l'instrument et endommager des parties en plastique. Placer le GeneReader en dehors de la lumière directe du soleil.</p>

<p>CAUTION</p> 	<p>Damage to the instrument (C7)</p> <p>Do not use bleach, solvents or reagents containing acids, alkalis or abrasives to clean the GeneReader.</p>
DE	<p>Beschädigung des Gerätes</p> <p>Zur Reinigung des GeneReader darf keine Chlorbleiche, Lösungsmittel, saure oder alkalische Reagenzien oder Scheuermittel verwendet werden.</p>
FR	<p>Détérioration de l'instrument</p> <p>Ne pas utiliser de l'eau de Javel, des produits solvants ou abrasifs, ou de réactifs contenant des acides ou des solutions alcaliques pour nettoyer le GeneReader.</p>

Index

- Cautions 11
- Cleaning the workstation hard disk 61
- Cooling air outlet 24
- Creating a service archive folder 62
- Environmental conditions 66
- External features 22
- File handling 41
- Flow cell door 22
- Flow cell setup 41
- Fluidic drawer 24
- General information 9
- GeneReader Flow Cell 25
- GeneReader principle 21
- GeneReader Sample to Insight workflow 20
- Getting started 31
- Hood 23
- Intended use 9
- Internal features 25
- Loading and running the GeneReader 43
- Long-term storage 58
- Maintenance equipment and reagents 48
- Maintenance wash 48
- Mechanical data and hardware features 66
- Monthly cleaning 59
- Operating procedures 34
- Ordering information 73, 74
- Power cord socket 24
- Power requirements 28
- Power switch 24
- Requirements 10
- Routine maintenance 48
- Run finished 43
- Safety 11
- Safety information
 - chemical 16
 - electrical safety 14
 - environment 15
 - maintenance safety 17
 - mechanical hazards 17
 - proper use 12
 - symbols on the GeneReader 18
 - waste disposal 16
- Site requirements 27
- Software license 69
- Software workflow 35
- Software-related errors 64
- Staggered loading of flow cells 43
- Status lights 22
- Status of Configuration 39
- Status of Flow Cells 37
- Status of Fluids 38
- System
 - installation 27
- System delivery and installation 27
- Technical assistance 9
- Technical data 66
- Troubleshooting 64
- Unloading reagents and flow cells 43
- USB ports 24
- User interface 36
- Using the GeneReader software 35
- Warnings 11
- Waste Electrical and Electronic Equipment (WEEE) 70
- Weekly maintenance 51
- Workflow procedures 41
- Workstation equipment 25
- Workstation requirements 30
- Workstation specifications 67

Trademarks: QIAGEN[®], Sample to Insight[®], QIAGEN GeneRead[®], QIAGEN Clinical Insight[™], QCI[™], GeneRead QIAcube[®], GeneRead[™] (QIAGEN Group); GeneReader[™], QIAGEN GeneReader[®] (Intelligent Bio-Systems, Inc.); Adobe[®], Reader[®] (Adobe Systems Incorporated); Microsoft[®], Windows[®] (Microsoft Corporation); Mikrozid[®] (Schülke & Mayr GmbH); Incidin[®] (Ecolab, Inc.); Intel[®], Xeon[®] (Intel Corporation); NVIDIA[®], Quadro[®] (NVIDIA Corporation).
Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. HB2503-001 12/2017

© 2017 QIAGEN, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com