

Automated purification of mycoplasma DNA from respiratory samples

Martin Vondracek, Kristina Fahlander, and Tonje Fossheim*

Karolinska University Hospital, Solna, Stockholm, Sweden

*QIAGEN AS, Oslo, Norway

Mycoplasma DNA was purified from 4 different kinds of respiratory research samples using the BioRobot® M48 workstation and the MagAttract® DNA Mini M48 Kit with the QIAGEN Supplementary Protocol "Purification of bacterial DNA from primary samples using the MagAttract DNA Mini M48 Kit". DNA purified from nasopharyngeal swabs, nasopharyngeal secretions, sputum, and tracheal secretions performed well in real-time PCR, and showed identical results when compared to DNA purified by an established manual method.

Mycoplasmas are small, atypical bacteria with flexible cell walls and can be very difficult to detect and identify using traditional culture-based techniques. *Mycoplasma pneumoniae* is the pathogen most commonly associated with primary atypical pneumonia in humans, and is a common cause of both community and hospital acquired lung infections in children and immunocompromised individuals. Infection is often asymptomatic but can result in life-threatening pneumonia or long-term sequelae. In order to ensure accurate life science and clinical research and effective development of diagnostic tests, a rapid and standardized method of processing a range of respiratory research sample types is required.

Table 1. Comparison of Manual and Automated DNA Purification for Quantitative Amplification of *M. pneumoniae* DNA by Real-Time PCR

Sample number	Sample type [†]	Real-time PCR analysis*		
		Manual DNA purification and in-house PCR: <i>M. pneumoniae</i>	Automated DNA purification and real-time PCR: <i>M. pneumoniae</i>	Automated DNA purification and real-time PCR: control amplicon
1	NG-swab	+	+	+
2	NG-swab	+	+	+
3	NG-swab	+	+	+
4	NG-swab	+	+	+
5	NG-swab	+	+	+
6	NG-swab	+	+	+
7	NG-fluid	+	+	+
8	Sputum	+	+	+
9	Sputum	-	-	+
10	ET-fluid	-	-	+
11	Positive control	+	+	+
12	Negative control	-	-	+

* + = amplicon detected; - = amplicon not detected.

† NG-swab = Nasopharyngeal swab, NG-fluid = Nasopharyngeal fluid, ET-fluid = Endotracheal fluid.



Material and methods

Four different respiratory sample types (Table 1) were processed. Nasopharyngeal secretions, sputum, and tracheal secretions were processed directly. Nasopharyngeal swabs were incubated in 1x PBS for 30 minutes at room temperature. Samples (200 μ l) were processed according to the Bact_200ul protocol using the MagAttract DNA Mini M48 Kit and the App. Package, M48, Infectious Disease. DNA was eluted in 100 μ l water.

M. pneumoniae DNA was amplified using 5 μ l purified DNA as template in 25 μ l real-time PCR spiked with internal positive control (IPC) DNA. A *M. pneumoniae*-specific probe (FAM™) and a probe specific for the internal amplification control (JOE™) were used to detect and identify amplified DNA. A true negative test result is confirmed by a positive control-amplification. Negative results in both test and control amplifications indicate the presence of PCR inhibitors.

Results

Real-time PCR results (Figures 1 and 2) show that automated DNA purification using the BioRobot M48 workstation results in high-performance DNA that is well-suited for sensitive real-time PCR. *M. pneumoniae*-specific amplicons were detected in all positive samples, but not in negative and control samples (Table 1), indicating both cross-contamination-free purification and that purified DNA was not degraded. Internal positive controls were detected in all samples and showed very similar C_T profiles (Figure 2), indicating that purified DNA was free of PCR inhibitors (Figure 1).

Reliable Amplification of Mycoplasma DNA from Respiratory Samples




Figure 1 Real-time PCR of *M. pneumoniae* DNA using purified DNA representative of purifications from 3 different types of respiratory sample containing *M. pneumoniae*. **Red:** 6 x nasopharyngeal swab; **Blue:** nasopharyngeal fluid; **Green:** sputum; **Yellow:** endotracheal fluid; **Broken lines:** 2 x positive control samples isolated with Bact_200ul protocol on the BioRobot M48 workstation; **Gray:** negative controls.

Reliable Real-Time PCR Controls Using Inhibitor-Free Purified DNA




Figure 2 Internal positive control DNA was amplified efficiently from all samples and no significant variation was observed in C_T values, indicating purified DNA was consistently free from PCR inhibitors. **Red:** nasopharyngeal swab samples; **Blue:** nasopharyngeal fluid; **Yellow:** endotracheal fluid; **Green:** sputum; **Gray:** negative controls; **Broken lines:** positive control isolated with Bact_200ul protocol on the BioRobot M48 workstation.

Conclusions

High performance *M. pneumoniae* DNA was purified from a range of respiratory research samples, using the Bact_200ul protocol on the BioRobot M48 workstation, in combination with the MagAttract DNA Mini M48 Kit and the App. Package, M48, Infectious Disease. Using this DNA in sensitive downstream applications, such as quantitative real-time PCR, allows:

- **Standardized amplification of *M. pneumoniae* DNA** — from respiratory samples and swabs
- **Rapid purification** — with no enrichment culture
- **Low CVs and reliable controls** — using inhibitor-free DNA

Ordering Information

Product	Contents	Cat. no.
MagAttract DNA Mini M48 Kit (192)	For 192 DNA preps: MagAttract Suspension B, Buffers, Proteinase K	953336
App. Package, M48, Infectious Disease	Software protocol package for infectious disease applications, v. 2.0, on the BioRobot M48 workstation	9016145
BioRobot M48	Robotic workstation for automated purification of nucleic acids using MagAttract M48 Kits; installation, 1-year warranty on parts and labor	9000708

Contact QIAGEN today and discover efficient purification of DNA from respiratory samples!

Trademarks: QIAGEN®, BioRobot®, MagAttract® (QIAGEN Group); FAM™, JOE™ (Applera Corp.). QIAGEN robotic systems are not available in all countries; please inquire. The BioRobot M48 workstation and MagAttract Kits are intended as general purpose devices. No claim or representation is intended for their use in identifying any specific organism or for a specific clinical use. It is the user's responsibility to validate the performance of the BioRobot M48 workstation and MagAttract Kits for any particular use, since performance characteristics have not been validated for any specific organism. The BioRobot M48 workstation and MagAttract Kits may be used in clinical diagnostic laboratory systems after the laboratory has validated their complete system as required by CLIA '88 regulations in the U.S. or equivalents in other countries. The PCR process is covered by U.S. Patents 4,683,195 and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG. 1026536 06/2004 © 2004 QIAGEN, all rights reserved.

Australia = QIAGEN Pty Ltd = Orders 03-9840-9800 = Fax 03-9840-9888 = Technical 1-800-243-066
Belgium = QIAGEN Benelux B.V. = Orders 0800-79612 = Fax 0800-79611 = Technical 0800-79556
Canada = QIAGEN Inc. = Orders 800-572-9613 = Fax 800-713-5951 = Technical 800-DNA-PREP (800-362-7737)
France = QIAGEN S.A. = Orders 01-60-920-920 = Fax 01-60-920-925 = Technical 01-60-920-930
Germany = QIAGEN GmbH = Orders 02103-29-12000 = Fax 02103-29-22000 = Technical 02103-29-12400
Italy = QIAGEN S.p.A. = Orders 02-33430411 = Fax 02-33430426 = Technical 02-33430414
Japan = QIAGEN K.K. = Telephone 03-5547-0811 = Fax 03-5547-0818 = Technical 03-5547-0811
Luxembourg = QIAGEN Benelux B.V. = Orders 8002-2076 = Fax 8002-2073 = Technical 8002-2067
The Netherlands = QIAGEN Benelux B.V. = Orders 0800-0229592 = Fax 0800-0229593 = Technical 0800-0229602
Switzerland = QIAGEN AG = Orders 061-319-30-30 = Fax 061-319-30-33 = Technical 061-319-30-31
UK and Ireland = QIAGEN Ltd. = Orders 01293-422-9111 = Fax 01293-422-922 = Technical 01293-422-999
USA = QIAGEN Inc. = Orders 800-426-8157 = Fax 800-718-2056 = Technical 800-DNA-PREP (800-362-7737)

