

October 2026

QIASymphony® Connect PAXgene® Blood RNA Instructions for Use (Labware List)

Sample and eluate tubes/racks that can be used with the QIASymphony PAXgene Blood RNA Kit and the QIASymphony Connect instrument (software version 1.0.1; labware package SOW-1048)

Trademarks: PAXgene®, PreAnalytiX® (PreAnalytiX GmbH)
QIAGEN®, Sample to Insight®, QIASymphony® (QIAGEN Group);
BD® (Becton, Dickinson and Company);
Corning® (Corning, Inc.);
Eppendorf®, LoBind® (Eppendorf AG);
Sarstedt® (Sarstedt AG and Co.);

PreAnalytiX GmbH, 8634 Hombrechtikon, CH.

HB-0139-L08-001 © 2026 PreAnalytiX GmbH, all rights reserved. Unless otherwise noted, PreAnalytiX, the PreAnalytiX Logo and all other trademarks are property of PreAnalytiX GmbH, Hombrechtikon, CH.

PreAnalytiX Distributors

PreAnalytiX products are manufactured and distributed by QIAGEN and BD for PreAnalytiX.

General information

For Research Use Only. Not for use in diagnostics procedures. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

Important: Before using this Labware List, confirm that it is consistent with the labware package installed on your QIASymphony system.

Legend



Recommended: These tubes are for use in combination with this protocol.



User-defined: User is responsible for defining and/or validating as the user deems appropriate.



Not recommended: These tubes are not for use in combination with this protocol.

x/■

Minimum sample volume (µL) required per sample per protocol (denoted by x)/clot detection possible.

x/□

Minimum sample volume (µL) required per sample per protocol (denoted by x)/clot detection not possible.

Important: Other tubes are available that can be used in combination with other QIASymphony Kits but are not recommended for any protocols that can be used with this kit.

“Sample” drawer, tube carrier

					Protocol
Supplier	Material	Example cat. no.	Name in touchscreen	Insert	PAXgene RNA
BD	PAXgene Blood RDNA Tube 16 x 100 mm	762165	BD#762165 PAXgeneRNA 16x100	No insert needed	

"Eluate" drawer

Supplier	Material	Example cat. no.	Category	Name in touchscreen	Adapter on Elution slots†	Protocol
						PAXgene RNA
QIAGEN	Elution Microtubes CL 96	19588	Deep Well	QIA#19588 *EMTR	Elution Microtube Rack QS, 96-wells, cat. no. 9020730 (cooling slot 1)	
Eppendorf®	1.5 mL Eppendorf LoBind® Snap Cap Safe-Lock tube	0030108.051	Tube 1.5 mL	EP#0030108.051* T1.5 Snap Cap	Micro Tube Snap Cap, 24-wells, cat. no. 9020731 (cooling slot 1)	
Sarstedt®	1.5 mL Microtube, PP, non-skirted	72.607	Tube, 1.5 mL Adapter V2	SAR#72.607* T1.5 Screw	Micro Tube Screw Cap QS, 24-wells, cat. no. 9020674 (cooling slot 1)	
Sarstedt	2.0 mL Microtube, PP, non-skirted	72.693	Tube 2.0 mL/Tube_2.0mL AdapterV1 (no BC)	SAR#72.693* T2.0 Screw	Micro Tube Screw Cap QS, 24-wells, cat. no. 9020674 (cooling slot 1)	
				SAR#72.693 T2.0 Screw	24-Tube 1.5/2.0 ml QS (non-cooling slots 2-4)	
Starlab®	1.5 mL Microtube, graduated conical tube, non-skirted	E1415-2231	Tube, 1.5 mL/Tube_1.5mL AdapterV1 (no BC)	SAR#72.693 T2.0 Screw	24-Tube 1.5/2.0 mL QS (non-cooling slots 2-4)	
				SL#E1415-2231* T1.5 Screw	Micro Tube Screw Cap QS, 24-wells, cat. no. 9020674 (cooling slot 1)	

* Indicates labware that can be cooled using a cooling adapter with barcode.

Visit www.qiagen.com/QIAsymphony for ordering information.

Document revision history

Date	Description
R1, October 2026	Initial release for labware package SOW-1048



For up-to-date licensing information and product-specific disclaimers, see the respective PreAnalytiX or QIAGEN kit handbook or user manual. PreAnalytiX and QIAGEN kit handbooks and user manuals are available at www.preanalytix.com and www.qiagen.com or can be requested from QIAGEN Technical Service or your local distributor.

**Better samples
More to explore**

Explore more at: www.preanalytix.com

HB-0139-L08-001 10/2026

 **PreAnalytiX**
A QIAGEN / BD Company