Novembar 2024.

Uputstvo za upotrebu testa PartoSure[®] Test

Verzija 1

Za korišćenje u in vitro dijagnostici

0197

TTDT-1-20-IVDR

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Reference

1. Antenatal Corticosteroids Revisited: Repeat Courses. NIH Consensus Statement: NIH, 2000.
2. ACOG Committee on Practice Bulletins- Obstetrics. ACOG Practice Bulletin no. 127: Management of Preterm Labor. Clinical management guidelines for obstetrician-gynecologists. *Obstet Gynecol* 2012; 119:1308-1316.
3. Berghella V, Hayes E, Visintine J, Baxter JK. Fetal fibronectin testing for reducing the risk of preterm birth. *Cochrane Database of Systematic Reviews* 2008, Issue 4. Art. No.: CD006843. (Systematic review).
4. Li HY, Chang TS. Vaginal fluid creatinine, human chorionic gonadotropin and alpha-fetoprotein levels for detecting premature rupture of membranes. *Zhonghua Yi Xue Za Zhi (Taipei)*. 2000 Sep;63(9):686-90.
5. Würz H, Geiger W, Künzig HJ, Jabs-Lehmann A, Bohn H, Lüben G. Radioimmunoassay of SP1 (pregnancy-specific beta 1-glycoprotein) in maternal blood and in amniotic fluid normal and pathologic pregnancies. *J Perinat Med*. 1981;9(2):67-78.
6. Kishida T, Hirao A, Matsuura T, Katamine T, Yamada H, Sagawa T, Fujimoto S. Diagnosis of premature rupture of membranes with an improved alpha- fetoprotein monoclonal antibody kit. *Clin Chem*. 1995 Oct;41(10):1500-3.
7. Fuks, Boris et al. "Devices and methods for detecting amniotic fluid in vaginal secretions." Patent 7,709,272. 04 May 2010.
8. Abbassi-Ghanavati M, Greer LG, Cunningham FG. Pregnancy and laboratory studies: a reference table for clinicians. *Obstet Gynecol*. 2009 Dec;114(6):1326-31.
9. CLSI: EP07 Interference Testing in Clinical Chemistry; Approved Guideline. 3rd Edition.

10. CLSI: EP12-A2 User Protocol for Evaluation of Qualitative Test Performance; Approved Guideline - Second Edition.
11. CLSI: EP37 Supplemental Tables for Interference Testing in Clinical Chemistry. 1st Edition.
12. Premarket Approval (PMA) of Adeza Biomedical corporation's Fetal Fibronectin Enzyme Immunoassay Kit (P90048). Food & Drug Administration (FDA), Department of Health & Human Services. 21 Sep 1995.
13. Hainer, BL, Gibson, MV. Vaginitis: Diagnosis and Treatment. Am Fam Physician. 2011 Apr 1;83(7):807-815.
14. Roy S, Sharma M, Ayyagari A, Malhotra S. A quantitative study of bacterial vaginosis. Indian J Med Res 1994; 100:172–6.
15. ACOG Committee Opinion. Use of progesterone to reduce preterm birth. American College of Obstetricians and Gynecologists. Obstet Gynecol. 2003 Nov;102 (5 Pt 1):1115-6.
16. ACOG Practice Bulletin, 107. Induction of Labor. Obstet Gynecol 2009;114:386–97.
17. Nikolova T., et al. Comparison of a novel test for placental alpha microglobulin-1 with fetal fibronectin and cervical length measurement for the prediction of imminent spontaneous preterm delivery in patients with threatened preterm labor. J Perinat Med. 2015. Jul;43(4):395-402.
18. Çekmez Y., et al. Use of cervicovaginal PAMG-1 protein as a predictor of delivery within seven days in pregnancies at risk of premature birth. BMC Pregnancy Childbirth. 2017. Jul 26;17(1):246-250.
19. Wing D. A., et al. Placental Alpha Microglobulin-1 Compared With Fetal Fibronectin to Predict Preterm Delivery in Symptomatic Women. Obstet Gynecol. 2017. Dec;130 (6):1183–91.
20. Melchor J. C., et al. Predictive performance of PAMG-1 vs fFN test for risk of spontaneous preterm birth in symptomatic women attending an emergency obstetric unit: retrospective cohort study. Ultrasound Obstet Gynecol. 2018. Mar;51: 644–649.

21. Ravi M., et al. Evaluation of the quantitative fetal fibronectin test and PAMG-1 test for the prediction of spontaneous preterm birth in patients with signs and symptoms suggestive of preterm labor. *J Matern Fetal Neonatal Med.* 2019. Dec;32(23): 3909-3914.
22. Dawes L. K., et al. The Biomarkers for Preterm Birth Study—A prospective observational study comparing the impact of vaginal biomarkers on clinical practice when used in women with symptoms of preterm labor. *Acta Obstet Gynecol Scand.* 2020. Feb;99(2): 249-258.
23. Nikolova T., et al. Evaluation of a novel placental alpha microglobulin-1 (PAMG-1) test to predict spontaneous preterm delivery. *J. Perinat. Med.* 2014; 42(4): 473–477.
24. Bolotskikh V. and Borisova, V. Combined value of placental alpha microglobulin-1 detection and cervical length via transvaginal ultrasound in the diagnosis of preterm labor in symptomatic patients. *J Obstet Gynaecol Res.* 2017. Aug;43(8): 1263-1269.
25. Santipap M. and Phupong V. Combination of three-dimensional ultrasound measurement of foetal adrenal gland enlargement and placental alpha microglobulin-1 for the prediction of the timing of delivery within seven days in women with threatened preterm labour and preterm labour. *J Obstet Gynaecol.* 2018 Nov;38(8): 1054-1059.
26. Nikolova T., et al. Prediction of spontaneous preterm delivery in women presenting with premature labor: a comparison of placenta alpha microglobulin-1, phosphorylated insulin-like growth factor binding protein-1, and cervical length. *Am J Obstet Gynecol.* 2018 Dec;219(6): 610.e1-610.e9.
27. Lotfi G., et al. Comparison of the effectiveness of a PAMG-1 test and standard clinical assessment in the prediction of preterm birth and reduction of unnecessary hospital admissions. *J Matern Fetal Neonatal Med.* 2019 Mar;32(5): 793-797.
28. Kashanian M., et al. Evaluation of placental alpha microglobulin-1(PAMG1) accuracy for prediction of preterm delivery in women with the symptoms of spontaneous preterm labor; a comparison with cervical length and number of contractions. *J Matern Fetal Neonatal Med.* 2022. Feb;35(3): 534-540.

29. Konoplyannikov A. G., et al. Combination of the placental alpha-1 microglobulin test and ultrasonic cervical length measurement to predict the time of preterm birth. J Matern Fetal Neonatal Med. 2022. Feb;35(3): 541-545.
30. Gokce A., et al. Role of cervical length and placental alpha microglobulin-1 to predict preterm birth. J Matern Fetal Neonatal Med. 2020 Nov; 22;1-5. (Online ahead of print).
31. Barinov S. V., et al. Clinical and biochemical markers of spontaneous preterm birth in singleton and multiple pregnancies. Matern Fetal Neonatal Med. 2021 Feb 24; 1-6. (Online ahead of print).
32. Chawanpaiboon S., et al. Placental α -microglobulin-1 in cervicovaginal fluid and cervical length to predict preterm birth by Thai women with symptoms of labor. Asian Biomed (Res Rev News) 2021; 15(3):119–127.
33. Kehl S., et al. The use of PAMG-1 testing in patients with preterm labor, intact membranes and a short sonographic cervix reduces the rate of unnecessary antenatal glucocorticoid administration. J Perinat Med. 2021 Jul 19;49(9): 1135-1140.
34. Di Renzo G.C., et al. Preterm Labor and Birth Management: Recommendations from the European Association of Perinatal Medicine. J Matern Fetal Neonatal Med. 2017.30(17): 2011-2030.

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Each pack contains:

Svako pakovanje sadrži:

The PartoSure Test is a rapid, non-instrumented, qualitative immuno-chromatographic test for the *in vitro* detection of placental alpha-microglobulin-1 (PAMG-1) in vaginal secretions of pregnant women using a sterile vaginal swab and is not intended for self-testing.

PartoSure Test je brzi, kvalitativni i imunohromatografski test bez instrumenata za *in vitro* otkrivanje alfa mikroglobulina-1 posteljice (PAMG-1) u vaginalnom iscetku trudnice pomoću sterilnog štapića za uzorkovanje iz vagine koji se isporučuje u kompletu.

Warning: Do not fold or bend the test pouch and/or the test strip.

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11/2024 HB-3629-001 1124336SR -IVDR © 2024 QIAGEN, sva prava zadržana.