

August 2013

## Important Note

Dear REPLI-g<sup>®</sup> Mitochondrial DNA Kit user,

The REPLI-g Mitochondrial DNA Kit can now be used for amplifying mitochondrial DNA from other species as well as humans. To enable this, the REPLI-g mt Reaction Buffer no longer contains the human primers. Instead, the primers are provided in a separate tube labeled REPLI-g Human mt Primer Mix. Information on correct use of the REPLI-g Human mt Primer Mix is provided in the REPLI-g Mitochondrial DNA Kit Handbook. Please note that if the kit is used for species other than human, the appropriate primer mix for that species (not included in the kit) should be used. The other kit contents remain unchanged.

If you have any questions, please contact one of the QIAGEN Technical Service Departments or distributors (see back cover of handbooks or visit [www.qiagen.com](http://www.qiagen.com)).

Best regards,

QIAGEN

Trademarks: QIAGEN<sup>®</sup>, REPLI-g<sup>®</sup> (QIAGEN Group).

