

Scoring Sheet — The Pre-Screen Assay Suite

Date:	Protein:	Protein vol.:	μ l
Operator:	Buffer:	Solution vol.:	μ l
Plate ID:	Additives:	Additive vol.:	μ l

Date of observation

Location	Crystallization condition					
A1	4.3 M Sodium chloride					
A2	1.5 M Sodium chloride					
A3	2.0 M Ammonium sulfate, 30% (v/v) Glycerol					
A4	2.0 M Ammonium sulfate, 15% (v/v) Glycerol					
A5	30% (v/v) PEG 400					
A6	12% (v/v) PEG 400					
B1	3.2 M Ammonium sulfate					
B2	1.1 M Ammonium sulfate					
B3	40% (v/v) Isopropanol					
B4	15% (v/v) Isopropanol					
B5	30% (w/v) PEG 4000					
B6	12% (w/v) PEG 4000					
C1	3.0 M Sodium/Potassium phosphate pH 7.0					
C2	1.0 M Sodium/Potassium phosphate pH 7.0					
C3	65% (v/v) MPD					
C4	23% (v/v) MPD					
C5	30% (w/v) PEG 8000					
C6	12% (w/v) PEG 8000					
D1	0.1 M Sodium acetate pH 4.6, 2.0 M Ammonium sulfate					
D2	0.1 M Tris pH 8.5, 2.0 M Ammonium sulfate					
D3	0.1 M Sodium acetate pH 4.6, 50% (v/v) MPD					
D4	0.1 M Tris pH 8.5, 50% (v/v) MPD					
D5	0.1 M Sodium acetate pH 4.6, 20% (w/v) PEG 4000					
D6	0.1 M Tris pH 8.5, 20% (w/v) PEG 4000					

Order EasyXtal and NeXtal products online at www.qiagen.com/crystallization

Trademarks: QIAGEN®, EasyXtal®, NeXtal® (QIAGEN Group) 09/2008 © 2006–2008 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Canada ■ 800-572-9613

China ■ 0086 21 3865 3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800 787980

Japan ■ 03-5547-0811

Korea (South) ■ 1544 7145

Luxembourg ■ 8002 2076

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 65-67775366

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



Sample & Assay Technologies