
Important Note

March 2018

Dear valued QIAxcel Advanced customer,

In addition to the “Electrical safety” information listed on pages 2–3 of the *QIAxcel® Advanced User Manual*, please be aware of the following:

The instrument must be placed in a location that allows clear access to the power cord for easy connection/disconnection.

Best regards,

Your QIAxcel team

Wichtiger Hinweis

März 2018

Sehr geehrte Kundin/sehr geehrter Kunde von QIAxcel Advanced,

zusätzlich zu den Hinweisen zur elektrischen Sicherheit auf den Seiten 2–3 des *QIAxcel® Advanced-Benutzerhandbuchs* sollten Sie Folgendes beachten:

Bei der Aufstellung des Geräts müssen Sie darauf achten, dass das Stromkabel leicht zugänglich bleibt und einfach angeschlossen bzw. getrennt werden kann.

Mit freundlichen Grüßen,

Ihr QIAxcel Team

Note importante

Mars 2018

Cher client QIAxcel Advanced,

En plus des informations de « Sécurité électrique » figurant aux pages 2-3 du *Manuel d'utilisation QIAxcel® Advanced*, veuillez prendre connaissance des éléments suivants :

L'instrument doit être placé dans un endroit qui permet un accès dégagé au cordon d'alimentation pour un branchement / débranchement facile.

Cordialement,

Votre équipe QIAxcel

