



March 2024

Important Note

Changed Pall Acro® Filter

Dear valued QIAGEN® customer,

The #1027058 Pall Acro Filter PTFE 0.2 µm 50 mm has been replaced with 17805 Sartorius Midisart® 2000 Filter PTFE 0.2 µm 64 mm for the QIAvac Connecting System, cat. no. 19419.

Impact on function or performance was tested and our equivalency studies have shown that the updated product can be used equivalently.

Rest assured that this change does not have an impact on quality, safety, function, and performance of the system.

We apologize for any inconvenience this change may cause.

If you have any questions, you can reach us by phone or get online support. Visit www.qiagen.com/support for localized contact details.

Kind regards,

QIAGEN

www.qiagen.com

Trademarks: QIAGEN®, Sample to Insight® (QIAGEN Group); Acro® (Pall Corporation); Midisart® (Sigma-Aldrich). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

03-2024 1134864 HB-3564-001 © 2024 QIAGEN, all rights reserved.