



August 2025

Release Note: QIAquant® Software Version 1.0.4

At QIAGEN®, we are continuously making improvements to provide the highest level of performance with our instruments. Consequently, we are pleased to inform you that an upgrade to your QIAquant instrument operating software is now available. Version 1.0.4 includes the following improvement:

Bug fix

- **Group Name Import:** Fixed an issue where sample group assignments were misaligned after importing a sample layout from an Excel file. Group information now imports correctly.

Note: The QIAquant Software is compatible with QIAquant instruments with serial numbers that start with 3107 (produced before November 2021).

If you have any questions or need assistance, please contact your local QIAGEN representative or contact our Technical Support Center at www.qiagen.com/support/technical-support