

Important Note

Expanded Access: Key PyroMark® Q24 assays now available for the PyroMark Q48 Autoprep

Dear valued PyroMark customer,

We are pleased to inform you that six established PyroMark Q24 kits, previously available exclusively for the PyroMark Q24 platform, are now also available in versions compatible with the PyroMark Q48 Autoprep:

Kit	Cat. No.	Corresponding PyroMark Q24 Kit	Cat. No.
BRAF Pyro Kit Q48	1065142	BRAF Pyro Kit	970470
KRAS Pyro Kit Q48	1065145	KRAS Pyro Kit	970460
EGFR Pyro Kit Q48	1065143	EGFR Pyro Kit	970480
NRAS Pyro Kit Q48	1065144	NRAS Pyro Kit	970530
RAS Extension Pyro Kit Q48	1088848	RAS Extension Pyro Kit	970590
MGMT Pyro Kit Q48	1067482	MGMT Pyro Kit	970061

The main difference between the PyroMark Q24 and PyroMark Q48 formats is the kit composition. PyroMark Q24 kits consist of two component boxes: Box 1 contains the primer and PCR components, and Box 2 contains PyroMark Q24-specific reagents and buffers. Going forward, Box 1 of the PyroMark Q24 kits will include one additional tube for each sequencing primer to ensure compatibility with the automatic primer loading of the PyroMark Q48 Autoprep.

PyroMark Q48 kits consist only of the components included in Box 1 of the PyroMark Q24 kits. Reagents and buffers required for the pyrosequencing reaction are not included in the PyroMark Q48 kits and must be purchased separately as either PyroMark Q48 Advanced Reagents (cat. no. 974002) or PyroMark Q48 Advanced CpG Reagents (cat. no. 974022).

A handbook containing the protocol for running assays on the PyroMark Q48 Autoprep can be downloaded from the Resources section of the PyroMark Q48 website.

We hope these new kit options provide you with greater flexibility and support a smooth transition from the PyroMark Q24 to the PyroMark Q48 Autoprep for customers who wish to continue with pyrosequencing: (www.qiagen.com/PyroMark-Q48-Autoprep).

If you have any questions, please contact QIAGEN Technical Support or your QIAGEN sales representative.

Thank you for your confidence in our Sample to Insight solutions.

Best regards,

Your QIAGEN team

Trademarks: QIAGEN®, Sample to Insight®, PyroMark® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

06/2026 HB-3850-001 © 2026 QIAGEN, all rights reserved.