

QIASymphony® SP-compatible kits at a glance



QIASymphony Kits	Catalog no.	Sample type	Target nucleic acid
QIASymphony DSP DNA Mini Kit (192)*	937236	Human whole blood, buffy coat, tissues and FFPE tissues, as well as viral DNA from human whole blood	Genomic and viral DNA
QIASymphony DSP DNA Midi Kit (96)*	937255		
QIASymphony DSP Circulating DNA Kit (192)*	937556	Plasma, urine	ccfDNA
QIASymphony DSP Circulating DNA Kit (96)*	937555		
QIASymphony DSP Circulating DNA Maxi Kit (192)*	937566		
QIASymphony DSP Virus/Pathogen Mini Kit (192)*	937036	Serum, plasma or CSF Respiratory samples such as swabs, as well as urine and urogenital swabs	Viral DNA, viral RNA and bacterial DNA
QIASymphony DSP Virus/Pathogen Midi Kit (96)*	937055		
QIASymphony DNA Mini Kit (192)	931236	Human whole blood, buffy coat, tissue lysates FFPE tissue lysates, cultured cell lysates, bacterial culture lysates, viral nucleic acids from human whole blood	Genomic DNA
QIASymphony DNA Midi Kit (96)	931255		
QIASymphony DNA Maxi Kit (96)	937266	Human whole blood	
QIASymphony Virus/Bacteria Mini Kit (192)	931036	Plasma, serum, CSF, respiratory samples (BAL, sputum, aspirates, dried swabs, transport media) and urogenital samples (urine, urogenital swabs, transport media)	Viral nucleic acid and bacterial DNA
QIASymphony Virus/Bacteria Midi Kit (96)	931055		
QIASymphony PowerFecal® Pro DNA Kit (192)	938036	Stool, gut, soil samples	Microbial DNA
QIASymphony DNA Investigator Kit (192)	931436	Forensic samples e.g., trace samples, blood, semen, swabs, saliva, bones and hair	Genomic and mitochondrial DNA
QIASymphony RNA Kit (192)	931636	Cultured cells, tissue, FFPE tissues	Total RNA
QIASymphony PAXgene® Blood RNA Kit (96)	762635	Human whole blood collected in PAXgene Blood RNA Tubes	Total RNA
QIASymphony mericon® Bacteria Kit (360)	931156	Primary bacterial enrichment cultures, bacteria from food or environment sample enrichment cultures, bacteria concentrated by filtration	Bacterial DNA extraction
QIASymphony Certal Residual DNA Kit (96)	931855	Bioprocess purification buffer samples Cell culture supernatant samples, biologics purified using chromatography procedures	Residual DNA and viral nucleic acid extraction
QIASymphony Certal Vaccine NA Kit (96)	931955	Vaccine samples; cell culture supernatants, fermentation broths	Residual DNA and viral nucleic acid extraction

* These products are only available in selected countries. Confirm availability with your local representative.

The following products are intended for in vitro diagnostic use: QIASymphony DSP DNA Kits, QIASymphony DSP Circulating DNA Kits and QIASymphony DSP Virus/Pathogen Kits.

The listed products are intended for non-clinical applications. These products are not intended for the diagnosis, prevention or treatment of a disease:

QIASymphony DNA Kits, QIASymphony Virus/Bacteria Kits, QIASymphony PowerFecal Pro DNA Kit, QIASymphony DNA Investigator Kit, QIASymphony RNA Kit, QIASymphony PAXgene Blood RNA Kit, QIASymphony mericon Bacteria Kit, QIASymphony Certal Vaccine NA Kit and QIASymphony Certal Residual DNA Kit.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIASymphony®, mericon®, PowerFecal® (QIAGEN Group); PAXgene® (PreAnalytiX GmbH). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. PROM-24124-001 QPRO-14806 01/2026 © 2026 QIAGEN, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support www.support.qiagen.com | Website www.qiagen.com