

January 2026

PyroMark® Q48 Autoprep – Network Connectivity Notice

Dear valued customer,

We would like to share an important note regarding the network connectivity of the PyroMark® Q48 Autoprep. For security reasons, we strongly advise against connecting the instrument to any network. This recommendation applies both to connections with a local network and to direct connections to the internet.

If you have any questions, you can reach us by phone or get online support. Visit support.qiagen.com for localized contact details.

Kind regards,

QIAGEN

www.qiagen.com

Trademarks: QIAGEN®, Sample to Insight®, PyroMark® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

HB-1024-001 © 2026 QIAGEN, all rights reserved.