

June 2026

Important Note

Use of QX Alignment Marker 1.5 bp/5 kb with the QX DNA High Sensitivity Kit and QIAxcel® ScreenGel®

Dear valued QIAxcel Connect customer,

For QIAxcel ScreenGel 2.1 and later, we are offering an Alignment Marker file to analyze DNA fragments up to 5000 bp with the QIAxcel DNA High Sensitivity Kit (available for download from www.qiagen.com/qiaxcel-connect-system). The QX Alignment Marker 1.5 bp/5 kb (cat. no. 929524) is to be used in a 1:50 dilution with QX DNA HS Dilution Buffer, replacing the QX DNA HS Alignment Marker in the QX DNA High Sensitivity Kit.

If you have further questions, please contact QIAGEN Technical Support. Visit support.qiagen.com for additional resources and contact information.

Best regards,

QIAGEN

Trademarks: QIAGEN®, Sample to Insight®, QIAxcel®, ScreenGel® (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

06/2026 HB-3854-001 © 2026 QIAGEN, all rights reserved.