



June 2026

Important Note

Pre-inserted Filter-Tips & Holders in EZ[®] Kits

Dear valued customer,

At QIAGEN, we are committed to streamline workflows, improving process safety, and reducing environmental impact. That's why we are pleased to share that Filter-Tips are now pre-inserted into their Holders, significantly reducing worktable preparation time before each run. Furthermore, the Filter-Tips design has been optimized to further improve the tip handling by the instrument.

Applicable kits	Old kit components	New kit components
DSP kits*	1 x Disposable Tip Holders (50), mat. no. 1044261	2x Filter-Tips & Holders (25), mat. no. 1136384
	1 x Disposable Filter-Tips (50), mat. no. 1044260	
Other kits**	1 x Disposable Tip Holders (50), mat. no. 1029837	2x Filter-Tips & Holders (25), mat. no. 1136383
	1 x Disposable Filter-Tips (50), mat. no. 1023728	

These changes do not affect kit performance or require the installation of new protocols.

For your further convenience and process safety, all Filter-Tips & Holders are oriented with their upper edges facing the bag opening, enabling faster and more reliable handling during the setup.

If you have further questions, please contact your local QIAGEN representative or visit www.qiagen.com/support for Technical Support.

Sincerely,
QIAGEN

* DSP kits are intended for in vitro diagnostic use: EZ1 DSP Virus Kit, EZ1 DSP DNA Blood Kit
** The listed products are intended for non-clinical applications. These products are not intended for the diagnosis, prevention, or treatment of a disease: EZ1&2 DNA Blood 200 µL Kit, EZ1&2 DNA Blood 350 µL Kit, EZ1&2 DNA Tissue Kit, EZ1&2 Virus Mini Kit v2.0, EZ1&2 RNA Tissue Mini Kit, EZ1&2 DNA FFPE Kit, EZ1&2 DNA FFPE UNG Kit, EZ1&2 ccdDNA Kit, EZ2 RNA FFPE Kit, EZ2 All Prep DNA/RNA FFPE Kit, EZ2 RNA/miRNA Tissue/Cells Kit, EZ2 PowerFecal Pro DNA/RNA Kit, EZ1&2 DNA Investigator Kit, EZ2 DNA Investigator Sep&Prep Kit

Trademarks: QIAGEN[®], Sample to Insight[®], EZ[®], EZ2[®], EZ1&2[®], AllPrep[®], Investigator[®], PowerFecal[®] (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.
1140290 06/2026 HB-3851-01 © 2026 QIAGEN, all rights reserved.