

# Perceive saliva as the next horizon in -omic sciences

With the PAXgene® Saliva Collector



Searching for a better way to collect samples?



Non-invasive and painless approach



User and lab personnel safety



All-in-one device from collection, stabilization to transportation



Reliable and reproducible testing results



Easy to use, at-home collection

Explore more at [www.preanalytix.com](http://www.preanalytix.com)

 **PreAnalytiX**


A QIAGEN / BD Company

Designed to be convenient and cost friendly

Nucleic Acid	Stabilization Period	Downstream Applications	Unique Features
DNA	At least 24 months at up to 25°C	dPCR, real-time qPCR, NGS	<ul style="list-style-type: none"> <li>Bacteriostatic properties preventing bacterial growth</li> <li>Minimized DNA degradation</li> </ul>
SARS-CoV-2 derived RNA	At least 96 hours at up to 20°C	RT-qPCR and more	<ul style="list-style-type: none"> <li>Enhanced lab personnel safety through viral inactivation, preventing infectiousness and replication*</li> <li>Preservation of SARS-CoV-2 derived RNA copy numbers</li> </ul>

\*in cell culture model

PAXgene Saliva Collector completes the comprehensive and verified saliva analysis workflow



#### Workflow Features & Benefits

- Direct sample processing
- Efficient nucleic acid detection & quantification
- High sample reproducibility and reliability
- High quality DNA with no qPCR inhibition, high purity and integrity
- Efficient SARS-CoV-2 derived RNA detection and quantification
- Standardized workflow from sample to insight

Product		Catalog No.
PAXgene Saliva Collectors (25)**	For 25 individual samples: PAXgene Saliva Collector in a single blister box for 2 ml human saliva collection	769040

\*\*For molecular biology applications. This product is not intended for the diagnosis, prevention, or treatment of a disease.



For up-to-date licensing information and product-specific disclaimers, see the respective PreAnalytiX or QIAGEN kit handbook or user manual. PreAnalytiX and QIAGEN kit handbooks and user manuals are available at [www.preanalytix.com](http://www.preanalytix.com) and [www.qiagen.com](http://www.qiagen.com) or can be requested from PreanalytiX Technical Services.



**BLOOD · TISSUE · BONE MARROW**  
The better the source, the more to eXplore.

Explore more at: [www.preanalytix.com](http://www.preanalytix.com)  
PROM-18937-001 1125442 08/2021

**PreAnalytiX**  
A QIAGEN / BD Company