

EZ2[®] Connect-compatible kits at a glance



EZ2 Connect/EZ2 Connect MDx/ EZ2 Connect Fx Kits	Catalog no.	Sample type	Target nucleic acid
EZ1 [®] DSP DNA Blood Kit ¹	62124	Whole blood	Genomic DNA
EZ1 DSP Virus Kit ¹	62724	Serum, plasma, CSF, stool, nasopharyngeal swabs, Gram-positive bacteria	Viral DNA, viral RNA and bacterial DNA
EZ1&2 [®] DNA Blood 200 µl Kit ²	951034	Whole blood	Genomic DNA
EZ1&2 DNA Blood 350 µl Kit ²	951054	Whole blood, buffy coat	Genomic DNA
EZ1&2 DNA Tissue Kit ²	953034	Tissue, buccal swabs, cultured cells, dried blood spots, bacteria in primary samples, bacterial culture or pellets	Genomic DNA, bacterial DNA
EZ1&2 Virus Mini Kit v2.0 ²	955134	Serum, plasma, CSF, urine, whole blood, stool, respiratory samples, transport media, dried swabs, mycoplasma	Viral DNA, viral RNA and bacterial DNA
EZ1&2 RNA Tissue Mini Kit ²	959034	Cultured animal or human cells, human white blood cells, frozen or stabilized easy-to-lyse or high cell density tissue	Total RNA
EZ1&2 DNA FFPE Kit ²	954404	FFPE tissue	Genomic DNA
EZ1&2 DNA FFPE UNG Kit ²	954414	FFPE tissue	Genomic DNA
EZ1&2 ccfDNA Kit ²	954854	Serum, plasma, urine	ccfDNA
EZ2 RNA FFPE Kit ²	959734	FFPE tissue	RNA
EZ2 AllPrep [®] DNA/RNA FFPE Kit ²	954734	FFPE tissue	DNA and RNA
EZ2 RNA/miRNA Tissue/Cells Kit ²	959035	Cultured animal or human cells, human white blood cells, frozen or stabilized tissue	Total RNA, including miRNA and other small RNA
EZ2 PowerFecal [®] Pro DNA/RNA Kit ²	954634	Stool	Microbial DNA, RNA or total nucleic acids
EZ1&2 DNA Investigator [®] Kit ³	952034	Forensic samples: Trace samples, blood, semen, saliva, bones and hair	DNA
EZ2 DNA Investigator Sep&Prep Kit ³	952134	Sexual assault samples	Sperm-derived DNA

¹ Kit can be used with EZ2 Connect MDx only

² Kit can be used with EZ2 Connect, EZ2 Connect Fx, EZ2 Connect MDx (in research mode)

³ Kit can be used with EZ2 Connect Fx only

The following products are intended for in vitro diagnostic use: EZ1 DSP Virus Kit, EZ1 DSP DNA Blood Kit and EZ2 Connect MDx.

The listed products are intended for non-clinical applications. These products are not intended for the diagnosis, prevention or treatment of a disease:

EZ1&2 DNA Blood 200 µl Kit, EZ1&2 DNA Blood 350 µl Kit, EZ1&2 DNA Tissue Kit, EZ1&2 Virus Mini Kit v2.0, EZ1&2 RNA Tissue Mini Kit, EZ1&2 DNA FFPE Kit, EZ1&2 DNA FFPE UNG Kit, EZ1&2 ccfDNA Kit, EZ2 RNA FFPE Kit, EZ2 AllPrep DNA/RNA FFPE Kit, EZ2 RNA/miRNA Tissue/Cells Kit, EZ2 PowerFecal Pro DNA/RNA Kit, EZ1&2 DNA Investigator Kit, EZ2 DNA Investigator Sep&Prep Kit, EZ2 Connect and EZ2 Connect Fx.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], Sample to Insight[®], EZ1[®], EZ2[®], EZ1&2[®], AllPrep[®], Investigator[®], PowerFecal[®] (QIAGEN Group). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. PROM-23806-003 10/2025 © 2025 QIAGEN, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support www.support.qiagen.com | Website www.qiagen.com