

QIAxpert[®] Potassium Dichromate (3)

For testing the photometric accuracy of the QIAxpert instrument

QIAxpert Potassium Dichromate is intended for use in testing the photometric accuracy and precision of the QIAxpert instrument in combination with QIAxpert slide consumables.

Storage conditions

QIAxpert Potassium Dichromate should be stored at room temperature (15–25°C) upon arrival and protected from light. Opened vials should be used within 1 hour. Please check the expiration date on the product label.

Quality control

In accordance with QIAGEN's ISO-certified Quality Management System, each lot of QIAxpert Potassium Dichromate is tested against predetermined specifications to ensure consistent product quality.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view, and print the SDS for each QIAGEN kit and kit component.

Product use limitations

QIAxpert Potassium Dichromate is intended to be used only in combination with the QIAxpert instrument.

Procedure

For complete instructions for use, please refer to the *QIAxpert User Manual*, which can be downloaded from the resources tab at www.qiagen.com/p/QIAxpert.



Ordering Information

Product	Contents	Cat. no.
QIAxpert Potassium Dichromate (3)	Solution for testing the photometric accuracy of the QIAxpert instrument; 3 x 1 ml vials	990701

QIAxpert Potassium Dichromate is intended for molecular biology applications. This product is not intended for the diagnosis, prevention, or treatment of a disease.

Trademarks: QIAGEN®, QIAxpert® (QIAGEN Group).
1081164 02/2014 © 2014 QIAGEN, all rights reserved.

www.qiagen.com
Australia ■ 1-800-243-800
Austria ■ 0800-281011
Belgium ■ 0800-79612
Brazil ■ 0800-557779
Canada ■ 800-572-9613
China ■ 800-988-0325

Denmark ■ 80-885945
Finland ■ 0800-914416
France ■ 01-60-920-930
Germany ■ 02103-29-12000
Hong Kong ■ 800 933 965
India ■ 1-800-102-4114
Ireland ■ 1800 555 049

Italy ■ 800-787980
Japan ■ 03-6890-7300
Korea (South) ■ 080-000-7145
Luxembourg ■ 8002 2076
Mexico ■ 01-800-7742-436
The Netherlands ■ 0800-0229592
Norway ■ 800-18859

Singapore ■ 1800-742-4368
Spain ■ 91-630-7050
Sweden ■ 020-790282
Switzerland ■ 055-254-22-11
Taiwan ■ 0080-665-1947
UK ■ 0808-2343665
USA ■ 800-426-8157

