

# Fungi ITS2 Primer 1 For Sequencing

Analyte Specific Reagent. Analytical and performance characteristics are not established.

## Contents

| Fungi ITS2 Primer 1 For Sequencing |        |
|------------------------------------|--------|
| Catalog no.                        | 979013 |
| Seq 1 Fungi (90 µl; 20 µM)         | 1 vial |
| Product Sheet                      | 1      |

## Shipping and Storage

Fungi ITS2 Primer 1 For Sequencing is shipped at ambient temperature and should be stored at –20°C upon receipt.

## Product Description

Fungi ITS2 Primer 1 For Sequencing consists of Seq 1 Fungi for sequencing of a hypervariable region of ITS2 of fungal species. The purity of the product was determined by anion-exchange HPLC and found to be >97%.

## Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs). These are available online in convenient and compact PDF format at [www.qiagen.com/Support/MSDS.aspx](http://www.qiagen.com/Support/MSDS.aspx) where you can find, view, and print the MSDS for each QIAGEN kit and kit component.

## 24-hour emergency information

Emergency medical information in English, French, and German can be obtained 24 hours a day from: Poison Information Center Mainz, Germany, Tel: +49-6131-19240



---

Sample & Assay Technologies

Trademarks: QIAGEN® (QIAGEN Group).

1056951 03/2009 © 2009 QIAGEN, all rights reserved.

**www.qiagen.com**

**Canada** ■ 800-572-9613

**China** ■ 0086-21-3865-3865

**Denmark** ■ 80-885945

**Finland** ■ 0800-914416

**France** ■ 01-60-920-930

**Germany** ■ 02103-29-12000

**Hong Kong** ■ 800 933 965

**Ireland** ■ 1800 555 049

**Italy** ■ 800-787980

**Japan** ■ 03-6890-7300

**Korea (South)** ■ 1544 7145

**Luxembourg** ■ 8002 2076

**Mexico** ■ 01-800-7742-639

**The Netherlands** ■ 0800 0229592

**Norway** ■ 800-18859

**Singapore** ■ 65-67775366

**Spain** ■ 91-630-7050

**Sweden** ■ 020-790282

**Switzerland** ■ 055-254-22-11

**UK** ■ 01293-422-911

**USA** ■ 800-426-8157

**Australia** ■ 1-800-243-800

**Austria** ■ 0800/281010

**Belgium** ■ 0800-79612

**Brazil** ■ 0800-557779



---

**Sample & Assay Technologies**