

March 2017

GeneRead™ Sequencing Q Wash Buffers (2): Product Sheet

Catalog number 1106081

For next-generation sequencing applications using the
QIAGEN GeneReader™ instrument

For research use only

Not for use in diagnostic procedures

GeneRead Sequencing Q Wash Buffers are included in the GeneRead Fast Sequencing Q Kit (cat. no. 185221) to perform sequencing runs on the QIAGEN GeneReader instrument. The wash buffers support the sequencing chemistry process and maintenance procedures.

IMPORTANT: When used in combination with the GeneRead Fast Sequencing Q Kit (cat. no. 185221), the GeneRead Sequencing Q Wash Buffers will allow the user to process 2 flow cells.

For details about reagent preparation and use, refer to the *QIAGEN GeneRead® Fast Sequencing Q Kit Handbook* and the *QIAGEN GeneReader User Manual for Advanced Process Flow (APF) Instrument Configuration*.

Shipping and Storage

GeneRead Sequencing Q Wash Buffers (cat. no. 1106081) are shipped on cool packs and should be stored immediately upon receipt at 2°C to 8°C. Check buffers for precipitates before use and re-dissolve at 37°C, if necessary. Buffers should be protected from light exposure. We recommend storing all reagents in the dark.

Ordering Information

Product	Contents	Cat. no.
GeneRead Sequencing Q Wash Buffers (2)	Includes 1 x 267.9 ml Wash Buffer 9, 2 x 249.6 ml of Wash Buffer 11 Premix, 2 x 27.7 ml Wash Buffer 11 Add-on, and 1 X 170 ml of Maintenance Wash Buffer	1106081
Related products		
GeneRead Fast Sequencing Q Kit (2)	Includes Reagents and Buffers, Add-Ons, and 2 flow cells supplied for up to 2 single flow cell runs on the GeneReader	185221
GeneReader Platform	Next-generation sequencing instrument: includes installation and training, 1 year warranty on parts and labor	9002312

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIAGEN GeneRead®, GeneRead™ (QIAGEN Group); GeneReader™ (Intelligent Bio-Systems, Inc.).

Limited License Agreement for GeneRead Sequencing Q Wash Buffers

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. QIAGEN grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at www.qiagen.com. Some of these additional protocols have been provided by QIAGEN users for QIAGEN users. These protocols have not been thoroughly tested or optimized by QIAGEN. QIAGEN neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, QIAGEN makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. QIAGEN specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. QIAGEN may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see www.qiagen.com.

HB-2363-001 1106876 03/2017 © 2017 QIAGEN all rights reserved.

