


轻松掌控您的 NGS 工作流程

想要一个让您高枕无忧的 GeneReader™ NGS 系统？

使用 Horizon® 多重定量 FFPE 参考标准品即可对 NGS 工作流程进行质控



无缝衔接安装、培训和设备调试

- 确保整个 NGS 工作流程成功安装和调试
- 采用标准品进行安装培训，更保证培训效果
- 利用验证过的数据集合和确定的流程指导 QC 工作

从核酸提取到数据解读，整个 GeneReader 工作流程无懈可击

- Horizon FFPE 参比样本中, 利用 GeneReader NGS 系统能准确地检出 GeneRead QIAact Actionable Panel 覆盖的所有突变位点, 等位基因频率低至 5%
- 不同样本, 试剂批次, 操作人员, 仪器和实验场地, 均能确保检测结果一致性




图 1. 在不同场地, 对不同批次 Horizon 多重定量 FFPE 参考标准品进行变异频率检测的可重复性。

Learn more about the GeneReader NGS System and Horizon Reference Standards at www.GeneReaderNGS.com.


关于最新的许可信息和产品特定的免责声明, 请阅读相关的 QIAGEN 试剂盒手册或操作指南。QIAGEN 试剂盒手册和操作指南可在 www.qiagen.com 下载, 或向 QIAGEN 技术服务或当地的经销商索取。

Trademarks: QIAGEN®, Sample to Insight®, GeneRead® (QIAGEN Group), GeneReader™ (Intelligent Bio-Systems, Inc.), Horizon® (Horizon Discovery Group plc.).
Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.
© 2016 QIAGEN, all rights reserved.
The GeneReader NGS System is for Research Use Only. Not intended for diagnostic procedures

凯杰企业管理 (上海) 有限公司

电话: 021-3865 3865 ■ 技术支持热线: 800-988-0325 400-880-0325

TechService-CN@qiagen.com ■ www.qiagen.com



LS201606004



凯杰诊断