

QIAGEN® Fecal Kits Guide

Isolate high-quality DNA and RNA from fecal samples

QIAGEN fecal DNA and RNA isolation kits isolate high-quality DNA or RNA from even the most challenging fecal samples. With these kits, purification of nucleic acids from difficult samples, including those with PCR inhibitors, is made simple and protocols can be automated on several instruments which saves time and standardizes methods.

QIAamp® PowerFecal® DNA Kit



- Best for DNA microbiome applications
- Use for stool, gut and biosolid samples
- Based on the same technology as the DNeasy® PowerSoil® Kit
- Protocol can be automated on the QIAcube® instrument

Cat. no. 12830-50

QIAamp 96 PowerFecal QIAcube HT Kit



- Best for high-throughput pathogen detection
- New product for stool, gut and biosolids
- This kit has a combination of Inhibitor Removal Technology® and optimized extraction chemistry
- Contains 96-well bead beating plates and 96-well silica-spin filter plates
- Automated protocol on the QIAcube HT

Cat. no. 51531

RNeasy® PowerMicrobiome® Kit



- Best for total RNA isolation from stool, gut and biosolid samples
- Isolates DNA and RNA if DNase step is omitted
- Protocol can be automated on the QIAcube

Cat. no. 26000-50

AllPrep® PowerViral® DNA/RNA Kit



- Best for viral and bacterial nucleic acid isolation from waste water and stool samples

Cat. no. 28000-50

AllPrep PowerFecal DNA/RNA Kit



- Best for simultaneous DNA and RNA isolation from stool, gut or biosolid samples
- This new kit combines Inhibitor Removal Technology with optimized extraction chemistry
- Contains new optimized beads for increased yields

Cat. no. 80244

MagAttract® PowerMicrobiome DNA/RNA EP/KF Kits



- Best for high-throughput isolation using automated liquid handling systems
- Uses magnetic particle technology
- Designed for isolation of nucleic acids from stool, gut and biosolid samples
- Automated protocol on the KingFisher® Flex or Duo and other liquid handling systems

Cat. no. 27500-4-EP

Cat. no. 27600-4-KF

Visit www.qiagen.com or contact QIAGEN today to discover more about QIAGEN fecal kits!

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QIAamp®, QIAcube®, AllPrep®, DNeasy®, Inhibitor Removal Technology®, MagAttract®, PowerFecal®, PowerMicrobiome®, PowerSoil®, PowerViral®, RNeasy® (QIAGEN Group); KingFisher® (Thermo Fisher Scientific, Inc.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.
1107569 PROM-11140-002 07/2017 © 2017 QIAGEN, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com