


pQE-100 DoubleTag™ Vector

Positions of elements in bases	
Vector size (bp)	3512
Start of numbering at <i>Xho</i> I (CTCGAG)	1–6
T5 promoter/lac operator element	7–87
T5 transcription start	61
6xHis-tag coding sequence	127–144
Multiple cloning site	145–192
Tag-100 epitope coding sequence	193–229
Lambda <i>t</i> ₀ transcriptional termination region	259–353
<i>rrnB</i> T1 transcriptional termination region	1115–1212
ColE1 origin of replication	1689
β-lactamase coding sequence	3307–2447



pQE-100 DoubleTag™

